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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for the coloring of peripheral blood cells, bone marrow or for cytological study of cellular elements. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Dye second May Grünwald - Giemsa
May Grünwald dye - Giemsa used to stain cells is a mixture of dyes with neutral characteristics, which stain the nuclear components and cytoplasmic cells.

REAGENTS

Dye - May Grünwald - Store between 2 and 8°C. Contains: Eosin and Methylene Blue 180 mg/dL.

Dye - Giemsa - Store between 15 and 30°C. Contains: Azur II 80 mg/dL and Azur Eosinato II 100 mg/dL.

PRESENTATION

Reagent	Volume
C004 Dye - May Grünwald	500 mL
C002 Dye - Giemsa	500 mL

Reagents are sold separately.

EQUIPMENTS AND OPERATIONAL INPUTS

Support for coloring, watch or stopwatch, pipettes, slides. Can be found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transportation temperature should be between 15 to 30°C. Protect from light and avoid moisture.

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed.

3- Water used in material cleaning must be recent and free of contaminants.

4- Handle reagent with care, avoid contact with skin and mucous.

5- Slides for the fabrication of films should be clean and free of fat.

6- A drop of blood used for making the smear should not be too large because the larger the droplet the more thicker the smear.

7- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

8- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

9- Do not use the product in case of damaged packaging.
10- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Peripheral blood collected by fingerstick or veinpuncture, preferably without the use of anticoagulant or **EDTA Bioclin**. The slides made from blood collected with Anticoagulants should be made up to 30 minutes after collection to avoid cellular deformations.

PROCESS DESCRIPTION

TECHNIQUE

Place slide on staining rack.

1- Cover the smear with 20 drops of May Grünwald dye. Leave for 3 minutes.

2- After this time, add 20 drops of tap water. Mix the solution in different parts of the blade. Prevent the leakage of the dye. Wait 2 minutes.

3- Discard the mixture that coats the slide and, without washing, cover it with 20 drops of Giemsa diluted solution, prepared at the time of coloring (a drop of Giemsa for each milliliter of distilled water - on average it is used 3 mL for each slide). Wait for 12 to 15 minutes.

4- Discard the dye, rinse the slide under running water.

5- Allow it to dry upright.

6- Check with immersion objective.

RESULTS

Erythrocytes: Pinkish coloring.

Platelets: Bluish coloring.

Neutrophiles: Dark blue nucleus and pale pink cytoplasm with granules ranging from pink to blue.

Lymphocytes: Violet blue nucleus and blue cytoplasm.

Basophiles: Purple to dark blue nucleus and cytoplasm with bulky dark blue granules.

Monocytes: Violet blue nucleus and pale blue cytoplasm.

Eosinophiles: Blue nucleus and pale pink cytoplasm with bulky granules ranging from red to orange.

PROCEDURE LIMITATIONS

Smears that present a blue or green color are alkaline due to the prolonged action of the dye, thick smears or insufficient washing of the slide.

Smears that are present red color are overly acidic due to weak staining, washing time or excessive acidity of the dye.

Precipitation may appear when the slides are unwashed, dust or drying of the dye occurs during the staining process.

Eventually the water used by the laboratory does not provide a good coloring. In this case use a Phosphate Tampon (pH 7,2), using the following formula:

1 - Solution A	
KH ₂ PO ₄	9,1 g
Distilled H ₂ O q.s.p.	1.000 mL

2 - Solution B	
Na ₂ HPO ₄	9,5 g
Distilled H ₂ O q.s.p.	1.000 mL

USAGE TAMPON

Solution A	28 mL
Solution B	72 mL

PRODUCT PERFORMANCE

CONTROL QUALITY

Accuracy

COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY

May Grünwald - Giemsa kit was compared with other smears coloring methods commercially available. According to the coloring comparison of 10 clinical samples, the method shows an excellent correlation.

Precision

REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 5 different samples, obtaining the following results:

Sample	Nº of Repetitions	Expected Results	Found Results
01	20	Standard Coloring	Standard Coloring
02	20	Standard Coloring	Standard Coloring
03	20	Standard Coloring	Standard Coloring
04	20	Standard Coloring	Standard Coloring
05	20	Standard Coloring	Standard Coloring

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 5 different samples, obtaining the following results:

Sample	N° of Repetitions	Expected Results	Found Results
01	20	Standard Coloring	Standard Coloring
02	20	Standard Coloring	Standard Coloring
03	20	Standard Coloring	Standard Coloring
04	20	Standard Coloring	Standard Coloring
05	20	Standard Coloring	Standard Coloring

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QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca

CEP 31565-130 - Belo Horizonte - MG - Brasil

Phone.: + 55 (31) 3439.5454 - Fax: + 55 (31) 3439.5455

E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 - Made in Brazil

OBELIS S.A.

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1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service


Phone.: 0800 0315454

E-mail: sac@bioclin.com.br

ANVISA registration for May Grünwald - Giemsa kit:
10269360112

Review: July/2013

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED