MICROALBUMINÚRIA

REF K078

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for quantitative determination of Microalbumin in human samples of urine by turbidimetry analysis, for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Immunoturbidimetry

Turbidimetric assay for quantification of Microalbumin in urine sample. The Microalbumin present in human urine sample reacts with a specific antibody. Turbidity is induced by the formation of immune complexes, which are directly proportional to its concentration in the sample.

REAGENTS

Number 1 - Buffer - Store between 2 and 8°C. Contains: PEG 6.000 < 15%, buffer and preservative.

Number 2 - Antiserum - Store between 2 and 8°C. Contains: Goat antiserum for Human Albumin < 2 g/L, buffer and preservative.

Number 3 - Calibrator - Store between 2 and 8°C. Contains: Serum, stabilizer and preservative. **Potentially Infectious.**

Warning: The concentration of the Calibrator varies with the lot - See the bottle label.

Number 4 - Control - Store between 2 and 8°C. Contains: Serum, stabilizer and preservative. **Potentially Infectious.**

Warning: The concentration of the Calibrator varies with the lot - See the bottle label.

PRESENTATION

Presentation	Reagent			
	N° 1	N° 2	N° 3	N° 4
1	25 mL	2 mL	1 mL	1 mL
2	2 x 25 mL	2 x 2 mL	1 mL	1 mL
3	3 x 25 mL	3 x 2 mL	1 mL	1 mL
4	4 x 25 mL	4 x 2 mL	1 mL	1 mL
5	40 mL	10 mL	1 mL	1 mL
7	20 mL	5 mL	1 mL	1 mL
8	2 x 20 mL	2 x 5 mL	1 mL	1 mL
9	3 x 20 mL	3 x 5 mL	1 mL	1 mL
10	4 x 20 mL	4 x 5 mL	1 mL	1 mL

EQUIPMENTS AND OPERATIONAL INPUTS

This kit is specific to automatic equipment.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be 2 to 8°C. Transportation, at temperatures up to 30°C should not exceed 5 days. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain accurate results.

3- The water used in the material rinsing must be recent and deionized.

4- Reagents n° 1, 2, 3 e 4 contain Sodium Azide, irritant to skin and mucous. Handle carefully.

5- Do not freeze the reagents, the freezing of antibody or buffer can affect their performance.

6- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

7- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

8- Do not use the product in case of damaged packaging.

9- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Turbid urine or with visible sediment must be centrifuged before use. Centrifuge for 10 minutes at 3000 rpm. Use the supernatant for testing.

The Microalbumin in urine is stable for 48 hours if refrigerated at 2 and 8°C. For longer periods, freeze the sample at -20°C for up to 3 months, and only once. All samples are considered potentially infectious, and therefore must be handled following the rules of biosafety.

WORKING REAGENT PREPARATION

Reagents are ready to use.

TECHNIQUE

The calibration stability of the Microalbumin kit installed on refrigerated equipment is at least 15 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

For calibration and control of the reaction, using calibrator and internal controls in the kit.

This kit is for use in automated biochemical analyzers. **Check the application sheet for this equipment on www.bioclin.com.br or through SAC.**

CALCULATIONS

24 Hour Urine

For determination of Microalbumin in 24-hours urine should be performed dosage Microalbumin in the sample following the original technique of kit. The result found by the technique is given in mg/L. For expression of the result in mg/24 hours using the formula below.

$$\text{Microalbumin (mg/24horas)} = \frac{\text{Microalbumin (mg/L)}}{\text{Urinary Volume (L)*}} \times 24 \text{ hours (L)*}$$

*To convert the urine volume of mL to L divide the value given in mL by 1000.

Example:

Results of Microalbumin found in urine samples of 24 hours = 56,0 mg/L

Urinary volume of 24 hours = 1,9 L

Microalbumin (mg/24 hours) = 56,0 x 1,9 = 106,4 mg/24 hours

Relationship Microalbumin/Creatinine (mg/g)

To determine the relationship between Microalbumin and urinary Creatinine due, firstly, to perform dosage Microalbumin and Creatinine in the same urine sample. A urine sample can be random or first morning.

The urine Creatinine values are given in mg/dL. To convert the obtained concentration in mg/dL to g/L, simply divide the result by 100.

The result of Microalbumin and Creatinine ratio is given as mg/g.

Example:

Creatinine value found in sample = 120 mg/dL

Microalbumin value found in the sample (mg/L) = 60,0 mg/L

$$\text{Urinary Creatinine (g/L) found} = \frac{\text{Value Creatinine}}{100}$$

$$\text{Urinary Creatinine (g/L)} = \frac{120 \text{ mg/dL}}{100} = 1,2 \text{ g/L}$$

$$\text{Microalbumin / Creatinine (mg/g)} = \frac{\text{Microalbumin (mg/L)}}{\text{Creatinine (g/L)}}$$

$$\text{Microalbumin / Creatinine (mg/g)} = \frac{60,0 \text{ mg/L}}{1,2 \text{ g/L}} = 50,0 \text{ mg/g}$$

PROCEDURE LIMITATIONS

Clinical diagnosis should be made not only with results of one test should be considered at the same time, clinical data of patients.

INTERFERENCES

Triglycerides up to 2500 mg/dL, Bilirubin up to 15 mg/dL and Hemoglobin up to 125 mg/L did not significantly interfere with results. Some drugs and substances affect the concentration of Microalbumin in urine.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's calibrator is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 927D.

REFERENCE VALUES

Reference values below should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

Random Urine or First Morning

Microalbumin: 0 - 25 mg/L (IFCC)

Relation Microalbumin and Urinary Creatinine (mg/g)

Reference Value	Men (mg/g)	Women (mg/g)
Normal	< 17	< 25
High	17 - 250	25 - 355
Very High	> 250	> 355

24 Hour Urine

Reference Value	Men and Women (mg/24 hours)
Normal	< 30
High	30 - 300
Very High	> 300

Urine 4, 6 or 12 Hours

Urines to 4, 6 or 12 hours perform the same calculations used for 24-hour urine, but with the total volume of urine collected during the time interval.

Collection Time	Reference Value Men e Women	
Urine of 4 Hours	< 5 (mg/4 hours)	Or < 20 µg/min
Urine of 6 Hours	< 7,5 (mg/6 hours)	
Urine of 12 Hours	< 15 (mg/12 hours)	

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy RECOVERY

The recovery analysis was performed with 05 determinations samples. Accuracies were calculated and met in good agreement with the reference values, obtaining a recovery between 91 and 102%.

COMPARISON OF METHODS AND SPECIFICITY METHODOLOGY

Microalbumin kit was compared with another method for measurement of Microalbumin commercially available. 40 analyzes were performed and the results were evaluated. The linear equation obtained was $Y = 0,9458X + 5,0277$ and coefficient correlation 0,9981. With these results we can conclude that the kit shows good methodological specificity.

Precision

REPEATABILITY

The repeatability was calculated from 10 successive determinations, using 2 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2
Average Concentration (mg/dL)	20,20	130,31
Standard Deviation (mg/L)	0,49	0,96
Coefficient of Variation (%)	2,40	0,74

REPRODUCIBILITY

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 2 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2
Average Concentration (mg/dL)	18,89	129,11
Standard Deviation (mg/L)	1,46	1,31
Coefficient of Variation (%)	7,75	1,02

Sensitivity

The sensitivity was calculated from 10 determinations of a sample free of Microalbumin. The average found was 0,251 mg/L, with standard deviation of 0,087 mg/L. The sensitivity, wich indicates the method detection limit, corresponds the average plus 3 times the standard deviation and is equal to 0,512 mg/L.

Linearity

The linearity depends on the concentration value of the calibrator used and can not be greater than 400 mg/L. For higher values, dilute sample 1:5 with distilled or deionized water and repeat the determination. Multiply the result by the factor dilution.

DIAGNOSTIC SIGNIFICANCE

Diabetic nephropathy, which is accompanied by irreversible damage kidneys and persistent proteinuria is the leading cause of death people with insulin dependent *diabetes mellitus*. A signal early diabetic nephropathy is the presence of Albumin in urine at low concentrations. Early detection of minor kidney damage is very important, because the rapid onset of treatment allows reversion of the condition.

BIBLIOGRAPHIC REFERENCES

1. Mount, JN, J. Clin. Pathology, 22, 12 (1986).
2. Schmidt, A., et al., Diabetic Medicine, 5, 126 (1988).
3. Killingsworth, LM; SAVORY, J: Automated Immunochemical Procedures for Measurement of Immunoglobulins IgG, IgA and IgM in human serum, Clin. Chem. 1971, 17, 936.
4. PESCE, AJ, KAPLAN, LA: Methods in Clinical Chemistry CV Mosby Company, 1987.
5. Hellsing, K., Profiles in the Biological Fluids, 1973, 23: 579.
6. BLOM, M. and HOJME, H., Clinical Chemistry, 1976, 22: 657.
7. HILLS, LP and TIFFANY, IT, Clinical Chemistry, 1980, 26: 1459.
8. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

 **OBELIS S.A.**

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for Microalbumin kit: 10269360206

Review: October/2020

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MANUFACTURED BY

BATCH CODE

CONTROL

DATE OF MANUFACTURE

POSITIVE CONTROL

USED BY
(last day of month)

NEGATIVE CONTROL

TEMPERATURE LIMITATION
(store at)

BIOLOGICAL RISK

CONTAINS SUFFICIENT
FOR <N> TESTS

INFLAMMABLE

CONSULT INSTRUCTIONS
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC DEVICE

POISON

EUROPEAN AUTHORIZED
REPRESENTATIVE

CE MARK

KEEP AWAY
FROM SUNLIGHT

DO NOT USE IF
PACKAGE IS
DAMAGED