MUCOPROTEÍNAS

REF K028

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for determination of Mucoproteins in samples of serum. Colorimetric test. For *in vitro* diagnostic only.

PRINCIPLE OF ACTION

Methodology: Winzler Modified

The proteins precipitate in a Perchloric Acid solution (Reagent Nº 1), resulting a glycoprotein fraction called Seromuroid or Mucoproteins. These are precipitate in filtered with Phosphotungstic Acid (Reagent Nº 2), later dissolved and dosed through the contend of Tyrosine with the Folin Reagent. The blue color developed is proportional to the concentration of Mucoproteins in the sample.

REAGENTS

Number 1 - Perchloric Acid - Store between 15 and 30°C. Contains: Perchloric Acid < 10 mol/L.

Number 2 - Phosphotungstic Acid - Store between 15 and 30°C. After use, close the vial to avoid contamination with CO₂ from the air. Contains: Phosphotungstic Acid < 35mmol/L and stabilizer.

Number 3 - Sodium Carbonate Stock - Store between 15 and 30°C. Contains: Sodium Carbonate < 4 mol/L.

Number 4 - Folin Reagent - Store between 15 and 30°C. Contains: Sodium tungstate < 600 mmol/L, Sodium Molybdate < 40 mmol/L, Phosphoric Acid < 2 mol/L, Lithium Sulphate < 3 mol/L and activators.

Number 5 - Standard - Store between 15 and 30°C. After handling it is advisable to keep between 2 and 8°C to avoid evaporation. Contains: Tyrosine 40 mg/dL (5 mg/dL Mucoproteins) and stabilizer.

PRESENTATION

Reagent	Presentation 1	Presentation 2
Nº 1	75 mL	300 mL
Nº 2	23 mL	90 mL
Nº 3	50 mL	200 mL
Nº 4	10 mL	40 mL
Nº 5	2 mL	8 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, water-bath at 37°C, watch and stopwatch, funnel, high retention glass paper (Whatman Nº 50, Toyo Nº 4, Green S807 or SS589 blue track), pipettes, test tubes, volumetric balloon, Biocontrol N and Biocontrol P Bioclin. These items are found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperatures shall be between 15 and 30°C. Protect from light and avoid moisture.

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must to be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- The water level in water bath must by higher than the reagent level in test tubes.

6- The filtered deproteinization must be clear. Use the recommended filter paper. In case of high concentration of Mucoproteins, if the filter present turbidity, filter it again.

7- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

8- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

9- Do not use the product in case of damaged packaging.

10- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLE

Serum obtained free of hemolysis. The analyte is stable for 7 days between 2 and 8°C.

PROCESS DESCRIPTION

STABILITY AFTER OPEN

The stability after preparation of the working reagent (Sodium Carbonate) of the Mucoproteínas kit is at least 6 months between 15 and 30°C. This stability may vary according to the test conditions and the environment.

Therefore, it is suggested to monitor the product's performance using control sera.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin kits.

PREPARATION OF USAGE SODIUM CARBONATE

Presentation 1

Transfer 50 mL of Reagent Nº 3 (Sodium Carbonate Stock) to a 250 mL volumetric balloon. Complete the volume with distilled or deionized water and homogenize it. Store it in a plastic flask. The reagent is stable for 06 months in temperatures between 15 and 30°C.

Presentation 2

Transfer 200 mL of Reagent Nº 3 (Sodium Carbonate Stock) to a 1000 mL volumetric balloon. Complete the volume with distilled or deionized water and homogenize it. Store it in a plastic flask. The reagent is stable for 06 months in temperatures between 15 and 30°C.

DEPROTEINIZATION

In a 18x180 test tube, marked with the letter S (Sample), add:

Sodium Chloride 0.85%	6.0 mL
Serum	1.0 mL
Reagent Nº 1	3.0 mL

Roughly agitate in a vortex for exact 2 minutes and filter it through the high retention filter paper (Whatman Nº 50, Toyo Nº 4, Green S807 or SS589 blue track).

MUCOPROTEINS PRECIPITATION

In centrifuge tube marked with the letter "5" (Sample), add:

Filtered	4.0 mL
Reagent Nº 2	0.8 mL

Agitate and let it rest during 10 minutes. Centrifuge at 2500 - 4000 rpm for 5 minutes. Disposal all supernatant. Pour the tube over an absorbent paper, let it rest for 10 minutes to drain all the liquid excess. This sample tube will be used in Colorimetry.

COLORIMETRY

	Blank	Standard	Sample
Reagent Nº 5	---	50 mL	---
Usage Carbonate	5.0 mL	5.0 mL	5.0 mL

Agitate the tube "S" (Sample) to dissolve the precipitate. Then, add:

Reagent Nº 4	200 mL	200 mL	200 mL
--------------	--------	--------	--------

Agitate the tubes immediately and place them in water-bath at 37°C for 15 minutes. Determinate the absorbance at 680 nm (640 - 700 nm), matching the zero with the tube "B" (Blank). The color reaction is stable for 2 hours.

CALCULATIONS

Mucoproteins (mg/dL) = $\frac{\text{Sample Absorbance}}{\text{Standard Abs.}}$ x 5
in Tyrosine fraction

The reaction follows the Lambert Beer Law, it can use the Calibration Factor:

Calibration Factor = $\frac{\text{Standard Concentration (5 mg/dL)}}{\text{Standard Absorbance}}$

Mucoproteins (mg/dL) = $\frac{\text{Sample Absorbance}}{\text{x Calibration Factor}}$
in Tyrosine fraction

To calculate the total Mucoprotein value, multiply the Tyrosine fraction obtained result by the 23.8.

PROCEDURE LIMITATIONS

The use of the high retention paper filter is indispensable for obtaining correct results.

INTERFERENT

No interference was observed for Bilirubin up to 5 mg/dL, Hemoglobin up to 30 mg/dL and Triglycerides up to 900 mg/dL.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

REFERENCE VALUES

The reference values in mg/dL, for this method, were obtained through the determination of Mucoproteins in healthy populations of male and female.

Mucoprotein Tyrosine fraction: 2.4 to 4.5 mg/dL.

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

METHODS COMPARISON AND METHODOLOGICAL SPECIFICITY

The Mucoproteins kit was compared with other method for determination of Mucoprotein dosage commercially available. 07 analyses were performed and the results were evaluated. The linear equation obtained was $Y = 1.006X - 0.020$ and the correlation coefficient 0.992. With these results, we can conclude that the kit shows good methodology specificity.

Precision

REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Concentration (mg/dL)	2.13	3.41	4.69
Standard Deviation (mg/dL)	0.07	0.09	0.08
Coefficient of Variation (%)	3.08	2.50	1.68

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Concentration (mg/dL)	2.11	3.41	4.71
Standard Deviation (mg/dL)	0.02	0.01	0.02
Coefficient of Variation (%)	0.99	0.29	0.32

Sensitivity

The sensitivity was calculated from 20 determinations of a sample free of the Mucoproteins. The average found was 0.015 mg/dL, with standard deviation of 0.001 mg/dL. The sensitivity, which indicates the method limit detection, correspond the average plus 3 times the standard deviation, and is equal to 0.017 mg/dL.

Linearity

The reaction is linear up to the concentration of 15 mg/dL. For higher values, dilute the filtered (1:2 or 1:3) with distilled or deionized water. Execute a new dosage and multiply the obtained value by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

The Mucoprotein are known as acute phase protein. Considerably increase of acute inflammatory process and are an important index of the rheumatoid activity, because if they kept elevated, when other proofs are already normalized. Increased values are observed in tuberculoses, collagen diseases, mellitus diabetes, rheumatic fever, gout and neoplasia.

Reduced values are observed in nephrosis, adrenocortical insufficiency, hypophysary and hepatic.

NUMBER OF TESTS

Presentation 1 - 25 Tests / 1 mL of Sample
50 Tests / 500 mL of Sample

Presentation 2 - 100 Tests / 1 mL of Sample
200 Tests / 500 mL of Sample

BIBLIOGRAPHIC REFERENCES

- 1 - WINZLER, R. J.; DEVOR, A. W.; MEHL, J. W. and SMYTH, I. M.: J. Clin. Invests. 27:609, 1948.
- 2 - Weimer, H. E.; MOSHIN, J. R., Am. Ver. Tuberc. Pulmonary Diseases, 68:594, 1952.
- 3 - WINZLER, R. J., Methods of Biochemical Analysis, Vol. 2 nd. D. Glick, Interscience, N. Y., p.279, 1955.
- 4 - TONKS, D. B., Chem, 9:217, 1963.
- 5 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil

Phone: +55 (31) 3439.5454

E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 - Made in Brazil

EC REP **OBELIS S.A.**

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1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service

Phone: 0800 0315454

E-mail: sac@bioclin.com.br

ANVISA registration for the Mucoproteins kit: 10269360100

Review: October/2020

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MANUFACTURED BY

BATCH CODE

CONTROL

DATE OF MANUFACTURE

POSITIVE CONTROL

USED BY
(last day of month)

NEGATIVE CONTROL

TEMPERATURE LIMITATION
(store at)

BIOLOGICAL RISK

CONTAINS SUFFICIENT
FOR <N> TESTS

INFLAMMABLE

CONSULT INSTRUCTIONS
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC DEVICE

POISON

EUROPEAN AUTHORIZED
REPRESENTATIVE

CE MARK

KEEP AWAY
FROM SUNLIGHT

DO NOT USE IF
PACKAGE IS
DAMAGED