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1 - KILLINGSWORTH, L. M; SAVORY, J.: Automated Immunochemical Procedures for Measurement of Immunoglobulins IgG, IgA and IgM in human serum, Clin. Chem., 1971; 17: 936.

2 - PESCE, A. J.; KAPLAN, L. A. : Methods in Clinical Chemistry C. V. Mosby Company, 1987.

3 - HELLSING, K., Profiles in the Biological Fluids, 1973; 23: 579.

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USAGE INSTRUCTIONS

FUNCTION

Reagent used for calibration of turbidimetrics tests, for *in vitro* diagnostic use only.

REAGENTS

Number 1 - Multicalibrator - Store between 2 and 8°C. Contains: Serum Proteins and Sodium Azide 0,9%. **Potentially infectious.**

Warning: The concentration of the analytes varies according to the batch.

See table that accompanies the product.

PRESENTATION

Reagent N° 1	1 mL
--------------	------

EQUIPMENTS AND OPERATIONAL INPUTS

Automated or semi-automated biochemical devices, watch or stopwatch, pipettes, test tubes, kits and water bath at 37°C. They can be found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport, at temperatures between 15 and 30°C, should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- Reagent 1 and 2 should be handled cautiously since they are subject to biological contamination.

6- Handle with care, as the reagent contains Sodium Azide, which is irritant to skin and mucous.

7- The reagent has been tested for antibodies anti-HIV and HBs Antigen using state of the art methods and found negative results. The risk of infection cannot be excluded and the reagent should be handled as cautiously as observed for the patient's serum. **Potentially infectious.**

8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

9- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

10- Do not use the product in case of damaged packaging.

11- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

PROCESS DESCRIPTION

The utilization of the calibrator must comply with specific technique for the equipment to be used in calibrations.

The concentration of analytes varies each lot. See table that accompanies the product.

BIBLIOGRAPHIC REFERENCES

1 - KILLINGSWORTH, L. M; SAVORY, J. : Automated Immunochemical Procedures for Measurement of Immunoglobulins IgG, IgA and IgM in human serum, Clin. Chem., 1971; 17: 936.

2 - PESCE, A. J.; KAPLAN, L. A. : Methods in Clinical Chemistry C. V. Mosby Company, 1987.

3 - HELLSING, K., Profiles in the Biological Fluids, 1973; 23: 579.

4 - BLOM, M. and HJOME, H., Clinical Chemistry, 1976; 22: 657.

5 - HILLS, L. P. and TIFFANY, T. I., Clinical Chemistry, 1980; 26: 1459.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

EC	REP	OBELIS S.A.
----	-----	--------------------

Bd. Général Wahis, 53
1030 Brussels, Belgium

ANVISA registration for Multical kit: 10269360077

Review: November/2017

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED