## PROTEÍNA C REATIVA

REF K059

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## FINALIDADE

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### USAGE INSTRUCTIONS

#### FUNCTION

Method for quantitative determination of C Reactive Protein (CRP) in serum or plasma (Heparin or EDTA) samples. Immunoturbidimetric test, for *in vitro* diagnostic only.

#### PRINCIPLE OF ACTION

**Methodology:** Immunoturbidimetry

The reaction permits to quantify by a turbidimetric method, the CRP concentration in the sample. The polystyrene particles coated with anti-CRP mix with the sample forming aggregates in the presence of CRP. The agglutination process causes an increase of particle size and consequently an increase in absorbance, which is measured by comparison with the calibrator of known concentration.

#### REAGENTS

**Number 1 - Buffer** - Store between 2 and 8°C. Contains: Buffer < 200 mmol/L, stabilizers, surfactant and preservative.

**Number 2 - CRP Latex** - Store between 2 and 8°C. Contains: Suspension of Latex particles sensitized with anti-CRP antibodies < 10 µg/dL, stabilizers and preservative.

**Number 3 - Calibrator** - Store between 2 and 8°C. Contains: C Reactive Protein, stabilizers and preservative. **Potentially infectious.**

**Warning: Concentrations of C Reactive Protein vary according to the lot - See bottle label.**

#### PRESENTATION

Presentation	Reagent		
	N° 1	N° 2	N° 3
1	18 mL	2 mL	1,0 mL
2	45 mL	5 mL	1,0 mL
3	2 x 45 mL	2 x 5 mL	2 x 1,0 mL
4	20 mL	5 mL	1,0 mL
5	2 x 20 mL	2 X 5 mL	1,0 mL
6	3 x 20 mL	3 x 5 mL	2 x 1,0 mL
7	4 x 20 mL	4 x 5 mL	2 x 1,0 mL
8	40 mL	10 mL	1,0 mL
9	2 x 40 mL	2 x 10 mL	2 x 1,0 mL
10	3 x 40 mL	3 x 10 mL	3 x 1,0 mL
11	4 x 40 mL	4 x 10 mL	4 x 1,0 mL
12	60 mL	15 mL	2 x 1,0 mL
13	2 x 60 mL	2 x 15 mL	3 x 1,0 mL

#### EQUIPMENTS AND OPERATIONAL INPUTS

Biochemistry automated equipment, Multical and Multicontrol Bioclin. These items are found at markets specialized on Laboratories of Clinical Analysis.

#### TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

#### SPECIAL CARE

**1- For professional *in vitro* diagnostic use only.**

**2-** Strictly follow the methodology proposed to obtain exact results.

**3-** Water used in material cleaning must be recent and free of contaminants.

**4-** Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

**5-** Reagent N° 3 should be handled cautiously, since it is subject to biological contamination.

**6-** Handle with care, as the reagent contains Sodium Azide, which is irritant to skin and mucous.

**7-** The biological materials were tested for HIV and HBsAg using state of the art methods and found negative results. The risk of infection cannot be excluded and the reagent should be handled as cautiously as observed for the patient's serum. **Potentially Infectious.**

**8-** We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

**9-** To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.

**10-** Do not use the product in case of damaged packaging.

**11-** It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

#### SAMPLES

Serum or plasma (Heparin or EDTA) obtained free of hemolysis and intense lipemia. Samples are stable between 2 and 8°C for up to 8 days or 3 months at -20°C.

Samples with presence of fibrin should be centrifuged before testing.

#### PROCESS DESCRIPTION

The calibration stability of the C Reactive Protein kit installed on refrigerated equipment is at least 22 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

#### TECHNIQUE

The reagents are ready to use. Homogenize well the Reagent N° 2 before use.

For calibration of the reaction, utilize preferentially the kit's internal calibrator or the Bioclin Multical kit. For control use the Bioclin Multicontrol Kit.

This kit is for use in automated biochemical analyzers only. **Check the application sheet for this equipment on the website [www.bioclin.com.br](http://www.bioclin.com.br) or through SAC 0800 031 5454.**

#### PROCEDURE LIMITATIONS

Hemolysis, icterus and lipemia low or moderate do not interfere in test performance. Gross lipemia and turbid samples should be processed with a greater dilution. Values below 2 mg/L lead to few reproductive outcomes.

#### INTERFERENT

No interference was observed for Bilirubin up to 60 mg/dL, Lipids up to 900 mg/dL Hemoglobin up to 500 mg/dL and Rheumatoid Factor up to 300 IU/mL.

#### INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present an analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

#### TRACEABILITY

The kit's calibrator is traceable to the reference material NIBSC 85/506 (Human C Reactive Protein 1<sup>st</sup> International Standard - WHO International Standard).

#### REFERENCE VALUES

The reference values for this method were obtained through the determination of C-Reactive Protein in healthy populations of male and female.

Less than 8 mg/L.

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

#### PRODUCT PERFORMANCE CONTROL QUALITY

##### Accuracy

##### RECOVERY

The recovery analysis was performed with 05 determinations of samples. Accuracies were calculated and were found in good agreement with the reference values, obtaining a recovery between 97 and 106%.

#### COMPARISON OF METHODS AND METHODOLOGY SPECIFICITY

The C Reactive Protein kit was compared with other commercially available method for dosage of C Reactive Protein. 41 analyzes were performed and the results were evaluated. The linear equation obtained was  $Y = 0,9998X - 1,2981$  and correlation coefficient 0,9976. With these results we can conclude that the kit shows good methodological specificity.

**Precision****REPEATABILITY**

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
<b>Concentration (mg/L)</b>	9,47	17,22	77,82
<b>Standard Deviation (mg/L)</b>	0,08	0,10	0,51
<b>Coefficient of Variation (%)</b>	0,81	0,60	0,65

**REPRODUCIBILITY**

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
<b>Concentration (mg/L)</b>	9,47	16,73	77,15
<b>Standard Deviation (mg/L)</b>	0,09	0,40	0,70
<b>Coefficient of Variation (%)</b>	0,96	2,39	0,90

**Sensitivity**

Sensitivity was calculated from 40 determinations of a sample free of C Reactive Protein. The average found was 0,06 mg/L, with a standard deviation of 0,05 mg/L. The sensitivity, which indicates the Method Detection Limit, corresponds the average plus 3 times the standard deviation, and is equal to 0,21 mg/L.

**Linearity**

The reaction is linear to the highest point on the curve. For samples with higher values, it is recommended to dilute the sample with Sodium Chloride 0,85%, repeat the dosage and multiply the result by the dilution factor.

**High Dose Hook Effect**

Hook effect wasn't observed with a high concentration of PCR up to 400 mg/L.

**DIAGNOSTIC SIGNIFICANCE**

The C-reactive protein is a useful indicator of ongoing inflammatory process, either of infectious origin (pneumonia, tuberculosis) or noninfectious (active rheumatic fever, rheumatoid arthritis, lupus erythematosus).

Is also present in several other pathological conditions, such as acute myocardial infarction, neoplasm, intense trauma, viruses, burns. The determination of its plasma concentration is an effective test to monitor the therapy used and the prognosis of inflammation.

**NUMBER OF TESTS**

Variable according to the automated equipment used. Check the number of tests on schedule.

**BIBLIOGRAPHIC REFERENCES**

- 1 - PESCE, A. J.; KAPLAN, L. A.: Methods in Clinical Chemistry, C. V. Mosby Company, 1987.
- 2 - HELLSING, K., Profides in the Biological Fluids, 1973; 23, 579.
- 3 - BLOM, M. and HJOME, H., Clinical Chemistry, 1976; 22, 657.
- 4 - HILLS, L. P. and TIFFANY, T. I., Clinical Chemistry, 1980; 26, 1459.
- 5 - GALVIN, J. P. Et al Clin. Lab Assays 73, 4th, 1983.
- 6 - PASSING, H.; BABLOCK, W., J. Clin. Chem. Clin. Biochem 21, 709; 1983.
- 7 - NOKAMA, Y.; NAKAMURA, R. M., J. Clin. Lab. Anal., 1:15; 1987.
- 8 - OKOMURA, J. M.; MIYAGI, J. M.; TERADA, K. J., Clin. Lab. Anal. 4, 231; 1990.
- 9 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

**QUALITY ASSURANCE**

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaing, when stored and transported under appropriate conditions.

**QUIBASA QUÍMICA BÁSICA Ltda**  
 Rua Teles de Menezes, 92 - Santa Branca  
 CEP 31565-130 - Belo Horizonte - MG - Brasil  
 Phone.: +55 (31) 3439.5454  
 E-mail: bioclin@bioclin.com.br  
 CNPJ: 19.400.787/0001-07 - Made in Brazil

**OBELIS S.A.**

Bd. Général Wahis, 53  
 1030 Brussels, Belgium

**CUSTOMER SERVICE**

Customer Advisory Service  
 Phone.: 0800 0315454  
 E-mail: sac@bioclin.com.br

ANVISA registration number for C Reactive Protein kit:  
 10269360070

**Review:** July/2020

**UNIVERSAL SYMBOLOLOGY**

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED