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REF K079

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## FINALIDADE

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1. Thompson D, Milford-Ward A, Whicher JT. The value of acute phase protein measurements in clinical practice. Ann Clin Biochem 1992;29:123-31.
2. Gabay C, Kushner I. Acute-phase proteins and other systemic responses to inflammation. N Engl J Med 1999;340:448-54.
3. Hansson LO, Lindquist L. C-reactive protein: its role in the diagnosis and follow-up of infectious diseases. Curr Opin Infect Diseases 1997; 10:196-201.
4. Sipe JD. Acute-phase proteins in osteoarthritis. Semin Arthritis Rheum 1995;25:75-86.
5. Dati F, Schumann G, Thomas L, Aguzzi F, Baudner S, Bienvenu J et al. Consensus of a group of professional societies and diagnostic companies on guidelines for interim reference ranges for 14 proteins in serum based on the standardization against the IFCC/BCR/CAP reference material (CRM 470). Eur J Clin Chem Clin Biochem 1996;34:517-20.
6. Schlebusch H, Liappi N, Klein G. High sensitive CRP and creatinine: reference intervals from infancy to childhood. Poster presented at AACC/CSCC, July/August 2001, Chicago, Illinois.
7. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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1. Thompson D, Milford-Ward A, Whicher JT. The value of acute phase protein measurements in clinical practice. Ann Clin Biochem 1992;29:123-31.
2. Gabay C, Kushner I. Acute-phase proteins and other systemic responses to inflammation. N Engl J Med 1999;340:448-54.
3. Hansson LO, Lindquist L. C-reactive protein: its role in the diagnosis and follow-up of infectious diseases. Curr Opin Infect Diseases 1997; 10:196-201.
4. Sipe JD. Acute-phase proteins in osteoarthritis. Semin Arthritis Rheum 1995;25:75-86.
5. Dati F, Schumann G, Thomas L, Aguzzi F, Baudner S, Bienvenu J et al. Consensus of a group of professional societies and diagnostic companies on guidelines for interim reference ranges for 14 proteins in serum based on the standardization against the IFCC/BCR/CAP reference material (CRM 470). Eur J Clin Chem Clin Biochem 1996;34:517-20.
6. Schlebusch H, Liappis N, Klein G. High sensitive CRP and creatinine: reference intervals from infancy to childhood. Poster presented at AACC/CSCC, July/August 2001, Chicago, Illinois.
7. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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## C REACTIVE PROTEIN ULTRASENSITIVE

REF K079

### USAGE INSTRUCTIONS

#### FUNCTION

Method for quantitative determination of C-Reactive Protein (CRP) in very low concentrations. Immunoturbidimetric test, for *in vitro* diagnostic only.

#### PRINCIPLE OF ACTION

**Methodology:** Immunoturbidimetry

The reagent allows the quantitative determination of CRP in human serum by antigen-antibody reaction. In the presence of a polymer activator, which increases sensitivity and speed of the assay, C-Reactive Protein forms an insoluble complex with specific antibody, causing turbidity, whose intensity increases proportional to the concentration of CRP in the sample.

#### REAGENTS

**Number 1 - Buffer** - Store between 2 and 8°C. Contains: Buffer < 200 mmol/L, stabilizers, surfactant and preservative.

**Number 2 - Antiserum** - Store between 2 and 8°C. Contains: Suspension of Latex particles sensitized with anti-CRP antibodies < 10 µg/dL, stabilizers and preservative.

**Number 3 - Calibrator** - Store between 2 and 8°C. Contains: C Reactive Protein, stabilizers and preservative. **Potentially infectious.**

**Warning: Concentrations of C Reactive Protein vary according to the lot - See bottle label.**

#### PRESENTATION

| Presentation | Reagent   |           |          |
|--------------|-----------|-----------|----------|
|              | N°1       | N°2       | N°3      |
| 1            | 10 mL     | 10 mL     | 1 mL     |
| 2            | 25 mL     | 25 mL     | 1 mL     |
| 3            | 2 x 25 mL | 2 x 25 mL | 2 x 1 mL |
| 4            | 40 mL     | 10 mL     | 1 mL     |
| 5            | 2 x 40 mL | 2 x 10 mL | 2 x 1 mL |
| 6            | 3 x 40 mL | 3 x 10 mL | 2 x 1 mL |
| 7            | 4 x 40 mL | 4 x 10 mL | 2 x 1 mL |
| 8            | 1 x 20 mL | 1 x 5 mL  | 1 x 1 mL |
| 9            | 2 x 20 mL | 2 x 5 mL  | 2 x 1 mL |
| 10           | 3 x 20 mL | 3 x 5 mL  | 2 x 1 mL |
| 11           | 4 x 20 mL | 4 x 5 mL  | 2 x 1 mL |

#### EQUIPMENTS AND OPERATIONAL INPUTS

Biochemistry automated equipment, Multical and Multicontrol Bioclin. They can be found at markets specialized on Clinical Analysis Laboratories.

#### TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

#### SPECIAL CARE

**1- For professional *in vitro* diagnostic use only.**

**2- Strictly follow the methodology proposed to obtain exact results.**

**3- Water used in material cleaning must be recent and free of contaminants.**

**4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers, that may alter the results significantly.**

**5- Do not use particulate matter and turbidity reagents.**

**6- Handle with care, all the reagent that contains Sodium Azide, which is irritant to skin and mucous.**

**7- All the raw material of product is tested and should be non-reactive for HBsAg, HCV and HIV 1 & 2. However, these tests do not provide total assurance of the absence of infectious agents. The manual manipulation of any product containing serum is potentially capable of transmitting diseases. Therefore, we must take due care in handling the biosafety of these products.**

**8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.**

**9- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.**

**10- Do not use the product in case of damaged packaging.**

**11- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.**

#### SAMPLES

Serum or plasma (Heparin or EDTA) obtained free of hemolysis and intense lipemia. Samples are stable between 2 and 8°C for up to 8 days or 3 months at -20°C.

Samples with presence of fibrin should be centrifuged before testing.

#### PROCESS DESCRIPTION

The calibration stability of the C Reactive Protein Ultrasensitive kit installed on refrigerated equipment is at least 40 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

#### TECHNIQUE

The reagents are ready to use.

Homogenize well the Reagent N° 2 before use.

For calibration of the reaction, utilize preferentially the kit's internal calibrator or the Bioclin Multical kit. For control use the Bioclin Multicontrol Kit.

This kit is for use in automated biochemical analyzers only. **Check the application sheet for this equipment on the website [www.bioclin.com.br](http://www.bioclin.com.br) or through SAC 0800 031 5454.**

#### INTERFERENCE

No interference was observed for Bilirubin up to 24 mg/dL, Lipids up to 360 mg/dL, Hemoglobin up to 200 mg/dL and Rheumatoid Factor up to 120 UI/mL.

#### INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present an analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

#### TRACEABILITY

The kit's calibrator is traceable to the reference material NIBSC 85/506 (Human C-Reactive Protein 1<sup>st</sup> International Standard - WHO International Standard).

#### REFERENCE VALUES

The reference values for this method were obtained through the determination of C-Reactive Protein in healthy populations of male and female. Values are indicators of cardiovascular risk:

- High Risk: greater than 3,0 mg/L
- Medium Risk: 1,0 to 3,0 mg/L
- Low Risk: less than 1,0 mg/L

#### IMPORTANT

- The interpretation of values for cardiovascular risk can only be done in the absence of systemic inflammatory conditions, since their evaluation is also a parameter for definition of inflammatory conditions.

- CRP should not be used in isolation but must be associated with clinical test or complementary exams.

- These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

- The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

#### PRODUCT PERFORMANCE QUALITY CONTROL

##### Accuracy RECOVERY

The recovery analysis was performed with 05 determinations of samples. Accuracies were calculated and were found in good agreement with the reference values, obtaining a recovery between 94 and 107%.

##### COMPARISON OF METHODS AND METHODOLOGY SPECIFICITY

The C Reactive Protein Ultrasensitive kit was compared with other commercially available methods for determination of CRP dosage in extremely low concentrations. 41 analyses were performed and the results were evaluated. The linear equation obtained was  $Y = 0,9982X + 0,0374$  with a correlation coefficient equal to 0,9988. With these results we can conclude that the kit has good methodological specificity.

## Precision

### REPEATABILITY

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

|                              | Sample 1 | Sample 2 | Sample 3 |
|------------------------------|----------|----------|----------|
| Concentration (mg/L)         | 0,99     | 2,02     | 10,71    |
| Standard Deviation (mg/L)    | 0,03     | 0,03     | 0,05     |
| Coefficient of Variation (%) | 2,98     | 1,45     | 0,45     |

### REPRODUCIBILITY

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

|                              | Sample 1 | Sample 2 | Sample 3 |
|------------------------------|----------|----------|----------|
| Concentration (mg/L)         | 1,00     | 2,04     | 10,66    |
| Standard Deviation (mg/L)    | 0,04     | 0,04     | 0,08     |
| Coefficient of Variation (%) | 4,51     | 2,04     | 0,79     |

## Sensitivity

Sensitivity was calculated from 40 determinations of a sample free of C-Reactive Protein. The average found was 0,031 mg/L with standard deviation of 0,015 mg/L. The sensitivity, which indicates the method detection limit, corresponds the average plus 3 times the standard deviation, and is equal to 0,076 mg/L.

## Linearity

The reaction is linear up to the highest point on the curve. For samples with higher values, it is recommended to dilute the sample with NaCl 0,85%, repeat the dosage and multiply the result by the dilution factor.

## High Dose Hook Effect

Hook effect wasn't observed with a high concentration of PCR up to 250 mg/L.

## DIAGNOSTIC SIGNIFICANCE

The levels of C-Reactive Protein (CRP) have been recommended for cardiac risk assessment. CRP levels above the reference value, in apparently healthy patients, may signal coronary heart disease, and thus is an important indicator for AMI and stroke. Importantly, the assessment of CRP as an indicator of cardiac risk can be employed only in the absence of inflammatory conditions, since in those situations, their values are greatly increased.

## NUMBER OF TESTS

Variable according to the automated equipment used. Check the number of tests on schedule.

## BIBLIOGRAPHIC REFERENCES

1. Thompson D, Milford-Ward A, Whicher JT. The value of acute phase protein measurements in clinical practice. *Ann Clin Biochem* 1992;29:123-31.
2. Gabay C, Kushner I. Acute-phase proteins and other systemic responses to inflammation. *N Engl J Med* 1999;340:448-54.
3. Hansson LO, Lindquist L. C-reactive protein: its role in the diagnosis and follow-up of infectious diseases. *Curr Opin Infect Diseases* 1997; 10:196-201.
4. Sipe JD. Acute-phase proteins in osteoarthritis. *Semin Arthritis Rheum* 1995;25:75-86.
5. Dati F, Schumann G, Thomas L, Aguzzi F, Baudner S, Bienvendu J et al. Consensus of a group of professional societies and diagnostic companies on guidelines for interim reference ranges for 14 proteins in serum based on the standardization against the IFCC/BCR/CAP reference material (CRM 470). *Eur J Clin Chem Clin Biochem* 1996;34:517-20.
6. Schiebisch H, Liappis N, Klein G. High sensitive CRP and creatinine: reference intervals from infancy to childhood. Poster presented at AACC/SCCC, July/August 2001, Chicago, Illinois.
7. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

## QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

### QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31585-130 - Belo Horizonte - MG - Brasil  
Phone.: +55 (31) 3439.5454  
E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

EC REP

OBELIS S.A.

Bd. Général Wahis, 53  
1030 Brussels, Belgium

## CUSTOMER SERVICE

Customer Advisory Service  
Phone.: 0800 0315454  
E-mail: sac@bioclin.com.br

ANVISA registration for C Reactive Protein Ultrasensitive kit: 10269360207

Review: January/2020

## UNIVERSAL SYMBOLOGY

|  |                                       |  |                                        |
|--|---------------------------------------|--|----------------------------------------|
|  | CATALOG NUMBER                        |  | MANUFACTURED BY                        |
|  | BATCH CODE                            |  | CONTROL                                |
|  | DATE OF MANUFACTURE                   |  | POSITIVE CONTROL                       |
|  | USED BY<br>(last day of month)        |  | NEGATIVE CONTROL                       |
|  | TEMPERATURE LIMITATION<br>(store at)  |  | BIOLOGICAL RISK                        |
|  | CONTAINS SUFFICIENT<br>FOR <N> TESTS  |  | INFLAMMABLE                            |
|  | CONSULT INSTRUCTIONS<br>FOR USE       |  | CORROSIVE                              |
|  | IN VITRO DIAGNOSTIC DEVICE            |  | POISON                                 |
|  | EUROPEAN AUTHORIZED<br>REPRESENTATIVE |  | CE MARK                                |
|  | KEEP AWAY<br>FROM SUNLIGHT            |  | DO NOT USE IF<br>PACKAGE IS<br>DAMAGED |