FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Immunochromatographic method for rapid qualitative determination of Occult Blood on feces. For *in vitro* diagnostic only.

PRINCIPLE OF ACTION

Methodology: Immunochromatography

The method uses pre-covered membranes with monoclonal antibody to identify selectively the human hemoglobin in feces samples. After the collection and preparation of the sample, the buffered extracted sample is added to the membrane. If there is human hemoglobin present in the sample it will bind to a conjugate antibody with colloidal gold and migrates through the membrane by capillarity until it finds a second fixed antibody in the region if the test line membrane. In that case, a red line will appear in the test line region. In the absence of human hemoglobin, no line will appear in the test zone. The conjugate not captured continues its way until the control line, confirming the correct function of the test.

REAGENTS

Number 1 - Cassette: Store between 15 and 30°C. **Do not freeze.** Containing:

- Control Region: Polyclonal goat Antibody and mouse Anti-IgG o immobilized in membrane.
- Test Region: Monoclonal mouse Antibody Human Anti-Hemoglobin.
- Conjugate: Monoclonal mouse Antibody Anti-Hemoglobin conjugate with Coloidal Gold.

Number 2 - Diluent: Contained in **analysis tube**; Store between 15 and 30°C. Contains: Buffer PBS.

PRESENTATION

Presentation	Packing Number of Cassettes	Diluent
1	1	1 x 2 mL
2	5	5 x 2 mL
3	10	10 x 2 mL
4	20	20 x 2 mL
5	25	25 x 2 mL
6	30	30 x 2 mL
7	40	40 x 2 mL
8	50	50 x 2 mL
9	100	100 x 2 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Watch or stopwatch. These itens are found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The transport and storage temperatures must be between 15 and 30°C. Avoid moisture.

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain accurate results.

3- The stripes do not show contaminants, but because the samples are biological, safety routines should be followed accordingly and the use of disposable gloves is essential.

4- Use clean and dry containers for storage of samples.

5- Open the aluminum envelope, to remove the cassette, only by the time of use.

6- The kit should not be frozen or exposed to high temperatures, as this will cause deterioration.

7- The reagents and the sample must be equilibrated with the environment before starting the test.

8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

9- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

10- Do not use the product in case of damaged packaging.

11- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Only for tests with human feces.

The feces samples should be collected in clean and dry recipients. Avoid contamination from water and other substances from the bathroom.

The samples shall be kept in a fresh place under control temperatures. Samples must be tested within 6 hours after collection. To store it for a longer period, freeze at -20°C for 1 week or 3 days at 2 - 8°C.

The stool sample inside the analysis tube (R2 - Diluent) is stable for 3 days at room temperature or for 7 days at 2 to 8°C.

The samples must always be in environment temperature to carry out the test.

Samples should not be collected during the menstrual period.

Samples of patients suffering from hemorrhoids, bleeding or presenting blood in the urine may produce inconsistent results.

The use of medicines irritating to the gastric mucosa (anti-inflammatories, corticoids or iron and vitamin C based) should be stopped at least 2 days before collecting the sample.

Alcohol in excess can also cause gastrointestinal irritation. Therefore, the use of such substances should not be performed at least 72 hours before testing.

Dietary restriction is not necessary for testing.

PROCESS DESCRIPTION

The stability of the Hidden Blood cassette (after opening the aluminized sachet) is 2 hours. This stability may vary according to the conditions of the environment after opening the aluminum sachet.

TECHNIQUE

The sample should be at room temperature between 15 and 30°C (room temperature) before starting the test.

1- Open the analysis tube containing the diluent solution and insert the dose handle in three distinct points of the sample.

2- Close the analysis tube with the dose handle infused by the feces, vigorously agitate from 3 to 5 minutes to emulsify the feces, and guarantee the disperssion and dissolution (a vortex is allowed).

3- Remove the cassette from its envelope and appropriately identify it so.


4- Open the lower extremity of the analysis tube, so it will function as a dropper.

5- Cassette: Place 2 drops in the circular cassette window allowing the first drop to be absorbed so that it put the next one.


6- Read the result between 10 and 15 minutes. Do not read the result after 15 minutes.

INTERPRETATION OF RESULTS



Reagent Test: Formation of two red lines, one in the control region (C) and one in the test region (T) after 10 minutes. Do not interpret after 15 minutes.

Non-Reagent Test: Formation of only one red line, after 10 minutes, in the Control region. Do not interpret after 15 minutes.

Inadequate Testing: The absence of line formation in the Control region (C) indicates an error in the procedure or deterioration of the cassette. In this case, repeat the test using a new cassette.

PROCEDURE LIMITATIONS

1- The test is qualitative although there is a proportionality in the intensity of the color obtained in the test line.

2- The test results must be interpreted along with the patient's clinical evaluation available and other diagnostics procedures.

3- It is important to use the correct volume of sample, since lower or higher volumes may determine erroneous results.

4- The reading time of the reaction must be followed according to established technique in order to avoid misinterpretations results.

5- The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

INTERFERENT

No interference was observed for Ascorbic Acid up to 20 mg/dL and Bilirubin up to 20 mg/dL.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal program of quality control, where procedures, norms, limits and tolerance for variations are clearly established. It is important to note that all measuring systems characteristic analytical variability, which should be monitored by the laboratories themselves. Therefore, it is the use of controls is precision and accuracy of dosages.

PRODUCT PERFORMANCE
QUALITY CONTROL

Accuracy
COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY
The Occult Blood kit was compared with another method commercially available to measure Ocult Blood. 648 analysis were conducted. The following results were founded:

Comparison of Methods		Test Kit Results		Subtotal
		+	-	
Reference Kit Results	+	28	0	28
	-	2	618	620
Subtotal		30	618	648

The results lead to the following values:
Sensitivity: > 99,9% (28/28)
Specificity: 99,68% (618/620)
Concordance: 99,69% (646/648)

Precision
REPEATABILITY
The repeatability was calculated from 20 successive determinations, using 5 samples with different concentrations, obtaining the following results:

Samples	Nº of Repeatability	Expected Result	Obtained Result
Feces 01	20	100% Reagent	100% Reagent
Feces 02	20	100% Non-Reagent	100% Non-Reagent
Feces 03	20	100% Reagent	100% Reagent
Feces 04	20	100% Non-Reagent	100% Non-Reagent
Feces 05	20	100% Non-Reagent	100% Non-Reagent

REPRODUCIBILITY
The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 5 samples with different concentrations, obtaining the following results:

Samples	Nº of Days	Nº of Repeatability per Day	Expected Result	Obtained Result
Feces 01	03	20	100% Reagent	100% Reagent
Feces 02	03	20	100% Non-Reagent	100% Non-Reagent
Feces 03	03	20	100% Reagent	100% Reagent
Feces 04	03	20	100% Non-Reagent	100% Non-Reagent
Feces 05	03	20	100% Non-Reagent	100% Non-Reagent

Sensitivity
The test is capable of detecting concentrations of human hemoglobin equals or higher than 50 ng/mL.

Specificity
The Occult Blood kit, was evaluated for possible cross reaction with suine, chicken, bovine, mouse, duck, fish and goat blood. No significant interference was found in 20 samples with positive or negative results on the feces.
DIAGNOSTIC SIGNIFICANCE
The presence of fecal Occult Blood is directly associated with gastrointestinal disorders such as Hemorrhoid, Anal Fissure, Diverticulitis, Polyps, Crohn's Disease, Colorectal Cancer, Gastric or Duodenal Ulcers, Intestinal Infections, among others. Small amounts of blood in the stool or bleeding detectable only after cleaning the anus with toilet paper are the most common forms of rectal bleeding. In 90% of cases, the etiology is benign and corresponds mainly hemorrhoids and anal fissures. When the amount of blood in the stool is moderate to large, or when there melaena (bloody stools digested), the source of bleeding is usually more internal, usually colon or stomach. Colorectal cancer is a major cause. Immunologic tests used as screening and developed to detect human hemoglobin are more accurate and do not require special diets for patients. Although the screening test is not specific to determine what is causing the disease bleeding, it serves as an important early alert to the clinical research. Early diagnosis and prompt treatment showed significant reduction of complications and mortality from colorectal cancer.

BIBLIOGRAPHIC REFERENCES
1- NaKama H. et al., A comparative study of immunochemical fecal testes for detection of colorectal adenomatous polyps. Hepato-Gastroenterology. 2000, 47:386-389.
2- Uchida K. et al., Immunochemical detection of human blood in feces. Clinica Chimica Acta. 1990, 189:267-274.
3- QUIBASA: dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE
Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda**
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

EC	REP
----	-----

OBELIS S.A.
Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE
Customer Advisory Service
Phone: 0800 0315454
Email: sac@bioclin.com.br

ANVISA registration for the Occult Blood kit: 10269360202

Review: July/2020

UNIVERSAL SYMBOLOGY

 CATALOG NUMBER

 BATCH CODE

 DATE OF MANUFACTURE

 USED BY
(last day of month)

 TEMPERATURE LIMITATION
(store at)

 CONTAINS SUFFICIENT
FOR <N> TESTS

 CONSULT INSTRUCTIONS
FOR USE

 IN VITRO DIAGNOSTIC DEVICE

 EUROPEAN AUTHORIZED
REPRESENTATIVE

 KEEP AWAY
FROM SUNLIGHT

 MANUFACTURED BY

 CONTROL

 POSITIVE CONTROL

 NEGATIVE CONTROL

 BIOLOGICAL RISK

 INFLAMMABLE

 CORROSIVE

 POISON

 CE MARK

 DO NOT USE IF
PACKAGE IS
DAMAGED