FINALIDADE

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REFERÊNCIAS BIBLIOGRÁFICAS
1 - Hirsh J , Deykin D, Poller L: Therapeutic Range For Oral Anticoagulant Therapy. Chest 89/12 (Suppl) 115-155, 1986.
2 - International Committee For Standardization In Haematology And International Committee On Trombosis And Haemostasis. Amer. J. Clin. Path . 88:779 1985.
3 - National Committee For Clinical Laboratory Standards: Colleccion, Transport And Processing Of Blood Specimens For Coagulation Testing And Performance Of Coagulation Assays, 1991. Nccls.
4 - Douketis Jd, Lane A , Milne J , Ginsberg Js. Trombosis Research , 1998 ; 92 : 11-17.
5 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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REFERENCIAS BIBLIOGRÁFICAS
1 - Hirsh J , Deykin D, Poller L: Therapeutic Range For Oral Anticoagulant Therapy. Chest 89/12 (Suppl) 115-155, 1986.
2 - International Committee For Standardization In Haematology And International Committee On Trombosis And Haemostasis. Amer. J. Clin. Path . 88:779 1985.
3 - National Committee For Clinical Laboratory Standards: Colleccion, Transport And Processing Of Blood Specimens For Coagulation Testing And Performance Of Coagulation Assays, 1991. Nccls.
4 - Douketis Jd, Lane A, Milne J, Ginsberg Js. Trombosis Research, 1998; 92: 11-17.
5 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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USAGE INSTRUCTIONS

FUNCTION

Method for determination of manual or automated of Prothrombin Time in citrated plasma. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Coagulometry

Consists of measuring the Clotting Time of plasma after addition a source of Thromboplastin (factor III, tissue factor) and Calcium. The factor VII is then activated to form a complex (complex FT-FVIIa) that will activate factors IX and factor Xa XO with Phospholipids factor tissue, Calcium and factor Va form the Prothrombin activator complex which transforms Prothrombin (factor II) into Thrombin (factor IIa). Thrombin acts on Fibrinogen, transforming it into Fibrin. This test assesses the factors of Prothrombin Complex (factors II, V, VII and X).

REAGENT

Number 1 - Clot Former - Store between 2 and 8°C. Contains: Rabbit brain extract 2.6%, Calcium Chloride, buffer and preservative.
Presentation ready for use.

PRESENTATION

Presentation	Reagent 1
1	10 x 5 mL
2	20 x 2 mL
3	10 x 2 mL
4	20 x 5 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Log-log paper, pipette and tips, tubes, coagulometer or stopwatch, Biocontrol N and P Coagulation Bioclin. These itens are found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature is between 2 and 8°C and the transport temperature shall be up to 30°C, not exceeding 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- **For professional *in vitro* diagnostic use only.**
- 2- Follow strictly the proposed methodology to obtain results accurate.
- 3- Handle with care reagents because they contain Sodium Azide, irritating to skin and mucous membranes.
- 4- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 5- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Safety Data Sheet Safety of Chemicals) available on the website www.bioclin.com.br or upon request by the SAC (Advisory Service Customer) of Quibasa.
- 6- Do not use the product in case of damaged packaging.
- 7- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Material: Plasma

Anticoagulant: Use the kit as anticoagulant Citrate Bioclin A001 as written in usage instructions of the same, or 3.8% Sodium Citrate (0.130 M) at a ratio of 9 parts blood to 1 part anticoagulant (Example: 4.5 mL blood + 0.5 mL of anticoagulant).

Collection: Blood should be obtained by venipuncture, avoiding hemolysis, prolonged tourniquet, blistering and aspiration tissue fluid.

Preparation: Blood should be mixed with the anticoagulant soon after collection. Centrifuge immediately at 3000 rpm for 15 minutes. Remove the plasma without cellular pipette red or yellow layer. The samples should be tested in less than three hours. If the test can not be done in this period the plasma should be frozen for at maximum 1 week at -20°C.

PROCESS DESCRIPTION

The calibration stability of the PT Bioclin kit installed on refrigerated equipment is at least 14 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Coagulation Bioclin kits.

- 1- Homogenize Reagent N°1 before use.
- 2- Preheat Reagent N°1 to 37°C for 4 to 5 minutes. This procedure should be performed only with the volume that will be used to perform the tests. The reagent must remain warm for the entire procedure.
- 3- Separate the reaction tubes that will be used to perform the tests.
- 4- Pipette 100 µL of the Sample/Control into the respective reaction tube and incubate at 37°C for 2 to 3 minutes in a water bath or thermal block.
- 5- Pipette 200 µL of Reagent N°1 into the Sample/Control tube (preheated to 37°C) in the Sample/Control tube, simultaneously starting the timer.
- 6- Leave the tube in the water bath at 37°C, shaking gently. Keep in the bath for about 09 seconds.
- 7- Remove the tube from the bath, gently tilt once per second and stop the timer at the time of the clot formation. Note the clotting time.

Note: After opened, reagent N° 1 is stable for 14 days at 2 to 8°C.

CALCULATIONS AND RESULTS

The results of Prothrombin Time may be expressed as:

Quick Time: In seconds.

Percentage of activity: Related to a normal plasma (considered 100% activity), and for this you must draw a curve activity, using a normal plasma pool and freshly harvested and diluted with saline (0.9%).

Curve

Dilution	1:1	1:2	1:3	1:4	1:8	1:10
Pool (mL)	0.2	0.2	0.2	0.2	0.2	0.2
Physiologic Solution (mL)	---	0.2	0.4	0.6	1.4	1.8
Activity	100%	50%	33%	25%	12.5%	10%

In a log-log paper, plot the graph by putting the seconds value on the “y” axis and the percentage value on the “x” axis. This procedure should be repeated on every batch.

INR (International Normalized Ratio)

To obtain INR must be related to the prothrombin time ratio the patient by the time of normal plasma with the ISI (Sensitivity Index International), as shown below. The value of ISI is sent with the kit.

INR = (PT patient / PT normal)^{ISI}
Patient PT: 13.5 sec
Normal PT (reference): 12 sec
ISI: 1.20

Example:
INR = (13.5 / 12)^{1.20} = 1.15

INTERFERENCES

Excessive hemolysis and lipemia. Use of plasma with Heparin and EDTA.
The PT can be increased by substances such as Anticoagulants Heparin and Warfarin, Oral Contraceptives, Ethanol, Tetracycline, Corticosteroids and EDTA.
The reduction of PT is observed in patients taking: Vitamin K, Anti-Histamines, Caffeine and Phenobarbital.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The PT Bioclin kit is traceable to the reference material NIBSC RBT/05 (International Standard for Thromboplastin - WHO International Standard).

REFERENCE VALUES

Prothrombin time (quick) - 10 to 14 seconds
Activity (%) - 70 to 100%

INR - 1.0 to 1.08 in healthy people
- 2.0 to 3.5 in patients taking oral anticoagulants

Each laboratory should establish its reference values, and the results vary with respect to the reagent used. The results should be expressed in INR to normalize the variation of results between laboratories.

ISI values are determined by comparison with a reference material for thromboplastin primary. As low as the ISI is, the more sensitive is the reagent.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE
QUALITY CONTROL

Accuracy
COMPARISON OF METHODS AND SPECIFICITY
METHODODOLOGY
PT Bioclin kit was compared with another method for evaluating Prothrombin Time commercially available.
Were performed 07 analysis and the results were evaluated. The obtained linear equation was $Y = 1.0804X - 1.2695$ and correlation coefficient 0.9954. With these results, we can conclude that the kit shows good methodological specificity.

Specificity
The specificity of the kit is provided by the addition of an exclusive source thromboplastin, specific activator of factor VII, which eliminates the possibility of direct activation of the cascade of other factors coagulation.

Accuracy
REPEATABILITY
The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Time (sec)	12.65	24.85	17.45
Standard Deviation (%)	0.49	0.75	0.51
Coefficient of Variation (%)	3.87	3.00	2.93

REPRODUCIBILITY
The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Time (sec)	12.57	24.77	17.40
Standard Deviation (%)	0.19	0.08	0.05
Coefficient of Variation (%)	1.51	0.31	0.29

Sensitivity
The analytical sensitivity of the kit is given by the value of Sensitivity Index International (ISI) and ranges from 1.0 to 1.2.

DIAGNOSTIC SIGNIFICANCE
The test measures the Prothrombin Time, this increased time is due to congenital or acquired deficiencies of factors I, II, V, VII and X. Being thus, the test of choice for assessing the reduction of the factors vitamin K dependent, and indicated in the control patients using oral anticoagulants, hemorrhagic disease of the newborn, hepatic failure, intravascular coagulation, disordered eating and absorption of vitamin K. Evaluates both the extrinsic pathway as the common pathway of coagulation.

NUMBER OF TESTS
Presentation 1 - 250 tests
Presentation 2 - 200 tests
Presentation 3 - 100 tests
Presentation 4 - 500 tests

BIBLIOGRAPHIC REFERENCES
1 - Hirsh J, Deykin D, Poller L: Therapeutic Range For Oral Anticoagulant Therapy. 89/12 Chest (Suppl) 115-155, 1986.
2 - International Committee For Standardization In Haematology International Committee On Thrombosis And Haemostasis. Amer. J. Clin. Path. 88:779 1985.
3 - National Committee For Clinical Laboratory Standards: Colleccion, Transport And Processing Of Blood Coagulation Specimens For Testing And Performance Of Coagulation Assays, 1991. NCCLS.
4 - Douketis JD, Lane A, Milne J, Ginsberg JS. Thrombosis Research, 1998, 92: 11-17.
5 - Quibasa: Data from the Department of Research and Development.

QUALITY ASSURANCE
Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda**
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

 **OBELIS S.A.**
Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE
Customer Advisory Service
Phone: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA resgistration for PT Bioclin kit: 10269360151

Review: June/2020

UNIVERSAL SYMBOLOGY

 CATALOG NUMBER	 MANUFACTURED BY
 BATCH CODE	 CONTROL
 DATE OF MANUFACTURE	 POSITIVE CONTROL
 USED BY (last day of month)	 NEGATIVE CONTROL
 TEMPERATURE LIMITATION (store at)	 BIOLOGICAL RISK
 CONTAINS SUFFICIENT FOR <N> TESTS	 INFLAMMABLE
 CONSULT INSTRUCTIONS FOR USE	 CORROSIVE
 IN VITRO DIAGNOSTIC DEVICE	 POISON
 EUROPEAN AUTHORIZED REPRESENTATIVE	 CE MARK
 KEEP AWAY FROM SUNLIGHT	 DO NOT USE IF PACKAGE IS DAMAGED