FINALIDADE

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- 1 - The Committee on Enzymes of the Scandinavian society for Clinical Chemistry and Clinical Physiology Scand J. Clin Lab Invest., 1974, 33, 291 - 306.
- 2 - BERGMEYER, Bowers and cols., Clin. Chim., Acta., 1.976, 70, 19 - 42, 1977, 21 - 22.
- 3 - BERGMEYER, HV.; SCHEIBE, P.; WAHLEFELD, A.W., Clin., Chem., 1.978, 24, 58 - 73.
- 4 - Expert Panel on Enzymes of the International Federation of Clinical Chemistry : Part 3. Revised IFCC Method for Aspartate Aminotransferase, Clin., Chem., 1.978, 24, 720.
- 5 - Scandinavian Committee on Enzymes Scand J. Clin Lab Invest, 1.981, 41, 107 - 116.
- 6 - BURTIS; CARL, A.;ASHWOOD; EDWARD, R., Clin. Chem.,Tietz Text Book of; 2a ed., 1.986, 788 - 797.
- 7 - PESCE, A., J.; KAPLAN, L., A., Methods in Clin. Chem., 1.987.
- 8 - IFCC Reference Procedure for the Measurement of Catalytic Concentration of Alanine Aminotransferase. Clin Chem Lab Med 2002; 40(7):718-24.
- 9 - Soldin SJ, Brugnara C, Wong EC: Pediatric Reference Intervals, 5.ed.Washington: AACCC Press, 2005.p.3-4.
- 10 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- 11 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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REFERENCIAS BIBLIOGRÁFICAS
1 - The Committee on Enzymes of the Scandinavian society for Clinical Chemistry and Clinical Physiology Scand J. Clin Lab Invest., 1974, 33, 291 - 306.
2 - BERGMAYER, Bowers and cols., Clin. Chim., Acta., 1.976, 70, 19 - 42, 1977, 21 - 22.
3 - BERGMAYER, HV.; SCHEIBE, P.; WAHLEFELD, A.W., Clin., Chem., 1.978, 24, 58 - 73.
4 - Expert Panel on Enzymes of the International Federation of Clinical Chemistry : Part 3. Revised IFCC Method for Aspartate Aminotransferase, Clin., Chem., 1.978, 24, 720. 5 - Scandinavian Committee on Enzymes Scand J. Clin Lab Invest, 1.981, 41, 107 - 116.
6 - BURTIS; CARL, A.;ASHWOOD; EDWARD, R., Clin. Chem.,Tietz Text Book of; 2a ed., 1.986, 788 - 797.
7 - PESCE, A., J.; KAPLAN, L., A., Methods in Clin. Chem., 1.987.
8 - IFCC Reference Procedure for the Measurement of Catalytic Concentration of Alanine Aminotransferase. Clin Chem Lab Med 2002; 40(7):718-24.
9 - Soldin SJ, Brugnara C, Wong EC: Pediatric Reference Intervals, 5.ed.Washington: AACC Press, 2005.p.3-4.
10 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
11 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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USAGE INSTRUCTIONS


FUNCTION

Method for determination of Alanine Amino Transferase (ALT) in samples of serum or plasma (EDTA or Heparin). Kinetic test, for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Kinetic (UV)

Kinetic determination (UV) of ALT, according to the reaction:

ALT catalyzes the transfer of Amino groups of Alanine for α -Ketoglutarate, leading to formation of Pyruvate and Glutamate. Pyruvate in the presence of LDH reacts with NADH, reducing the Lactate and oxidizes NADH to NAD⁺. The rate of oxidation is proportional to ALT activity in the sample.

REAGENTS

Number 1 - Substrate - Store between 2 and 8°C. Contains: Tris Buffer < 200 mmol/L (pH 7.8), LDH 2400 U/L, L-Alaninine < 500 mmol/L, stabilizer and preservative.

Number 2 - Coenzyme - Store between 2 and 8°C. Contains: Tris Buffer < 200 mmol/L (pH 10), Alpha-Ketoglutarate < 100 mmol/L, NADH < 5 mmol/L and preservative.

PRESENTATION

Presentation	Reagent N° 1	Reagent N° 2
1	54 mL	6 mL
2	27 mL	3 mL
3	2 x 27 mL	2 x 3 mL
4	4 x 27 mL	4 x 3 mL
5	1 x 40 mL	1 x 10 mL
6	2 x 40 mL	2 x 10 mL
7	4 x 40 mL	4 x 10 mL
8	2 x 40 mL	1 x 20 mL
9	4 x 40 mL	2 x 20 mL
10	4 x 27 mL	1 x 12 mL
11	6 x 36 mL	1 x 24 mL
12	6 x 36 mL	2 x 12 mL
13	3 x 36 mL	1 x 12 mL
34	4 x 80 mL	4 x 20 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Thermostated spectrophotometer, pipettes, watch or stopwatches, test tubes, Biocontrol N, Biocontrol P Bioclin and Bical Bioclin. These items are found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures between 15 and 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning and preparation of the working reagent must be recent and free of contaminants.

4- It is important, for the good development of the test, a rigorous control of time, temperature and pH.

5- Reagent N° 1 contains Sodium Azide, and should be handled with care.

6- Lipemic and icteric samples increase the absorbance at 340 nm. In this case, you should dilute the sample 1:2 with Saline Solution. Multiply the results by 2.

7- Readings of absorbance less than 0.800 of Reagent Working indicate loss of it. In this case, do not use the reagent.

8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

9- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Advisory Service Customer) of Quibasa.

10- Do not use the product in case of damaged packaging.

11- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or plasma with EDTA or heparin, obtained free of hemolysis. The sample is stable for 3 days between 2 and 8°C, and 3 months at -20°C ¹⁰.

PROCESS DESCRIPTION

The calibration stability of the Transaminase ALT Kinetic kit installed on refrigerated equipment is at least 12 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

PREPARATION OF WORKING REAGENT

Mix 4 parts of Reagent N° 1 with 1 part of Reagent N° 2. Working Reagent is stable during 72 hours between 15 and 30°C and 14 days between 2 and 8°C.

REACTION CONDITIONS

Is indispensable condition the use of thermostated cuvette at 37°C, 1 cm optical path and reading at 340 nm (334 - 365 nm).

TECHNIQUE

Bioclin recommends using the kit Biocal Bioclin as calibrator and as control serum, Biocontrol N and P Bioclin Kits.

Add 100 μ L of Sample to 1.0 mL of Working Reagent, mix and transfer to a thermostated cuvette at 37°C and wait for 1 minute. Make the initial reading, simultaneously starting the timer. Repeat readings after 1, 2 and 3 minutes. Calculate the average differences in absorbance per minute ($\Delta A/\text{min.}$) and use it to calculate the result.

CALCULATIONS

ALP (U/L) 340 nm = $\Delta A/\text{min.}$ x 1746
 334 nm = $\Delta A/\text{min.}$ x 1780
 365 nm = $\Delta A/\text{min.}$ x 3235

Results are expressed in U/L.

For an average variation in absorbance ≥ 0.15 in 340 nm and 334 nm or ≥ 0.080 in 365 nm, repeat the determination diluting the sample with NaCl 0.85%. Multiply the results obtained by the dilution factor.

PROCESS LIMITATIONS

The kinetic method is based on the molar absorptive. By this reason, the readings must be conducted in a spectrophotometer that satisfies the following conditions: Wavelength 340 nm
Semi trajectory of the pass band 10 nm
Stray light less than 0.5%
1 cm thermostated cuvette

INTERFERENT

No interference was observed by Bilirubin until 19 mg/dL, Hemoglobin until 180 mg/dL and Triglycerides until 650 mg/dL.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit calibration can be made using the theoretical calibration factor, based on the molar absorptivity of the NADH, or through the BIOCAL calibrator. Bioclin recommends the usage of the BIOCAL calibrator which is traceable to the reference material ERM-AD454.

REFERENCE VALUES

The reference values, in U/L, for this method were obtained through the determination of ALT in healthy populations of male and female.

	Male (U/L)	Female (U/L)
1 - 30 days	20 - 54	21 - 54
1 - 6 month	26 - 55	26 - 61
7 - 12 month	26 - 59	26 - 55
1 - 3 years	19 - 59	24 - 59
4 - 11 years	24 - 49	24 - 49
12 - 15 years	24 - 59	19 - 44
Adults	≤ 45	≤ 34

These values should be use as guidance, being that every laboratory must create their own range of values for reference, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE
QUALITY CONTROL

Accuracy
RECOVERY

The recovery analysis was performed with 05 determinations of sample. Accuracies were calculated and were found in good agreement with the reference values, obtaining a recovery between 95% and 102%.

COMPARISON OF METHODS AND METHODOLOGICAL
SPECIFICITY

The Transaminase ALT Kinetic kit was compared with another method commercially available to measure Alanine Amino Transferase. 42 analyzes were conducted and the results were evaluated. The linear equation obtained was Y = 0.9693X + 1.33, with correlation coefficient 0.9995. With these results we can conclude the kit shows good methodological specificity.

Precision
REPEATABILITY

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (U/L)	34.05	146.95	90.55
Standard Deviation (U/L)	0.22	0.64	0.50
Coefficient of Variation (%)	0.65	0.43	0.56

REPRODUCIBILITY

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (U/L)	34.11	144.88	89.83
Standard Deviation (U/L)	0.47	1.97	1.11
Coefficient of Variation (%)	1.36	1.36	1.24

Sensitivity

Sensitivity of Transaminase ALT Kinetic kit was calculated from 40 determinations of a sample free of Alanine Amino Transferase. The average found was of 0.050 U/L with a standard deviation of 0.316 U/L. The sensitivity, which indicates the method detection limit, corresponds the average plus 3 times the standard deviation and is equal to 0.998 U/L.

Linearity

This reaction is linear up to the concentration of 400 U/L. For samples with values higher than 400 U/L it is recommended to dilute the sample with Sodium Chloride 0.85%, repeat the dosage and multiply the results obtained by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

The increased activity of the Aspartate Amino Transferase enzymes - AST (location cytho-mitochondrial) and Alanine Amino Transferase - ALT (both cytoplasmic) reflects changes in various tissues.

The highest activity of ALT is located in the liver tissue. Minor activities occur in skeletal muscle, heart, kidney and pancreas. Its activity is increased in viral hepatitis and toxic (30 - 50 or 100 times the reference values - RV), as well in other liver diseases (DH), associated with necrosis liver. In the DH-associated chronic cell necrosis, due to increased release of AST - mitochondria, can be inversion ratio of ALT/AST. Occurs, further increasing their levels in infectious mononucleosis and the intra and extra liver cholestasis.

BIBLIOGRAPHIC REFERENCES

- 1 - The Committee on Enzymes of the Scandinavian society for Clinical Chemistry and Clinical Physiology Scand J. Clin Lab Invest., 1974, 33, 291 - 306.
- 2 - BERGMEYER, Bowers and cols., Clin. Chim., Acta., 1.976, 70, 19 - 42, 1977, 21 - 22.
- 3 - BERGMEYER, HV.; SCHEIBE, P.; WAHLEFELD, A.W., Clin., Chem., 1.978, 24, 58 - 73.
- 4 - Expert Panel on Enzymes of the International Federation of Clinical Chemistry : Part 3. Revised IFCC Method for Aspartate Aminotransferase, Clin., Chem., 1.978, 24, 720.
- 5 - Scandinavian Committee on Enzymes Scand J. Clin Lab Invest, 1.981, 41, 107 - 116.
- 6 - BURTIS; CARL, A.;ASHWOOD; EDWARD, R., Clin. Chem.,Tietz Text Book of; 2a ed., 1.986, 788 - 797.
- 7 - PESCE, A., J.; KAPLAN, L., A., Methods in Clin. Chem., 1.987.
- 8 - IFCC Reference Procedure for the Measurement of Catalytic Concentration of Alanine Aminotransferase. Clin Chem Lab Med 2002; 40(7):718-24.
- 9 - Soldin SJ, Brugnara C, Wong EC: Pediatric Reference Intervals, 5.ed.Washington: AACC Press, 2005.p.3-4.
- 10 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- 11 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

EC REP **OBELIS S.A.**

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE
Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for Transaminase ALT Kinetic kit:
10269360145

Review: July/2020

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

BATCH CODE

DATE OF MANUFACTURE

USED BY
(last day of month)

TEMPERATURE LIMITATION
(store at)

CONTAINS SUFFICIENT
FOR <N> TESTS

CONSULT INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC DEVICE

EUROPEAN AUTHORIZED
REPRESENTATIVE

KEEP AWAY
FROM SUNLIGHT

MANUFACTURED BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

INFLAMMABLE

CORROSIVE

POISON

CE MARK

DO NOT USE IF
PACKAGE IS
DAMAGED