TRANSAMINASE AST (TGO) CINÉTICA

REF K048

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FINALIDADE

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- 1 - The Committee on Enzymes of the Scandinavian Society for Clinical Chemistry and Clinical Physiology Scand J. Clin Lab Invest., 1974, 33, 291 - 306.
- 2 - BERGMEYER, Bowers and cols., Clin. Chim., Acta., 1.976, 70, 19 - 42, 1977, 21 - 22.
- 3 - BERGMEYER, HV.; SCHEIBE, P.; WAHLEFELD, A.W., Clin., Chem., 1.978, 24, 58 - 73.
- 4 - Expert Panel on Enzymes of the International Federation of Clinical Chemistry: Part 3. Revised IFCC Method for Aspartate Aminotransferase, Clin., Chem., 1.978, 24, 720.
- 5 - Scandinavian Committee on Enzymes Scand J. Clin Lab Invest, 1.981, 41, 107 - 116.
- 6 - BURTIS; CARL, A.;ASHWOOD; EDWARD, R., Clin. Chem.,Tietz TEXT BOOK of; 2a ed., 1.986, 788 - 797.
- 7 - PESCE, A., J.; KAPLAN, L., A., Methods in Clin. Chem., 1.987.
- 8 - SOLDIN, S.J., BRUGNARA, C., WONG, E.C.: Pediatric Reference Intervals, 5.ed.Washington: AACC Press, 2005.p.3-4.
- 9 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31

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- 1 - The Committee on Enzymes of the Scandinavian Society for Clinical Chemistry and Clinical Physiology Scand J. Clin Lab Invest., 1974, 33, 291 - 306.
- 2 - BERGMAYER, Bowers and cols., Clin. Chim., Acta., 1.976, 70, 19 - 42, 1977, 21 - 22.
- 3 - BERGMAYER, HV.; SCHEIBE, P.; WAHLEFELD, A.W., Clin., Chem., 1.978, 24, 58 - 73.
- 4 - Expert Panel on Enzymes of the International Federation of Clinical Chemistry: Part 3. Revised IFCC Method for Aspartate Aminotransferase, Clin., Chem., 1.978, 24, 720.
- 5 - Scandinavian Committee on Enzymes Scand J. Clin Lab Invest, 1.981, 41, 107 - 116.
- 6 - BURTIS; CARL, A.;ASHWOOD; EDWARD, R., Clin. Chem.,Tietz TEXT BOOK of; 2a ed., 1.986, 788 - 797.
- 7 - PESCE, A., J.; KAPLAN, L., A., Methods in Clin. Chem., 1.987.
- 8 - SOLDIN, S.J., BRUGNARA, C., WONG, E.C.: Pediatric Reference Intervals, 5.ed.Washington: AACC Press, 2005.p.3-4.
- 9 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31

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USAGE INSTRUCTIONS



FUNCTION

Method for determination of Aspartate Amino Transferase (AST). Kinetic test, for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Kinetic (UV)

Kinetic determination (UV) of AST, according to the reaction:

AST catalyzes the transfer of Amino groups of Aspartate for α - Ketoglutarate, leading to formation of Glutamate and Oxalacetate. Oxalacetate in the presence of MDH reacts with NADH, reducing itself into Malate and the NADH oxidates to NAD⁺. The rate of oxidation is proportional to AST activity in the sample.

REAGENTS

Reagent N° 1 - Substrate - Store between 2 e 8°C. Contains: Tris Buffer < 200 mmol/L, Sodium Aspartate < 450 mmol/L, D-Lactate Dehydrogenase < 5 KU, Malate Dehydrogenase < 5 KU, chelator, surfactant and preservative.

Reagent N° 2 - Coenzyme - Store between 2 e 8°C. Contém: Tris Buffer < 40 mmol/L, Alpha Ketoglutaric Acid < 180 mmol/L, NADH < 5 mmol/L, surfactant and preservative.

PRESENTATION

Presentation	Reagent N° 1	Reagent N° 2
1	54 mL	6 mL
2	27 mL	3 mL
3	2 x 27 mL	2 x 3 mL
4	4 x 27 mL	4 x 3 mL
5	1 x 40 mL	1 x 10 mL
6	2 x 40 mL	2 x 10 mL
7	4 x 40 mL	4 x 10 mL
8	2 x 40 mL	1 x 20 mL
9	4 x 40 mL	2 x 20 mL
10	4 x 27 mL	1 x 12 mL
11	6 x 36 mL	1 x 24 mL
12	6 x 36 mL	2 x 12 mL
13	3 x 36 mL	1 x 12 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Thermostated spectrophotometer, pipettes, watch or stopwatch, test tubes, Biocontrol N, Biocontrol P and Biocal Bioclin. They can be found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport, at temperatures up to 30°C, should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must to be recent and free of contaminants.

4- It is important, for the good development of the test, a rigorous control of time and temperature.

5- Reagents N° 1 and N° 2 contains Sodium Azide, and should be handled with care.

6- Lipemic and icteric samples increase the absorbance at 340 nm. In this case, you should dilute the sample 1:2 with Saline Solution. Multiply the results by 2.

7- Readings of less than 0,800 absorbance from Working Reagent indicate loss of itself. In this case, do not use the reagent.

8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

9- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Advisory Service Customer) of Quibasa.

10- Do not use the product in case of damaged packaging.

11- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or plasma with EDTA or Heparin, obtained free of hemolysis. The enzyme is stable for 3 days between 2 and 8°C and 3 months at -20°C °.

PROCESS DESCRIPTION

The calibration stability of the Transaminase AST Kinetic kit installed on refrigerated equipment is at least 30 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

PREPARATION OF WORKING REAGENT

Mix 4 parts of Reagent N°1 with 1 part of Reagent N°2. Working Reagent is stable during 72 hours between 15 and 30°C and 14 days between 2 and 8°C.

REACTION CONDITIONS

It is indispensable to use a thermostated cuvette at 37°C, 1 cm optical path and reading at 340 nm (334 - 365 nm).

TECHNIQUE

Bioclin recommends using the kit Biocal as calibrator and as control serum, Biocontrol N and P Bioclin kits.

Add 100 µL of Sample to 1,0 mL of Working Reagent, mix and transfer to a thermostated cuvette at 37°C and wait for 1 minute. Make the initial reading, simultaneously starting the stopwatch. Repeat readings after 1, 2 and 3 minutes. Calculate the average differences in absorbance per minute (ΔA/min.) and use it to calculate the result.

DESCRIPTION OF THE CALCULATIONS

AST (U/L) 340 nm = ΔA/min. x 1746

334 nm = ΔA/min. x 1780

365 nm = ΔA/min. x 3235

Results are expressed in U/L.

For an average variation in absorbance ≥ 0,15 in 340 nm and 334 nm or ≥ 0,080 in 365 nm, repeat the determination diluting the sample with NaCl 0,85%. Multiply the results obtained by the dilution factor.

PROCESS LIMITATIONS

The kinetic method is based on the molar absorptiveness, by this reason the readings must be conducted in a spectrophotometer that satisfies the following conditions:

Wavelength 340 nm

Semi trajectory of the pass band 10 nm

Stray light less than 0,5%

1 cm thermostated cuvette

INTERFERENT

No interference was observed by Bilirrubin until 19 mg/dL, Hemoglobin until 45 mg/dL, Triglycerides until 650 mg/dL and Pyruvate until 0,2 mmol/L.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit calibration can be made using the theoretical calibration factor based on the molar absorptivity of the NADH, or through the BIOCAL calibrator. Bioclin recommends the usage of the BIOCAL calibrator which is traceable to the reference material ERM-AD457.

REFERENCE VALUES

The reference values in U/L, for this method were obtained through the determination of AST in healthy populations of male and female.

Age	Male (U/L)	Female (U/L)
1 - 7 days	26 - 98	20 - 93
8 - 30 days	16 - 67	20 - 69
1 - 6 months	16 - 62	16 - 61
7 - 12 months	16 - 52	16 - 60
1 - 3 years	16 - 57	16 - 57
4 - 6 years	10 - 47	10 - 47
7 - 15 years	10 - 41	5 - 36
Adults	11 - 39	10 - 37

These values should be use as guideline, each laboratory should create their own range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and / or treatment of the patient.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

RECOVERY

The recovery analysis was performed with samples of 05 determinations. Accuracies were calculated and were found in good agreement with the reference values, obtaining a recovery between 94% and 103%.

COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY

The Transaminase AST Kinetic kit was compared with another kit commercially available to measure Aspartate Amino Transferase. 42 analyzes were conducted and the results were evaluated. Was obtained the linear equation $Y = 1,002X - 1,2973$, with correlation coefficient equal to 0,9994. With these results we can conclude the kit shows good methodological specificity.

Precision

REPEATABILITY

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration(U/L)	166,29	71,08	20,76
Standard Deviation (U/L)	1,14	0,58	0,30
Coefficient of Variation (%)	0,69	0,82	1,44

REPRODUCIBILITY

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration(U/L)	166,41	71,05	20,11
Standard Deviation (U/L)	0,12	0,02	0,82
Coefficient of Variation (%)	0,07	0,03	4,07

Sensitivity

Sensitivity was calculated from 40 determinations of a sample free of Aspartate Amino Transferase. The average found was of 2,10 U/L with a standard deviation of 0,29 U/L. The sensitivity, which indicates the method detection limit, corresponds the average plus 3 times the standard deviation, and is equal to 2,97U/L.

Linearity

This reaction is linear up to the concentration of 400 U/L. For samples with values higher than 400 U/L it is recommended to dilute the sample with Sodium Chloride 0,85%, repeat the dosage and multiply the results obtained by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

The increased activity of enzymes Aspartate Amino Transferase - AST (location citomitocondrial) reflects changes in several tissues. This enzyme is found in high concentrations in heart, liver, skeletal muscle, kidney and pancreas. Their action on plasma increases after 6 to 8 hours of myocardial infarction, reaching a peak in 24 to 48 hours after the onset. Considerable increases occur in viral hepatitis, toxic, necrotic liver diseases - 3 to 50 or 100 times the reference values (RV), mononucleosis (20 times the RV), intrahepatic cholestasis - liver (20 times the VR) and muscular dystrophies (8 times the VR). In chronic liver diseases associated with cell necrosis due to increase in AST, can occur inversion of ALT (TGP) /AST.

BIBLIOGRAPHIC REFERENCES

- 1 - The Committee on Enzymes of the Scandinavian society for Clinical Chemistry and Clinical Physiology Scand J. Clin Lab Invest., 1974, 33, 291 - 306.
- 2 - BERGMEYER, Bowers and cols., Clin. Chim., Acta., 1.976, 70, 19 - 42, 1977, 21 - 22.
- 3 - BERGMEYER, HV.; SCHEIBE, P.; WAHLEFELD, A.W., Clin., Chem., 1.978, 24, 58 - 73.
- 4 - Expert Panel on Enzymes of the International Federation of Clinical Chemistry: Part 3. Revised IFCC Method for Aspartate Aminotransferase, Clin., Chem., 1.978, 24, 720.
- 5 - Scandinavian Committee on Enzymes Scand J. Clin Lab Invest, 1.981, 41, 107 - 116.
- 6 - BURTIS; CARL, A.;ASHWOOD; EDWARD, R., Clin. Chem.,Tietz TEXT BOOK of; 2a ed., 1.986, 788 - 797.
- 7 - PESCE, A., J.; KAPLAN, L., A., Methods in Clin. Chem., 1.987.
- 8 - SOLDIN, S.J., BRUGNARA, C., WONG, E.C.: Pediatric Reference Intervals, 5.ed.Washington: AACC Press, 2005.p.3-4.
- 9 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

OBELIS S.A.

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for Transaminase AST Kinetic kit:
10269360119

Review: January/2020

UNIVERSAL SYMBOLGY

CATALOG NUMBER

MANUFACTURED BY

BATCH CODE

CONTROL

DATE OF MANUFACTURE

POSITIVE CONTROL

USED BY
(last day of month)

NEGATIVE CONTROL

TEMPERATURE LIMITATION
(store at)

BIOLOGICAL RISK

CONTAINS SUFFICIENT
FOR <N> TESTS

INFLAMMABLE

CONSULT INSTRUCTIONS
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC DEVICE

POISON

EUROPEAN AUTHORIZED
REPRESENTATIVE

CE MARK

KEEP AWAY
FROM SUNLIGHT

DO NOT USE IF
PACKAGE IS
DAMAGED