TRANSAMINASE TGO

REF K034

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FINALIDADE

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USAGE INSTRUCTIONS


FUNCTION

Method for determination of Oxalacetic Transaminase (AST). Colorimetric test. For *in vitro* diagnostic only.

PRINCIPLE OF ACTION

Methodology: Reitmam and Frankel

The Transaminase catalyzes the transfer of the Amino group from an Alpha-Amino Acid to a Alpha-Keto Acid.

The Oxaloacetate reacts with the Dinitrophenylhydrazine. The intensity of the color of the formed Hydrazone, in na alkaline environment, is directly proportional to the quantity of the Oxaloacetate, in a determinate time, is function of the enzymatic activity.

REAGENTS

Number 1 - AST Substrate - Store between 2 and 8°C. Contains: Phosphate Buffer 0,1 mol/L (pH 7,4), Aspartic L Acid 0,153 mol/L, Alpha - Ketoglutarate Acid 2 mmol/L and Sodium Azide 7,7 mmol/L.

Number 2 - Color Reagent - Store between 2 and 8°C. Contains: 2-4 Dinitrophenylhydrazine 1 mmol/L, Chloridric Acid 1,0 mol/L.

Number 3 - Concentrated Sodium Hydroxide - Store between 2 and 8°C. Contains: Sodium Hidroxide 5 mol/L.

Number 4 - Standard - Store between 2 and 8°C. Homogenized well before use. Contains: Sodium Pyruvate 2,53 mmol/L (22 mg/dL Pyruvate), Phosphate Buffer 0,1 mol/L (pH 7,4) and Sodium Azide 7,7 mmol/L.

PRESENTATION

REAGENT	VOLUME
Nº 1	50 mL
Nº 2	50 mL
Nº 3	40 mL
Nº 4	4 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, water-bath (37°C), watch and stopwatch, pipettes and test tubes, monolog or graph paper, Biocontrol N and Biocontrol P Bioclin. They can be found at the market specialized on Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature is between 2 and 8°C. The transport, at temperature between 15 and 30°C, should not exceeding 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- Water used in material cleaning must to be recent and free of contaminants.
- 4- Although the Substrate contain a preservative, it should avoid bacterial contamination.

- 5- Due to the large concentration of AST in the erythrocytes, the hemolysis prejudice the test, supplying results wrongly elevated.
- 6- Handle the Reagents Nº 1 and Nº 4 carefully, they contain Sodium Azide , irritating to skin and mucous membrane.
- 7- The Reagent Nº 3 is caustic; avoid the skin contact. The Reagent Nº 2 contain acid and shall be handle with care, because it is irritating to skin and mucous membrane, causing minor burns.
- 8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 9- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 10- Do not use the product in case of damaged packaging.
- 11- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLE

Serum or plasma (collected with heparin or EDTA) The analytical is stable for 4 days between 2 and 8°C and 2 weeks at 10°C negative.

PROCESS DESCRIPTION

USAGE SODIUM HYDROXIDE PREPARATION (0,4 N)

Transfer quantitatively the contend of flask Nº 3 (Concentrated Sodium Hydroxide - 40 mL) for a volumetric balloon of 500 mL and complete until the mark with distilled or deionized water free of CO². Homogenize well and store in plastic flask. Stable in temperatures between 15 and 30°C.

CALIBRATION CURVE

Do not use the factor method for AST dosages, because the absorbencies are not directly proportional to the enzymatic activity. Prepare with great accuracy the calibration curve described bellow.

Tube	Standard	AST Substrate	Distilled Water	AST U/mL
1	---	1,0 mL	200 µL	0
2	100 µL	0,9 mL	200 µL	24
3	200 µL	0,8 mL	200 µL	61
4	300 µL	0,7 mL	200 µL	114
5	400 µL	0,6 mL	200 µL	190

Add to each tube 1,0 mL of Reagent Nº 2 (Color Reagent), homogenize and let it rest for 20 minutes at environment temperature. Optionally, this time can be reduced to 5 minutes, incubating in water-bath at 56°C. Add 10,0 mL of Usage Sodium Hydroxide (0,4 N). Homogenize and let it rest for 5 minutes. The color remains stable for 60 minutes. Determine the absorbencies or transmittances at 505 nm (490 - 540 nm), matching the zero with distilled or deionized water. Trace the calibration curve, placing the absorbencies values in order (graph paper) or transmittance (monolog paper), and in abscissa the values in Units/mL.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin kits.

Mark the test tubes to dosage AST and proceed as follows:

	AST
Reagent Nº 1	250 µL
Incubate at 37°C per 3 minutes	
Sample	100 µL
Homogenize and incubate at 37°C per 30 minutes	
Reagent Nº 2	250 µL
Homogenize and let it rest per 20 minutes in environment temperature or per 5 minutes at 56°C	
Usage Sodium Hidroxide	2,5 mL

Homogenized by inversion and let it rest per 5 minutes in environment temperature. Transfer to the cuvettes and measure the absorbance or transmittance in 505 nm (490 - 540 nm), matching the zero with distilled or deionized water. The color is stable per 30 minutes. Calculate the AST values through the calibration curve.

PROCEDURE LIMITATIONS

Test periodically the calibration to detect any change in response of the spectrophotometer or from the colorimeter.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material ERM-AD457.

REFERENCE VALUES

The reference values in U/mL for this method were obtained through the determination of AST in healthy populations of male and female.

AST.....4 to 36 units/mL
To convert the values of units/mL to international units (UI), multiply by 0,482.
These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

PRODUCT PERFORMANCE
QUALITY CONTROL

Accuracy
RECOVERY

The Transaminase AST kit recovery analysis was made with 05 samples determinations. The accuracies were calculated and if found in good concordance with the reference values, obtaining a recovery between 95 and 103%.

COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY

The Transaminase AST kit was compared with other method commercially available for dosage of AST. 07 analyzes were performed and the results were evaluated. The linear equation obtained was $Y = 1,013X - 0,013$ and correlation coefficient 0,997. With these results, we can conclude that the kit shows good methodological specificity.

Precision

REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (U/mL)	24,05	31,65	56,90
Standard Deviation (U/mL)	0,76	0,81	0,79
Coefficient of Variation (%)	3,16	2,57	1,39

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (U/mL)	23,95	31,75	56,67
Standard Deviation (U/mL)	0,10	0,10	0,36
Coefficient of Variation (%)	0,42	0,31	0,64

Sensitivity

The sensitivity was calculated from 20 determinations of a sample free of AST. The average found was 3,377 U/mL, with standard deviation of 0,006 U/mL. The sensitivity, wich indicates the method detection limit, corresponds the average plus 3 times the standard deviation, and is equal to 3,394 U/mL.

Linearity

This reaction is linear up to the concentration of 180 U/L. For samples with values higher than 180 U/L it is recommended to dilute the sample with Sodium Chloride 0,85%, repeat the dosage and multiply the results obtained by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

The increase on the Aspartate Amine Transferase - AST enzymes activity, of cytomitochondrial location, reflects the changes of several tissues.

This enzyme is found in high concentration in the heart, liver, skeletal muscle, kidney and pancreas, and its activity in the plasma increases 6 to 8 hours, after the myocardium stroke, reaching the apex in 24 hours to 48 hours, after the assail.

Considerable increases occur in viral hepatitis, toxics, necrotic hepatic diseases (3 to 50 times the values references - RV), mononucleosis (20 times the RV), intra-hepatic cholestasis (20 times the RV) and muscle dystrophies (8 times the RV). In the chronic HD associated to cell necrosis, due the increase of the AST liberation (mitochondrial), can occur the inverted relation ALT/AST. Still, the increase of its levels in infecting mononucleosis, and in intra cholestasis and extra-hepatics.

NUMBER OF TESTS

200 tests/100 µL of sample/3 mL of reagent

BIBLIOGRAPHIC REFERENCES

1. CABAL, P., LEEPER, R., WROBLEWISKI, F.: Am. J. Clin. Path. 26:1101, 1956;
2. LA DUE, J. S., WROBLEWISKI, F., KARMEN, A.: Science 120:497,01954;
3. MOHUN, A. F., COOK, I. J. Y.: Am. J. Clin. Path. 27:235, 1957;
4. WROBLEWISKI, F., CABAUD, P., Am. J. Clin. Path. 27:235, 1957;
5. TONKS, D. B.: Clin. Chem. 9: 217, 1983.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

OBELIS S.A.

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for the Transaminase AST kit: 10269360139

Review: September/2016

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MANUFACTURED BY

BATCH CODE

CONTROL

DATE OF MANUFACTURE

POSITIVE CONTROL

USED BY
(last day of month)

NEGATIVE CONTROL

TEMPERATURE LIMITATION
(store at)

BIOLOGICAL RISK

CONTAINS SUFFICIENT
FOR <N> TESTS

INFLAMMABLE

CONSULT INSTRUCTIONS
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC DEVICE

POISON

EUROPEAN AUTHORIZED
REPRESENTATIVE

CE MARK

KEEP AWAY
FROM SUNLIGHT

DO NOT USE IF
PACKAGE IS
DAMAGED