TRANSAMINASE TGP

REF K035

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FINALIDADE

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USAGE INSTRUCTIONS


FUNCTION

Method for determination of Pyruvic Transaminase (ALT). Colorimetric test. For *in vitro* diagnostic only.

PRINCIPLE OF ACTION

**Methodology:** Reitmamm and Frankel  
The Transaminase catalyzes the transfer of the Amino group from an Alpha-Amino Acid to a Alpha-Keto Acid.

The pyruvate reacts with the dinitrophenylhydrazine. The intensity of the color of the formed hydrazone, in an alkaline environment, is directly proportional to the quantity of the Pyruvate, in a determinate time, in function of the enzymatic activity.

REAGENTS

**Number 1 - ALT Substrate** - Store between 2 and 8°C. Contains: Phosphate Buffer 0,1 mol/L (pH 7,4), L-Alanine 0,2 mol/L, Alpha-Ketoglutarate Acid 2 mmol/L, and Sodium Azide 7,7 mmol/L.

**Number 2 - Color Reagent** - Store between 2 and 8°C. Contains: 2,4-Dinitrophenylhydrazine 1,0 mmol/L and Chloridric Acid 1,0 mol/L.

**Number 3 - Concentrated Sodium Hydroxide** - Store between 2 and 8°C. Contains: Sodium Hydroxide 5 mol/L.

**Number 4 - Standard** - Store between 2 and 8°C. Homogenized well before use. Contains: Sodium Pyruvate 2,53 mmol/L (22 mg/dL Pyruvate), Phosphate Buffer 0,1 mol/L (pH 7,4) and Sodium Azide 7,7 mmol/L

PRESENTATION

| Reagent | Volume |
|---------|--------|
| N° 1    | 50 mL  |
| N° 2    | 50 mL  |
| N° 3    | 40 mL  |
| N° 4    | 4 mL   |

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, water-bath at 37°C, watch and stopwatch, pipettes, test tubes, monolog or millimeter paper, Biocontrol N and Biocontrol P Bioclin. The instruments and glassware can be found in the specialized market of clinical analysis laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 and 8°C. The transport at temperature between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- Water used in material cleaning must be recent and free of contaminants.

4- Although the Substrate contain preservative, should avoid bacterial contamination.

5- Handle the Reagents N° 1 and N° 4 carefully, they contain Sodium Azide, irritating to skin and mucous membrane.

6- The Reagent N° 3 is caustic; avoid the skin contact. The Reagent N° 2 contain acid and shall be handle with care, because it is irritating to skin and mucous membrane, causing minor burns.

7- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

8- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Advisory Service Customer) of Quibasa.

9- Do not use the product in case of damaged packaging.

10- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLE

Serum or plasma with EDTA or heparin, obtained free of hemolysis. The sample is stable for 3 days between 2 and 8°C, and 3 months at -20°C °.

PROCESS DESCRIPTION

USAGE SODIUM HYDROXIDE PREPARATION (0,4 N)

Transfer quantitatively the contend of flask N° 3 (40 mL) for a volumetric balloon of 500 mL and complete until the mark with distilled or deionized water free of CO<sub>2</sub>. Homogenize well and store in plastic flask. Stable in temperatures between 15 and 30°C, until the validity of the product printed on the label.

CALIBRATION CURVE

Do not use the factor method for ALT dosages, because the absorbencies are not directly proportional to the enzymatic activity. Prepare with great accuracy the calibration curve described below.

| Tube | Standard | ALT Substrate | Distilled Water | ALT U/mL |
|------|----------|---------------|-----------------|----------|
| 1    | ---      | 1,0 mL        | 200 µL          | 0        |
| 2    | 100 µL   | 0,9 mL        | 200 µL          | 28       |
| 3    | 200 µL   | 0,8 mL        | 200 µL          | 57       |
| 4    | 300 µL   | 0,7 mL        | 200 µL          | 97       |
| 5    | 400 µL   | 0,6 mL        | 200 µL          | 150      |

Add to each tube 1,0 mL of Color Reagent N° 2, homogenize and let it rest for 20 minutes at environment temperature. Optionally, this time can be reduced to 5 minutes, incubating in water-bath at 56°C. Add 10,0 mL of Usage Sodium Hydroxide (0,4 N). Homogenize and let it rest for 5 minutes. The color remains stable for 60 minutes. Determine the absorbencies or transmittances at 505 nm (490 - 540 nm), matching the zero with distilled or deionized water. Trace the calibration curve, placing in orderly absorbance values (millimeter paper) or transmittance (monolog paper), and in abscissa, the values in U/mL.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin kits.

Mark the test tubes for the ALT dosage and proceed as follows:

|                                                                                                |        |
|------------------------------------------------------------------------------------------------|--------|
|                                                                                                | ALT    |
| Reagent N° 1                                                                                   | 250 µL |
| Incubate at 37°C per 3 minutes                                                                 |        |
| Sample                                                                                         | 50 µL  |
| Homogeneize and incubate at 37°C per 30 minutes                                                |        |
| Reagent N° 2                                                                                   | 250 µL |
| Homogeneize and let it rest per 20 minutes in environment temperature or per 5 minutes at 56°C |        |
| Usage Sodium Hydroxide                                                                         | 2,5 mL |

Homogenized by inversion and let it rest per 5 minutes in environment temperature. Transfer to the cuvettes and measure the absorbance or transmittance in 505 nm (490 - 540 nm), matching the zero with distilled or deionized water. The color is stable per 30 minutes. Calculate the ALT values through the calibration curve.

PROCEDURE LIMITATIONS

Test periodically the calibration to detect any change in response of the spectrophotometer or from the colorimeter.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material ERM-AD454.

REFERENCE VALUES

The reference values in U/mL for this method were obtained through the determination of ALT in healthy populations of male and female.

ALT.....4 to 32 U/mL

To convert the values of U/mL to international units (UI), multiply by 0,482.

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE  
QUALITY CONTROL

Accuracy  
RECOVERY

The Transaminase ALT kit recovery analysis was made with 5 samples determinations. The accuracies were calculated and if found in good concordance with the reference values, obtaining a recovery between 95 and 105%.

COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY

The Transaminase ALT kit was compared with other commercially available methods for measurement of LT. 07 analyzes were performed and the results were evaluated. The linear equation obtained was  $Y = 0,993X + 0,707$  and the correlation coefficient 0,999. With these results, we can conclude that the kit presents good methodological specificity.

Precision  
REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

|                              | Sample 1 | Sample 2 | Sample 3 |
|------------------------------|----------|----------|----------|
| Average Concentration (U/mL) | 14,45    | 23,10    | 35,40    |
| Standard Deviation (U/mL)    | 0,51     | 0,79     | 0,94     |
| Coefficient of Variation (%) | 3,53     | 3,41     | 2,66     |

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

|                              | Sample 1 | Sample 2 | Sample 3 |
|------------------------------|----------|----------|----------|
| Average Concentration (U/mL) | 14,50    | 23,35    | 35,17    |
| Standard Deviation (U/mL)    | 0,05     | 0,23     | 0,21     |
| Coefficient of Variation (%) | 0,34     | 0,98     | 0,59     |

Sensitivity

The sensitivity was calculated from 20 determinations of a sample free of ALT. The average found was 3,389 U/mL, with standard deviation of 0,012 U/mL. The sensitivity, which indicates the method detection limit, corresponds the average plus 3 times the standard deviation, and is equal to 3,424 U/mL.

Linearity

This reaction is linear up to the concentration of 140 U/L. For samples with values higher than 140 U/L it is recommended to dilute the sample with Sodium Chloride 0,85%, repeat the dosage and multiply the results obtained by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

The increase activity of the Alanine Amino Transferase - ALT (of cytoplasmic origin) reflect the changes of several tissues. The higher activity of this enzyme is located in the hepatic tissue. Lower activities occur in the skeletal muscle, heart, kidney and pancreas. Its activity is increased in viral and toxic hepatitis

(30 - 50 or 100 times the RV) as well as other hepatic diseases (HD), associated to hepatic necrosis. In the chronic HD, associated to cell necrosis, due the increase of the AST liberation (mitochondrial), can occur the inverted relation ALT/AST. Still, the increase of its levels in infecting mononucleosis, in intra cholestasis and extra-hepatics.

NUMBER OF TESTS

200 Tests / 50 µL of Sample / 3 mL of Reagent.

BIBLIOGRAPHIC REFERENCES

- 1 - CABAL, P., LEEPER, R., WROBLEWISKI, F.: Am. J. Clin. Path. 26:1101, 1956;
- 2 - LA DUE, J. S., WROBLEWISKI, F., KARMEN, A.: Science 120:497,01954;
- 3 - MOHUN, A. F., COOK, I. J. Y.: Am. J. Clin. Path. 27:235, 1957;
- 4 - WROBLEWISKI, F., CABAUD, P., Am. J. Clin. Path. 27:235, 1957;
- 5 - TONKS, D. B.: Clin. Chem. 9: 217, 1983.
- 6 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455  
E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

OBELIS S.A.

Bd. Général Wahis, 53  
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service  
Phone: 0800 0315454  
E-mail: sac@bioclin.com.br

ANVISA registration for the Transaminase ALT kit: 10269360140

Review: March/2019

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MANUFACTURED BY

BATCH CODE

CONTROL

DATE OF MANUFACTURE

POSITIVE CONTROL

USED BY  
(last day of month)

NEGATIVE CONTROL

TEMPERATURE LIMITATION  
(store at)

BIOLOGICAL RISK

CONTAINS SUFFICIENT  
FOR <N> TESTS

INFLAMMABLE

CONSULT INSTRUCTIONS  
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC DEVICE

POISON

EUROPEAN AUTHORIZED  
REPRESENTATIVE

CE MARK

KEEP AWAY  
FROM SUNLIGHT

DO NOT USE IF  
PACKAGE IS  
DAMAGED