# Bioclin

## TTPA BIOCLIN

REF **K090**

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#### FINALIDADE

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FUNCTION

Method for measuring manual or automatic time of Activated Partial Thromboplastin Time in citrated plasma. For in vitro diagnostics use only.

PRINCIPLE OF ACTION

**Methodology:** Coagulometry

Consist on the determination of time of citrated plasma after adding reagents that contain a Plasmatic Activator (Ellagic Acid) and Phospholipids, that will act as substitutes to platelets. Initially plasma is incubated with the Contact Activator, at 37°C for 3 minutes. Following, add the Clot Forming (Calcium Chlorite) for recalcification, triggering the intrinsic coagulation mechanism. This step should be timed until the formation of the clot.

REAGENTS

**Number 1 - Contact Activator (Activated Cephalin)** - Store between 2 and 8°C. Contains: Rabbit brain extract, Ellagic Acid and preservative.

**Number 2 - Clot Former (Calcium Chlorite Solution)** - Store between 2 and 8°C. Contains: Calcium Chloride and preservative.

**Presentation ready for use.**

PRESENTATION

| Presentation | Reagent 1   | Reagent 2 |
|--------------|-------------|-----------|
| 1            | 6 x 2.5 mL  | 4 x 4 mL  |
| 2            | 6 x 5 mL    | 8 x 4 mL  |
| 3            | 12 x 2.5 mL | 8 x 4 mL  |
| 4            | 12 x 5 mL   | 16 x 4 mL |

EQUIPMENTS AND OPERATIONAL INPUTS

Coagulometer, water bath at 37°C, watch or stopwatch, pipettes, test tubes, Biocontrol N Coagulation and Biocontrol P Coagulation. These itens are found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

Storage temperature should be between 2 and 8°C. Transportation under temperatures between 15 and 30°C should exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

**1- For professional *in vitro* diagnostic use only.**

**2-** Strictly follow the methodology proposed to obtain exact results.

**3-** The level of water in water bath must be superior to the level of reagents in test tubes.

**4-**Reagent N° 1 contains Sodium Azide, irritant to skin and mucous, handle with care.

**5-** Avoid the warming of reagent that will not be used. Already heated reagent should not be returned to the original bottle.

**6-** We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

**7-** To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.

**8-** Do not use the product in case of damaged packaging.

**9-** It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

**Material:** Plasma

**Anticoagulant:** Use as coagulant the Bioclin Citrate kit, according to instructions of usage, or Sodium Citrate 3.8% (0.130 M) or 3.2% (0.109 M), proportionally in 9 parts of Blood to 1 of Anticoagulant (for example: 4.5 mL of Blood + 0.5 mL of Anticoagulant).

**Collection:** Blood must be obtained through vein puncture, avoiding hemolysis, prolonged tourniquet, bubble formation and tissue liquid aspiration. When evaluating only the coagulation test, always collect 2 tubes of citrate and discard the first one

**Preparation:** Blood should be mixed with anticoagulant right after collection. Centrifuge immediately at 3000 rpm, during 15 minutes.

Remove plasma without pipetting red cells and the yellow layer.

Samples must be tested in less than 3 hours. If tests can not be made during this period plasma should be frozen for at maximum 1 week at -20°C.

TECHNIQUE

Bioclin recommends, for use of the kit, to use the Biocontrol N Coagulation and Biocontrol P Coagulation kits as the control serum.

**1-** Preheat the Reagent N° 2 to 37°C. This procedure should be performed only with the volume that will be used to perform the tests. The reagent must remain warm for the entire procedure.

**2-** Separate the reaction tubes that will be used to perform the tests.

**3-** Pipette 100 µL of the Sample/Control into the respective reaction tube and preheat to 37°C for 2 minutes in a water bath.

**4-** Pipette 100 µL of Reagent N° 1 into the Sample/Control tube and incubate at 37°C, for 2 to 3 minutes.

**5-** Pipette 100 µL of Reagent N° 2 (preheated to 37°C) in the Sample/Control tube, simultaneously triggering the timer.

**6-** Leave the tube in the water bath at 37°C, shaking gently. Keep in the bath for about 20 seconds.

**7-** Remove the tube from the bath, gently tilt it once a second and stop the timer at the time of the clot formation. Note the clotting time.

CALCULATIONS

Calculate the average between coagulation time of each samples/controls.

Release all results in seconds.

CURVE OF SENSIBILITY TO HEPARIN

This procedure is useful to monitor therapies with Heparin.

A curve of sensibility to Heparin must be built adding the quantity of Heparin known to a mixture of normal plasmas (pool).

As a suggestion the following values of Heparin can be used:

**0.0 / 0.1 / 0.2 / 0.4 / 0.6**

Determine the dilution time dosages, plotting in a graphic the times and concentrations.

Each laboratory must establish its own calibration curve.

PROCEDURE LIMITATIONS

The concentration of citrate may interfere with coagulation time because it directly affects the concentration of calcium. Thus, APTT values are higher in tubes of 3.8% than in tubes of 3.2%.<sup>4</sup>

Citrate at 3.8% concentration will release results higher than 3.2%. Therefore, it is necessary that the laboratory standardize so that the monitoring of the patient's diagnosis is always performed with the same anticoagulant concentration.

The concentration of sodium citrate should be adjusted in patients with hematocrit values above 55%. In samples with high hematocrit, the blood/anticoagulant ratio (9: 1) is not maintained and may cause excess citrate to the volume of blood present in the tube. To calculate the citrate volume, use the formula described in item 4 of the Laboratory Diagnostic Manual for Hereditary Coagulopathies and Plaquetopathies of the Ministry of Health.<sup>5</sup>

INTERFERENCES

Icteric, lipemic or hemolyzed samples.

Usage of plasma containing Heparin, EDTA or Sodium Oxalate.

Conditions and substances that may alter TTPA: Oral contraceptives, Estrogen, pregnancy, drugs Coumarin, Heparin, Asparaginase, Acetylsalicylic Acid, Metronidazole, Digitalis and Tetracycline.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

REFERENCE VALUES

Interval should be used only as guidance. It is recommended that each laboratory establish, in population served, its own reference range, using a pool of normal plasmas or reference plasmas.

| Age                                     | Interval (seconds) |
|-----------------------------------------|--------------------|
| 2 months                                | 26 – 47            |
| 5 months                                | 26 – 46            |
| Children older than 6 months and Adults | 27 - 38            |

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

COMPARISON OF METHODS, METHODOLOGICAL SPECIFICITY AND SENSITIVITY

APTT Bioclin kit was compared with other method commercially available for dosage of Activated Partial Thromboplastin Time. 07 analyzes were performed and the results were evaluated. The linear equation obtained was Y = 1.0094X - 1.5788 and the correlation coefficient 0.9974. With these results, we can conclude that the kit shows good methodological specificity.

**Precision**

**REPEATABILITY**

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

|                              | Sample 1 | Sample 2 | Sample 3 |
|------------------------------|----------|----------|----------|
| Time (s)                     | 32.15    | 64.90    | 44.55    |
| Standard Deviation (s)       | 1.14     | 0.97     | 0.94     |
| Coefficient of Variation (%) | 3.54     | 1.49     | 2.12     |

**REPRODUCIBILITY**

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

|                              | Sample 1 | Sample 2 | Sample 3 |
|------------------------------|----------|----------|----------|
| Time (s)                     | 32.67    | 64.32    | 44.65    |
| Standard Deviation (s)       | 0.48     | 0.93     | 0.22     |
| Coefficient of Variation (%) | 1.47     | 1.44     | 0.49     |

**Sensitivity**

The kit has been evaluated for its sensitivity to the factors of deficiency in coagulation and obtained the following results:

| Factor          | VIII |      | IX   |      | XI    |      | XII   |      | Pre-kallikr. |
|-----------------|------|------|------|------|-------|------|-------|------|--------------|
| % Activ.        | < 1  | 20   | < 1  | 20   | < 1   | 20   | < 1   | 20   | < 1          |
| Result. in sec. | 82.3 | 44.7 | 83.7 | 41.0 | 134.3 | 47.6 | < 200 | 36.1 | 69.4         |

**Specificity**

The specificity of the kit is guaranteed by the addition of one single activation source of coagulant factors, Ellagic Acid.

**DIAGNOSTIC SIGNIFICANCE**

The Activated Partial Thromboplastin Time (APTT) is a test that monitors both the intrinsic pathway and a common pathway of coagulation and is useful in monitoring of Heparin therapy and to verify the deficiency all coagulation factors except factor VII. Are widely used in deficiencies identified in stage 1 of the coagulation factors VIII, IX, XI, XII, Fletcher factor.

| Presentation | Reagent 1   | Reagent 2 | N° of tests |
|--------------|-------------|-----------|-------------|
| 1            | 6 x 2.5 mL  | 4 x 4 mL  | 150         |
| 2            | 6 x 5 mL    | 8 x 4 mL  | 300         |
| 3            | 12 x 2.5 mL | 8 x 4 mL  | 300         |
| 4            | 12 x 5 mL   | 16 x 4 mL | 600         |

**BIBLIOGRAPHIC REFERENCES**

1- International Committee For Standardization in Hematology Trombosis and Haemostasis. Amer. J. Clin. Path.

2- National Committee For Clinical Laboratory Standards: Coleccion, Transport And Processing Of Blood Specimens For Coagulation Test And Performance Of Coagulation Assays. Approved Guideline Fourth Edition. NCCLS document H12-A4, 2003

3- Banez,el, D.a., Koepek, Laboratory Monitoring Of Heparin Therapy. Effect Of Different Salts of Heparin On The Activated Partial Thromboplastin Time. Amer. J. Clin. Path.

4- ADCOCK, D. M.; KRESSIN, D. C. Effect of 3.2% vs. 3.8% sodium citrate concentration on routine coagulation testing. **AM. J. Clin. Pathol.** 107(1): 105-10, 1997.

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**QUALITY ASSURANCE**

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda**  
Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone.: +55 (31) 3439.5454  
E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

|    |     |
|----|-----|
| EC | REP |
|----|-----|

**OBELIS S.A.**

Bd. Général Wahis, 53  
1030 Brussels, Belgium

**CUSTOMER SERVICE**  
Customer Advisory Service  
Phone.: 0800 0315454  
E-mail: sac@bioclin.com.br

ANVISA registration number for APTT Bioclin kit: 10269360152

**Review:** June/2020

UNIVERSAL SYMBOLOGY

 CATALOG NUMBER

 BATCH CODE

 DATE OF MANUFACTURE

 USED BY  
(last day of month)

 TEMPERATURE LIMITATION  
(store at)

 CONTAINS SUFFICIENT  
FOR <N> TESTS

 CONSULT INSTRUCTIONS  
FOR USE

 IN VITRO DIAGNOSTIC DEVICE

|    |     |
|----|-----|
| EC | REP |
|----|-----|

 EUROPEAN AUTHORIZED  
REPRESENTATIVE

 KEEP AWAY  
FROM SUNLIGHT

 MANUFACTURED BY

|         |
|---------|
| CONTROL |
|---------|

 CONTROL

|           |
|-----------|
| CONTROL + |
|-----------|

 POSITIVE CONTROL

|           |
|-----------|
| CONTROL - |
|-----------|

 NEGATIVE CONTROL

 BIOLOGICAL RISK

 INFLAMMABLE

 CORROSIVE

 POISON

 CE MARK

 DO NOT USE IF  
PACKAGE IS  
DAMAGED