# Bioclin

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REF K211

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#### FINALIDADE

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### USAGE INSTRUCTIONS


#### FUNCTION

Method for quantitative determination of Latent Binding Capacity of the Iron in serum or heparinized plasma samples. Colorimetric test, for *in vitro* diagnostic use only.

#### PRINCIPLE OF ACTION

**Methodology:** Colorimetric test/Ferrozine

A known concentration of Iron ions is incubated with serum and will bind specifically with the Transferrin in insaturated Iron-binding sites. The exceeding Iron ions are measured through the reaction with Ferrozine. The difference between the exceeding Iron and the total quantity of iron initially added to the serum is equivalent to the quantity bounded to the Transferrin, which will be the Latent Binding Capacity of the Iron – CLLF.

#### REAGENTS

**Number 1 - Buffer** - Store between 2 and 8°C. Contains: Tris Buffer < 500 mmol/L, Iron III Chloride < 50 mmol/L, Sodium Bicarbonate < 500 mmol/L, surfactants and stabilizers.

**Number 2 - Color Reagent** - Store between 2 and 8°C. Contains: Ferrozine < 25 mmol/L, reducing agents and stabilizer.

**Number 3 - Calibrator** - Store between 2 and 8°C. Contains: Serum lyophilized and preservative. **Potentially Infectant.**

**Attention:** Concentration varies by lot.

#### PRESENTATION

Reagent 1	Reagent 2	Reagent 3
1 x 20 mL	1 x 5 mL	5 mL
2 x 20 mL	2 x 5 mL	5 mL
3 x 20 mL	3 x 5 mL	5 mL
4 x 20 mL	4 x 5 mL	5 mL
5 x 20 mL	5 x 5 mL	5 mL
1 x 40 mL	1 x 10 mL	5 mL
2 x 40 mL	2 x 10 mL	5 mL
4 x 40 mL	4 x 10 mL	5 mL
2 x 40 mL	1 x 20 mL	5 mL
4 x 40 mL	2 x 20 mL	5 mL
1 x 20 mL	1 x 5 mL	3 mL
2 x 20 mL	2 x 5 mL	3 mL
3 x 20 mL	3 x 5 mL	3 mL
4 x 20 mL	4 x 5 mL	3 mL
5 x 20 mL	5 x 5 mL	3 mL
1 x 40 mL	1 x 10 mL	3 mL
2 x 40 mL	2 x 10 mL	3 mL

4 x 40 mL	4 x 10 mL	3 mL
2 x 40 mL	1 x 20 mL	3 mL
4 x 40 mL	2 x 20 mL	3 mL
1 x 20 mL	1 x 5 mL	-
2 x 20 mL	2 x 5 mL	-
3 x 20 mL	3 x 5 mL	-
4 x 20 mL	4 x 5 mL	-
5 x 20 mL	5 x 5 mL	-
1 x 40 mL	1 x 10 mL	-
2 x 40 mL	2 x 10 mL	-
4 x 40 mL	4 x 10 mL	-
2 x 40 mL	1 x 20 mL	-
4 x 40 mL	2 x 20 mL	-

#### EQUIPMENTS AND OPERATIONAL INPUTS

Biochemistry auto-chemistry analyzers, Biocontrol N, Biocontrol P Bioclin and Biocal Bioclin. These items are found at markets specialized on Laboratories of Clinical Analysis.

#### TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperature should be between 2 and 8°C. The transport temperature up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

#### SPECIAL CARE

- For professional *in vitro* diagnostic use only.
- Reagent N° 2 must be kept protected from light.
- Water used in material cleaning must be recent and free of contaminants.
- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- Calibration should be repeated periodically to verify any changes in the response of the colorimeter or spectrophotometer.
- Hemolysis, even discrete, interferes in the dosage.
- To avoid possible analytical errors, it is advisable that the distilled water used for reconstitution of the calibrator is of adequate quality. For this, it is recommended that water used in clinical laboratories follow the specifications of the Clinical Laboratory Standards Institute (CLSI).
- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.
- Do not use the product in case of damaged packaging.
- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

#### SAMPLE

Serum or heparinized plasma free of hemolysis.

Serum samples and heparinized Plasma are stable for 7 days at 2 to 8°C, or 1 month at -20°C.

The sample for therapeutic control should always be collected at the same time.

The sample should not be frozen more than once.

#### PROCESS DESCRIPTION

#### REAGENT PREPARATION

The reagents n° 1 and 2 are ready to use.

The calibration stability of the UIBC Automated kit installed on refrigerated equipment is at least 5 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control sera.

#### RECONSTITUTION OF CALIBRATOR

**Before reconstitute the Calibrator should be noted the volume indicated on the bottle and proceed as described below:**

**Presentation with Calibrator 3 mL:** Open carefully the bottle and add exactly **3 mL of distilled water**, as volume described in the label of the bottle.

**Presentation with Calibrator 5 mL:** Open carefully the bottle and add exactly **5 mL of distilled water**, as volume described in the label of the bottle.

Mix. Leave the bottle resting for 30 minutes. Every 10 minutes, agitate the bottle carefully using circular motion to prevent material sticking to the wall of the bottle. Avoid formation of foam.

After reconstitution, observe the stability according to the table below.

-20°C	+4°C	+25°C
30 days	2 days	8 hours

#### Warning:

- The reconstituted Calibrator should not be thawed more than once. Therefore, for minor routines, fractionate the material into aliquots before freezing.

#### TECHNIQUE

For the use of the kit, use as calibrator the Biocal Bioclin kit and as control serum the Biocontrol P and Biocontrol N Bioclin kits.

This kit is for use in automated biochemical analyzers only. **Check the application sheet for this equipment on the website [www.bioclin.com.br](http://www.bioclin.com.br) or through SAC 0800 031 5454.**

#### CALCULATIONS

The calculations are also based on the concentrations obtained with the Fe AUTOMATED kit (K209).

The results obtained with the UIBC AUTOMATED kit correspond to the Latent Iron-Binding Capacity (LIBC).

To determine the Total Iron-Binding Capacity (TIBC) and the Transferrin Saturation Index (TSI) follow the calculations as described below:

$$\text{TIBC (mcg/dL)} = \text{LIBC} + \text{Serum Iron}$$

$$\text{TSI (\%)} = \frac{\text{Serum Iron}}{\text{TIBC (mcg/dL)}} \times 100$$

The results will be expressed in mcg/dL.

**INTERFERENCES**

There was no sign of interference of Ascorbic Acid up to 30 mg/dL, conjugated and free Bilirubin up to 20 mg/dL, lipemia up to 2000 mg/dL of Triglycerides, RF up to 350 IU/mL, Copper up to 15 mg/dL, Zinc up to 7.5 mg/dL and Hemoglobin < 200 mg/dL.

In rare cases, samples from patients with gammopathies may show false results.

**INTERNAL QUALITY CONTROL**

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present an analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

**TRACEABILITY**

The Iron-Binding Capacity values of the kit's internal calibrator and the Biocal were defined by iron and transferrin dosages. Therefore, they are traceable to the primary iron standards (NIST SRM 937) and transferrin (ERM DA470k / IFCC).

**REFERENCE VALUES**

The reference values for this method were obtained through the determination of the Iron-Binding Capacity in healthy populations of male and female.

Latent Capacity (CLLF) ..... 140 - 280 µg/dL  
Total Binding Capacity (CTLF) ..... 250 - 410 µg/dL  
Transferrin Saturation (IST) ..... 20 - 50%

To convert the values of µg/dL to µmol/L (SI), multiply the obtained results by 0.179.

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

**PRODUCT PERFORMANCE****QUALITY CONTROL****Accuracy****COMPARISON OF METHODS**

UIBC Automated kit was compared with other methods for determination of Iron-Binding Capacity dosage commercially available measurement. 42 analyzes were performed and the results were evaluated. The equation linear obtained was  $Y = 1.0107X - 7.0284$  and the correlation coefficient 0.984. With these results, we can conclude that the kit shows good methodological specificity.

**Precision****REPEATABILITY**

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
<b>Average Concentration (µg/dL)</b>	177.15	180.88	188.23
<b>Standard Deviation (µg/dL)</b>	8.66	7.10	8.65
<b>Coefficient of Variation (%)</b>	4.89	3.93	4.60

**REPRODUCIBILITY**

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
<b>Average Concentration (µg/dL)</b>	168.52	173.52	176.20
<b>Standard Deviation (µg/dL)</b>	7.68	6.39	10.44
<b>Coefficient of Variation (%)</b>	4.56	3.68	5.93

**Sensitivity**

The sensitivity was calculated from 40 determinations of a sample free of Iron-Binding Capacity. The average found was -1.350 µg/dL, with deviation standard of 2.860 µg/dL. The sensitivity, which indicates the method limit detection, corresponds the average plus 3 times the standard deviation, and is equal to 7.231 µg/dL.

**Linearity**

The reaction is linear until the concentration of 600 µg/dL or 108 µmol/L.

For higher values, dilute the sample with NaCl solution at 0.85% and repeat the dosage. Multiply the results by the dilution factor applied.

**DIAGNOSTIC SIGNIFICANCE**

Several diseases involving Iron in the organism may be evaluated from the determination of the Iron Binding Capacity combined with the determination of Iron Serum. The sum of those two values represents the maximum iron concentration that the proteins are capable of binding. In a iron deficiency, the iron binding capacity is increased, and in severe or malignant inflammatory diseases this binding capacity is reduced.

**BIBLIOGRAPHIC REFERENCES**

- FAIRBANKS, V.F.; KLEE, G.G. Biochemical aspects of hematology. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 1642-1710.
- WICH, M.; PINGERRA, W.; LEHMANN, P. Clinical aspects and laboratory. Iron metabolism, anemias. 5th ed. Wien, New York: Springer; 2003.
- GUDER, W.G.; ZAWTA, B. et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001. p. 46.
- DATI, F.; SHUMANN, G.; THOMAS, L.; AGUZZI, F.; BAUDNER, S.; BIENVENU, J. et al. Consensus of a group of professional societies and diagnostic companies on guidelines for interim reference ranges for 14 proteins in serum based on the standardization against the IFCC/BCR/CAP reference material (CRM 470). Eur J Clin Chem Clin Biochem 1996; 34:517-20.
- THOMAS L ed. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft, 1998:273-5.
- BAKKER AJ, MÜCKE M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007; 45(9):1240-1243.
- Bioclin – Dados de arquivo
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

**QUALITY ASSURANCE**

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

**QUIBASA QUÍMICA BÁSICA Ltda**  
Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone.: +55 (31) 3439.5454  
E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

**CUSTOMER SERVICE**

Customer Advisory Service  
Phone.: 0800 0315454  
E-mail: sac@bioclin.com.br

ANVISA registration for the UIBC Automation kit: 10269360321

Review: September/2020

**UNIVERSAL SYMBOLOGY**

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED