Bioclin

URÉIA ENZIMÁTICA

REF

K047

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FINALIDADE

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USAGE INSTRUCTIONS


FUNCTION

Method for determination of Urea in samples of serum, plasma (EDTA, Heparin or Citrate) and Urine. Enzymatic colorimetric test, only for *in vitro* diagnostic.

PRINCIPLE OF ACTION

Methodology: Enzymatic Colorimetric

The Urea is hydrolyzed by Ammonium ions and CO₂ by Urease.

In alkaline pH and in presence of Salicylate and Sodium Hypochloride, the Ammonia originates a greenish compound, which color intensity is proportional to the Urea concentration in the analyzed samples.

REAGENTS

Number 1 - Standard - Store between 2 and 8°C. Contains: Urea 70.0 mg/dL and acidifying.

Number 2 - Enzyme - Store between 2 and 8°C. Contains: Phosphate Buffer < 100 mmol/L, Mercaptosuccinic Acid < 0.18 mmol/L, Albumin < 60 mmol/L, EDTA Dissodate Salt < 5 mmol/L, Urease < 500 KU, stabilizers and preservative.

Number 3 - Buffer Stock - Store between 2 and 8°C. Contains: Phosphate Buffer < 100 mmol/L, Sodium Nitroprusside < 38 mmol/L, stabilizers and preservative.

Number 4 - Oxidizing Stock - Store between 2 and 8°C. Contains: Sodium Hydroxide < 3.0 mol/L and Sodium Hypochloride < 2%.

PRESENTATION

REAGENT	VOLUME
Nº 1	3 mL
Nº 2	25 mL
Nº 3	100 mL
Nº 4	50 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, water-bath at 37°C, watch and stopwatch, pipettes, test tubes, volumetric flask and Biocontrol N e Biocontrol P Bioclin. These itens are found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperatures shall be between 2 and 8°C. The transport temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must to be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- Heavy Metal and fluoride (above 2 mg/dL) are Urease inhibitors.

6- Do not smoke next to the dosage area, because the smoke contain Ammonia steam who contaminate the sample, leading to fake increased values.

7- The Reagent Nº 4 is caustic, handle carefully. The Reagent Nº 2 contain Sodium Azide, irritating to skin and mucous membrane.

8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

9- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

10- Do not use the product in case of damaged packaging.

11- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLE

Serum, plasma (EDTA, Heparin or Citrate) and Urine.

Do not use anticoagulants containing ammonia.

The Serum or Plasma sample is stable for 7 days at 2 to 8°C⁵, and 3 months at -20°C.

Urine sample is stable for 7 days at 2 to 8°C and 4 weeks at -20°C.⁵

To dose Urea in urine, collect the 24 hours sample in a flask containing 2.0 mL of HCl at 50% (V/V). Centrifuge the sample before technique begins and proceed the analysis in a few hours, because the urine excrement urea is easily decomposed by bacterial action.

PROCESS DESCRIPTION

PREPARATION OF WORKING BUFFER

Add the contend of the Reagent Nº 3 (Buffer Stock - 100 mL) in 400 mL of distilled or deionized water and homogenize well. Store in amber flask well sealed. The working buffer is stable for 12 months between 2 and 8°C.

PREPARATION OF WORKING REAGENT

Add 1 share of Reagent Nº 2 (Enzyme) in 20 parts of Working Buffer and homogenized it well. Example: Add 1.0 mL of Reagent Nº 2 in 20 mL of Working Buffer and homogenize it. Store in an amber flask. It is stable for 20 days between 2 and 8°C.

PREPARATION OF WORKING OXIDIZING

Dilute the Reagent Nº 4 (Oxidizing Stock) in 450 mL of distilled or deionized water. Keep it in an amber flask. Stable for 12 months between 2 and 8°C.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and Biocontrol P Bioclin Kits.

The calibration stability of the Urea Enzymatic kit installed on refrigerated equipment is at least 20 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

Mark 3 test tubes with letters B (Blank), A (Sample) and P (Standard), and proceed as follows:

	Blank	Standard	Sample
Sample	---	---	10 µL
Reagent Nº 1	---	10 µL	---
Working Reagent	1.0 mL	1.0 mL	1.0 mL
Homogenize and place in water-bath at 37°C, for 5 minutes.			
Working Oxidizing	1.0 mL	1.0 mL	1.0 mL
Homogenize and place in water-bath at 37°C, for 5 minutes.			

Read sample and standard absorbance at 600 nm (570 - 610 nm), matching the zero with Blank. The color is stable for 60 minutes. For the dosage of urea in urine, follow the described technique, using diluted sample 1:50 (100 µL of urine + 4.9 mL of distilled or deionized water). Multiply the obtained result by 50.

CALCULATIONS

$$\text{Urea (mg/dL)} = \frac{\text{Sample Absorbance}}{\text{Standard Absorbance}} \times 70$$

The reaction follows the Lambert Beer Law, it can use the Calibration Factor:

$$\text{Calibration Factor} = \frac{\text{Standard Concentration (70mg/dL)}}{\text{Standard Absorbance}}$$

$$\text{Urea (mg/dL)} = \text{Sample Absorbance} \times \text{Calibration Factor}$$

Urine

$$\text{Urea (g/24h)} = \frac{\text{Test Absorbance} \times \text{Factor X 50 X Volume (L)}}{100}$$

The results for serum or plasma will be express in mg/dL and for urine in g/24h.

PROCEDURE LIMITATIONS

The water, contamination of glassware, environment and/or from sample with Ammonia can lead to fake elevated results.

INTERFERENT

Some drugs, such as Aminoglycosides, Cephalosporins, Allopurinol, Methyl dopa, Furosemide and Propranolol, can interfere on the results, elevating the urea values.

Moderate hemolysis and Bilirubin up to 20 mg/dL do not interfere significantly in the results.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 912A.

REFERENCE VALUES

The reference values for this method were obtained through the determination of urea in healthy populations of male and female.

Serum or Plasma	
1 day - 12 months	2 - 34 mg/dL
1 - 13 years	8 - 36 mg/dL
Adults	15 – 40 mg/dL
Urine	
All ages	26 - 43 g/24h

To convert the mg/dL values to mmol/L (SI) multiply by 0.166. These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served. The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY

The Urea Enzymatic kit was compared with other commercially available methods for measurement of Urea Enzymatic. 07 analyzes were performed and the results were evaluated. The linear equation obtained was $Y = 0.973X + 1.466$ and the correlation coefficient 0.999. With these results, we can conclude that the kit shows good methodological specificity.

Precision

REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration(mg/dL)	24.90	37.85	46.00
Standard Deviation (mg/dL)	0.79	0.75	0.79
Coefficient of Variation (%)	3.16	1.97	1.73

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration(mg/dL)	24.95	37.88	45.98
Standard Deviation (mg/dL)	0.18	0.15	0.18
Coefficient of Variation (%)	0.72	0.40	0.38

Sensitivity

The sensitivity was calculated from 20 determinations of Urea. The average found was 0.0225 mg/dL, with standard deviation of 0.0005 mg/dL. The sensitivity, which indicates the detection limit of the method corresponds the average plus 3 times the standard deviation, and is equal to 0.024 mg/dL.

Linearity

The reaction is linear up to 250 mg/dL. For sample with value above 250 mg/dL, we recommend the sample dilution with deionized or distilled water, repeat the dosage and multiply the obtained value by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

The Urea, main product of proteins and amino acids catabolism, have its serum concentration affected by diet and by hydration states, constituting an rough indication of the kidney function condition.

Increased values of plasmatic urea may be originated from pre-kidney, kidney and post kidney causes. Pre Kidney cause: result of excrement malfunction, observed in the decompensate heart, hemorrhagic shock, acute dehydration, elevated protein catabolism (burns, fever). Kidney cause: consequence of the acute or chronic kidney disease (glomerulonephritis, pyelonephritis, tubular necroses) with reduced glomerular filtration, where we could observe plasmatic levels of urea 300 mg/dL or more. Post kidney: generally resultant of an urinary tract obstruction, being common in kidney lithiasis, and in tumors by bladder compression.

The decrease of serum urea occurs only in a few situations, such as acute hepatic insufficiency, in inanition, and in last trimester pregnancy.

NUMBER OF TESTS

500 Tests / 10 μ L of Sample / 2 mL of Reagent

BIBLIOGRAPHIC REFERENCES

- 1 - BERTHELOT, M. P. E., Clin. Appl., 1968, 284.
- 2 - BALLETER, W. D. E.; BUSITMAN, C. S.; TIDWELL, P. W., Anal. Chim., 33-59.
- 3 - WINDMANN, F. K.; TURNER, K., Clin. Chem, 1987, 21:1754-1770.
- 4 - BERGMAYER, HU., Methods of Enzymatic Analysis, vol. 9, VCH Publishers, 1985, 449-453.
- 5 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- 6 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

EC REP **OBELIS S.A.**

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for the Urea Enzymatic kit: 10269360101

Review: July/2020

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MANUFACTURED BY

BATCH CODE

CONTROL

DATE OF MANUFACTURE

POSITIVE CONTROL

USED BY
(last day of month)

NEGATIVE CONTROL

TEMPERATURE LIMITATION
(store at)

BIOLOGICAL RISK

CONTAINS SUFFICIENT
FOR <N> TESTS

INFLAMMABLE

CONSULT INSTRUCTIONS
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC DEVICE

POISON

EUROPEAN AUTHORIZED
REPRESENTATIVE

CE MARK

KEEP AWAY
FROM SUNLIGHT

DO NOT USE IF
PACKAGE IS
DAMAGED