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600 Pruebas (LCR)

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1 - GRADWOHL: Clinical Laboratory Methods and Diagnosis - Sonnen Writh and Jerret, Mosby Cd., 1980.

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USAGE INSTRUCTIONS

FUNCTION

Screening method to detect syphilis reagin. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Reaction of Flocculation

The combination of lecithin, cholesterol and cardioliolipin possui immunological similarity with antigens of *Treponema pallidum*, consisting of a non-treponemal antigen. The interaction of reagin in the sample with this antigen produces flocculation that can be detected by microscope optical.

REAGENTS

Reagent N° 1 - VDRL antigen in suspension - Store between 2 and 8°C. Do not freeze. Homogenize it well before use. Contains: Cardioliolipin 0,44 µmol/L, Lecithin 3,12 µmol/L and Cholesterol 23,2 µmol/L in Phosphate Buffer (pH 6,0) 10 mmol/L.

PRESENTATION

Reagent	Volume
N° 1	6 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Excavated slide, microscope, adjust the rotary shaker level at 180 rpm, pipettes, watch or stopwatch. They can be found at markets specialized on Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperature should be between 2 and 8°C. The transport temperature between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Antigen (Reagent N° 1) and the sample must be at room temperature just before use.
- 3- Do not freeze the reagent.
- 4- Do not use plasma.
- 5- The reagent should not come into contact with rubber materials.
- 6- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 7- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Advisory Service Customer) of Quibasa.
- 8- Do not use the product in case of damaged packaging.
- 9- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or cerebrospinal fluid (clear and free of clot fragments). Do not use hemolyzed, lipemic or suspected microbial contamination samples, as they may produce nonspecific agglutination. Turbid samples after thawing should not be used. **No need to inactivate the sample.** Fasting is recommended 8 hours before blood collection.

Serum: stable for 5 days at 2 to 8 ° C or 30 days at -20 ° C.

Cerebrospinal fluid: stable for 1 day between 2 and 8 ° C. Do not freeze.

DESCRIPTION OF PROCESS

Reagent N° 1 is ready for use. Shake it before performing the test.

TECHNIQUE

TESTING WITH SERUM

In a hollowed plate pipette:

Sample	50 µL
Reagent N° 1	20 µL

Agitate it manually or in a circular motion or in a mixer for 4 minutes at 180 rpm. After that, examine the microscope at 100 X.

Besides the sample without dilution it is recommended a test using a 1:8 dilution with 0,85% NaCl (for all samples) in order to avoid the prozone effect. See Process Limitations.

SEMI QUANTITATIVE EVIDENCE

Cavity N°	1	2	3	4	5	6
Dilution	1:2	1:4	1:8	1:16	1:32	1:64

Pipette 50 µL of NaCl 0,85% into each cavity from the carved slide. Cavity N° 1, pipette 50 µL of sample to be titrated, homogenize. Transfer 50 µL from the cavity N° 1 to N° 2 and so on until the desired dilution. Discard the excess 50 µL from last cavity. Then add to each 20 µL of Reagent N°1 to each cavity. Agitate manually using circular motions or rotary mixer for 4 minutes at 180 rpm. Look under the microscope at 100 X increase soon after stirring. The title comes from the highest dilution sample in which flocculation occurred.

TESTING WITH CEREBROSPINAL FLUID

- 1- Dilute Reagent N° 1 1:2 with NaCl 10%. Homogenize by inversion or rotation for 5 minutes. This solution is stable for 2 hours.
- 2- Proceed as described below, using centrifuged cerebrospinal fluid. It is not necessary to inactivate.

In a hollowed plate pipette:

Sample	50 µL
Reagent N° 1 (diluted)	10 µL

Agitate it manually or in a circular motion or in a mixer for 4 minutes at 180 rpm. After that, examine the microscope at 100 X.

Besides the sample without dilution it is recommended a test using a 1:8 dilution with 0,85% NaCl (for all samples) in order to avoid the prozone effect. See Process Limitations.

RESULTS

Positive - Reactive: Flocculation occurs with formation of clumps of varying sizes. Suspension with a heterogeneous aspect. In this case, proceed to dilution of the sample and performing the test semi-quantitative.

Negative - Non-Reactive: Absence of flocculation, suspension of homogeneous aspect.

PROCEDURE LIMITATIONS

It is estimated that 1 to 2% of patients with secondary syphilis have negative or weakly positive VDRL when using the undiluted serum. Positivation of the test occurs in higher dilutions due to the prozone effect. Recommendation was therefore to test all undiluted samples and diluted 1:8 with NaCl 0,85%.

There may yield false-positive in certain pathologies such as: tuberculosis, infectious mononucleosis, leprosy, hepatitis, collagen disease, autoimmune diseases, malaria and cancer.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

**PRODUCT PERFORMANCE
QUALITY CONTROL**

Sensitivity

In 163 truly positive samples, there were no false negative results found, thus, the sensitivity was 100%.

Specificity

In 80 truly negative samples, there was only 1 false positive result found, thus, the specificity was 98,75%.

DIAGNOSTIC SIGNIFICANCE

The method is suitable for a screening test. Final diagnosis should not be based solely on laboratory data. Results should be correlated with the signs and symptoms from clinical patient.

Treated patients infected with syphilis may have positive results in low titers (serological scar).

VDRL can be used to monitor therapy, represented by the rapid response by the fall in securities by negative changes and slower or even nonexistent for treponemal tests.

NUMBER OF TESTS

300 Tests (Serum)
600 Tests (CSF)

BIBLIOGRAPHIC REFERENCES

- 1 - GRADWOHL'S Clinical Laboratory Methods and Diagnosis - Sonnen Writhe and Jerret, Cd Mosby, 1980.
- 2 - Manual of Diagnosis reaction to el de la Syphilis - Nº 311 - World Health Organization
- 3 - FERREIRA, A.W.; Laboratory Diagnosis of Major Infectious and autoimmune diseases, 2nd ed., Publisher Guanabara Koogan SA, 2001.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda**
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

 **OBELIS S.A.**

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for VDRL Ready for Use kit: 10269360120

Review: June/2018

English 2/2

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED