FINALIDADE

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USE INSTRUCTIONS

PURPOSE

Rapid immunochromatographic test for qualitative determination of canine distemper antigen in samples of ocular or nasal discharge.

PRINCIPLE OF ACTION

Methodology: Immunochromatography.

Distemper Ag VET FAST is an *in vitro* diagnostic test for rapid and qualitative detection of distemper virus antigens in ocular or nasal secretions. The cassette consists of a nitrocellulose membrane, impregnated in the test line with antibody against canine distemper virus. The conjugate consists of antibody against distemper virus, conjugated with colloidal gold particles. When the sample is added to the cassette, it migrates by capillary action until it reaches the conjugate. In this region, the antigens present in the sample interact with the conjugate, forming the complex: antigen - conjugate.

The sample continues to migrate through the immunochromatographic strip until reaching the test line, where the antigens also interact with the antibody against canine distemper virus immobilized in the test line, forming the complex: immobilized antibody - antigen - conjugate. This shows the formation of a pink-red line on the test line; indicating positive result.

Sample excess continues to migrate through the test strip until it reaches the control line, which is composed of impregnated anti-conjugate antibody. Thus, the formation of a second line (control line) is observed; indicating that the test yielded a valid result.

REAGENTS

Reagent Number 1: Cassette with pipette in aluminized sachet - Store between 15 and 30°C. **Do not freeze.** It contains:

01 plastic cassette with 01 test strip containing:

- Control region (C): Antibody anti-mouse antibody immobilized.
- Test region (T): Immobilized canine distemper virus monoclonal antibody.
- Conjugate: Monoclonal mouse antibody against canine distemper virus, colloidal gold conjugate.

Reagent Number 2: Diluent - Store between 15 and 30°C. **Do not freeze.** Contains: Buffer solution and preservative.

Reagent 3: Sterile flexible cotton swab.

PRESENTATIONS

Reagent	Presentation	Number of Cassettes for Package	Diluent
Cassette + Diluent	VET005-1	10	10 x 0,5 mL
	VET005-2	25	25 x 0,5 mL
	VET005-3	10	10 x 0,5 mL
	VET005-4	25	25 x 0,5 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Clock or stopwatch. They are in the specialized market of articles for Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperature of the product must be between 15 and 30°C. Keep away from direct sunlight.

Do not freeze.

The cassette is sensitive to air humidity. Therefore, it should always be kept in its inviolate primary packaging (aluminized sachet) until the time of use.

SPECIAL CARE

1- Product for *in vitro* diagnosis, for veterinary use only.

2- Strictly follow the methodology proposed for obtaining results.

3- Cassettes (test devices) do not contain contaminants; but since the samples are biological, the professional must strictly follow the safety rules and routines. The use of disposable gloves and other personal protective equipment is essential.

4- Use clean, dry containers to package biological samples.

5- Open the aluminum envelop to remove the cassette only at the time of use.

6- The kit cannot be frozen or exposed to high temperatures, since this will cause deterioration of its components.

7- The reagents, as well as the samples, should reach room temperature before starting the test.

8- Do not use reagents from different lots or kits.

9- Do not use kit beyond its expiration date. The expiry date is stated on the product label.

10- The cassette (test device), the sterile flexible cotton swab and the plastic pipette are disposable and cannot be reused.

11- The kit materials that had contact with the sample of the patient must have a special destination. The disposal of sterile flexible cotton swab, the plastic pipettes, the cassettes (test devices) and diluent vials should be disposed of in organic waste.

12- Secondary packaging (box). Instruction for use, aluminized sachet and silica sachet can be disposed of in trash common.

13- We recommend to apply the local, state and federal environmental protection so that the disposal of reagents and biological material is done in accordance with the current legislation.

14- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or through request by the Customer Service Department (SAC) of the Quibasa.

SAMPLES

Perform the test preferably with ocular (conjunctive) or nasal secretion samples.

Test the sample preferably immediately after collection.

If necessary, samples can be stored between 2 to 30°C for up to 2 days or at -20°C for longer periods.

DESCRIPTION OF PROCESS

1- The sample should be at room temperature before starting the test.

2- Remove the cassette (test device) from the protective packaging and place it on a clean, dry and level surface

3- Collect the eye discharge or nasal discharge samples using the sample collection swab.

Note: Use a gentle twisting motion during sample collection to completely wet the swab with the secretion.


4- Insert the swab into the diluent tube, dissolving the sample.

Note: Use rotary movements of the swab inside the tube and rub it against the tube wall to maximize dissolution of the secretion in the diluent.

5- Use the pipette that comes with the kit to collect the dissolved sample.

6- Add 4 drops of sample (80 µL), slowly, into the sample application well of the cassette.

7- Wait for the formation of the lines. Interpret the results between 5 and 10 minutes. Do not interpret after 10 minutes.

INTERPRETATION OF RESULTS

Reactive Test: Formation of two lines: one in the control region (C) and one in the test region (T). Do not interpret after 10 minutes.	
Not Reactive Test: Line formation only in control region (C). Do not interpret after 10 minutes.	
Invalid Test: The absence of line formation in the control region (C), with or without line formation in the test region (T), indicates procedure error or deterioration of the cassette. In this case, repeat the test using a new cassette.	

TEST LIMITATIONS

This product is intended to detect canine distemper virus infections. A negative result can be observed if the amount of analyte present in the sample is below the detection limit of the test. This can happen, depending on disease period in which the sample was collected, the type of sample used in the test and its storage conditions. A negative result may also occur in case of inefficiency in the sample extraction step in the diluent.

Therefore, for the diagnosis of the disease, the professional must take into account not only the result of this test; but also, clinical signs, animal health conditions, and other diagnostic results that may be available.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal program of quality control, where procedures, standards, limits and tolerance for variations are clearly established. It is important to note that all measuring systems characteristic analytical variability, which should be monitored by the laboratories themselves. Therefore, it is the use of controls, which allow the assessment of the precision and accuracy of dosages.

PRODUCT PERFORMANCE

QUALITY CONTROL

CLINICAL SENSITIVITY AND SPECIFICITY

This product has been tested comparing with other methods. A total of 108 samples characterized by canine distemper were analyzed; 41 samples were diagnosed as positive and 67 were diagnosed as negative. The results show that the clinical sensitivity of the product is 97.56% and the clinical specificity is 97.01%.

METHOD		REFERENCE		Total	
Cassette Distemper Ag	Result	Reactive	Not reactive		
	Reactive	40	2		42
	Not reactive	1	65		66
Total result		41	67	108	

Clinical sensitivity: (40/41): 97,56%

Clinical Specificity: (65/67): 97,01%

Precision: (40+59/41+61): 97,22%

ANALYTICAL SENSITIVITY

The product was tested with samples of tear secretion, which contained a determined virus load. The product was able to identify with positive result all the concentrations tested with a value equal to or greater than 103 TCID₅₀ / mL. Therefore, the analytical sensitivity (or detection limit) of 103 TCID₅₀ / mL was established.

Precision

REPEATABILITY

Repeatability was calculated from 10 successive determinations, using 2 different samples, obtaining the following results:

Samples	Nº of Repetitions	Expected Result	Found Result
Sample 1	10	Reactive	Reactive
Sample 2	10	Reactive	Reactive

REPRODUCIBILITY

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 2 different samples, obtaining the following results:

Samples	Nº of Days	Nº of Repetitions	Expected Result	Found Result
Sample 1	03	10	Reactive	Reactive
Sample 2	03	10	Reactive	Reactive

BIBLIOGRAPHIC REFERENCES

1- QUIBASA: Data from the Research and Development Department.

QUALITY GUARANTEE

Before being released for consumption, all Quibasa reagents are tested by the Quality Control Department. The quality of the reagents is assured until the validity date mentioned on the packaged, as long as they are stored and transported under the appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Tel.: (31) 3439.5454
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Industria Brasileira

CONSUMER SERVICE

Customer Support Service

Tel.: 0800 0315454

E-mail: sac@bioclin.com.br

Product licensed in the Ministry of Agriculture, Livestock and Supply since 11/11/2019 with the number 10.282/2019.

Technical manager: Dr. Hugo Vieira Fróes (CRMV/MG 8963).

Review: June/2020

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED