FIV Ac / FeLV Ag COMBO VET FAST

REF VET013

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FINALIDADE

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USE INSTRUCTIONS

PURPOSE

Fast combined immunochromatographic test for the qualitative determination of Antibodies against the virus of feline immunodeficiency (FIV) and feline leukemia virus Antigen (FeLV) in serum, plasma or whole blood samples.

PRINCIPLE OF ACTION

Methodology: Immunochromatography.

The FIV Ac / FeLV Ag Combo VET FAST is an in vitro diagnostic test for the rapid and qualitative detection of total Antibodies against feline immunodeficiency virus (FIV Ac) and feline leukemia virus Antigen (FeLV Ag) in feline blood, serum or plasma.

The test consists of two independent cassettes, linked by a support (combo cassette): one detects Anti-FIV Antibodies and the other detects the FeLV Antigen. Each cassette consists of a nitrocellulose membrane, impregnated in the test line with the feline immunodeficiency virus Antigen (FIV cassette) and Antibodies against feline leukemia virus (FeLV cassette). The conjugate consists of FIV Antigen (FIV cassette) and Anti-FeLV Antibody (FeLV cassette); both conjugated with colloidal gold particles.

When the sample is added to the cassettes, it migrates by capillary action through the test strip until it reaches the conjugate. The following occurs in this region:

In the FIV cassette: Anti-FIV Antibodies, if present in the sample, interact with the conjugate, formed of colloidal gold conjugated FIV Antigen. Therefore, formation of the complex occurs: Anti-FIV Antibody – FIV Antigen with colloidal gold. The sample continues to migrate through the test strip until it reaches the test line, where the Anti-FIV Antibody also interacts with the immobilized FIV Antigen in the test line, forming the complex: conjugate - Anti-FIV Antibody - Antigen immobilized. This shows the formation of a pink-red line on the test line; indicating positive result.

In the FeLV cassette: the FeLV Antigen, possibly present in the sample, interacts with the conjugate, formed of colloidal gold conjugated Anti-FeLV Antibodies. This is how complex formation occurs: FeLV Antigen - colloidal gold conjugated Anti-FeLV Antibody. The sample continues to migrate through the test strip until it reaches the test line, where the FeLV Antigen also interacts with the Anti-FeLV Antibody immobilized in this area, forming the complex: conjugate - FeLV Antigen - immobilized Anti-FeLV Antibody. This shows the formation of a pink-red line on the test line; indicating positive result.

In both cassettes, sample excess continues to migrate through the test strip until it reaches the control line, which is composed of impregnated Anti-conjugate Antibody. Thus, the formation of a second line (control line) is observed; indicating that the test yielded a valid result.

REAGENTS

Reagent Number 1 (R1): Cassette - Store between 15 and 30°C. **Do not freeze.** It contains:

A) FIV cassette:

- 01 plastic cassette with 01 test strip containing:
- Control Region (C): Goat immobilized Anti-mouse Antibody.
- Test Region (T): Immobilized FIV Antigen.
- Conjugate: FIV Antigen colloidal gold conjugate and mouse Antibody colloidal gold conjugate.

B) FeLV cassette:

- 01 plastic cassette with 01 test strip containing:
- Control Region (C): Immobilized Anti-mouse Antibody.
- Test Region (T): Immobilized monoclonal Antibody against feline leukemia virus.
- Conjugate: Anti-FeLV colloidal gold conjugate and mouse Antibody colloidal gold conjugate.

C) 01 Sachet of Silica

Reagent Number 2 (R2): Diluent - Store between 15 and 30°C. **Do not freeze.** Contains: Buffer solution and preservative.

Collection Accessories:

- Dispesable Plastic Pipette

PRESENTATIONS

Reagent	Presentation	Number of Cassettes for Package	Diluent
Cassette + Diluent	VET013-1	10	1 x 2,5 mL
	VET013-2	25	1 x 5,0 mL
	VET013-3*	10	10 x 0,5 mL
	VET013-4*	25	25 x 0,5 mL
	VET013-5	5	1X1,25 mL

*Individual presentations

EQUIPMENTS AND OPERATIONAL INPUTS

Clock or stopwatch. These itens are found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperature of the product must be between 15 and 30°C. Keep away from direct sunlight. **Do not freeze.** The combo cassette is sensitive to air humidity. Therefore, it should always be kept in its inviolate primary packaging (aluminized sachet) until the time of use.

SPECIAL CARE

- Product for *in vitro* diagnosis, for veterinary use only.**
- Strictly follow the methodology proposed for obtaining results.
- Cassettes (test devices) do not contain contaminants; but since the samples are biological, the professional must strictly follow the safety rules and routines. The use of disposable gloves and other personal protective equipment is essential.
- Use clean, dry containers to package biological samples.
- Open the aluminum envelop to remove the combo cassette only at the time of use.
- The kit cannot be frozen or exposed to high temperatures, since this will cause deterioration of its components.
- The samples should reach room temperature before starting the test.
- Do not use reagents from different lots or kits.
- Do not use kit beyond its expiration date. The expiry date is stated on the product label.
- The combo cassette (test devices) and the capillary tube are disposable and cannot be reused.
- Blood samples with clots, microcoagulants and fibrin mesh formation can cause the filter mesh to clog, making it impossible for the sample to pass through the test membrane correctly.
- The kit materials that had contact with the sample of the patient must have a special destination. The disposal of capillary tubes, the cassettes (test devices) and diluent vials should be disposed of in organic waste.
- Secondary packaging (box), Instruction for use, envelope aluminum and silica sachet can be disposed of in trash common.
- We recommend to apply the local, state and federal environmental protection so that the disposal of reagents and biological material is done in accordance with the current legislation.
- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or through request by the Customer Service Department (SAC) of the Quibasa.

SAMPLES

Perform the test with serum, plasma or whole blood.

Serum: Allow blood to clot and then centrifuge to separate serum.

Plasma: Collect whole blood in a tube containing Heparin, Citrate, or EDTA. Centrifuge the blood and separate the plasma.

Whole Blood: Collect the blood sample in a tube containing EDTA, Heparin, or Citrate using standard laboratory procedures.


Whole blood samples should be analyzed within 24 hours, can be stored at 2 and 8°C and never be frozen.

Serum and plasma samples should be analyzed immediately after collection, but if necessary, they can be stored at 2 and 8°C for up to 7 day sor at -20°C for longer periods.

VETERINARY USE

DESCRIPTION OF PROCESS

- The sample should be at room temperature before starting the test.
- Remove the combo cassette (test devices) from the protective packaging and place it on a clean, dry and level surface.
- Collect the sample using the disposable capillary tube that accompanies the kit.
- Tighten the capillary tube above the marked trace and touch the sample. Relieve the pressure so the sample is aspirated. Aspirate until the trace marked to obtain the volume of 20 µL.
- Dispense the sample by pressing the pipette again, touching the test membrane at the sample well.
- Add the samples in both cassettes (FIV Ac and FeLV Ag). **Note:** Alternatively, an automatic pipette can be used to collect and add the sample.
- Hold the Diluent bottle (Reagent N° 2) vertically and apply 2 drops (70 µL) in sample well. Add the Diluent in both cassettes.
- Wait for the lines to form. Interpret the results between 15 and 20 minutes. Do not interpret after 20 minutes.

INTERPRETATION OF RESULTS BY CASSETTE

Reagent Test: Formation of two lines: one in the control region (C) and one in the test region (T). Do not interpret after 20 minutes.	
No Reactive Test: Line formation only in control region (C). Do not interpret after 20 minutes.	
Invalid Test: The absence of line formation in the control region (C), with or without line formation in the test region (T), indicates procedure error or deterioration of the cassette. In this case, repeat the test using a new cassette.	

Note: Each cassette must be interpreted individually.

TEST LIMITATIONS

This product is intended for detection of infection by feline immunodeficiency virus (FIV - detection of Antivirus Antibodies) and by feline leukemia virus (FeLV - detection of viral Antigen). A no reactive result can be observed if the amount of analyte present in the sample is below the detection limit of the test. This may happen depending on the period of the disease in which the sample was collected, the type of sample used in the test or the storage conditions. In these cases, consult the responsible veterinarian and consider repeating the test after a few weeks. For the diagnosis of the disease, the professional must take into account not only the result of this test; but also, clinical signs, animal health conditions, and other diagnostic results that may be available.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal program of quality control, where procedures, standards, limits and tolerance for variations are clearly established. It is important to note that all measuring systems characteristic analytical variability, which should be monitored by the laboratories themselves. Therefore, it is the use of controls, which allow the assessment of the precision and accuracy of dosages.

PRODUCT PERFORMANCE
QUALITY CONTROL

CLINICAL SENSITIVITY AND SPECIFICITY

This product has been tested comparing with other methods. A total of 102 samples of FIV Ac were analyzed, which 41 samples were diagnosed as reagent and 61 were diagnosed as no reactive. The results show that the clinical sensitivity of the product is 97.56% and the clinical specificity is 96.72%.

METHOD		REFERENCE		Total
FIV Ac Cassette	Result	Reagent	No Reactive	
	Reagent	40	2	
	No Reactive	1	59	
Total Result		41	61	102

Clinical Sensitivity: (40/41): 97.56%
Clinical Specificity: (59/61): 96.72%
Precision: (99/102): 97.05%

A total of 104 samples of FeLV Ag were analyzed as well, which 46 samples were diagnosed as reactive and 58 were diagnosed as no reactive. The results show that the clinical sensitivity of the product is >99% and the clinical specificity is 98.27%.

METHOD		REFERENCE		Total
FeLV Ag Cassette	Result	Reagent	No Reactive	
	Reagent	46	1	
	No Reactive	0	57	
Total Result		46	58	104

Clinical Sensitivity: (46/46): >99%
Clinical Specificity: (57/58): 98.27%
Precision: (103/104): 99.04%

Precision
REPEATABILITY
Repeatability was calculated from 10 successive determinations, using 2 different samples of FIV Ac and 2 different samples of FeLV Ag. The following results was obtained:

Samples	N° of Repetitions	Expected Result	Found Result
FIV Ac Sample 1	10	Reagent	Reagent
FIV Ac Sample 2	10	Reagent	Reagent
FeLV Ag Sample 1	10	Reagent	Reagent
FeLV Ag Sample 2	10	Reagent	Reagent

REPRODUCIBILITY
Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 2 different samples of FIV Ac and 2 different samples of FeLV Ag. The following results were selected:

Samples	N° of Days	N° de Repetitions	Expected Result	Found Result
FIV Ac Sample 1	03	10	Reagent	Reagent
FIV Ac Sample 2	03	10	Reagent	Reagent
FeLV Ag Sample 1	03	10	Reagent	Reagent
FeLV Ag Sample 2	03	10	Reagent	Reagent

CLINICAL SIGNIFICANCE

The Feline Immunodeficiency Virus (FIV) and the Feline Leukemia Virus (FeLV) are the cause of important diseases that reduce the life expectancy of cats. FIV transmission occurs mainly through bites and scratches, while FeLV transmission occurs through contact with secretions such as saliva and sharing objects, and the incubation period of these diseases can vary from months to years.

In FIV infection, the main clinical signs observed are fever and generalized lymphadenopathy, chronic opportunistic infections with a high rate of recurrence, neurological alterations and tumors. After infection, there is a significant and persistent increase in antibodies against the virus, which makes antibody detection a good diagnostic tool. In FeLV infection, the appearance of tumors is common and animals often become susceptible to secondary infections with low efficacy in its treatment. The production of antibodies against FeLV is deficient and nonspecific, while the circulation of viral antigen is high and persistent, making detection of viral antigen a safe and effective tool for diagnosis.
It is common for asymptomatic cats to occur, so the diagnosis is essential for the application of the appropriate behavior with the animal and for the control of the spread of these diseases.

BIBLIOGRAPHIC REFERENCES

- 1- QUIBASA: Data from the Research and Development Department.
- 2- MARTINS ES, TOGNOLI GK, ILHA PHO, et al. Prevalência de imunodeficiência viral felina e leucemia viral felina no Distrito Federal. Archives of Veterinary Science, v.17, n.S1, p.274-276, 2012.
- 3- HARBOUR DA, CANEY SMA, SPARKES AH. Feline Immunodeficiency Virus Infection. p. 607-622. Cap.24. IN: Chandler EA, Gaskeç CJ, Kaskel RM. Feline Medicine and Therapeutics. Blackwell: Wiley. 2004.
- 4- BIEZUS G, MACHADO G, FERIAN PE, COSTA UM, et al. Prevalence of and factors associated with feline leukemia virus (FeLV) and feline immunodeficiency virus (FIV) in cats of the state of Santa Catarina, Brazil. Comparative Immunology, Microbiology and Infectious Disease. v.63, p.17-21, 2019.

QUALITY GUARANTEE

Before being released for consumption, all Quibasa reagents are tested by the Quality Control Department. The quality of the reagents is assured until the validity date mentioned on the packaged, as long as they are stored and transported under the appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda**
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Tel.: (31) 3439.5454
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Industria Brasileira

CONSUMER SERVICE
Customer Advisory Service
Tel.: 0800 0315454
E-mail: sac@bioclin.com.br

Product licensed in the Ministry of Agriculture, Livestock and Supply since 11/11/2019 with the number 10.283/2019.

Technical manager: Dra. Camila Eckstein (CRMV/MG 20611)

Review: September/2021

UNIVERSAL SYMBOLOGY

REF

CATALOG NUMBER

LOT

BATCH CODE

DATE OF MANUFACTURE

USED BY
(last day of month)

TEMPERATURE LIMITATION
(store at)

CONTAINS SUFFICIENT
FOR <N> TESTS

CONSULT INSTRUCTIONS
FOR USE

IVD

IN VITRO DIAGNOSTIC DEVICE

EC

REP

EUROPEAN AUTHORIZED
REPRESENTATIVE

KEEP AWAY
FROM SUNLIGHT

MANUFACTURED BY

CONTROL

CONTROL

CONTROL

CONTROL +

CONTROL +

POSITIVE CONTROL

CONTROL -

CONTROL -

NEGATIVE CONTROL

BIOLOGICAL RISK

INFLAMMABLE

CORROSIVE

POISON

CE

CE MARK

DO NOT USE IF
PACKAGE IS
DAMAGED