WRIGHT CORANTE HEMATOLÓGICO

REF: C006

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Staining method for peripheral blood cells, bone marrow or cytological study of cellular elements. Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Dye according to Wright
Wright Dye is a mixture of dye with neutral characteristics, which stain the nuclear and cytoplasmic components of cells.

REAGENTS

Number 1 - Wright Dye - Store between 15 and 30°C. Contains: Eosin Methylene Blue and solubilizer.

PRESENTATION

Reagent	Volume
Nº1	1000 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Support for staining, watch or stopwatch, microscope, pipettes and blades. They can be found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transportation temperature should be between 15 to 30°C. Protect from light and avoid moisture.

SPECIAL CARE

- For professional *in vitro* diagnostic use only.**
- Strictly follow the methodology proposed.
- Water used in material cleaning must to be recent and free of contaminants.
- Handle reagent with care, avoid contact with skin and mucous.
- The slides for the fabrication of smears should be clean and fat-free.
- A drop of blood used for making the smear should not be too large because the larger the droplet the more thicker the smear.
- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Advisory Service Customer) of Quibasa.

- Do not utilize the product in case of damages on the packaging.
- It is essential that the instruments and equipments utilized stay properly calibrated and subjected to periodic maintenances.

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Precipitation may appear when the slides are unwashed, dust or drying of the dye occurs during the staining process. Eventually the water used by the laboratory does not provide a good coloring. In this case use a Phosphate Tampon (pH 7,2), using the following formula:

1 - Solution A	
KH ₂ PO ₄	9,1 g
Distilled H ₂ O q.s.p.	1.000 mL

2 - Solution B	
Na ₂ HPO ₄	9,5 g
Distilled H ₂ O q.s.p.	1.000 mL

Usage Tampon

Solution A	28 mL
Solution B	72 mL

PRODUCT PERFORMANCE

CONTROL QUALITY

Accuracy

COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY

Wright Hemotologic Dye kit was compared with other smears coloring methods commercially available.

According to the coloring comparison of 20 clinical samples, the method shows an excellent correlation.

Precision

REPEATABILITY

The repeatability was calculated from 10 successive colorations, using 1 sample, obtaining the following results:

Number	Sample 1
1	Standard Coloring
2	Standard Coloring
3	Standard Coloring
4	Standard Coloring
5	Standard Coloring
6	Standard Coloring
7	Standard Coloring
8	Standard Coloring
9	Standard Coloring
10	Standard Coloring

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4 - NARASIMHA, A., KUMAR, H., PRASAD, CSBR. Anticoagulant induced artefacts in peripheral blood smears. Indian Journal of Hematology & Blood Transfusion. 2008;24(2):43-48.

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REPRODUCIBILITY

The reproducibility was calculated from 10 successive colorations for 3 consecutive days, using 1 sample, obtaining the following results:

Days	Nº of Repetitions	Expected Results	Found Results
1º Day	10	Standard Coloring	Standard Coloring
2º Day	10	Standard Coloring	Standard Coloring
3º Day	10	Standard Coloring	Standard Coloring

BIBLIOGRAPHIC REFERENCES

1 - WALF., Pratical Clinical Hematology Interpretation and Techniques, 375-376,1973.

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QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax.: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

EC

REP

OBELIS S.A.
Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for Wright Hematologic Dye kit: 10269360095

Review: October/2019

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED