BIOCAL

REF **K072**

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1. Dati F. Reference materials and guidelines for standardization of methods in laboratory medicine. In: Thomas L, editor. Clinical laboratory diagnostics. 1st ed. Frankfurt:TH-Books Verlagsgesellschaft; 1998. p. 1404-26.

2. Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).

3. Moss DW, Henderson AR. Enzymes. In: Burtis CA, Ashwood ER, editors. Tietz textbook of Clinical Chemistry. 2nd ed.philadelphia: W.B Saunders Company;1994, p.735-896.

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- 1. Dati F. Reference materials and guidelines for standardization of methods in laboratory medicine. In: Thomas L, editor. Clinical laboratory diagnostics. 1st ed. Frankfurt:TH-Books Verlagsgesellschaft; 1998. p. 1404-26.
- 2. Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).
- 3. Moss DW, Henderson AR. Enzymes. In: Burtis CA, Ashwood ER, editors. Tietz textbook of Clinical Chemistry. 2nd ed.philadelphia: W.B Saunders Company;1994. p.735-896.
- 4. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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BIOCAL

REF **K072**

USAGE INSTRUCTIONS CALBRATOR FOR BIOCHEMISTRY

FUNCTION

Reagent used for calibrating automated analytical testing methodologies. For *in vitro* diagnostic use only.

REAGENTS

Number 1 - Calibrator - Store between 2 and 8°C. Contains: Calibrator containing approximately 30 analytes embedded in a lyophilized human matrix and preservative **Potentially Infectious**.

Warning: The concentration of analytes varies each lot. See table that accompanies the product.

PRESENTATION

Presentation	Volume
1	1 x 3 mL
2	5 x 3 mL
3	1 x 5 mL
4	5 x 5 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Automated biochemistry analyzers, watch or stopwatch, pipettes, test tubes and Bioclin kits. These items are found at markets specialized on Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture.

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.**
- Strictly follow the methodology proposed to obtain exact results.
- Water used in material cleaning must be recent and free of contaminants.
- The reagent should be handled cautiously since it is subject to biological contamination.
- To avoid possible analytical errors, it is advisable that the distilled water used for reconstitution of the product is of adequate quality. For this, it is recommended that water used in clinical laboratories follow the specifications of the Clinical Laboratory Standards Institute (CLSI).
- Handle with care, as the reagent contains Sodium Azide, which is irritant to skin and mucous.
- The reagent has been tested for anti-HIV antibodies, anti-HCV and HBs antigen using state of the art methods and found negative results. The risk of infection cannot be excluded and the reagent should be handled as cautiously as observed for the patient's serum. **Potentially Infectious**.
- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

9- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

10- Do not use the product in case of damaged packaging.

11- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

PROCESS DESCRIPTION

RECONSTITUTION OF CALBRATOR

Before reconstitute the Calibrator should be noted the volume indicated on the bottle and proceed as described below.

Presentation 1 and 2: Open carefully the bottle and add exactly **3 mL of distilled water**, as volume described in the label of the bottle.

Presentation 3 and 4: Open carefully the bottle and add exactly **5 mL of distilled water**, as volume described in the label of the bottle.

Mix. Leave the bottle resting for 30 minutes. Every 10 minutes, agitate the bottle carefully using circular motion to prevent material sticking to the wall of the bottle. Avoid formation of foam.

After reconstitution, observe the stability according to the table below.

	-20°C	+4°C	+25°C
Bilirubine (protected from the light)	14 days	8 hours	4 hours
CK and CK-MB	30 days	6 hours	2 hours
Other parameters	30 days	2 days	8 hours

Warning:

- The reconstituted Calibrator should not be thawed more than once. Therefore, for minor routines, fractionate the material into aliquots before freezing.

TRACEABILITY

The concentrations of the analytes present in the BIOCAL calibrator are traceable to the standard reference material of NIST (National Institute of Standards and Technology), IFCC (International Federation of Clinical Chemistry) and samples with consensus values from nationals and internationals proficiency programs.

See Bioclin kits Usage Instructions.

REFERENCE VALUES

Concentration of analytes vary with each batch.
See table that accompanies the product.

BIBLIOGRAPHIC REFERENCES

- 1. Dati F. Reference materials and guidelines for standardization of methods in laboratory medicine. In: Thomas L, editor. Clinical laboratory diagnostics. 1st ed. Frankfurt:TH-Books Verlagsgesellschaft; 1998. p. 1404-26.
- 2. Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).
- 3. Moss DW, Henderson AR. Enzymes. In: Burtis CA, Ashwood ER, editors. Tietz textbook of Clinical Chemistry. 2nd ed.philadelphia: W.B Saunders Company;1994, p.735-896.
- 4. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda**
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for Biocal kit: 10269360123

Review: March/2022

UNIVERSAL SYMBOLOGY

	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER
	MADE BY		