BIOCLIN FAST COVID-19 Ag

[REF] K244

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FINALIDADE

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
USAGE INSTRUCTIONS

FUNCTION

Test for the qualitative determination of antigens for SARS-CoV-2, the virus that causes COVID-19 (severe acute respiratory syndrome coronavirus 2) in biological **nasopharyngeal swab samples**. Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Immunochromatography
 The BIOCLIN FAST COVID-19 Ag kit is an immunochromatographic test that detects antigens of the SARS-CoV-2 Virus employing mouse monoclonal antibodies anti N protein (nucleocapsid) and RBD (receptor binding domain) present in the sample. The test line contains Anti-N and Anti-RBD monoclonal antibodies. After adding the sample, colloidal gold particles conjugated Anti-N and Anti-RBD bind to viral antigens and migrate by capillarity across the membrane. The colloidal gold complex and viral antigen are captured in the test line (T), leading to the formation of a reddish line indicating the presence of SARS-CoV-2 in the sample. The non-formation of a reddish line in the T region indicates the absence of a viral particle in the sample or an amount of Virus below the detection limit. The sample-antigen mixture continues running through the membrane to the control line, where a reddish line is formed, both for positive and negative samples. The formation of the reddish line in the Control region (C) demonstrates that the test is suitable for use.

REAGENTS

Reagent Number 1 (R1) Cassette – Store at 2 to 30°C. **Do not freeze**. Contains:

- A) 01 plastic device with 01 test strip containing:
- Control region (C): Goat antibody Anti-mouse (> 0.06 µg).
- Test region (T): Anti-N (> 0.3 µg), Anti-RBD (> 0.3 µg)
- Conjugate: Colloidal gold particles conjugated with mouse monoclonal antibody Anti-N and Anti-RBD.

- B) 01 silica sachet

Reagent Number 2 (R2) Collection Tube - Store at 2 to 30°C. **Do not freeze**. Contains: Buffer <100 mmol/L pH 7.4 and preservative

Accessories for Sample Collection:

- Dropper cap (attached to the collection tube)
- Disposable Swab

PRESENTATION

Reagent	Nº of Cassettes per packaging	Nº of Collection tube per packaging
Cassette + Collection Tube (Included dropper cap and disposable swab) *	05	05
	10	10
	25	25
	50	50
	100	100
Reagent	Nº of tests per packaging	
Aluminized pouch, contains:	05	
1 Cassette	10	
1 Collecting Tube	25	
1 Dropper Cap	50	
1 Disposable Swab	100	

* The number of accessories is proportional to the number of cassettes and collecting tubes.

INSTRUMENTS AND OPERACIONAL INPUTS

Clock or stopwatch, PPE (glove, mask, goggles), and garbage for disposal. All the inputs listed are in the specialized market for articles for Clinical Analysis Laboratories.

CONDITIONS OF STORAGE AND TRANSPORTATION

The storage and transport temperature should be 2 to 30°C. Avoid moisture. **Do not freeze**.

WARNINGS

- 1- For professional *in vitro* diagnostic use only.**
2- Strictly follow the methodology proposed to obtain exact results.
3- The used device has potential infectious use. The safety routinemust be followed strictly and the use of disposable gloves is indispensable.
4- Use clean and dry recipients to settle the samples.
5- Open the aluminum envelop, to take the cassette off, only for its immediate use.
6- The kit must not be frozen or exposed to temperatures above 45°C this exposure may cause its deterioration.
7- The reagents, as well as the sample, must reach the environment temperature before the test begins.
8- Do not mix reagents between kits with the same batch or with different batch numbers.

9- Do not use the kit after the expiration date. The expiry date is indicated on the product label.

10- The cassette (test device), collecting tube, dropper cap, and swab are disposable and cannot be reused.

11- The kit materials that had contact with the patient's biological material must have a special destination. The collecting tubes and swabs, and Cassettes (test devices) must be discarded in biological waste.

12- Boxes, instructions for use, and aluminized pouches can be disposed of in regular trash.

13- We recommend applying local, state, and federal environmental protection standards so that the disposal of reagents and biological material is done following current legislation.

14- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Safety Information Sheet for Chemical Products) available on the website www.bioclin.com.br or upon request by Quibasa's Customer Service Department.

15- Do not use the product in case of damage to the packaging.

SAMPLES


Nasopharynx swab. The sample must be tested preferably after collection or should be stored between 2-8°C for 24 hours.

DESCRIPTION OF THE TECHNICAL PROCESS

The stability of the BIOCLIN FAST COVID-19 Ag cassette (after opening the aluminized pouch) is 2 hours. It is stability may vary according to the environmental conditions after opening the aluminized pouch.

TECHNIQUE

- A) Segregate the collection material: collecting tube (dropper cap attached), swab.
B) Nasopharynx sampling:
- Use the swab included in the kit to collect the nasopharynx sample.
- The sample should be tested preferably after collection.

- C) Make sure that the buffer of the collection tube is at the bottom of the tube and remove the dropper cap.


D) After collection, insert the swab into the collection tube.

E) Mix for 10 seconds to ensure good extraction through the buffer.

F) Remove excess buffer at the tip of the swab by pressing the walls of the collection tube against the head of the swab.

G) Insert the dropper cap into the collection tube.

H) Add three drops of the buffer (90-100 µL) to the sample well. Read the test between 15 and 20 minutes. Do not interpret after 20 minutes.

INTERPRETATION OF THE RESULT

Reagent Test: Formation of a reddish line in the test region (T) and another line in the control region (C). Do not interpret after 20 minutes.