Bioclin

BIOCLIN FAST COVID-19 Ag

[REF] K244

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Test for the qualitative determination of antigens for SARS-CoV-2, virus which causes COVID-19 (severe acute respiratory syndrome coronavirus 2) in biological **nasopharyngeal swab and nasal swab** samples. *Only for in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Immunochromatography
The BIOCLIN FAST COVID-19 Ag kit is an immunochromatographic test that detects antigens of the SARS-CoV-2 Virus employing mouse monoclonal antibodies anti N protein (nucleocapsid) and RBD (receptor binding domain) present in the sample. The test line contains Anti-N and Anti-RBD monoclonal antibodies. After adding the sample, colloidal gold particles conjugated Anti-N and Anti-RBD bind to viral antigens and migrate by capilarity across the membrane. The colloidal gold complex and viral antigen are captured in the test line (T), leading to the formation of a reddish line indicating the presence of SARS-CoV-2 in the sample. The non-formation of a reddish line in the T region indicates the absence of a viral particle in the sample or an amount of Virus below the detection limit. The sample-antigen mixture continues running through the membrane to the control line, where a reddish line is formed, both for positive and negative samples. The formation of the reddish line in the Control region (C) demonstrates that the test is suitable for use.

REAGENTS

Reagent Number 1 (R1) Cassette – Store at 2 to 30°C. **Do not freeze.** Contains:

- A) 01 plastic device with 01 test strip containing:
- Control region (C): Goat antibody Anti-mouse (> 0.06 µg).
 - Test region (T): Anti-N (> 0.3 µg), Anti-RBD (> 0.3 µg) - Conjugate: Colloidal gold particles conjugated with mouse monoclonal antibody Anti-N and Anti-RBD.

- B) 01 silica sachet

Reagent Number 2 (R2) Collection Tube - Store at 2 to 30°C. **Do not freeze.** Contains: Buffer < 100 mmol/L pH 7.4, and preservatives

Accessories for Sample Collection:

- **Dropper cap (attached to the collection tube)**
- **Disposable Swab.**

PRESENTATION

Reagent	Nº of Cassettes per packaging	Nº of Collection tube per packaging
Cassette + Collection Tube (Included dropper cap and disposable swab) *	05	05
	10	10
	25	25
	50	50
	100	100


DESCRIPTION OF THE TECHNICAL PROCESS

The BIOCLIN FAST COVID-19 Ag cassette's stability (after opening the aluminized pouch) is 2 hours. It is stability may vary according to the environmental conditions after opening the aluminized pouch.

TECHNIQUE

- A) Segregate the collection material: collecting tube (dropper cap attached), swab.
B) Sampling
- Use the swab included in the kit to collect the nasopharynx or nasal samples.
- The sample should be tested preferably after collection.


Nasopharynx Swab




- 1- Ask the patient to tilt his head slightly.
- 2- Carefully insert the swab into the patient's nostril until the swab head meets resistance against the nasopharynx surface.
- 3- Rotate the swab 5 times, in contact with the posterior nasopharynx, remove the swab.
- 4- Repeat steps 1, 2, and 3 in the second nostril.

Nasal Swab

- 1- Ask the patient to blow his nose.
- 2- Insert approximately 2 cm from the tip of the swab into one nostril.
- 3- Using circular movements, scrape the swab head 5 times into the patient's nostril.
- 4- Repeat steps 2 and 3 for the other nostril.

- C) Ensure that the collection tube's buffer is at the bottom of the tube and remove the dropper cap.
D) After collection, insert the swab into the collection tube.
E) Mix for 10 seconds to ensure good extraction through the buffer.
F) Remove excess buffer at the tip of the swab by pressing the walls of the collection tube against the head of the swab.
G) Insert the dropper cap into the collection tube.
H) Add three drops of the buffer (90-100 µL) to the sample well. Read the test between 15 and 20 minutes. Do not interpret after 20 minutes.

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- 1 - Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
- 2 - CENTERS FOR DISEASE CONTROL AND PREVENTION et al. Interim guidelines for collecting, handling, and testing clinical specimens for COVID-19. Retrieved June, v. 19, p. 2020, 2020. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13. A.Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
- 3 - Corman VM, Landt O, et al. Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
- 4 - Huang C, Wang Y, et al. Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb15;395(10223):497-506.
- 5 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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
INTERFERENT

For the assessment of interferents, substances were added to 2 samples of known results. No divergent results were found up to concentrations of: acetylsalicylic acid 65 mg/dL, bilirubin 30 mg/dL, hemoglobin 500 mg/dL, triglycerides 1000 mg/dL and mucin 100 µg/ mL and rheumatoid factor 1000 IU/mL.

PRODUCT PERFORMANCE

Accuracy

METHODS COMPARISON, METHODOLOGICAL SPECIFICITY, AND SENSITIVITY
To establish the sensitivity and specificity of the kit BIOCLIN FAST COVID-19 Ag were analyzed 292 known clinical samples. Clinical sensitivity and clinical specificity found were:

Sensitivity: (135/142) = 95,1% (IC 95%: 91,6% – 98,6%)
Specificity: (146/150) = 97,3% (IC 95%: 94,8% – 99,9%)

To establish the sensitivity and specificity for the BIOCLIN FAST COVID-19 Ag kit for nasal swab samples, 100 known clinical samples were analyzed. The clinical sensitivity and clinical specificity found for nasal swab samples were:

Sensitivity: (48/50) = 96,0% (IC 95%: 90,67% – 100%)
Specificity: (49/50) = 98,0% (IC 95%: 94,2% – 100%)

Precision

REPEATABILITY

The repeatability was calculated from 10 successive determinations, using 2 different samples, obtaining the following results:

Samples	Repetition Nº	Expected Results	Results Found
01	10	100% Reagent	100% Reagent
02	10	100% Non-Reagent	100% Non-Reagent

REPRODUCIBILITY

The reproducibility was calculated from 5 successive determinations for 3 consecutive days, using 2 different samples, obtaining the following results:

Samples	Nº fo days	Repetition Nº	Expected Results	Results Found
01	03	05	100% Reagent	100% Reagent
02	03	05	100% Non-Reagent	100% Non-Reagent

ANALYTICAL SENSITIVITY

Sensitivity 2 x 10^{3.1} TCID50/mL (Tissue Culture Infectious Dose 50%/mL).

DIAGNOSTIC SIGNIFICANCE

Coronaviruses (CoV) belong to the Coronaviridae family and are widely distributed infecting humans and other mammals. In humans, the most common symptoms presented are fever, cough, dyspnoea, myalgia, or fatigue. It can cause diseases ranging from the common cold to more serious illnesses, such as the Middle East Respiratory Syndrome (MERSCoV) and the Severe Acute Respiratory Syndrome (SARS-CoV). In December 2019, a series of pneumonia cases of unknown cause appeared in Wuhan in China, with clinical symptoms very similar to viral pneumonia, after complete sequencing of the respiratory samples it was shown to be Coronavirus infection, which was called the new coronavirus (2019-nCoV). In approximately 150 countries on all continents, cases of infection have already been reported.

BIBLIOGRAPHIC REFERENCES

- 1 - Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
- 2 - CENTERS FOR DISEASE CONTROL AND PREVENTION et al. Interim guidelines for collecting, handling, and testing clinical specimens for COVID-19. Retrieved June, v. 19, p. 2020, 2020. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation, and storage of specimens for molecular methods; approved guideline. CLSI document MM13. A.Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
- 3 - Corman VM, Landt O, et al. Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
- 4 - Huang C, Wang Y, et al. Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb15;395(10223):497-506.
- 5 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Department of Quality Control. The quality of reagents is assured until the expiration date stated on the presentation packaging when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte/MG - Brasil
Tel.: +55 31 3439 5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Indústria Brasileira

CUSTOMER SERVICE

Customer Advisory Service
Phone: +55 31 3439 5454
E-mail: customerservice@bioclin.com.br

ANVISA registration for Bioclin Fast COVID-19 Ag kit: 10269360335

Review: June/2021

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR 40- TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER