BIOLISA SÍFILIS Ac TOTAL

REF **K240**

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FINALIDADE

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2. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and sérum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.

3. Avelleira JCR, Bottino G. Sífilis: diagnóstico, tratamento e controle. An Bras Dermatol. 2006;81(2):111-26.

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7. HORVATH, A. Biology and natural history of syphilis. In: GROSS, G. Sexually transmitted infections and sexually transmitted diseases. [S.I.]: Springer, 2011. p.

8. TELELAB. Manual de Coleta de Sangue - Diagnóstico e monitoramento das DST, Aids e Hepatites Virais.

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REF K240

USE INSTRUCTIONS

FUNCTION
Test for qualitative determination of total antibodies (IgG, IgM and IgA) for *Treponema pallidum*, bacteria that causes Syphilis, in biological samples (serum, plasma or whole blood on filter paper) through enzyme immunoassay. Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION
Methodology: Enzyme Immunoassay or Immunoenzymatic
BIOLISA SÍFILIS Ac TOTAL is an immunoenzymatic solid-phase assay based on the "sandwich" principle for the detection of total antibodies (IgG, IgM and IgA) to *Treponema pallidum* in human samples of serum, plasma or whole blood on filter paper. IgG, IgM and IgA antibodies present in the sample bind to recombinant *Treponema pallidum* antigens coated on the microplate forming immunocomplexes. After the initial incubation, the microplate is washed to remove unbound materials. Recombinant antigens conjugated to Peroxidase are added to the microplate, which is then incubated. The enzyme-conjugated antigens bind to the IgG, IgM and IgA anti-*Treponema pallidum* antibodies present, attached to the antigen-coated plate. A new wash is performed to remove the surplus. After this stage, the Substrate is added and incubated, producing a blue color that indicates the amount of IgG, IgM and IgA anti-*Treponema pallidum* antibodies present in the sample. The Stop Solution is added to stop the reaction with a color change from blue to yellow, measured in a microplate reader.

- REAGENTS**
1 - Sensitized Plate - Store between 2 and 8 °C. Contains: Sensitized Plate with recombinant *Treponema pallidum* antigens and preservative.
2 - Conjugate - Store between 2 and 8 °C. Contains: Buffer Solution containing recombinant *Treponema pallidum* antigen linked to Peroxidase, surfactant, stabilizers, dye and preservative.
3 - Concentrated Washing - Store at 2 to 8 °C. Contains: Buffer Solution, surfactant and preservative.
4 - Sample Diluent - Store between 2 and 8 °C. Contains: Buffer Solution, stabilizers, surfactant and preservative.
5 - Substrate - Store between 2 and 8 °C. Contains: Buffer Solution containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.
6 - Stop Solution - Store between 2 and 8 °C. Contains: Hydrochloric Acid 1 M.
7 - Negative Control - Store between 2 and 8 °C. Contains: Buffer Solution, stabilizer, surfactant and preservative. **Potentially infective.**
8 - Positive Control - Store between 2 and 8 °C. Contains: Buffer Solution, IgG, IgM and IgA anti-*Treponema pallidum* antibodies, dye, stabilizers, surfactant and preservative. **Potentially infective.**
9 - Plate Sealers.
10 - Filter Paper Elution Buffer - Store between 2 and 8 °C. Contains: Buffer Solution, surfactant, stabilizer and preservative.
11 - Negative Extraction Control - Store between 2 and 8 °C. Contains: Non-reactive sample for IgG, IgM and IgA anti-*Treponema pallidum* antibodies impregnated in Filter Paper. **Potentially infective.**
12 - Positive Extraction Control - Store between 2 and 8 °C. Contains: Reactive sample for IgG, IgM and IgA anti-*Treponema pallidum* antibodies impregnated in Filter Paper. **Potentially infective.**

REAGENTS	1	2	3
	96 Cavities	192 Cavities	480 Cavities
1 - Sensitized Plate	1 Unit (96 cavities)	2 Units (192 cavities)	5 Units (480 cavities)
2 - Conjugate	1 Flask x 12 mL	2 Flasks x 12 mL	5 Flasks x 12 mL
3 - Concentrate Washing	1 Flask x 50 mL	2 Flasks x 50 mL	5 Flasks x 50 mL
4 - Sample Diluent	1 Flask x 42 mL	2 Flasks x 42 mL	5 Flasks x 42 mL
5 - Substrate	1 Flask x 12mL	2 Flasks x 12mL	5 Flasks x 12mL
9 - Stop Solution	1 Flask x 12mL	2 Flasks x 12mL	5 Flasks x 12mL
7 - Negative Control	1 Flask x 300 µL	2 Flasks x 300 µL	5 Flasks x 300 µL
8 - Positive Control	1 Flask x 300 µL	2 Flasks x 300 µL	5 Flasks x 300 µL
9 - Plate Sealers	3 Units	6 Units	15 Units
10 - Filter Paper Elution Buffer	1 Flask x 12 mL	2 Flasks x 12 mL	5 Flasks x 12 mL
11 - Negative Extraction Control	1 unit	2 units	5 units
12 - Positive Extraction Control	1 unit	2 units	5 units

EQUIPMENTS AND OPERATIONAL INPUTS
Materials in the kit:
- Reagents described in the above table
- Operating instructions (manual)

- Required materials not contained in the kit:**
1 - Pipette capable of dispensing volumes of 5 to 100 µL with lower coefficient of variation than 1.5%.
2 - Re-pipettor for repetitive pipetting volumes of 100 µL, with lower coefficient of variation than 1.5% or multichannel pipette (optional).
3 - Microplate washer (optional).
4 - ELISA reader capable of absorbance at 450 and 630 nm wavelength.
5 - Paper towel to dry cavities
6 - Stopwatch or watch.
7 - Flask to store the washing solution after diluted.
8 - Distilled or deionized water.
9 - Tools of Quality Control.
10 - Incubator 37°C ± 2°C.
11 - Paper shredder (3 mm diameter) for filter paper technique.

TRANSPORTATION AND STORAGE CONDITIONS
The storage temperature should be 2 to 8 °C. Transport at temperatures up to 30 °C should not exceed 5 days. Keep away from light and avoid humidity. **Do not freeze.**

- SPECIAL CARE**
1 - For professional *in vitro* diagnostic use only.
2 - Strictly follow the methodology proposed to obtain accurate results.
3 - The sachet containing the microplate should be opened only after it reaches room temperature. Place the strip with unused cavities in the sachet, seal and store be 2 to 8 °C.
4 - The water used in material cleaning must be recent and free of contaminants.
5 - Deionized and saturated columns release alkaline water, several ions and oxidizing and reducing agents that can significantly alter the results.
6 - Stop Solution contains Hydrochloric Acid which is a strong acid. Handle it with care.
7 - All the raw material of product is tested and should be nonreactive for HBsAg, Anti-HIV 1&2 and Anti-HCV. However, these tests do not provide total assurance of the absence of infectious agents. The manipulation of any product containing human serum is potentially capable of transmitting diseases. Therefore, we must take due care in handling the bio safety of these products.
8 - Always add reagents in the same order to minimize the difference in reaction time between the cavities.
9 - As a safety measure, you should cover the plate during the reaction.

- 10 - You must ensure that the bottom of the cavity is clean and dry, and there are no bubbles on the surface fluid before reading the plate. Do not let the cavities run dry during the test.
11 - Do not expose reagents, especially the Substrate, to strong light or Hypochlorite fumes during storage or incubation steps.
12 - We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
13 - To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
14 - Do not use the product in case of damaged packaging.
15 - It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES
Serum or Plasma (EDTA or Heparin)
Hemolysed or highly lipemic samples should not be used. The samples can be kept refrigerated, between 2 and 8°C, for a maximum period of 5 days. If samples cannot be analyzed within 5 days, they can be stored for up to 30 days at -20 °C.

Whole Blood on Filter Paper
Whole Blood (EDTA or Puncture)
Samples dried on filter paper can be stored at room temperature as long as they are protected from direct sunlight and in low humidity. For storage of up to 2 years, the samples must remain refrigerated, between 2 and 8 °C. Storage for more than 2 years must be carried out the temperature of -20 °C.

PROCESS DESCRIPTION
The results of the stability test show that the Biolisa kit SÍFILIS Ac TOTAL is stable after opening for at least 30 days. This stability may vary according to the conditions of the test and the environment. Therefore, it is suggested to monitor the product's performance using internal kit controls and technique validation criteria.

Washing Solution
Dilute the contents of the Reagent N° 3 (Concentrate Washing) in 1000 mL of distilled or deionized water. After preparation the solution may be stored at 2 to 30 °C until expiration date printed on the original bottle. In case of crystallization, heat it at 37 °C until dissolution.

Substrate
The Substrate is ready for use.

TECHNIQUE FOR USE ON EQUIPMENT
For use in automatic equipment, consult the SAC (Customer Service Department).
Samples of Serum or Plasma
Before starting the assay, place all reagents, Controls and Samples to stabilize at room temperature (15 – 30 °C) for at least 40 minutes.
1 - Select the cavities to be used considering: Controls and Samples (it is recommended to test in duplicate). Return the strips of the microplate will not be used, for the original sealed packaging.
2 - Select the first cavity for Blank (OPTIONAL).
3 - Pipette 100 µL of Sample Diluent, included in the cavity for Blank.
4 - Pipette 5 µL of Sample and Controls into the previously determined wells. Observe the color change of the diluent when adding the sample. A Color change indicates that the sample was successfully added to the well.
5 - Homogenize gently for ± 30 seconds, cover the wells with plate sealer.
6 - Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
7 - Remove the sealing of wells.

- 8 - Discard the contents of the wells by aspiration (Washer) or by decanting (manual). Use approximately 300 µL of Washing Solution, previously prepared to perform a total of five (5) washing cycles. To ensure drying of the plate, at the end of the wash, hit the board for a few seconds on absorbent paper. Note: Poor washing/drying can cause inadequate results.
9 - Pipette 100 µL of Conjugate into each well.
10 - Homogenize gently for ± 30 seconds, cover the wells with plate sealer.
11 - Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
12 - Remove the sealer from plate cavities.
13 - Repeat item 8.
14 - Pipette 100 µL of Substrate into all wells.
15 - Homogenize gently for ± 30 seconds, cover the wells with plate sealer.

- 16 - Incubate for 10 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
17 - Remove the sealer from plate cavities.
18 - Pipette 100 µL of Stop Solution into all wells.
19 - Mix gently for ± 30 seconds.
20 - Read using double filter: 450nm / 630nm within up to 15 minutes (maximum).

Whole Blood Sample on Filter Paper
Before starting the assay, place all reagents, Controls and Samples to stabilize at room temperature (15 – 30 °C) for at least 40 minutes.
1 - Select the wells to be used considering Positive Extraction Control, Negative Extraction Control and Whole Blood Samples on Filter Paper (it is recommended to test in duplicate). Return the strips of the microplate will not be used for the original sealed packaging.
2 - Select the first cavity for Blank (OPTIONAL).
3 - Add a 3 mm disk of Positive Extraction Control, Negative Extraction Control and Whole Blood Sample on Filter Paper, previously perforated, in the determined wells.
4 - Pipette 100 µL of Filter Paper Elution Buffer into each well. Included in the Blank cavity.
5 - Pipette 5 µL of Controls into previously determined wells.
6 - Homogenize gently for ± 30 seconds, cover the wells with plate sealer.
7 - Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
8 - Remove the sealing of wells.
9 - Discard the contents of the cavities by aspiration (Washing machine) or by decanting (Manual). Remove the paper discs, if necessary, with the aid of a needle or tip. Use approximately 300 µL of Washing Solution, previously prepared to perform a total of five (5) washing cycles. To ensure drying of the plate, at the end of the wash, hit the board for a few seconds on absorbent paper. Note: Poor washing/drying can cause inadequate results.
10 - Pipette 100 µL of Conjugate into all wells.
11 - Homogenize gently for ± 30 seconds, cover the wells with plate sealer.
12 - Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
13 - Remove the sealer from plate cavities.
14 - Repeat item 9.
15 - Pipette 100 µL of Substrate into each well.
16 - Homogenize gently for ± 30 seconds, cover the wells with plate sealer.
17 - Incubate for 10 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
18 - Remove the sealer from plate cavities.
19 - Pipette 100 µL of Stop Solution into all wells.
20 - Mix gently for ± 30 seconds.
21 - Read using double filter: 450nm / 630nm within up to 15 minutes (maximum).

TECHNIQUE VERIFICATION
Verify if the results obtained by the reading of the Blank and Controls are compatible to the values below:

ITEM	ABSORBANCE
Blank	< 0.200
Negative Control	< 0.200
Positive Control	> 1.000
Negative Extraction Control	< 0.200
Positive Extraction Control	> 0.500

If the values are out the expected values, you must repeat the technique.

CALCULATIONS

QUALITATIVE

Samples of Serum or Plasma

Calculate Cut Off according to the following formula:

Cut Off = (Average Absorbance Positive Control x 0,01) + 0,100

Example:

ITEM	ABSORBANCE
Positive Control (Reagent N° 8)	A1 = 1.944
	A2 = 1.945
Cut Off = (Average Absorbance Positive Control x 0,01) + 0,100	Cut Off = ((1.944 + 1.945) / 2 x 0,01) + 0,100
	Cut Off = 0,119

Calculate the Index by dividing the absorbance of the Sample by the Cut Off value.
Example:

ITEM	ABSORBANCE
Sample	1,41
Cut Off Value	0,119
Index = Sample / Cut Off Value	Index = 1.41 / 0.119
	Index = 11,84

Note: The data presented in the examples are for illustration only and cannot be used to calculate results.

Whole Blood Sample on Filter Paper

Calculate Cut Off according to the following formula:

Cut Off = (Average Absorbance of Positive Control x 0,01) + 0,100

Example:

ITEM	ABSORBANCE
Positive Control (Reagent N° 8)	A1 = 1.944
	A2 = 1.945
Cut Off = (Average Absorbance of Positive Control x 0,01) + 0,100	Cut Off = ((1.944 + 1.945) / 2 x 0,01) + 0,100
	Cut Off = 0,119

Note: Calculate the Index by dividing the absorbance of the sample by the Cut Off value. Example:

ITEM	ABSORBANCE
Sample	0,965
Cut Off Value	0,119
Index = Sample / Cut Off Value	Index = 0,965 / 0,119
	Index = 8,11

Note: The data presented in the examples are for illustration only and cannot be used to calculate results.

INTERPRETATION OF RESULTS

Samples of Serum, Plasma and Filter Paper Whole Blood

After calculating the index of the samples, consider the indices below to determine the results.

RESULTS	QUALITATIVE FOR SERUM SAMPLE, PLASMA AND TOTAL BLOOD ON FILTER PAPER
	INDEX
Negative	< 0.90
Positive	> 1.1
Undetermined	0.9 - 1.1

Note: In the case of an undetermined result, the sample must be reanalyzed. Samples that obtain results repeatedly Indeterminate conditions must be retested using an alternative method. If results remain indeterminate, a new sample in two weeks. The result of the last sample must prevail collected. The results provided by this kit must be interpreted responsible medical professional, and it is not the only criterion for the determining the diagnosis and/or treatment of the patient.

PROCEDURE LIMITATIONS

The interpretation of a diagnostic test should not be established with basis of a single trial. Other confirmatory tests should be included, before a sample is considered positive. A negative result does not exclude the possibility of exposure. All results must be interpreted in conjunction with other available clinical information before the definitive diagnosis of the disease.

INTERFERENT

No interference was observed for Triglycerides 1500 mg/dL, Acetylsalicylic Acid 20 mg/dL, Ascorbic Acid 2 g/dL, Creatine 200 mg/dL, Bilirubin 1 g/dL, Albumin 2 g/dL, Hemoglobin 1000 mg/dL, Oxalic Acid 60 mg/dL, Rheumatoid Factor 980 IU/ mL, C Reactive Protein 41.2 mg/dL and Anti Streptolysin O 1023 IU/ml.

CROSS REACTIVITY

Samples of Serum or Plasma

A cross-reactivity study was carried out, evaluating 100 samples negative for syphilis, but positive for other infections. Among them, 10 positive samples for each of the following infections: Cytomegalovirus, HIV, Zika, HBV, HCV, Dengue, Chikungunya, Toxoplasmosis, Malaria and Leptospirosis. No cross reactivity was observed with any of the samples. Despite the results found, the possibility of cross-reactivity cannot be completely ruled out. The final diagnosis should consider the patient's clinical data together with other laboratory data.

Whole Blood Sample on Filter Paper

A cross-reactivity study was carried out, evaluating 82 samples of whole blood on filter paper, which were negative for syphilis, but positive for other infections. Among them, 9 were positive for Zika, 11 for Toxoplasmosis, 13 for Rubella, 8 for Chagas, 5 for HIV, 10 for HCV, 7 for Cytomegalovirus, 9 for HBV, 7 for Malaria and 3 for Leptospirosis. No cross reactivity was observed with any of the samples. Despite the results found, the possibility of cross-reactivity cannot be completely ruled out. The final diagnosis should consider the patient's clinical data along with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

PRUDUCT PERFORMANCE

Accuracy

REPETIBILITY

The repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

REPETIBILITY	SAMPLE		
	1	2	3
Average	0.678	1.082	0.055
Standard Deviation	0.098	0.087	0.011
Coefficient of Variation (%)	14.475	8.041	19.861

REPRODUCIBILITY

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

REPRODUCIBILITY	SAMPLE		
	1	2	3
Average	0.519	1.125	0.060
Standard Deviation	0.082	0.110	0.010
Coefficient of Variation (%)	11.849	9.659	17.089

Clinical Sensitivity and Specificity

Serum and Plasma Sample

BIOLISA SÍFILIS Ac TOTAL Kit analyzed clinical samples in comparison with other methods of EIA. The results show that the clinical sensitivity of the BIOLISA SÍFILIS Ac TOTAL kit is > 99,9% and clinical specificity is > 99,9%.

	Expected Result	BIOLISA SÍFILIS Ac TOTAL
Positive Sample	165	165
Negative Sample	101	101
Total Tested Sample	266	

Clinical Sensitivity: >99,9% (165/165) - IC 95%: 97,8% - 100%

Clinical Specificity: >99,9% (101/101) - IC 95%: 96,4% - 100%

Whole Blood Samples on Filter Paper

BIOLISA SÍFILIS Ac TOTAL Kit analyzed clinical samples in comparison with other methods of EIA. The results show that the clinical sensitivity of the BIOLISA SÍFILIS Ac TOTAL kit is 98,21% and clinical specificity is > 99,9%.

	Expected Result	BIOLISA SÍFILIS Ac TOTAL
Positive Sample	56	55
Negative Sample	51	51
Total Tested Sample	107	

Clinical Sensitivity: 98,21% (56/55) IC 95%: 90,4% - 100%

Clinical Specificity: >99,9% (51/51) IC 95%: 93,0% - 100%

DIAGNOSTIC SIGNIFICANCE

Syphilis is a systemic infectious disease, of chronic evolution, caused by *Treponema pallidum*, which evolves slowly in three stages, characterized by lesions of the skin and mucous membranes. It can be transmitted by sexual contact, thus becoming a sexually transmitted disease and, more rarely, by fetal-placental contamination (congenital syphilis). The first manifestations occur after an average incubation period of 21 days, ranging from 3 to 90 days. The disease has three distinct phases, with characteristic manifestations in each one and a latency period (without symptoms) between the second and third phases. In the primary phase, there is a classic lesion on the genitals, the gateway to Treponema, called chancre. Chancretakes an average of 3 to 6 weeks to heal, and may leave no marks. The secondary phase occurs after 4 to 8 weeks of the chance onset, and this lesion may even be present. At this stage, the largest amount of treponemas occurs in the body. The tertiary phase is the phase of progressive and slow (chronic) inflammation with symptoms related to Organs predominantly compromised organs, it is destructive and disabling. The diagnosis can be made: by methods of detecting non-specific antibodies (using non-treponemic antigens) whose interpretation is visual, and by immunoenzymatic methods (ELISA) that detect the presence of specific antibodies against *Treponema pallidum* in samples of patients who are at different stages of the disease.

BIBLIOGRAPHIC REFERENCES

1. Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
2. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and sérum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
3. Avelleira JCR, Bottino G. Sífilis: diagnóstico, tratamento e controle. An Bras Dermatol. 2006;81(2):111-26.
4. Testes para diagnósticos da Sífilis. Secretaria de Ciência, Tecnologia e Insumos Estratégicos. Testes para diagnósticos da Sífilis – Conitec, N° 159 Maio 2015.
5. Montiel, Milagros; Arias, Julia; Pozo, Elieth y Mogollón, Adrian. Importancia de las pruebas específicas e inespecíficas para el diagnóstico de sífilis en donantes de sangre. Rev. Kasmera 36(2): 169 - 176, 2008
6. BRASIL. Ministerio da Saude. Secretaria de Vigilancia em Saude. Guia de Vigilancia em Saude/Ministerio da Saude, Secretaria de Vigilancia em Saude. Brasília: Ministerio da Saude, 2014a.812 p. Disponivel em: www. saude.gov.br/bvs.
7. HORVATH, A. Biology and natural history of syphilis. In: GROSS, G. Sexually transmitted infections and sexually transmitted diseases. [S.l]: Springer, 2011. p.
8. TELELAB. Manual de Coleta de Sangue - Diagnóstico e monitoramento das DST, Aids e Hepatites Virais.
9. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br
ANVISA registration for BIOLISA SÍFILIS Ac TOTAL kit: 10269360343

Review: October/2021

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MANUFACTURED BY

BATCH CODE

CONTROL

DATE OF MANUFACTURE

POSITIVE CONTROL

USED BY
(last day of month)

NEGATIVE CONTROL

TEMPERATURE LIMITATION
(store at)

BIOLOGICAL RISK

CONTAINS SUFFICIENT
FOR <N> TESTS

INFLAMMABLE

CONSULT INSTRUCTIONS
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC DEVICE

POISON

EUROPEAN AUTHORIZED
REPRESENTATIVE

CE MARK

KEEP AWAY
FROM SUNLIGHT

DO NOT USE IF
PACKAGE IS
DAMAGED