BIOLISA COVID-19 ANTICORPO NEUTRALIZANTE

[REF] K243

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FINALIDADE

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Bioclin

BIOLISA COVID-19 ANTIBODY NEUTRALIZER

[REF] K243

USAGE INSTRUCTIONS

FUNCTION

Test for qualitative determination of the neutralization capacity of total antibodies anti-SARS-CoV-2 in biological samples of serum or plasma (EDTA or Heparin) by inhibiting the formation of the S antigen and ACE2 protein complex (S-ACE2) by competition test. For use only *in vitro* diagnosis.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzyme
The BIOLISA COVID-19 ANTIBODY NEUTRALIZER kit is a solid phase immunoassay based on the principle of qualitative detection by competition of total antibodies anti-SARS-CoV-2 capable of inhibiting the S-ACE2 complex in serum and human plasma samples.

In the initial incubation, Conjugate 1, composed of the protein ACE2 conjugated to Biotin, binds to coated Streptavidin on the microplate. After the initial incubation, the microplate is washed to remove unbound materials. In the second incubation, serum and plasma samples containing total antibodies anti-SARS-CoV-2 compete with the ACE2 protein for binding with the Spike Antigen conjugated to Peroxidase present in Conjugate 2, preventing the formation of the S-ACE2 complex.

A New wash is performed to remove the total antibodies anti-SARS-CoV-2, linked or not linked to Spike Antigen.

After this step, the substrate is added and incubated. An blue color for samples indicates the absence of Total Antibodies anti-SARS-CoV-2 and complex formation between the Antigen Spike conjugated to Peroxidase and ACE2. On the other hand, a lower intensity of the blue color indicates a greater amount of Total anti-SARS-CoV-2 antibodies present in the sample prevented the formation of the S-ACE2 complex. The Stop Solution is added to stop the reaction if there is a color change from blue to yellow, measured on a reader microplate.

REAGENTS

- Sensitized Plate** - Store between 2 and 8°C. Contains: Plate sensitized with Streptavidin and preservative.
- Conjugate 1** - Store between 2 and 8°C. Contains: Buffer Solution containing Biotin-bound recombinant ACE2 protein, surfactant, stabilizers, dye and preservative.
- Concentrated Washing** - Store between 2 and 8 °C. Contains: Buffer Solution (Phosphate <0.5 mol/L, Potassium Chloride <100 mmol/L, Sodium Chloride <5 mol/L), surfactant and preservative.
- Sample Diluent** - Store at 2 to 8°C. Contains: Buffer Solution, stabilizers, surfactant, dye and preservative.
- Conjugate 2** - Store between 2 and 8°C. Contains: Buffer Solution containing Spike Antigen linked to Peroxidase, surfactant, stabilizers, dye and preservative.
- Substrate** - Store between 2 and 8°C. Contains: Solution Buffer containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.
- Stop Solution** - Store at 2 to 8°C. Contains: Acid Hydrochloride 1 M.
- Negative Control** - Store at 2 to 8°C. Contains: Buffer, stabilizer, surfactant and preservative solution. **Potentially infective.**
- Positive Control** - Store between 2 and 8°C. Contains: Buffer Solution, Total Anti-SARS-CoV-2 Antibodies, Dye, stabilizers, surfactant and preservative. **Potentially infecting.**
- Plate Sealers**

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	1	2	3
Reagentes	96 Cavidades	192 Cavidades	480 Cavidades
1- Placa Sensibilizada	1 Unidade (96 cavidades)	2 Unidades (192 cavidades)	5 Unidades (480 cavidades)
2- Conjugado 1	1 Frasco x 6 mL	2 Frascos x 6 mL	5 Frascos x 6 mL
3- Lavagem Concentrada	1 Frasco x 50 mL	2 Frascos x 50 mL	5 Frascos x 50 mL
4- Diluyente de Amostra	1 Frasco x 15 mL	2 Frascos x 15 mL	5 Frascos x 15 mL
5- Conjugado 2	1 Frasco x 12 mL	2 Frascos x 2 mL	5 Frascos x 12 mL
6- Substrato	1 Frasco x 12 mL	2 Frascos x 12 mL	5 Frascos x 12 mL
7- Solução de Parada	1 Frasco x 12 mL	2 Frascos x 12 mL	5 Frascos x 12 mL
8- Controle Negativo	1 Frasco x 300 µL	2 Frascos x 300 µL	5 Frascos x 300 µL
9- Controle Positivo	1 Frasco x 300 µL	2 Frascos x 300 µL	5 Frascos x 300 µL
10- Seladores de Placa	3 Unidades	6 Unidades	15 Unidades

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PRESENTATION

	1	2	3
REAGENTS	96 Cavities	192 Cavities	480 Cavities
1- Sensitized Plate	1 Unit (96 cavities)	2 Units (192 cavities)	5 Units (480 cavities)
2- Conjugate 1	1 Flask x 6 mL	2 Flasks x 6 mL	5 Flasks x 6 mL
3- Concentrated Washing	1 Flask x 50 mL	2 Flasks x 50 mL	5 Flasks x 50 mL
4- Sample Diluent	1 Flask x 15 mL	2 Flasks x 15 mL	5 Flasks x 15 mL
5- Conjugate 2	1 Flask x 12 mL	2 Flasks x 12 mL	5 Flasks x 12 mL
6- Substrate	1 Flask x t 12 mL	2 Flasks x 12 mL	5 Flasks x 12 mL
7- Stop Solution	1 Flask x 12 mL	2 Flasks x 12 mL	5 Flasks x 12 mL
8- Negative Control	1 Flask x 300 µL	2 Flasks x 300 µL	5 Flasks x 300 µL
9- Positive Control	1 Flask x 300 µL	2 Flasks x 300 µL	5 Flasks x 300 µL
10- Plate Sealers	3 Units	6 Units	15 Units

EQUIPMENTS AND OPERATIONAL INUTS

Materiais in the kit:

- Reagents described in the above table
- Operating instructions (manual)

Required materials not contained in the kit:

- Pipette capable of dispensing volumes of 10 to 100 µL with lower coefficient of variation than 1.5%.
- Re-pipettor for repetitive pipetting volumes of 10 to 100 µL, with lower coefficient of variation than 1.5% or multichannel pipette (optional).
- Microplate washer (optional).
- ELISA reader capable of absorbance at 450 and 630 nm wavelength.
- Paper towel to dry cavities
- Stopwatch or watch.
- Flask to store the Washing Solution after dilution.
- Distilled or deionized water.
- Tools of Quality Control.
- Incubator 37°C ± 2°C.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature and transport should be 2 to 8°C. Keep away from light and avoid humidity. **Do not freeze.**

SPECIAL CARE

- **For professional *in vitro* diagnostic use only.**
- Strictly follow the methodology proposed to obtain accurate results.
- The sachet containing the microplate should be opened only after it reaches room temperature. Place the strip with unused cavities in the sachet, seal and store be 2 to 8°C.
- The water used in material cleaning must to be recent and free of contaminants.
- Deionized and saturated columns release alkaline water, several ions and oxidizing and reducing agents that can signify cantly after the results.
- Stop Solution contains Hydrochloric Acid which is a strong acid. Handle it with care.
- All the raw material of product is tested and should be nonreactive for HBsAg, Anti-HIV 1&2 and Anti-HCV. However, these tests do not provide total assurance of the absence of infectious agents. The manipulation of any product containing human serum is potentially capable of transmitting diseases. Therefore, we must take due care in handling the bio safety of these products.

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PROCEDURE LIMITATIONS

The kit is a qualitative assay with the possibility of correlation of the inhibition capacity of the reaction by the analysis in percentage. There are still no clinical data that establish a relationship between the ability to inhibit the reaction presented in a sample and the clinical status of the patient. The results obtained with the BIOLISA COVID-19 ANTIBODY NEUTRALIZING kit demonstrate only that patients with results < 30% do not have neutralizing antibodies, and samples > 35% do have these antibodies. Results <30% do not exclude SARS-CoV-2 infection. Results > 35% may be due to infection by other coronaviruses, although no cross-reactions were evidenced in studies performed during development. This product simulates the immunodominant inhibition mechanism based on viral neutralization by neutralizing antibodies, blocking the formation of the S-ACE2 immune complex. Patients with positive virus may develop other neutralizing or even non-neutralizing antibodies against SARS-CoV-2, both undetectable by the proposed methodology. All results should be interpreted in conjunction with other clinical information available prior to definitive diagnosis of the disease.

INTERFERENT

No interference was observed for Acetilsalicylic Acid 20 mg/dL, Ascorbic Acid 2 g/dL, Creatine 200 mg/dL, Bilirubin 1 g/dL, Albumin 2 g/dL, Hemoglobin 1000 mg/dL, Oxalic Acid 60 mg/dL, Rematoid Factor 980 IU/mL, C Reactive Protein 41.2 mg/dL and Anti Streptolysin O 1023 IU/ml.

CROSS REACTIVITY

A study was carried out with 140 negative samples for COVID-19, but positive for other infections. Amongst them, 10 positive samples for influenza, 27 positive samples for Rhinovirus, 13 positive samples for Sycytial Virus Respiratory, 10 Zika positive samples, 10 samples positive for Dengue, 10 positive samples for Chikungunya, 10 positive samples for Toxoplasmosis, 10 positive samples for Cytomegalovirus, 10 positive samples for Rubella, 10 HIV positive samples, 10 positive samples for HBV and 10 HCV positive samples. It was not observed cross-reactivity with positive samples for Zika, Dengue, Chikungunya, Toxoplasmosis, Cytomegalovirus, Rubella, HIV, HBV and HCV. Despite the results found, it is not possible to completely rule out the possibility of cross-reactivity. The final diagnosis should consider the patient's clinical data along with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

PRUDUCT PERFORMANCE

Accuracy

REPETIBILITY

The repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

REPETIBILITY	Samples		
	1	2	3
Average	0.262	0.084	1.306
Standard Deviation	0.024	0.011	0.083
Coefficient of Variation (%)	9.141	12.586	6.357

REPRODUCIBILITY

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

REPRODUCIBILITY	Samples		
	1	2	3
Average	0.262	0.079	1.206
Standard Deviation	0.025	0.007	0.078
Coefficient of Variation (%)	10.125	8.936	6.446

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Clinical Sensitivity and Specificity

To carry out the sensitivity and specificity study 156 clinically known samples were analyzed, 96 positive and 60 negative. The results can be observed in the table below:

	Expected Result	BIOLISA COVID-19 ANTIBODY NEUTRALIZER
Positive Sample	96	93
Negative Sample	60	57
Total Samples Tested		156

Clinical Sensitivity: 96.87% (93/96) - 95% CI: 88.9% - 100% Clinical Specificity: 95% (57/60) - 95% CI: 89.62% - 100%

A new clinical sensitivity study was carried out with 70 positive samples for the PRNT gold standard, which are considered positive by the method when they present titer ≥ 1/20.

A new clinical specificity study was performed with 100 pre-pandemic samples.

Sensitivity (68/70) = 97.1% 95% CI: 93.3 - 100% Specificity (92/100) = 92.0% 95%CI: 86.9 - 97.1%

The study of detection of neutralizing antibodies in different groups of samples was also carried out. Below are the new data:

Study by sample group:

Grup*	Total evaluated samples	Positive Neutralizing Ac
Group of infected individuals with symptoms > 14 days	37	28
Group of uninfected individuals vaccinated with CORONAVAC (> 30 days after 2nd dose)	73	64

*The data presented above do not represent the sensitivity of the product, but the expression of neutralizing antibodies in different situations per group evaluated in the tested panel.

DIAGNOSTIC SIGNIFICANCE

Coronaviruses (CoV) belong to the Coronaviridae family, and are widely distributed infecting humans and other mammals. In humans, the most common symptoms presented are: fever, cough, dyspnoea, myalgia or fatigue. It can cause diseases ranging from the cold common to more serious diseases, such as the Eastern Respiratory Syndrome (MERS-CoV) and Severe Acute Respiratory Syndrome SARS-CoV. In December 2019, a series of cases of pneumonia caused disease appeared in Wuhan in China, with clinical symptoms very similar to viral pneumonia, after complete sequencing of the respiratory samples was found to be Coronavirus infection, which it was called the new coronavirus (2019-nCoV). In approximately 150 countries on all continents have already been reported cases of infection.

BIBLIOGRAPHIC REFERENCES

- Cancer Epidemiol Biomarkers Prev; 26(8): 1337–44.2017 AACR
- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
- Corman VM, Landt O, et al. Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
- Huang C, Wang Y, et al. Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb 15;395(10223):497- 506.
- Huo, J., Le Bas, A., Ruza, R.R. et al. Neutralizing nanobodies bind SARS-CoV-2 spike RBD and block interaction with ACE2. Nat Struct Mol Biol 27, 846–854 (2020).
- Quibasa: Dados do Departamento de Pesquisa e Desenvolvimento.

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- Cancer Epidemiol Biomarkers Prev; 26(8): 1337–44.2017 AACR
- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and sérum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
- Corman VM, Landt O, et al. Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
- Huang C, Wang Y, et al. Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb 15;395(10223):497- 506.
- Huo, J., Le Bas, A., Ruza, R.R. et al. Neutralizing nanobodies bind SARS-CoV-2 spike RBD and block interaction with ACE2. Nat Struct Mol Biol 27, 846–854 (2020).
- Quibasa: Dados do Departamento de Pesquisa e Desenvolvimento.

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QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte/MG - Brasil
Tel.: +55 31 3439 5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Indústria Brasileira

CUSTOMER SERVICE

Customer Advisory Service
Phone: +55 31 3439 5454
E-mail: customerservice@bioclin.com.br

ANVISA registration for Bioclin Fast COVID-19 Antibody Neutralizer kit: 10269360337

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