Toxoplasma PCR

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Review: May/2021

INDEX

Function	3
Principle of Action	3
Presentation	3
Reagents	4
Equipaments and Operational Inputs	4
Transportation and Storage Conditions	4
Special Care	4
Samples	5
Process Description	6
A . DNAExtraction	6
B . Preparation of the Reagents	6
C . Dilution of Quantitative Standard	7
D . PCR Preparation	7
E . Thermocycler settings for Real-Time PCR	8
F . Result Validation	9
G. Result Interpretation	10
Process Limitations	11
Product Performance / Quality Control	11
Comparison of Methods and Methodology Specificity	11
Repeatability	11
Reproducibility	12
Clinical Sensitivity	12
Analytical Sensitivity	12
Linearity	12
Diagnostic Significance	12
Bibliographic References	13
Customer Service	13
Universal Symbolology	14

FUNCTION

Quantitative detection of *Toxoplasma gondii* DNA using Real-Time Polymerase Chain Reaction (PCR) technology. Designed for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

The **Bio Gene Toxoplasma PCR** test is a *in vitro* assay based on the quantitative detection of Toxoplasma DNA by Real-Time PCR.

The Real-Time PCR method is used to amplify the pathogen DNA.

A thermocycler for Real-Time PCR is used to amplify and detect the fluorescent probe. Software calculates the Toxoplasma DNA concentration, expressed in copies/μL, using the standard curve generated from the quantitative standard, contained in the kit.

Reagent	Presentation	
	50 Tests	150 Tests
R1	1 x 55 μL *	1 x 165 μL *
R2	1 x 525 μL *	3 x 525 μL *
R3	1 x 525 μL	3 x 525 μL
R4	1 x 55 μL *	1 x 165 μL *
R5	1 x 600 μL *	1 x 600 μL *
R6	1 x 50 μL	1 x 150 μL
R7	1 x 500 μL *	1 x 500 μL *
R8	2 x 1,5 mL	2 x 1,5 mL
R9	1 x 1 mL	1 x 2 mL

*Lyophilized reagents. The volumes described above refer to the final volume after the reagents resuspension, as described in the item PROCESS DESCRIPTION, subitem B (Preparation of the Reagents).

REAGENTS

- R1. PCR Solution:** Primer, probe, TRIS-HCl.
R2. Mix Taq: Polymerase, dNTPs, MgCl₂.
R3. Mix Buffer: TRIS-HCl.
R4. IC PCR Solution: Primer, probe, TRIS-HCl.
R5. Internal Control: Plasmid, TRIS-HCl.
R6. Negative Control: TRIS-HCl.
R7. Standard A (2 x 10⁵ copies/μL): Plasmid, TRIS-HCl, EDTA.
R8. Diluent: TRIS-HCl, EDTA.
R9. Water: DNase/RNase free water.

EQUIPMENTS AND OPERATIONAL INPUTS**Materials in the kit:**

- Reagents described in the previous item.
- Usage instructions (manual).

Materials needed but not contained in the kit:

- 1- Real-time PCR thermocycler.
- 2- Laminar flow hood.
- 3- 1.5 ml microcentrifuge tubes.
- 4- PCR reaction tubes or plate.
- 5- Powder free latex gloves or similar material.
- 6- Microcentrifuge (3,000 - 12,000 rpm).
- 7- Vortex.
- 8- Micropipettes and sterile filter pipette tips (0.5-10μL, 10-100μL, 100-1000μL).
- 9- Nucleic acid extraction kit.

TRANSPORTATION AND STORAGE CONDITIONS

The required storage temperature is -20°C (-10 to -30°C). **After resuspension of lyophilized reagents, the product is stable for 6 months from the date of resuspension.** Avoid freezing and thawing.

The kit may be transported between 2 and 30°C for 7 days. Protect from light and avoid moisture.

SPECIAL CARE

- 1- **For professional *in vitro* diagnostic use only.**
- 2- Strictly follow the methodology proposed to obtain accurate results.
- 3- Handle and dispose of all biological samples, reagents and materials as it contain infectious agents. Avoid direct contact with biological samples

and reagents. Avoid spills and aerosol. Handle and dispose of wastes in accordance with appropriate security procedures.

4- It is require skilled personnel for the molecular biology procedures in order to minimize the risk of erroneous results, degradation of nucleic acids contained in the samples or even sample contamination by amplicons.

5- It is necessary to have separate areas for the extraction / preparation of amplification reactions and for the amplification / detection of amplified products. Never introduce an amplification product in the area intended for the extraction or preparation of amplification products.

6- All samples and reagents should be handled inside a laminar flow hood. It is necessary to use sterile and DNase/RNase free filter pipettes tips.

7- Avoid repeated thawing and freezing of the reagents.

8- Store the DNA samples at -20°C in case they will not be used immediately.

9- Do not use reagents after expiration date.

10- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

11- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

12- Do not use the product in case of damaged packaging.

13- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

This kit should be used with DNA samples extracted from whole blood, amniotic fluid and cerebrospinal fluid. Samples should be collected according to laboratory recommendations for molecular testing. They must be transported and stored as described below¹:

Amniotic fluid and cerebrospinal fluid: between 2 to 8° C for up to 1 day. For long-term storage, it is recommended to freeze samples at – 70° C.

Whole blood: between 2 to 8° C for up to 3 days.

DNA (extracted): Use DNA samples suitable for PCR amplification with adequate purity and concentration. For long-term storage, it is recommended to freeze samples at – 70° C. Avoid repeated freezing and thawing of the samples.

PROCESS DESCRIPTION

A. DNA Extraction

The nucleic acids (DNA) must be extracted in accord to the procedure of the chosen extraction kit.

For control of the extraction process, the **Internal Control (R5)** must be prepared (see item B,3) and added to the samples during the extraction as described below:

- 1- Add 4µL of the **Internal Control (R5)** to each tube containing the resuspended samples in extraction/lysis buffer.
- 2- Complete the extraction procedure according to the usage instructions of the extraction kit.

Note: Never add the Internal Control directly to pure biological sample. This may result in degradation of the sample.

B. Preparation of the Reagents

1- Before open, spin (*spin* pulse) each of the following reagents: **PCR Solution (R1)**, **IC PCR Solution (R4)**, **Internal Control (R5)** and **Standard A (R7)**.

2- Resuspend each vial of **Mix Taq (R2)*** reagent with 525µL of **Mix Buffer (R3)** reagent.

*For the presentation of 150 tests resuspend the **Mix Taq (R2)** vials as need of use.

*The **Mix Taq (R2)** reagent does not contain passive reference Dye (ROX).

3- Resuspend the following reagents: **PCR Solution (R1)**, **IC PCR Solution (R4)** and **Internal Control (R5)** with **Water (R9)** reagent, in according to the following table:

Reagent	Presentation	
	50 Tests	150 Tests
PCR Solution (R1)	55 µL	165 µL
IC PCR Solution (R4)	55 µL	165 µL
*Internal Control (R5)	600 µL	600 µL

*R5 and R7 reagents contain RNA/DNA template. They need to be handle (resuspended) in a separated area from the area designed to prepare other reagents in order to prevent contamination of other reagents.

C. Dilution of Quantitative Standard

- 1- Resuspend **Standard A (R7)** with 500 μL of **Diluent (R8)**.
- 2- Select 3 microtubes (not supplied in the kit) suitable for serial dilution of **Standard A (R7)**.
- 3- Pipette 90 μL of the **Diluent (R8)** in each microtube and name them as B, C and D respectively.
- 4- Pipette 10 μL of **Standard A (R7)** in microtube B and mix.
- 5- Replace the tip and pipette 10 μL of microtube B in microtube C and mix.
- 6- Replace the tip and pipette 10 μL of microtube C in microtube D and mix.
- 7- At the end of dilution procedure, it will be available 4 standards with the following concentrations:

Standard A - 2×10^5 copies/ μL

Standard B - 2×10^4 copies/ μL

Standard C - 2×10^3 copies/ μL

Standard D - 2×10^2 copies/ μL

D. PCR Preparation

- 1- Prior beginning the assay, select sufficient number of microtubes/wells to be used on the reaction, in accord to the number of samples, controls and quantitative standards to be used.
- 2- Prepare the final PCR solution in accord with reactions number required.

Reagents	1 Reaction	25 Reactions	50 Reactions	100 Reactions
Mix Taq (R2)	10 μL	250 μL	500 μL	1 mL
PCR Solution (R1)	1 μL	25 μL	50 μL	100 μL
CI PCR Solution (R4)	1 μL	25 μL	50 μL	100 μL
Water (R9)	3 μL	75 μL	150 μL	300 μL

To prepare different reaction number must be multiplied the volume of reagents for a reaction by the number of required reactions.

- 3- Pipette 15 μL of the final PCR solution in microtubes or wells previously determined.
- 4- Add 5 μL of extracted DNA sample or 5 μL of Quantitative Standards or 5 μL of **Negative Control (R6)** in previously determined tubes.
- 5- Mix well.

- 6-** The total volume of the reaction is 20 μ L. Each PCR run should include the relevant controls: Negative Control, Internal Control and Quantitative Standards.
- 7-** Transport the tubes/plate to the thermocycler and follow the instructions in the section E (Thermocycler Settings for Real-Time PCR).

E. Thermocycler Settings for Real-Time PCR

Check the operation manual of the Real-time PCR thermocycler to program the experiment.

1- Set the type of experiment:

Quantitative Test with Standard Curve.

2- Set the fluorescent probes:

Target	Detector	Quencher
Toxoplasma	FAM	NFQ-MGB
Internal Control	VIC	

Note: The Negative Control and Quantitative Standards do not present the Internal Control, because it is used for extraction and amplification control of the samples.

3- Set the quantitative standards:

Standard A - 2×10^5 copies/ μ L

Standard B - 2×10^4 copies/ μ L

Standard C - 2×10^3 copies/ μ L

Standard D - 2×10^2 copies/ μ L

4- Set the cycling conditions:

Stage	Temperature	Time	Cycles
1	95°C	2 minutos	1
2	95°C	10 segundos	50
	60°C	60 segundos	

Set "Data Collection" as "stage 2, step 2 (60@0:60)".

F. Result Validation

1- Standard Curve

Standard Curve	Permitted Range	Amplification / Detection
Correlation Coefficient (R^2)	$0,99 \leq R^2 \leq 1,00$	Valid

If the R^2 value does not fall within the limits of range, the result must be considered invalid and the test must be repeated.

2- Negative Control

CT Negative Control		Result	Amplification/ Detection
FAM	VIC		
Undetermined	Undetermined	Negative	Valid

3- Samples

Toxoplasma		Result	Detection
FAM	VIC		
Determined concentration	$25 \leq CT \leq 31$	Positive	Valid
	$CT < 25$ ou $CT > 31$	Positive	Invalid*
Undetermined concentration	$25 \leq CT \leq 31$	Negative	Valid
	$CT < 25$ ou $CT > 31$	Negative	Invalid*

*The CT values of Internal Control vary with the process conditions, as the extraction efficiency of the DNA/RNA, samples concentration and the thermocycler settings. Therefore, these conditions must be evaluated when the CT values are not appropriated and, in case it is relevant, the results can be validated.

Example: In some cases, the samples with high numbers of copies of DNA / RNA can inhibit the amplification of the Internal Control, resulting in CT value outside the optimal range, this does not invalidate the test results.

If the above requirements are not achieved, the assay must be considered invalid and must be repeated.

G. Result Interpretation

This kit is able to detect from 10 to 1.000.000 copies of DNA per reaction.

The thermocycler software calculates automatically the samples concentration.

Example: If the program shows a concentration of 2.00E+005, this means that the sample concentration is 2.0×10^5 copies/ μ L.

Sample Result in Copies/ μ L (FAM)	Copies per Reaction
$\geq 1 \times 10^6$	$> 1.000.000$
$2 \leq \text{Quantity} \leq 9 \times 10^5$	Quantity obtained
< 2	< 10

The non-detection of the pathogen DNA does not exclude the presence of infection when the pathogen title is below the detection limit of this kit.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only parameter for the determination of diagnosis and/or treatment of the patient.

The results obtained must be evaluated considering the clinical data and laboratory tests of the patient.

1- Conversion calculation for copies/mL

$$\text{Result in copies/mL} = \frac{\text{Copies}/\mu\text{L} \times \text{Elution volume } (\mu\text{L})^{**}}{\text{Extracted sample volume (mL)}^{**}}$$

*Result of sample quantification provided by the equipment in copies/ μ L.

**According to the protocol of the extraction kit used.

Example: Elution volume: 30µL
Sample volume: 200µL (0.2 mL)
Quantification of the sample: 5000 copies/µL

$$\text{Result in copies/mL} = \frac{(5000 \text{ copies/}\mu\text{L} \times 30\mu\text{L})}{0.2 \text{ mL}}$$

$$\text{Result in copies/mL} = 750,000 \text{ copies/mL}$$

PROCESS LIMITATIONS

Cross contamination that eventually occurs during the sample collection, processing, transportation and storage may give false results.

The heparin use in the samples collection for molecular tests is not recommended, due to the PCR inhibition which can generate false negative results.

The samples with presence of PCR inhibitors should be diluted to reduce their effect on the reaction.

PRODUCT PERFORMANCE QUALITY CONTROL

Accuracy

COMPARISON OF METHODS AND METHODOLOGY SPECIFICITY

The **Bio Gene Toxoplasma PCR** test was compared with another protocol for the quantification of the Toxoplasma DNA. 50 analyzes were performed and the results were evaluated. Comparing the two tests, the difference between the results obtained was ≤ 0.5 log. This show that the test has good methodological specificity.

Precision

REPEATABILITY

The repeatability was calculated from 10 successive determinations, using 2 positive samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2
Average Concentration (Log)	2.412	3.411
Standard Deviation (Log)	0.043	0.070
Coefficient of variation (%)	1.804	2.056

REPRODUCIBILITY

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 1 sample with different concentrations, obtaining the following results:

	Day 1	Day 2	Day 3
Average Concentration (Log)	3.417	3.377	3.317
Standard Deviation (Log)	0.064	0.069	0.068
Coefficient of variation (%)	1.888	2.057	2.069

Clinical Sensitivity

The **Bio Gene Toxoplasma PCR** showed clinical sensitivity of 99.9% and a clinical specificity of 99.9%.

Method		Clinical Data		Total
Bio Gene Toxoplasma PCR	Result	Positive	Negative	
	Positive	31	0	31
	Negative	0	25	25
TOTAL		31	25	56

Analytical Sensitivity

It was found in the performance study that the analytical sensitivity of the kit is 0.41 copies/ μ L (- 0.387 log). This means the assay is capable to detect approximately 2 target molecules in 5 μ L DNA extraction product, used on the amplification reaction.

Linearity

The linearity found in the performance studies was 5.230 log (2.03×10^5 copies/ μ L).

DIAGNOSTIC SIGNIFICANCE

Toxoplasma gondii is the etiological agent of toxoplasmosis can be transmitted through organ transplantation, blood transfusion, contact with feces from infected cats and ingestion of contaminated food. In adults the infection is usually asymptomatic, but there are cases in immunosuppressed with damage to the central nervous system. In pregnant, there is a risk

of fetal infection throughout pregnancy in the first trimester the risk of infection is lower, however, when the outcomes are serious occurs, such as miscarriage and hydrocephalus. Infections acquired later in pregnancy occur more often and tend to be less severe, but can lead to congenital disorders, such as cerebral calcifications and impaired learning. Early diagnosis in pregnant is essential to avoid or minimize consequences to the neonate. Therefore, the quantitative detection of toxoplasmosis by real time PCR is very important to be a sensitive method, precise and responsive.

BIBLIOGRAPHIC REFERENCES

1. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
2. Shen, L.J; Li, W; He, Y.H. Molecular diagnosis of toxoplasmosis. Chinese Journal of Schistosomiasis Control. 2011 Aug;23(4):465-7, 469. Review. Chinese.
3. Sterkers, Y; Varlet-Marie, E; Marty, P; Bastien, P; ANOFEL Toxoplasma-PCR Quality Control Group. Diversity and evolution of methods and practices for the molecular diagnosis of congenital toxoplasmosis in France: a 4-year survey. Clin Microbiol Infect. 2010 Oct;16(10):1594-602.
4. Wallon, M; Franck, J; Thulliez, P. et al. Accuracy of real-time polymerase chain reaction for *Toxoplasma gondii* in amniotic fluid. Obstet Gynecol. 2010 Apr;115(4):727-33.

CUSTOMER SERVICE

Customer Advisory Service

Phone: 0800 031 5454

E-mail: sac@bioclin.com.br

ANVISA registration for Bio Gene Toxoplasma PCR: 10269360314

UNIVERSAL SYMBOLOGY

	CUSTOM NUMBER		MANUFACTURED BY
	WATER CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	LOAD BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (temp °C)		BIOHAZARD RISK
	CHEMICAL SUPPRESSION FOR qPCR TESTS		FLAMMABLE
	COMPLETE INSTRUCTIONS FOR USE		CORROSIVE
	NO OTHER INFORMATION ENTRY		POISON
	NO OPEN FLAME SOURCES NEARBY DURING USE		CE MARK
	NO OPEN FLAME SOURCES NEARBY DURING USE		NO OPEN FLAME SOURCES NEARBY DURING USE

Bioclin · QUIBASA

Rua Teles de Menezes . 92 . Belo Horizonte . MG . Brasil . CEP: 31565-130

Tel +55 31 3439 5454 . www.bioclin.com.br

FARM. RESP. Sílvia Wandalsen Arndt - CRF MG 7422

C.N.P.J.: 19.400.787/0001-07 - Indústria Brasileira