FINALIDADE

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- Kaltwasser JP, Werner E. Diagnosis and clinical evaluation of iron overload. Baillieres Clin Haematol 1989; 2; 363-89.
- Wick M, Pingerra W, Lehmann P. Iron metabolism: diagnosis and therapy of anemias. 3rd ed. Vienna, New York: Springer Verlag, 1996.
- Lee MH, Means RT Jr. Extremely elevated serum ferritin levels in a university hospital: associated diseases and clinical significance. Am J Med 1996; 98:566-71.
- Worwood M. The laboratory assessment of iron status – an update. Clin Chim Acta 1997; 259:3-23.
- Baynes RD, Cook JD. Current issues in iron deficiency. Curr Opin Hematol 1996; 3:145-9.
- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento

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REFERENCIAS BIBLIOGRÁFICA
1. Kaltwasser JP, Werner E. Diagnosis and clinical evaluation of iron overload. Baillieres Clin Haematol 1989; 2; 363-89.
2. Wick M, Pingerra W, Lehmann P. Iron metabolism: diagnosis and therapy of anemias. 3rd ed. Vienna, New York: Springer Verlag, 1996.
3. Lee MH, Means RT Jr. Extremely elevated serum ferritin levels in a university hospital: associated diseases and clinical significance. Am J Med 1996; 98:566-71.
4. Worwood M. The laboratory assessment of iron status – an update. Clin Chim Acta 1997; 259:3-23.
5. Baynes RD, Cook JD. Current issues in iron deficiency. Curr Opin Hematol 1996; 3:145-9.
6. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
7. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento

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USAGE INSTRUCTIONS

FUNCTION

Method for quantitative determination of Ferritin, in photometric systems in samples of serum. Immunoturbidimetric test, for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Immunoturbimetry

Determination of ferritin concentration's ending point through photometric measuring of antigen-antibody reaction, between latex particles marked with antibodies anti-ferritin and ferritin present in sample.

REAGENTS

Number 1 - Buffer - Store between 2 and 8°C. Contains: Buffer < 300 mmol/L and stabilizers.

Number 2 - Ferritin Latex - Store between 2 and 8°C. Contains: Latex particles coated with Anti-Ferritin rabbit Antibodies < 10 g/L, buffer < 300 mmol/L and stabilizers.

Number 3 - Calibrator - Store between 2 and 8°C. Contains: Ferritin, stabilizers and preservative. **Potentially infectious.**

Number 4 - Control - Store between 2 and 8°C. Contains: Ferritin, stabilizers and preservative. **Potentially infectious.**

Warning: Concentration of Ferritin control and calibrator may vary according to the lot - See bottle label.

PRESENTATION

Reagent	Presentation				
	1	2	3	4	5
N°1	18 mL	45 mL	2 x 45 mL	40 mL	2 x 40 mL
N°2	2 mL	5 mL	2 x 5 mL	10 mL	2 x 10 mL
N°3	0.5 mL	1.0 mL	2 x 1.0 mL	1.0 mL	2 x 1.0 mL
N°4	0.5 mL	1.0 mL	2 x 1.0 mL	1.0 mL	2 x 1.0 mL

EQUIPMENTS AND OPERATIONAL INPUTS

This kit must be used in automated equipments only. These itens are found at markets specialized on Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must to be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- Reagents N° 3 and 4 should be handled cautiously although negative for HBsAg, anti-HCV and anti-HIV 1/2 as they are subject to biological contamination.

6- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

7- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Advisory Service Customer) of Quibasa.

8- Do not use the product in case of damaged packaging.

9- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum. Ferritin is stable for 07 days if kept between 2 and 8°C or up to 1 year if kept at -20°C ⁶. Sample must be frozen only once. Do not use highly lipemic or hemolysed samples. Sample with fibrin should be centrifuged before the test.

PROCESS DESCRIPTION

The calibration stability of the Ferritin kit installed on refrigerated equipment is at least 07 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

TECHNIQUE

For reaction control and calibration, use only calibrators and controls from this kit.

This kit is for use in automated biochemical analyzers only. **Check the application sheet for this equipment on the website www.bioclin.com.br or through SAC 0800 031 5454.**

PROCEDURE LIMITATIONS

The photometric scale of some instruments can be exceeded by highly lipemic samples, containing levels elevated ferritin.

The performance of the product depends directly on the performance of the automated equipment used. Therefore, to ensure best results of the test is very important to check periodically the optical and pipetting system of the equipment.

INTERFERENT

No interference was observed by Ascorbic Acid up to 30 mg/dL, Bilirubin up to 40 mg/dL, Hemoglobin up to 5000 mg/dL and Rheumatoid Factor up to 750 IU/mL. Lipids with values ≥ 2.5 g/dL may interfere.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

Do not compare the obtained results by different methodologies.

TRACEABILITY

The kit's calibrator is traceable to the reference material NIBSC 94/572 (Ferritin Human Recombinant - WHO International Standard).

REFERENCE VALUES

The reference values in µg/L or ng/mL, for this method were obtained through the determination of Ferritin in healthy populations of male and female.

Men: 30 - 220 µg/L or ng/mL

Women: 20 - 110 µg/L or ng/mL

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

RECOVERY

The recovery analysis was performed with samples of 05 determinations. Accuracies were calculated and were found in good agreement with the reference values, obtaining a recovery between 94% and 107%.

COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY

The Ferritin kit was compared with other commercially available methods for measurement of ferritin. 07 analyzes were performed and the results were evaluated. The linear equation obtained was Y = 0,995X + 0,676, and correlation coefficient 0,999. With these results we can conclude the kit shows good methodological specificity.

Precision
REPEATABILITY
The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Concentration (µg/L)	60.65	110.85	220.00
Standard Deviation (µg/L)	0.75	0.88	0.86
Coefficient of Variation (%)	1.23	0.79	0.39

REPRODUCIBILITY
The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Concentration (µg/L)	60.77	110.90	220.02
Standard Deviation (µg/L)	0.16	0.05	0.03
Coefficient of Variation (%)	0.26	0.05	0.01

High Dose Hook Effect
Hook effect wasn't observed with a high concentration of Ferritin sample up to 9000 µg/mL.

Sensitivity
The sensitivity was calculated from 20 determinations of a sample free of Ferritin. The average found was of 5,090 µg/dL with a standard deviation of 0,0013 µg/dL. The sensitivity, which indicates the method detection limit, corresponds the average plus 3 times the standard deviation and is equal to 5,094 µg/L.

Linearity
This kit is linear between 15 µg/L and the highest point in the calibration curve. For higher values, dilute the sample with a solution of NaCl 9 g/L. Mutiply the results by the dilution factor.

SIGNIFICANCE DIAGNOSTIC
Ferritin is the most important iron storage protein and is found in all cells, especially those involved in the synthesis of ferric compounds, metabolism and it's reserve.
The level of ferritin is the most accurate indicator of the amount of iron stored in the body.
His clinic is very useful in differential diagnosis of microcytic and hypochromic anemia iron deficiency, anemia of other etiologies.
In these cases, the ferritin decreases before changes of serum iron levels and red blood cell morphological changes.
However, to be part of the group of acute phase proteins, ferritin is increased in response to infection, trauma and acute inflammation.
The elevation occurs in 24 to 48 hours, with a peak in the third day, and it stays for a few weeks, which may hinder its interpretation.
Its levels can rise in excess of iron in transfused patients and cancer, especially in leukemias and lymphomas and breast cancer, liver, lung, colon and prostate.
Elevates also in megaloblastic and hemolytic anemias and in liver injury, especially alcoholic beverages. About 25% of patients with chronic hepatitis have increased ferritin.

BIBLIOGRAPHIC REFERENCES

1. Kaltwasser JP, Werner E. Diagnosis and clinical evaluation of iron overload. Baillieres Clin Haematol 1989; 2; 363-89.
2. Wick M, Pingerra W, Lehmann P. Iron metabolism: diagnosis and therapy of anemias. 3rd ed. Vienna, New York: Springer Verlag, 1996.
3. Lee MH, Means RT Jr. Extremely elevated serum ferritin levels in a university hospital: associated diseases and clinical significance. Am J Med 1996; 98:566-71.
4. Worwood M. The laboratory assessment of iron status – an update. Clin Chim Acta 1997; 259:3-23.
5. Baynes RD, Cook JD. Current issues in iron deficiency. Curr Opin Hematol 1996; 3:145-9.
6. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
7. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento

QUALITY ASSURANCE
Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda**
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

EC

REP

OBELIS S.A.
Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE
Customer Advisory Service
Phone: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for Ferritin kit: 10269360135

Review: November/2021

UNIVERSAL SYMBOLOGY

 CATALOG NUMBER

 BATCH CODE

 DATE OF MANUFACTURE

 USED BY
(last day of month)

 TEMPERATURE LIMITATION
(store at)

 CONTAINS SUFFICIENT
FOR <N> TESTS

 CONSULT INSTRUCTIONS
FOR USE

IVD

IN VITRO DIAGNOSTIC DEVICE

EC

REP

EUROPEAN AUTHORIZED
REPRESENTATIVE

 KEEP AWAY
FROM SUNLIGHT

 MANUFACTURED BY

CONTROL

CONTROL

CONTROL +

POSITIVE CONTROL

CONTROL -

NEGATIVE CONTROL

 BIOLOGICAL RISK

 INFLAMMABLE

 CORROSIVE

 POISON

CE

CE MARK

 DO NOT USE IF
PACKAGE IS
DAMAGED