BIOLISA SARS-CoV-2 IgG

 K234-1

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FINALIDADE

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USE INSTRUCTIONS

FUNCTION

Test for qualitative determination of IgG antibodies to SARS- CoV-2 (virus that causes COVID-19 disease) in biological samples (serum or plasma) through enzyme immunoassay test. Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzymatic The BIOLISA SARS-CoV-2 IgG kit is an immunoenzymatic assay in phase based on the principle of indirect qualitative detection of antibodies IgG for SARS-CoV-2 in serum, plasma and whole blood on filter paper human. IgG antibodies present in the sample bind to antigens recombinants of SARS-CoV-2 coated on the microplate forming immunocomplexes. After the initial incubation, the microplate is washed to remove unbound materials. Anti-IgG antibodies conjugated to Peroxidase is added to the microplate which is then incubated. The enzyme-conjugated antibodies bind to IgG Anti- SARS-CoV-2 present, connected to the antigen-coated plate. New wash is performed to remove surplus. After this stage, the Substrate is added and incubated, producing a blue color that indicates the amount of IgG Anti-SARS-CoV-2 Antibodies present in the sample. The Solution Stop is added to stop the reaction if there is a change in color from blue to yellow, measured in a microplate reader.

REAGENTS

1 - Sensitized Plate - Store between 2 to 8°C. Contains: Sensitized Plate with recombinant SARS-CoV-2 antigen and preservative.

2 - Conjugate - Store between 2 to 8°C. Contains: Buffer Solution containing Peroxidase-linked human anti-IgG antibody, surfactant, stabilizers, dye and preservative.

3 - Concentrated Washing - Store between 2 to 8°C. Contains: Buffer Solution, surfactant and preservative.

4 - Sample Diluent - Store between 2 to 8°C. Contains: Buffer Solution, stabilizers, surfactant and preservative.

5 - Substrate - Store between 2 to 8°C. Contains: Buffer Solution containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.

6 - Stop Solution - Store at 2 to 8°C. Contains: Hydrochloric Acid 1 M.

7 - Negative Control - Store between 2 to 8°C. Contains: Buffer Solution, stabilizer, surfactant and preservative. **Potentially infective.**

8 - Positive Control - Store between 2 to 8°C. Contains: Buffer Solution, IgG Anti - SARS-CoV-2 antibodies, dye, stabilizers, surfactant and preservative. **Potentially infective.**

9 - Plate Sealers

PRESENTATION

REAGENTS	1	2	3
	96 Cavities	192 Cavities	480 Cavities
1 - Sensitized Plate	1 Unit (96 cavities)	2 Unit (192 cavities)	5 Unit (480 cavities)
2 - Conjugate	1 Flask x 12 mL	2 Flask x 12 mL	5 Flask x 12 mL
3 - Concentrate Washing	1 Flask x 50 mL	2 Flask x 50 mL	5 Flask x 50 mL
4 - Sample Diluent	1 Flask x 100 mL	2 Flask x 100 mL	5 Flask x 100 mL
5 - Substrate	1 Flask x 12 mL	2 Flask x 12 mL	5 Flask x 12 mL
6 - Stop Solution	1 Flask x 12 mL	2 Flask x 12 mL	5 Flask x 12 mL
7 - Negative Control	1 Flask x 300 µL	2 Flask x 300 µL	5 Flask x 300 µL
8 - Positive Control	1 Flask x 300 µL	2 Flask x 300 µL	5 Flask x 300 µL
9 - Plate Sealers	3 Unit	6 Unit	15 Unit

EQUIPMENTS AND OPERATIONAL INPUTS

Materials in the kit:

- Reagents described in the above table
- Operating instructions (manual)

Required materials not contained in the kit:

- 1- Pipette capable of dispensing volumes of 5 to 500 µL with lower coefficient of variation than 1.5%.
- 2- Re-pipettor for repetitive pipetting volumes of 500 µL, with lower coefficient of variation than 1.5% or multichannel pipette (optional).
- 3- Microplate washer (optional).
- 4- ELISA reader capable of absorbance at 450 and 630 nm wavelength.
- 5- Paper towel to dry cavities.
- 6- Stopwatch or watch.
- 7- Flask to store the washing solution after dilution.
- 8- Distilled or deionized water.
- 9- Tools of Quality Control.
- 10- Incubator 37°C ± 2°C.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature and transport should be 2 to 8°C. Keep away from light and avoid humidity. **Do not freeze.**

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain accurate results.
- 3- The sachet containing the microplate should be opened only after it reaches room temperature. Place the strip with unused cavities in the sachet, seal and store be 2 to 8°C.
- 4- The water used in material cleaning must to be recent and free of contaminants.
- 5- Deionized and saturated columns release alkaline water, several ions and oxidizing and reducing agents that can significantly alter the results.
- 6- Stop Solution contains Hydrochloric Acid which is a strong acid. Handle it with care.

7- All the raw material of product is tested and should be nonreactive for HBsAg, Anti-HIV 1&2 and Anti-HCV. However, these tests do not provide total assurance of the absence of infectious agents. The manipulation of any product containing human serum is potentially capable of transmitting diseases. Therefore, we must take due care in handling the bio safety of these products.

8- Always add reagents in the same order to minimize the difference in reaction time between the cavities.

9- As a safety measure, you should cover the plate during the reaction.

10- You must ensure that the bottom of the cavity is clean and dry, and there are no bubbles on the surface fluid before reading the plate. Do not let the cavities run dry during the test.

11- Do not expose reagents, especially the Substrate, to strong light or Hypochlorite fumes during storage or incubation steps.

12- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

13- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

14- Do not use the product in case of damaged packaging.

15- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or Plasma (EDTA or Heparin).

Hemolysed or highly lipemic samples should not be used. The samples can be kept refrigerated, between 2 and 8°C, for a maximum period of 5 days. If samples cannot be analyzed within 5 days, they can be stored for up to 30 days at -20°C.

PROCESS DESCRIPTION

PREPARATION OF WORKING REAGENTS

Washing Solution

Dilute the contents of the Reagent N° 3 (Concentrate Washing) in 1000 mL of distilled or deionized water. After preparation the solution may be stored at 2 to 30°C until expiration date printed on the original bottle. In case of crystallization, heat it at 37°C until dissolution.

Substrate

The Substrate is ready for use.

STABILITY AFTER OPEN

The results of the stability test show that the Biolisa kit SARS- CoV-2 IgG is stable after opening for at least 30 days. This stability may vary according to the conditions of the test and the environment. Therefore, it is suggested to monitor the product's performance using internal kit controls and technique validation criteria.

TECHNIQUE

For use in automatic equipment, consult the SAC (Customer Service Department).

Before starting the assay, place all reagents, Controls and Samples to stabilize at room temperature (15 – 30°C) for at least 40 minutes.

1- Select the cavities to be used considering: Controls and Samples (it is recommended to test in duplicate). Return the strips of the microplate will not be used, for the original sealed packaging.

2- Select the first cavity for Blank (OPTIONAL).

3- Prepare a 1:101 dilution of Samples and Controls in microtube by adding 5 µL plus 500 µL of sample diluent. Homogenize

Observe the color change of the diluent when adding the sample. A Color change indicates that the sample was successfully added to the microtube.

4- Pipette 100 µL of Diluted Sample and controls into the previously determined wells. In the Blank cavity, pipette 100µL of Sample Diluent

5- Homogenize gently for ± 30 seconds, cover the wells with plate sealer.

6- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.

7- Remove the sealing of wells.

8- Discard the contents of the wells by aspiration (Washer) or by decanting (Manual). Use approximately 300 µL of Washing Solution, previously prepared to perform a total of five (5) washing cycles, with shaking for 5 seconds. To ensure drying of the plate, at the end of the wash, hit the board for a few seconds on absorbent paper.

Note: Poor washing/drying can cause inadequate results.

9- Pipette 100 µL of Conjugate into all wells. Except in the Blank cavity.

10- Homogenize gently for ± 30 seconds, cover the wells with plate sealer.

11- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.

12- Remove the sealer from plate cavities.

13- Repeat item 8.

14- Pipette 100 µL of Substrate into all wells. Included in the cavity for Blank.

15- Homogenize gently for ± 30 seconds, cover the wells with plate sealer.

16- Incubate for 10 minutes ± 2 minutes in an incubator at 37°C ± 2°C.

17- Remove the sealer from plate cavities.

18- Pipette 100 µL of Stop Solution into all wells. Included in the cavity for Blank.

19- Mix gently for ± 30 seconds.

20- Read using double filter: 450nm / 630nm within up to 15 minutes (maximum).

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NÚMERO DE TESTES

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- Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
- Corman VM, Landt O, et al. Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
- Huang C, Wang Y, et al. Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb 15;395(10223):497-506.
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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1. Cancer Epidemiol Biomarkers Prev; 26(8): 1337–44.2017 AACR
2. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and sérum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
3. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
4. Corman VM, Landt O, et al. Detection of 2019 novel coronavirus (2019- nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
5. Huang C, Wang Y, et al. Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb 15;395(10223):497-506.
6. QUIBASA: Dados do Departamento de Pesquisa Desenvolvimento.

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TECHNIQUE VERIFICATION

Verify if the results obtained by the reading of the Blank and Controls are compatible to the values below:

ITEM	ABSORBANCE
Blank	< 0.200
Negative Control	< 0.200
Positive Control	> 1.200

If the values are out the expected values, you must repeat the technique.

CALCULATIONS QUALITATIVE

Calculate Cut Off according to the following formula:

Cut Off = (Average Absorbance Positive Control x 0.02) + 0.360
Example:

ITEM	ABSORBANCE
Positive Control (Reagent N° 8)	A1 = 1.990
	A2 = 2.020
Cut Off = (Average Absorbance Positive Control x 0.02) + 0.360	$((1.990 + 2.020) / 2) \times 0.02 + 0.360 = 0.400$

Calculate the Index by dividing the absorbance of the sample by the Cut Off value.

Example:

ITEM	ABSORBANCIA
Sample	1.456
Cut Off Value	0.400
Index = Sample / Cut Off Value	$1.456 / 0.400 = 3.064$

INTERPRETATION OF RESULTS

After calculating the index of the samples, consider the indices below to determine the results.

RESULTS	QUALITATIVE INDEX
Negative	< 0.90
Positive	> 1.1
Undetermined	0.9 – 1.1

Note: In the case of an undetermined result, the sample must be reanalyzed. Samples that obtain results repeatedly indeterminate conditions must be retested using an alternative method. If results remain indeterminate, a new sample in two weeks. The result of the last sample must prevail collected. For a more complete assessment of the diagnosis and clinical correlation, it is suggested that each sample be tested for IgG and IgM. The results provided by this kit must be interpreted responsible medical professional, and it is not the only criterion for the determining the diagnosis and / or treatment of the patient. Note: The data presented in the examples are for illustration only and cannot be used to calculate results.

PROCEDURE LIMITATIONS

The interpretation of a diagnostic test should not be established with basis of a single trial. Other confirmatory tests should be included, before a sample is considered positive. A negative result does not exclude the possibility of exposure. All results must be interpreted in conjunction with other available clinical information before the definitive diagnosis of the disease.

INTERFERENT

No interference was observed for Acetylsalicylic Acid 20 mg/ dL, Ascorbic Acid 2 g/dL, Creatine 200 mg/dL, Bilirubin 1 g/dL, Albumin 2 g/dL, Hemoglobin 1000 mg/dL, Oxalic Acid 60 mg/dL, Rematoid Factor 980 IU/mL, C Reactive Protein 41.2 mg/dL and Anti Streptolysin O 1023 IU/ml. Plasma samples stored for periods longer than those recommended by this instruction for use (30 days) may show fibrin and fibronectin precipitation that may interfere with the test. These samples must not be used. Samples collected with anticoagulants other than those mentioned in this instruction for use (EDTA and Heparin) should not be used.

CROSS REACTIVITY

A cross-reactivity study was carried out, evaluating 132 samples negative for SARS-CoV-2, but positive for other infections. Among they 12 samples positive for Influenza, 27 positive for Rhinovirus, 13 positive samples for Respiratory Syncytial Virus, 10 positive samples for Zika, 10 positive samples for Dengue, 10 positive samples for Chikungunya, 10 samples positive for Toxoplasmosis, 10 samples positive for Rubella, 10 HIV positive samples, 10 samples positive for HBV and 10 positive samples for HCV. It was not observed cross-reactivity with positive samples for Influenza, Rhinovirus, Respiratory syncytial virus, Zika, Dengue, Chikungunya, Toxoplasmosis, Rubella, HIV, HBV and HCV. Despite the results found, it is not can completely rule out the possibility of cross-reactivity. O final diagnosis should consider the patient's clinical data together with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

PRUDUCT PERFORMANCE QUALITY CONTROL**Accuracy REPETIBILITY**

The repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

REPETIBILITY	SAMPLE		
	1	2	3
Average	2.039	1.164	0.233
Standard Deviation	0.032	0.029	0.015
Coefficient of Variation (%)	1.59	2.49	6.40

REPRODUCIBILITY

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

REPRODUCIBILITY	SAMPLE		
	1	2	3
Average	2.050	1.148	0.230
Standard Deviation	0.032	0.030	0.014
Coefficient of Variation (%)	1.54	2.61	6.21

Clinical Sensitivity and Specificity

BIOLISA SARS-CoV-2 IgG Kit analyzed clinical samples in comparison with other methods of EIA. The results show that the clinical sensitivity of the BIOLISA SARS-CoV-2 IgG kit is 92.6% and clinical specificity is 92.0%.

	Expected Result	BIOLISA SARS-CoV-2 IgG
Positive Sample	54	50
Negative Sample	50	46

Clinical Sensitivity: (50/54)= 92.6% (IC 95%: 85.8% - 99.3%)

Clinical Specificity: (46/50)= 92.0% (IC 95%: 84.7% - 99.2%)

DIAGNOSTIC SIGNIFICANCE

Coronaviruses (CoV) belong to the Coronaviridae family, and are widely distributed infecting humans and other mammals. In humans, the most common symptoms presented are: fever, cough, dyspnoea, myalgia or fatigue. It can cause diseases ranging from the cold common to more serious diseases, such as the Eastern Respiratory Syndrome (MERS-CoV) and Severe Acute Respiratory Syndrome SARS-CoV. In December 2019, a series of cases of pneumonia caused disease appeared in Wuhan in China, with clinical symptoms very similar to viral pneumonia, after complete sequencing of the respiratory samples was found to be Coronavirus infection, which it was called the new coronavirus (2019-nCoV). In approximately 150 countries on all continents have already been reported cases of infection.

NUMBER OF TESTS

Presentation 1 - 96 tests

Presentation 2 - 192 tests

Presentation 3 - 480 tests

BIBLIOGRAPHIC REFERENCES

1. Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
2. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/ LAB/99.1 rev. 2, 2002:31.
3. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
4. Corman VM, Landt O, et al. Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
5. Huang C, Wang Y, et al. Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb 15;395(10223):497-506.
6. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca

CEP 31565-130 - Belo Horizonte - MG - Brasil

Phone: +55 (31) 3439.5454

E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service

Phone.: 0800 0315454

E-mail: sac@bioclin.com.br

ANVISA registration for BIOLISA SARS-CoV-2 IgG kit: 10269360330

Review: February/2022

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

BATCH CODE

DATE OF MANUFACTURE

USED BY
(last day of month)

TEMPERATURE LIMITATION
(store at)

CONTAINS SUFFICIENT
FOR <N> TESTS

CONSULT INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC DEVICE

EUROPEAN AUTHORIZED
REPRESENTATIVE

KEEP AWAY
FROM SUNLIGHT

MANUFACTURED BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

INFLAMMABLE

CORROSIVE

POISON

CE MARK

DO NOT USE IF
PACKAGE IS
DAMAGED