FIV Ac / FeLV Ag COMBO VET FAST

[REF] VET013

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FINALIDADE

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USE INSTRUCTIONS

PURPOSE

Fast combined immunochromatographic test for the qualitative determination of Antibodies against the virus of feline immunodeficiency (FIV) and feline leukemia virus Antigen (FeLV) in serum, plasma or whole blood samples.

PRINCIPLE OF ACTION

Methodology: Immunochromatography.

The FIV Ac / FeLV Ag Combo VET FAST is an in vitro diagnostic test for the rapid and qualitative detection of total Antibodies against feline immunodeficiency virus (FIV Ac) and feline leukemia virus Antigen (FeLV Ag) in feline blood, serum or plasma.

The test consists of two independent cassettes, linked by a support (combo cassette): one detects Anti-FIV Antibodies and the other detects the FeLV Antigen. Each cassette consists of a nitrocellulose membrane, impregnated in the test line with the feline immunodeficiency virus Antigen (FIV cassette) and Antibodies against feline leukemia virus (FeLV cassette). The conjugate consists of FIV Antigen (FIV cassette) and Anti-FeLV Antibody (FeLV cassette); both conjugated with colloidal gold particles.

When the sample is added to the cassettes, it migrates by capillary action through the test strip until it reaches the conjugate. The following occurs in this region:

In the FIV cassette: Anti-FIV Antibodies, if present in the sample, interact with the conjugate, formed of colloidal gold conjugated FIV Antigen. Therefore, formation of the complex occurs: Anti-FIV Antibody – FIV Antigen with colloidal gold. The sample continues to migrate through the test strip until it reaches the test line, where the Anti-FIV Antibody also interacts with the immobilized FIV Antigen in the test line, forming the complex: conjugate - Anti-FIV Antibody - Antigen immobilized. This shows the formation of a pink-red line on the test line; indicating positive result.

In the FeLV cassette: the FeLV Antigen, possibly present in the sample, interacts with the conjugate, formed of colloidal gold conjugated Anti-FeLV Antibodies. This is how complex formation occurs: FeLV Antigen - colloidal gold conjugated Anti-FeLV Antibody. The sample continues to migrate through the test strip until it reaches the test line, where the FeLV Antigen also interacts with the Anti-FeLV Antibody immobilized in this area, forming the complex: conjugate - FeLV Antigen - immobilized Anti-FeLV Antibody. This shows the formation of a pink-red line on the test line; indicating positive result. In both cassettes, sample excess continues to migrate through the test strip until it reaches the control line, which is composed of impregnated Anti-conjugate Antibody. Thus, the formation of a second line (control line) is observed; indicating that the test yielded a valid result.

REAGENTS

Reagent Number 1 (R1): Cassette - Store between 15 and 30°C. **Do not freeze.** It contains:

A) FIV cassette:

01 plastic cassette with 01 test strip containing:

- Control Region (C): Goat immobilized Anti-mouse Antibody.
- Test Region (T): Immobilized FIV Antigen.
- Conjugate: FIV Antigen colloidal gold conjugate and mouse Antibody colloidal gold conjugate.

B) FeLV cassette:

01 plastic cassette with 01 test strip containing:

- Control Region (C): Immobilized Anti-mouse Antibody.
- Test Region (T): Immobilized monoclonal Antibody against feline leukemia virus.
- Conjugate: Anti-FeLV colloidal gold conjugate and mouse Antibody colloidal gold conjugate.

C) 01 Sachet of Silica

Reagent Number 2 (R2): Diluent - Store between 15 and 30°C. **Do not freeze.** Contains: Buffer solution and preservative.

Collection Accessories:

- Disposable Plastic Pipette

PRESENTATIONS

Reagent	Presentation	Number of Cassettes for Package	Diluent
Cassette + Diluent	VET013-1	10	1 x 2.5 mL
	VET013-2	25	1 x 5.0 mL
	VET013-3*	10	10 x 0.5 mL
	VET013-4*	25	25 x 0.5 mL
	VET013-5	5	1x1,25 mL

*Individual presentations

EQUIPMENTS AND OPERATIONAL INPUTS


Clock or stopwatch. These items are found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperature of the product must be between 15 and 30°C. Keep away from direct sunlight. **Do not freeze.** The combo cassette is sensitive to air humidity. Therefore, it should always be kept in its inviolate primary packaging (aluminized sachet) until the time of use.

SPECIAL CARE

- Product for *in vitro* diagnosis, for veterinary use only.**
- Strictly follow the methodology proposed for obtaining results.
- Cassettes (test devices) do not contain contaminants; but since the samples are biological, the professional must strictly follow the safety rules and routines. The use of disposable gloves and other personal protective equipment is essential.
- Use clean, dry containers to package biological samples.
- Open the aluminum envelop to remove the combo cassette only at the time of use.
- The kit must not be frozen or exposed to high temperatures, since this will cause deterioration of its components.
- The samples should reach room temperature before starting the test.

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QUALITY GUARANTEE

Before being released for consumption, all Quibasa reagents are tested by the Quality Control Department. The quality of the reagents is assured until the validity date mentioned on the packaged, as long as they are stored and transported under the appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Tel.: (31) 3439.5454
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Industria Brasileira

CONSUMER SERVICE

Customer Advisory Service
Tel.: 0800 0315454
E-mail: sac@bioclin.com.br

Product licensed in the Ministry of Agriculture, Livestock and Supply since 11/11/2019 with the number 10.283/2019.

Technical manager: Dr. Hugo Vieira Fróes (CRMV/MG 8963).

Review: July/2021

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal program of quality control, where procedures, standards, limits and tolerance for variations are clearly established. It is important to note that all measuring systems characteristic analytical variability, which should be monitored by the laboratories themselves. Therefore, it is the use of controls, which allow the assessment of the precision and accuracy of dosages.

PRODUCT PERFORMANCE

QUALITY CONTROL

CLINICAL SENSITIVITY AND SPECIFICITY

This product has been tested comparing with other methods. A total of 102 samples of FIV Ac were analyzed, which 41 samples were diagnosed as reagent and 61 were diagnosed as no reactive. The results show that the clinical sensitivity of the product is 97.56% and the clinical specificity is 96.72%.

BIBLIOGRAPHIC REFERENCES

1- QUIBASA: Data from the R&D Research and Development Department.