## VETLISA BABESIA IgM

REF **VET065**

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#### FINALIDADE

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- TABOADA J, MERCHANT SR. Babesiosis of Companion Animals and Man. Veterinary Clinics of North America: Small Animal Practice. n.1. v. 21, p. 103-123. Janeiro, 1991.
- COMAZZI S, PALTRINIERI S, MANFREDI MT, AGNES F. Diagnosis of canine babesiosis by Percoll gradient separation of parasitized erythrocytes. Journal of Veterinary Diagnostic and Investigation. v.11, p. 102-104. 1999.
- BOOZER AL, MACINTIRE DK. CANINE BABESIOSIS. Vet Clin Small Anim, v.33, p.885-904, 2003.
- BRANDÃO LP, HAGIWARA MK, MYIASHIRO SI. Humoral immunity and reinfection resistance in dogs experimentally inoculated with *Babesia canis* and either treated or untreated with imidocarb dipropionate. Vet Parasitol., v.114, p.253-265, 2003.
- VERCAMMEN F, DE DEKEN R, MAES L. Duration of protective immunity in experimental canine babesiosis after homologous and heterologous challenge. Veterinary Parasitology, v.68, p. 51-55, 1997.
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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## VETERINARY USE

#### FUNCTION

Test for qualitative determination of IgM antibodies against *Babesia vogeli* in serum or plasma of dogs, through immunoenzymatic assay, in microplate. For *in vitro* diagnostics only.

#### PRINCIPLE OF ACTION

**Methodology:** Enzyme immunoassay or immunoenzymatic assay.

VETLISA BABESIA IgM is a solid phase immunoenzymatic assay, based on the immunocapture principle for the qualitative detection of IgM antibodies against *Babesia vogeli* in serum or plasma of dogs. Antibodies against *Babesia vogeli* present in the sample bind to the recombinant antigens of *Babesia vogeli* coated on the microplate, forming immobilized antigen-antibody immunocomplexes. After the initial incubation, the microplate is washed to remove unbound materials. The Conjugate, formed by anti-dog IgM antibody linked to peroxidase, is added and binds to IgM antibodies linked to the antigens immobilized on the plate. After a second incubation, the microplate is washed to remove unbound materials. Then, the Substrate is added and incubated. The Stop Solution is added to stop the reaction, promoting a change of color from blue to yellow. The result should be evaluated by reading the absorbance in a microplate reader.

#### REAGENTS

**1- Sensitized Plate** - Store between 2 and 8°C. Contains: Polystyrene plate, divided into 12 strips of 8 wells each, impregnated with *Babesia vogeli* recombinant antigen.

**2- Concentrated Conjugate (100X)** - Store at 2 to 8°C. Contains: Dog anti-IgM antibody conjugated to peroxidase and protein stabilizer.

**3- Concentrated Washing (20X)** - Store at 2 to 8°C. Contains: Buffer solution, surfactant and preservative.

**4- Diluent** - Store at 2 to 8°C. Contains: Buffer solution, stabilizers, dye, surfactant and preservative.

**5- Substrate TMB** - Store between 2 and 8°C. Contains: Solution containing 3, 3', 5, 5' - tetramethylbenzidine (TMB), urea peroxide and preservative.

**6- Stop Solution** - Store at 2 to 8°C. Contains: Sulfuric Acid Solution 1 mol/L.

**7- Negative Control** - Store between 2 and 8°C. Contains: Serum containing negative serum for IgM antibodies against *Babesia vogeli* and protein stabilizer. **Potentially infective.**

**8- Positive Control** - Store between 2 and 8°C. Contains serum containing positive serum for IgM antibodies against *Babesia vogeli* and protein stabilizer. **Potentially infective.**

**9- Plate Sealers**

#### PRESENTATIONS

COMPONENTS	PRESENTATION			
	1	2	3	4
	96 cavities	192 cavities	480 cavities	960 cavities
<b>1- Sensitized Plate</b>	1 unit	2 units	5 units	10 units
<b>2- Concentrated Conjugate</b>	1 x 300 µL	2 x 300 µL	5 x 300 µL	10 x 300 µL
<b>3- Concentrated Washing</b>	1 x 50 mL	2 x 50 mL	5 x 50 mL	10 x 50 mL
<b>4- Diluent</b>	1 x 60 mL	2 x 60 mL	5 x 60 mL	10 x 60 mL
<b>5- Substrate TMB</b>	1 x 15 mL	2 x 15 mL	5 x 15 mL	10 x 15 mL
<b>6- Stop Solution</b>	1 x 15 mL	2 x 15 mL	5 x 15 mL	10 x 15 mL
<b>7- Negative Control</b>	1 x 1 mL	2 x 1 mL	5 x 1 mL	10 x 1 mL
<b>8- Positive Control</b>	1 x 1 mL	2 x 1 mL	5 x 1 mL	10 x 1 mL
<b>9- Plate Sealers</b>	3 units	6 units	15 units	30 units

#### EQUIPMENT AND OPERATIONAL INPUTS

**Materials contained in the kit:**

- Reagents described in the previous table
- Instructions for use (manual)

#### Necessary materials not contained in the kit:

- 1- Pipettes capable of dispensing volumes of 5 to 500 µL with a variation coefficient of less than 1.5%.
- 2- Repipettor for repetitive pipetting of volumes of 500 µL with coefficient of variation less than 1.5% or multichannel pipette (optional).
- 3- Microplate washer (optional).
- 4- ELISA reader with absorbance capacity at 450 and 630 nm wavelength.
- 5- Absorbent paper to dry the wells.
- 6- Stopwatch or watch.
- 7- Bottle to store the Washing Solution after dilution.
- 8- Distilled or deionized water.
- 9- Quality Control Tools.
- 10- Incubator at 37 ° C ± 2 ° C.

#### STORAGE AND TRANSPORT CONDITIONS

The product storage temperature is 2 to 8 °C. Transport can be done at room temperature (up to 30 °C) for up to 5 days. Keep away from light and avoid humidity. **Do not freeze.**

#### SPECIAL CARE

- 1- For professional in vitro diagnostic veterinary use only.
- 2- Strictly follow the proposed methodology to obtain accurate results.
- 3- The sachet containing the microplate must be opened only after reaching room temperature. Replace the unused microcavity strips in the sachet, seal and store at 2 to 8 °C.
- 4- The water used to clean the material must be fresh and free from contaminants.
- 5- Saturated deionizing columns release alkaline water, several ions and oxidizing and reducing agents that can significantly alter the results.
- 6- Do not mix reagents from different kits or lots. Do not use expired kit components.
- 7- Do not eat, drink or smoke at the test site.

**8-** Do not pipette reagents or sample(s) using the mouth. Do not use the same tip to pipette different samples.

**9-** The Positive and Negative Controls must be retested for each new test performed.

**10-** The Stop Solution contains Sulfuric Acid which is a strong acid. Therefore, handle it with due care.

**11-** The professional must strictly follow the safety rules and routines when handling biological samples. The use of disposable gloves and other personal protective equipment is essential.

**12-** Pipet the reagents always in the same order to minimize the difference in the reaction time between the wells.

**13-** As a safety measure, cover the plate during the reaction.

**14-** Make sure that the bottom of the cavity is clean and dry and that there are no bubbles on the liquid surface before reading the plate. Do not allow the wells to dry out during the test.

**15-** Handle the kit components with due care, in order to avoid contamination. Use new or properly cleaned tips and containers for handling and do not expose reagents, especially Substrate, to strong light or Hypochlorite vapors during storage or incubation steps.

**16-** We recommend applying local, state and federal environmental protection standards so that reagents and biological material are disposed of in accordance with current legislation.

**17-** To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Service of Customer Service) of Quibasa.

**18-** Do not use the product in case of damage to the packaging.

**19-** It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

#### SAMPLES

**Serum or Plasma (EDTA or Heparin).**

Hemolyzed or highly lipemic samples should not be used.

The samples can be kept refrigerated, between 2 and 8 °C, for a maximum period of 5 days. If the samples cannot be analyzed within 5 days, they can be stored for up to 30 days at -20 °C.

#### PROCESS DESCRIPTOR

##### Stability After Open

The results of the stability test show that the VETLISA BABESIA IgM kit is stable after being opened for at least 30 days. This stability may vary according to the test conditions and the environment. Therefore, it is suggested to monitor the product's performance using internal kit controls and the technique validation criteria.

##### Preparation of Work Reagents

###### WASHING SOLUTION

Dilute the contents of the Concentrated Washing (Reagent N° 3) in 1000 mL of distilled or deionized water. After preparation, the solution can be stored at 2 to 30 °C until the expiration date printed on the original bottle. It can be stored at room temperature. If crystallization occurs, heat to 37°C until dissolved.

###### CONJUGATE SOLUTION

Dilute the Concentrated Conjugate (Reagent N° 2) in the proportion of 1: 101 in Diluent (Reagent N° 4). Prepare the solution when performing the test.

To perform a test using all the wells of the kit, mix 110 µL of the Concentrated Conjugate in 11 mL of Diluent.

To perform an assay using 8 wells (1 strip), mix 10 µL of the Concentrated Conjugate in 1 mL of Diluent.

**IMPORTANT:** The diluted conjugate solution cannot be stored. Therefore, prepare only the amount necessary to perform the test.

##### DILUTION OF SAMPLES

In a test tube, dilute 15 µL of the sample in 300 µL of Diluent, if the test is going to be performed in duplicate. Cap the tube and vortex gently or mix manually by inversion. Dilutions cannot be stored.

**ATENÇÃO:** The Positive and Negative Controls are ready to use.

##### SUBSTRATE

The Substrate is ready to use.

##### Technique

Before starting the test, place all reagents, samples and controls to stabilize at room temperature (15 - 30 °C) for at least 40 minutes.

**1-** Separate the cavities to be used considering: Controls and Samples (it is recommended to test in duplicate). Return the unused strips of the microplate to the original sealed package.

**2-** Separate the first cavity for Blank (OPTIONAL).

**3-** Pipette 100 µL of the Negative Control and Positive Control in the previously determined wells. **Note: The controls are ready for use, there is no need to dilute them.**

**4-** Pipette 100 µL of the samples previously diluted in the previously determined wells. In the Blank cavity (OPTIONAL), if you have made the option, pipette only 100 µL of the diluent.

**5-** Gently mix for ± 10 seconds. Cover the wells with plate sealer.

**6-** Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.

**7-** Remove the sealer from the cavities.

**8-** Discard the contents of the cavities by aspiration (Washer) or by decantation (Manual).

Use approximately 300 µL of Washing Solution, previously prepared, and perform a total of five (5) wash cycles. Shake for 3 seconds with each wash. To guarantee the drying of the plate, at the end of the washing, beat the plate for a few seconds on absorbent paper.

**Note:** Poor washing / drying can cause inadequate results

**9-** Pipette 100 µL of Conjugate previously diluted in all wells.

**10-** Gently homogenize for ± 10 seconds. Cover the cavities with the plate sealer.

**11-** Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.

**12-** Remove the plate sealer from the cavities.

**13-** Repeat item 8.

**14-** Pipette 100 µL of Substrate into all wells.

**15-** Gently homogenize for ± 10 seconds. Cover the cavities with the plate sealer.

**16-** Incubate for 10 minutes ± 2 minutes in an incubator at 37°C ± 2°C.

**17-** Remove the plate sealer from the cavities.

**18-** Pipette 100 µL of Stop Solution into all wells.

**19-** Gently mix for ± 10 seconds.

**20-** Read using a double filter: 450 nm / 630 nm in up to 15 minutes (maximum).

**TECHNIQUE VERIFICATION**

Verify if the results obtained for reading the Blank and the Controls are compatible with the values presented below:

ITEM	ABSORBANCE (DOUBLE FILTER)
Blank	< 0,100
Negative Control	0,100 a 0,300
Positive Control	> 1.000

If the values are outside the expected values, the technique should be repeated.

**Calculations****QUALITATIVE**

Calculate the Cut Off according to the following formula:

$\text{Cut Off} = \text{Average Negative Control Absorbance} + 0.450$
---

Example:

ITEM	ABSORBANCE
Negative Control	0.212
	0.208
<b>Cut-Off = Average Negative Control Absorbance + 0.450</b>	$((0.212 + 0.208) / 2) + 0.450 = 0.660$

Calculate the Index by dividing the absorbance of the sample by the Cut Off value. Example:

ITEM	ABSORBANCE
Sample	1.680
Cut Off value	0.660
<b>Index = Sample / Cut Off value</b>	$1.680 / 0.660 = 2.55$

**INTERPRETATION OF RESULT**

After calculating the sample index, consider the indices below to determine the results.

RESULTS	INDEX
Negative	< 0.9
Positive	> 1.1
Indeterminate	≥0,9 - ≤1,1

**Note:** Indeterminate result indicates a test response within the indeterminate range (gray zone, between cut-off points) and, by itself, does not confirm or rule out infection. It may reflect an immunological window, early seroconversion, antibody decline, or pre-analytical interferences. Interpretation should take into account the patient's history and clinical presentation, with repetition of the test (preferably using a new sample), re-sampling in 7–14 days, and, depending on the suspicion and stage of the disease, performance of confirmatory tests. The diagnostic conclusion is the responsibility of the Veterinarian, integrating clinical signs and complementary examinations.

**Note:** The data presented in the examples are for illustration only and cannot be used to calculate the results.

**PROCESS LIMITATIONS**

The presence of IgM antibodies against *Babesia vogeli* should be evaluated along with clinical information and other data from complementary tests. IgM production usually begins between 7 and 14 days after infection and persists for approximately 10–14 days. The presence of IgM indicates that the dog has demonstrated an immune response to exposure to *Babesia vogeli*; however, the absence of IgM antibodies does not rule out infection. Factors such as the stage of the disease, the individual's immune response, the amount of antibodies available in the sample, and the presence of interfering factors should be considered when interpreting the result.

**INTERFERENCES**

No interference was observed for Triglycerides 1500 mg/dL, Ascorbic Acid 2 g/dL, Creatine 200 mg/dL, Bilirubin 1 g/dL, Albumin 2 g/dL, Hemoglobin 1000 mg/dL, Oxalic Acid 60 mg/dL, Factor Rheumatoid 980 IU/mL, C Reactive Protein 41.2 mg/dL and Anti-Streptolysin O 1023 IU/MI.

**CROSS REACTIVITY**

A cross-reactivity study was carried out, evaluating 45 negative samples for IgM antibodies against *Babesia vogeli*, but positive for other infections. Among then, it was tested 12 positive samples for Brucellosis, 2 positive samples for Anaplasmosis, 19 positive samples for Distemper, 3 positive samples for Parvovirus and 9 positive samples for Ehrlichiosis. No cross reactivity was observed with any of the samples. Despite the results found, the possibility of cross-reactivity cannot be completely ruled out. The final diagnosis should consider the patient's clinical data along with other laboratory data.

**INTERNAL QUALITY CONTROL**

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to note that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. For that, it is recommended to use controls, which allow to evaluate the precision and accuracy of the measurements.

**PRODUCT PERFORMANCE****Accuracy****REPEATABILITY**

Repeatability was calculated from 10 successive determinations, using 3 samples with different values. The following results were obtained:

REPEATABILITY	SAMPLE		
	1	2	3
<b>Average</b>	2.903	0.438	0.147
<b>Standard Deviation</b>	0.090	0.022	0.005
<b>Coefficient of Variation (%)</b>	3.07	5.01	4.03

**REPRODUCIBILITY**

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 3 samples with different values. The following results were obtained:

REPRODUCIBILITY	SAMPLE		
	1	2	3
<b>Average</b>	2.864	0.413	0.147
<b>Standard Deviation</b>	0.136	0.026	0.008
<b>Coefficient of Variation (%)</b>	4.73	6.36	5.72

**Clinical Sensitivity and Specificity**

In the VETLISA BABESIA IgM, clinical samples were analyzed in comparison with another method of Enzyme immunoassay. The results show that the clinical sensitivity of the kit is 95.7% and the clinical specificity is 96.1%.

Method	Reference		Total	
	Result	Positive		Negative
VETLISA BABESIA IgM	Positive	22	2	24
	Negative	1	51	52
<b>Total Result</b>		23	53	76

Clinical Sensibility: 95.7% (22/23) - IC 95%: 87.5% - 100%

Clinical Specificity: 96.1% (51/53) - IC 95%: 90.9% - 100%

Accuracy: 96.0% [(22+51) / (23+53)]

**CLINICAL MEANING**

Babesiosis is a disease caused by protozoa of the genus *Babesia*, which can affect domestic dogs, wild animals, and humans. Transmission occurs primarily through the bite of

infected ticks during a blood meal. The incubation period between infection and the onset of clinical signs ranges from 7 to 30 days, and the disease can manifest as peracute, acute, chronic, or asymptomatic. The most frequently observed clinical signs include anemia, thrombocytopenia, hemoglobinuria, fever, pale mucous membranes, tachycardia, tachypnea, depression, anorexia, jaundice, and splenomegaly. In chronic cases, progressive weight loss and persistent anorexia are common. Due to the clinical similarities between different hemoparasitoses and the possibility of coinfections, since they share the same vector, the use of complementary diagnostic tools, such as serological tests, is essential. These tests should be used to support clinical evaluation, never in isolation. Laboratory results, combined with the history and clinical examination, should be interpreted carefully, with the diagnostic conclusion always being the veterinarian's responsibility.

**BIBLIOGRAPHIC REFERENCES**

- 1- TABOADA J, MERCHANT SR. Babesiosis of Companion Animals and Man. Veterinary Clinics of North America: Small Animal Practice. n.1. v. 21, p. 103-123. Janeiro, 1991.
- 2- COMAZZI S, PALTRINIERI S, MANFREDI MT, AGNES F. Diagnosis of canine babesiosis by Percoll gradient separation of parasitized erythrocytes. Journal of Veterinary Diagnostic and Investigation. v.11, p. 102-104. 1999.
- 3- BOOZER AL, MACINTIRE DK. CANINE BABESIOSIS. Vet Clin Small Anim, v.33, p.885-904, 2003.
- 4- BRANDÃO LP, HAGIWARA MK, MYIASHIRO SI. Humoral immunity and reinfection resistance in dogs experimentally inoculated with *Babesia canis* and either treated or untreated with imidocarb dipropionate. Vet Parasitol., v.114, p.253-265, 2003.
- 5- VERCAMMEN F, DE DEKEN R, MAES L. Duration of protective immunity in experimental canine babesiosis after homologous and heterologous challenge. Veterinary Parasitology, v.68, p. 51-55, 1997.
- 6- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

**QUALITY ASSURANCE**

Before being released for consumption, all reagents produced by Quibasa Química Básica Ltda are tested by the Quality Control Department. The quality of the reagents is ensured until the expiration date mentioned on the packaged, when stored and transported under the appropriate conditions.

**QUIBASA QUÍMICA BÁSICA Ltda**

Rua Teles de Menezes, 92 – Santa Branca  
CEP 31.565 -130 - Belo Horizonte - MG - Brazil

Phone: (31) 3439.5454

E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 – Made in Brazil

**CUSTOMER SERVICE**

Customer Advisory Service

Phone: 0800 0315454

E-mail: sac@bioclin.com.br

**Product Licensed at the Ministry of Agriculture, Livestock and Supply from 03/03/2022 under the number 10.490/2022.**

**Technical manager:** Dra. Camila Eckstein (CRMV/MG 20611)

**Review:** January/2026

**UNIVERSAL SYMBOLOLOGY**

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER