VETLISA BABESIA IgG

REF VET038

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FINALIDADE

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PURPOSE

Test for qualitative determination of IgG antibodies against *Babesia canis* using dog serum or plasma, by solid phase enzyme immunoassay, in microplate. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzymatic The VETLISA Babesia IgG is an immunoenzymatic assay in solid phase, based on the principle of immunocapture for the qualitative detection of IgG antibodies against *Babesia canis* in serum or plasma. Antibodies against *Babesia canis* present in the sample, bind to the recombinant antigen of *Babesia canis* coated on the microplate, forming immobilized complexes: *Babesia canis* antigen – anti-*Babesia canis* antibodies. After the initial incubation, the microplate is washed to eliminate unbound materials. The conjugate, formed by the anti-dog IgG antibody conjugated to peroxidase, binds the immobilized antigen-antibody complex on the plate. After the second incubation, the microplate is washed and incubated with substrate. The intensity of the blue color produced by the addition of the substrate is proportional to the amount of antibodies against *Babesia canis* present in the sample. The stop solution is added to stop the reaction, promoting a color change to yellow, measured in a microplate reader.

KIT COMPONENTS

- 1- Sensitized Plate-** Store between 2 and 8°C. Contains: Polystyrene plate, divided into 12 strips of 8 wells each, impregnated with *Babesia canis* antigen.
- 2- Concentrated Conjugate (100x)** - Store between 2 and 8°C. Contains: Peroxidase-linked anti-dog IgG antibody solution.
- 3- Concentrated Wash Solution (20x)** - Store between 2 and 8°C. Contains: Phosphate Buffer Solution, surfactant and preservative.
- 4- Diluent** - Store between 2 and 8°C. Contains: Buffer solution, surfactant, stabilizer and preservative.
- 5- Substrate TMB** - Store between 2 and 8°C. Contains: Solution containing tetramethylbenzidine (TMB <1.0 mg / mL), Citric Acid Solution <5% and Urea Peroxide <1%.
- 6- Stop Solution** - Store between 2 and 8°C. Contains: 1M Sulfuric Acid.
- 7- Negative Control** - Store between 2 and 8°C. Contains: Non-reactive serum for canine Babesia, with preservative. **Potentially infectious.**
- 8- Positive Control** - Store between 2 and 8°C. Contains: Reactive serum for canine Babesia, with preservative. **Potentially infectious.**
- 9- Plate Sealers**

PRESENTATIONS

COMPONENTS	PRESENTATION		
	1	2	3
	96 wells	192 wells	480 wells
1- Sensitized Plate	1 unit (96 wells)	2 units (2 x 96 wells)	5 units (5 x 96 wells)
2- Concentrated Conjugate (100x)	1 x 350 µL	2 x 350 µL	5 x 350 µL
3- Concentrated Wash Solution (20x)	1 x 50 mL	2 x 50 mL	5 x 50 mL
4- Diluent	1 x 60 mL	2 x 60 mL	5 x 60 mL
5- Substrate TMB	1 x 15 mL	2 x 15 mL	5 x 15 mL
6- Stop Solution	1 x 15 mL	2 x 15 mL	5 x 15 mL
7- Negative Control	1 x 1,0 mL	2 x 1,0 mL	5 x 1,0 mL
8- Positive Control	1 x 1,0 mL	2 x 1,0 mL	5 x 1,0 mL
9- Plate Sealers	3 units	6 units	15 units

EQUIPMENTS AND OPERATIONAL INPUTS

- Materials contained in the kit:**
- Reagents described in the previous table.
 - Instructions for use (manual).

Materials needed, but not contained in the kit:

- 1-** Calibrated pipettes capable of dispensing volumes between 5 and 100 µL, with an accuracy greater than 1.5%.
- 2-** Repipetizers for repetitive pipetting of 100µL and 300µL volumes, with an accuracy greater than 1.5% (optional) or multichannel pipette.
- 3-** Microplate washing machine (optional).
- 4-** ELISA reader with absorbance capacity at 450 and 630 nm wavelength.
- 5-** Absorbent paper to dry the microcavities.
- 6-** Stopwatch or clock.
- 7-** Bottle to store the Washing Solution, after dilution.
- 8-** Distilled or deionized water.
- 9-** Quality Control Tools.
- 10-** Incubator with a capacity of 37°C ± 2°C.

STORAGE AND TRANSPORT CONDITIONS

The product storage temperature should be 2 to 8 °C. Transport can be carried out at room temperature (30°C) for up to 5 days. Keep away from light and avoid humidity. **Do not freeze.**

SPECIAL CARE

- 1-** Product for professional use.
- 2-** Product for in vitro diagnostic, for veterinary use only.
- 3-** Strictly follow the methodology proposed to obtain accurate results.
- 4-** The professional must strictly follow safety rules and routines when handling biological samples. The use of disposable gloves and other personal protective equipment is essential.
- 5- IMPORTANT:** Before starting the test, allow all kit components to reach room temperature. Open the envelope containing the wells only after reaching room temperature.
- 6-** Do not mix reagents from different kits or lots. Do not use expired kit components.
- 7-** Do not eat, drink, or smoke at the test site.
- 8-** Do not pipette reagents or samples by mouth. Do not use the same tip to pipette different samples.
- 9-** Positive and negative controls should be retested for each new test performed.
- 10-** The Stop Solution contains Sulfuric Acid, which is a strong acid. Therefore, handle it with care.
- 11-** Handle the kit components with care to avoid contamination. Take special care with the substrate, which is a light sensitive solution. Use new or appropriately clean tips and containers for handling and do not allow them to be exposed to strong light during storage or during the incubation periods of the assay.
- 12-** Use contaminant-free freshly distilled or deionized water when preparing solutions.
- 13-** Always add the reagents in the same order to minimize the difference in reaction time between the microcavities.
- 14-** Before reading the plate, ensure that the bottom of the microcavities are clean and dry. Make sure there are no bubbles on the surface of the liquid.
- 15-** We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 16-** To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa
- 17-** Do not use the product in case of damaged packaging.
- 18-** It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Use serum or plasma (EDTA or Heparin). Hemolyzed or highly lipemic samples should not be used. Samples may be refrigerated at between 2 and 8°C for a maximum of 5 days. If samples cannot be analyzed within 5 days, they can be stored for up to 30 days at -20°C (freezer).

PROCESS DESCRIPTION

PREPARATION OF WORK REAGENTS AND SAMPLES

1- Washing Solution

Dilute the content of the Concentrated Wash Solution (Reagent N° 3) in 1000 mL of freshly distilled or deionized water. Store at 2 to 8°C until the validity date printed on the original bottle. It can be stored at room temperature. If crystallization occurs, heat to 37°C until dissolution.

2- Conjugate Solution

Dilute the concentrated conjugate (Reagent N° 2) in the ratio of 1/100 in Diluent (Reagent N° 4). Prepare the solution at the moment of performing the test. To perform an assay using all wells in the kit, mix 110 µL of concentrated conjugate in 10,9 mL of Diluent. To perform an assay using 8 wells (1 strip), mix 10 µL of the Concentrated Conjugate in 990 µL of Diluent. **IMPORTANT:** The diluted conjugate solution cannot be stored. Therefore, prepare only the amount necessary to perform the test.

3- Sample and Controls Dilution

In a test tube, dilute 5 µL of the Sample, Negative Control (Reagent N° 7) and Positive Control (Reagent N° 8) in 400 µL of diluent. Cap the tube and gently vortex or manually mix by inversion. Dilutions cannot be stored.

TECHNIQUE

Before starting the assay, place all reagents, Controls and Samples to stabilize at room temperature (15 - 30°C) for at least 30 minutes. Return the unused strips to the original sealed packaging.

- 1-** Separate the wells to be used considering: Controls and Samples ((it is recommended to test in duplicate).
- 2-** Separate the first well for the Blank (OPTIONAL).
- 3-** Pipette 100 µL of the Samples, Negative Control and the Positive Control previously diluted into the previously determined wells.
- 4-** In the Blank cavity (OPTIONAL), if you have chosen, pipette only 100 µL of the diluent.
- 5-** Homogenize gently for ± 10 seconds. Cover the cavities with plate sealer.
- 6-** Incubate for 1 hour at 37°C.

7- Discard the content of the cavities by aspiration (washing machine) or by decantation (manual). Use approximately 300µL/well of previously diluted Wash Solution, for a total of five (5) wash cycles. Shake for three seconds after each wash cycle. To guarantee the drying of the plate, at the end of the wash, tap the plate for a few seconds on absorbent paper.

IMPORTANT: Poor washing/drying can cause inadequate results.

8- Pipette 100 µL of the Conjugate previously diluted into each well, even in the Blank cavity.

9- Homogenize gently for ± 10 seconds. Cover the cavities with plate sealer.

10- Incubate for 1 hour at 37°C.

11- Repeat item 7.

12- Add 100 µL of Substrate TMB into each well.

13- Homogenize gently for ± 3 seconds. Cover the cavities with plate sealer.

14- Incubate for 10 minutes at 37°C, protected from light.

15- Remove the plate sealer from the cavities.

16- Pipette 100 µL of Stop Solution into all wells.

17- Homogenize gently for ± 3 seconds

18- Read the absorbance in an ELISA reader on a dual filter 450 nm (primary filter) / 630nm (secondary filter) within a maximum of 10 minutes after adding the Stop Solution.

TECHNIQUE VERIFICATION

Verify if the results obtained from reading the Controls and the Blank are compatible with the specification above:

ITEM	ABSORBANCE (DUAL FILTER)
Blank	< 0,050
Negative Control	0,100 a 0,300
Positive Control	> 0,800

If the values are out the expected ones, the test should be repeated.

CUT OFF CALCULATION

If the results of the controls are valid, calculate the Cut off value with the following formula:

Cut Off = Negative Control Average Absorbance + 0,200

CALCULATION OF THE SAMPLE INDEX

Calculate the index by dividing the absorbance of the sample by the cut off value. Example:

ITEM	ABSORBANCE
Sample	0,992
Negative Control	0,210
Cut Off Value	0,210 + 0,200 = 0,410
Index: Sample Abs. / Cut Off Value	0,992 / 0,410 = 2,42

INTERPRETATION OF THE RESULT - QUALITATIVE

Samples are considered positive if the index value is greater than 1,1. They are considered negative if the index value is less than 0,9. They are considered indeterminate if the value of the index is greater than 0,9 and less than 1,1.

RESULTS	INDEX
Negative	< 0,9
Positive	> 1,1
Undetermined	≥ 0,9 y ≤ 1,1

IMPORTANT: Samples with indeterminate results must be retested. Samples that repeatedly show indeterminate results should be analyzed using an alternative method. If this sample shows a positive result, then it should be considered positive.

PRODUCT PERFORMANCE

REPEATABILITY

Repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following results:

REPEATABILITY	ABSORBANCE		
	1	2	3
Mean	0,780	0,980	0,090
Standard Deviation	0,037	0,077	0,002
Coefficient of Variation (%)	4,730	7,800	2,604

REPRODUCIBILITY

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 3 samples with different values, obtaining the following results:

REPRODUCIBILITY	ABSORBANCE		
	1	2	3
Mean	0,780	1,000	0,093
Standard Deviation	0,050	0,049	0,005
Coefficient of Variation (%)	6,414	4,879	5,072

CLINICAL SENSITIVITY AND SPECIFICITY

The kit VETLISA BABESIA IgG analyzed clinical samples compared to other EIA method. The results show that the kit VETLISA BABESIA IgG presents a clinical sensitivity of 97,4% and a clinical specificity of 97,8%.

Method		Reference		Total
VETLISA Babesia IgG	Result	Positive	Negative	
	Positive	88	1	89
	Negative	2	37	39
Total Result		90	38	128

Clinical sensitivity: 97,4% (37/38)

Clinical specificity: 97,8% (88/90)

Precision: 97,7% [(88+37)/(90+38)]

PROCEDURE LIMITATIONS

The interpretation of a diagnostic test should not be based on a single trial. Other confirmatory tests must be included before a sample is considered positive. A negative result does not exclude the possibility of exposure. However, all results should be interpreted in conjunction with vaccination history, available clinical and laboratory information.

INTERNAL QUALIT CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

BIBLIOGRAPHIC REFERENCES

1- QUIBASA: Data from the Research and Development Department.

QUALITY GUARANTEE

Before being released for consumption, all Quibasa reagents are tested by the Quality Control Department. The quality of the reagents is assured until the validity date mentioned on the packaged, as long as they are stored and transported under the appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda.

Rua Teles de Menezes, 92 - Santa Branca

CEP 31565-130 - Belo Horizonte - MG - Brasil

Tel.: (31) 3439.5454

E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 - Industria Brasileira

CONSUMER SERVICE

Customer Support Service

Tel.: 0800 0315454

E-mail: sac@bioclin.com.br

Product licensed in the Ministry of Agriculture, Livestock and Supply since 10/04/2019 with the number 10.235/2019.

Technical manager:

Dra. Camila Eckstein (CRMV/MG 20611)

Review: September/2021

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

BATCH CODE

DATE OF MANUFACTURE

USED BY
(last day of month)

TEMPERATURE LIMITATION
(store at)

CONTAINS SUFFICIENT
FOR <N> TESTS

CONSULT INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC DEVICE

EUROPEAN AUTHORIZED
REPRESENTATIVE

KEEP AWAY
FROM SUNLIGHT

MANUFACTURED BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

INFLAMMABLE

CORROSIVE

POISON

CE MARK

DO NOT USE IF
PACKAGE IS
DAMAGED