VETLISA ERLQUIOSE IgG

REF VET040

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Test for qualitative determination of IgG antibodies against *Ehrlichia canis* in dog serum or plasma, by immunoenzymatic assay, in microplate. For *in vitro* diagnostics only.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzymatic.
This product is a solid-phase enzyme immunoassay based on the immunocapture principle for the qualitative detection of IgG antibodies against *Ehrlichia canis* in dog serum or plasma. Antibodies against *Ehrlichia canis* present in the sample bind to the recombinant *Ehrlichia canis* antigen that coats the microplate, forming immobilized antigen-antibody complexes. After the initial incubation, the microplate is washed to remove unbound materials. The conjugate, formed by anti-dog IgG antibody linked to peroxidase, binds to the antigen-antibody complex immobilized on the plate. Then the microplate is washed and incubated with Substrate. The intensity of the blue color produced by the addition of the Substrate is proportional to the amount of antibodies against *Ehrlichia canis* present in the sample. Stopping Solution is added to terminate the reaction, causing a color change to yellow as measured on a microplate reader.

REAGENTS

- 1- Sensitized Plate** - Store between 2 and 8 °C. Contains: Polystyrene plate, divided into 12 strips of 8 wells each, impregnated with *E. canis* recombinant antigen.
- 2- Concentrated Conjugate (100X)** - Store between 2 and 8 °C. Contains: Peroxidase-linked anti-dog IgG antibody solution.
- 3- Concentrated Wash (20X)** - Store between 2 and 8 °C. Contains: Buffered solution, surfactant and preservative.
- 4- Diluent** - Store between 2 and 8 °C. Contains: Buffered solution, surfactant, stabilizer and preservative.
- 5- Substrate TMB** - Store between 2 and 8 °C. Contains: Solution containing tetramethylbenzidine (TMB <1.0 mg/mL), citric acid solution <5% and urea peroxide <1%.
- 6- Stop Solution** - Store between 2 and 8 °C. Contains: 1N sulfuric acid solution.
- 7- Negative Control** - Store between 2 and 8 °C. Contains: Negative serum for canine ehrlichiosis and preservative. **Potentially infectious.**
- 8- Positive Control** - Store between 2 and 8 °C. Contains: Positive serum for canine ehrlichiosis and preservative. **Potentially infectious.**
- 9- Plate Sealers**

PRESENTATIONS

Reagents	1	2	3
	96 Cavities	192 Cavities	480 Cavities
1- Sensitized Plate	1 Unit (96 cavities)	2 Units (192 cavities)	5 Units (480 cavities)
2- Concentrated Conjugate (100x)	1 x 350 µL	2 x 350 µL	5 x 350 µL
3- Concentrated Wash (20x)	1 x 50 mL	2 x 50 mL	5 x 50 mL
4- Diluent	1 x 60 mL	2 x 60 mL	5 x 60 mL
5- Substrate TMB	1 x 15 mL	2 x 15 mL	5 x 15 mL
6- Stop Solution	1 x 15 mL	2 x 15 mL	5 x 15 mL
7- Negative Control	1 x 1 mL	2 x 1 mL	5 x 1 mL
8- Positive Control	1 x 1 mL	2 x 1 mL	5 x 1 mL
9- Plate Sealers	3 Units	6 Units	15 Units

EQUIPMENT AND OPERATIONAL INPUTS

Materials contained in the kit:
- Reagents described in the previous table.
- Instructions for use (manual).

Necessary materials not contained in the kit:

- 1-** Pipettes capable of dispensing volumes of 5 to 500 µL with a variation coefficient of less than 1.5%.
- 2-** Repipettor for repetitive pipetting of volumes of 500 µL with coefficient of variation less than 1.5% or multichannel pipette (optional).
- 3-** Microplate washer (optional).
- 4-** ELISA reader with absorbance capacity at 450 and 630 nm wavelength.
- 5-** Absorbent paper to dry the wells.
- 6-** Stopwatch or watch.
- 7-** Bottle to store the Washing Solution after dilution.
- 8-** Distilled or deionized water.
- 9-** Quality Control Tools.
- 10-** Incubator at 37 °C ± 2 °C.

TRANSPORTATION AND STORAGE CONDITIONS

The product storage temperature should be 2 to 8 °C. Transport can be carried out at room temperature 30 °C for up to 5 days. Keep away from light and avoid humidity. **Do not freeze.**

SPECIAL CARE

- 1-** For professional veterinary *in vitro* diagnostic use only.
- 2-** Strictly follow the proposed methodology to obtain accurate results.
- 3-** The sachet containing the microplate must be opened only after reaching room temperature. Replace strips from unused microwells in the sachet, seal and store at 2 to 8 °C.
- 4-** The water used to clean the material must be fresh and free of contaminants.
- 5-** Saturated deionizer columns release alkaline water, different ions, and oxidizing and reducing agents that can significantly alter the results.
- 6-** The Stop Solution contains Sulfuric Acid which is a strong acid. So handle it with due care.

- 7-** The handling of any product that contains serum is potentially capable of transmitting diseases. Therefore, it is necessary to take proper biosafety precautions when handling these products.
- 8-** Pipette the reagents always in the same order to minimize the difference in reaction time between the microwells.
- 9-** For protection measure, the plate must be covered during the reaction.
- 10-** Ensure that the bottom of the cavity is clean and dry and that there are no bubbles on the surface of the liquid before reading the plate. Do not allow wells to dry during the assay.
- 11-** Do not expose reagents, especially Substrate, to strong light or Hypochlorite vapors during storage or incubation steps.
- 12-** We recommend applying local, state and federal regulations for environmental protection so that the disposal of reagents and biological material is done in accordance with current legislation.
- 13-** To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Safety Information Sheet for Chemical Products) available on the website www.bioclin.com.br or upon request by the SAC (Service of Customer Service) by Quibasa.
- 14-** Do not use the product in case of damage to the packaging.
- 15-** It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or Plasma (EDTA or Heparin)
Hemolyzed or highly lipemic samples should not be used. The samples can be kept refrigerated, between 2 and 8 °C, for a maximum period of 5 days. If the samples cannot be analyzed within 5 days, they can be stored for up to 30 days at a temperature of -20 °C.

PROCESS DESCRIPTION

Stability After Open
The results of the stability test show that the VETLISA Ehrlichiosis IgG kit is stable after being opened for at least 30 days. This stability may vary according to the test conditions and the environment. Therefore, it is suggested to monitor the performance of the product using internal controls of the kit and the criteria for validating the technique.

ATTENTION: The Positive and Negative Controls are ready to use.

Preparation of Works Reagents

WASHING SOLUTION
Dilute the contents of vial No. 3 (Concentrated Washing Solution) in 1000 mL of distilled or deionized water. After preparation, the solution can be stored at 2 to 30 °C for at least 30 days. It can be stored at room temperature. If crystallization occurs, heat to 37 °C until dissolved.

SUBSTRATE

The Substrate is ready for use.

CONJUGATE SOLUTION

Dilute the Concentrated Conjugate (Reagent No. 2) in the ratio of 1:101 in Diluent (Reagent No. 4). Prepare the solution at the moment of performing the test.

VETERINARY USE

To perform an assay using all cavities in the kit, mix 110 µL of Concentrated Conjugate in 11 mL of Diluent.
To perform an assay using 8 cavities (1 strip), mix 10 µL of the Concentrated Conjugate in 1 mL of Diluent.
IMPORTANT: The diluted conjugate solution cannot be stored. Therefore, prepare only the amount necessary to perform the test.

SAMPLE DILUTION

In a test tube, dilute 5 µL of the sample in 500 µL of Diluent (Reagent No. 4), if performing the assay in duplicate. Cap tube and vortex gently or mix manually by inversion. Dilutions cannot be stored.

TECHNIQUE

For use in automatic equipment, consult the Customer Service Department (SAC)
Before starting the assay, place all reagents, Controls and Samples to stabilize at room temperature (15 - 30 °C) for at least 30 minutes. Return the unused strips to the original sealed packaging.
1- Select the cavities to be used considering: Controls and Samples (it is recommended to test in duplicate).
2- Select the first cavity for the Blank (OPTIONAL).
3- Pipette 100 µL of the Negative Control and the Positive Control into the previously determined cavities. **Note: The controls are ready to use, there is no need to dilute them.**
4- Pipette 100 µL of the previously diluted samples into the previously determined cavities. In the Blank cavity, if you have chosen, pipette only 100 µL of the diluent.
5- Homogenize gently for ± 10 seconds. Cover the cavities with plate sealer.
6- Incubate for 30 minutes at 37 °C ± 2 °C.
7- Discard the content of the cavities by aspiration (washing machine) or by decantation (manual). Use approximately 300 µL/ cavity of previously diluted Wash Solution, for a total of five (5) wash cycles. Shake for three seconds after each wash cycle. To guarantee the drying of the plate, at the end of the wash, tap the plate for a few seconds on absorbent paper.
IMPORTANT: Poor washing/drying can cause inadequate results.
8- Pipette 100 µL of the Conjugate previously diluted into each cavity, even in the Blank cavity.
9- Homogenize gently for ± 10 seconds. Cover the cavities with plate sealer.
10- Incubate for 30 minutes at a 37 °C ± 2 °C.
11- Repeat item 7.
12- Add 100 µL of Substrate TMB into each cavity.
13- Homogenize gently for ± 3 seconds. Cover the cavities with plate sealer.
14- Incubate for 10 minutes at a 37 °C ± 2 °C, protected from light.
15- Remove the plate sealer from the cavities.
16- Pipette 100 µL of Stop Solution into all cavities.
17- Homogenize gently for ± 3 seconds
18- Read the absorbance in an ELISA reader on a dual filter 450 nm (primary filter) / 630nm (secondary filter) within a maximum of 10 minutes after adding the Stop Solution.

TECHNIQUE VERIFICATION

Check if the results obtained for reading the Blank and the Controls are compatible with the values presented below:

ITEM	ABSORBANCE (DUAL FILTER)
Blank	< 0.100
Negative Control	0.050 to 0.250
Positive Control	> 1.000

If the values are outside the expected values, the technique should be repeated.

Calculations

Calculate Cut Off according to the following formula:

Cut Off = Average Negative Control Absorbance + 0.400

Example:

ITEM	ABSORBANCE
	A1 = 0.151
Negative Control	A2 = 0.148
Cut Off = Average Negative Control Absorbance + 0,400	$(0.151 + 0.148) / 2 + 0.400 = 0.550$

Calculate the Index by dividing the Sample's absorbance by the Cut Off value.

Example:

ITEM	ABSORBANCIA
Sample	1.236
Cut Off Value	0.550
Índice = Samples / Cut Off Value	$1.236 / 0.550 = 2.25$

Note: The data presented in the examples are for illustration only and cannot be used to calculate the results.

INTERPRETATION OF THE RESULT

After calculating the sample index, consider the indices below to determine the results.

RESULTS	ÍNDICE
Negative	< 0.9
Positive	> 1.1
Undetermined	Between 0.9 e ≤ 1.1

Note: In case of an indeterminate result, the sample must be reanalyzed. Samples that repeatedly obtain indeterminate results should be retested using an alternative method. If results remain indeterminate, a new sample should be collected within two weeks. The result of the last sample collected must prevail. The results provided by this kit must be interpreted by the responsible veterinarian, and is not the only criterion for determining the diagnosis and/or treatment of the patient.

PROCESS LIMITATIONS

The interpretation of a diagnostic test should not be established on the basis of a single test. Other confirmation tests must be included before a sample is considered positive. A negative result does not exclude the possibility of exposure. Finally, all results should be interpreted together with vaccination history, clinical and laboratory information available.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, norms, limits and tolerance for variations are clearly established. It is important to emphasize that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. For that, it is recommended the use of controls, which allow evaluating the precision and accuracy of the measurements.

PRODUCT PERFORMANCE

Accuracy

REPEATABILITY

Repeatability was calculated from 10 successive determinations using 3 samples with different values. The following results were obtained:

REPEATABILITY	Sample		
	1	2	3
Mean	0.550	1.320	0.080
Standard Deviation	0.042	0.053	0.005
Coefficient of Variation (%)	7.623	3.987	6.447

REPRODUTIBILITY

Reproducibility was calculated from 10 successive determinations over 3 consecutive days using 3 samples with different values. The following results were obtained:

REPEATABILITY	Sample		
	1	2	3
Mean	0.529	1.303	0.077
Standard Deviation	0.047	0.063	0.004
Coefficient of Variation (%)	8.937	4.848	5.715

CLINICAL SENSITIVITY AND SPECIFICITY

The VETLISA EHRlichiosis IgG was used to analyze clinical samples in comparison with another method of indirect immunofluorescence. The results show that the clinical sensitivity of the kit is >99.9% and the clinical specificity is 97.3%.

METHOD		REFERENCE		Total
VETLISA EHRlichiosis IgG	Result	Negative	Positive	
	Negative	42	3	45
	Positive	0	107	107
Total		42	110	152

Clinical Sensitivity: >99.9% (42/42)

Clinical Specificity: 97.3% (107/110)

Accuracy: 98% $[(107 + 42) / (110 + 42)]$

CLINICAL SIGNIFICANCE

Ehrlichiosis is a disease caused by bacteria of the *Ehrlichia* genus, which parasitize blood cells such as leukocytes, erythrocytes, endothelial cells and platelets. Transmission occurs through the bite of the brown tick (*Rhipicephalus sanguineus*) that becomes contaminated with the bacteria after the blood meal in dogs with high parasitemia. The incubation period for ehrlichiosis in dogs varies between 8 to 20 days, and the cynical signs observed include petechiae in the skin and mucous membranes, lymph node enlargement, hepatomegaly and splenomegaly, apathy, weight loss, inappetence, anorexia, fever, runny nose and eye. When treatment does not result in total elimination from the body, further relapses may occur. The dog's immune response against *Ehrlichia* is essential to fight the infection and also allows for diagnosis by detecting antibodies in serological tests.

BIBLIOGRAPHIC REFERENCES

- 1- COUTO CG. Doenças Rickettsiais In: BIRCHARD, SHERDING, Manual Saunders: Clínica de pequenos animais. Ed. Roca: 139-42, 1998
- 2- ETtinger SJ, FElDman EC. Tratado de Medicina Interna Veterinária. 4 ed. São Paulo: Editora Manole. 1992.
- 3- ROSEZ KV, ALVES FR, BLEICH I. Erliquiose canina. Revista Cães & Gatos, n.96, p. 25-28. 2001.
- 4- WOLDEHIWET Z, RISTIC M. Rickettsial and chlamydial diseases of domestic animals, 1ed, Oxford, Pergamon Press, 1993, 427 p.
- 5- WOODY BJ, HOSKINS JD. Erlichial Diseases of dogs. Veterinary Clinics of North America: Small Animal Practice. n.1. v. 21, p.129-140. Janeiro 1991.

QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Quality Control Department. The quality of the reagents is ensured until the expiration date mentioned on the presentation packaging, as long as they are stored and transported under the appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes. 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Tel.: (31) 3439.5454

E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 - Industria Brasileira

CONSUMER SERVICE

Customer Support Service

Tel.: 0800 0315454

E-mail: sac@bioclin.com.br

Product licensed in the Ministry of Agriculture. Livestock and Supply since 06/09/2019 with the number 10.269/2019.

Technical manager: Dra. Camila Eckstein (CRMV/MG 20611)

Review: November/2021

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR «>» TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER