VETLISA LEISHMANIOSE FELINA IgG

REF VET047

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Test for qualitative determination of IgG antibodies against *Leishmania infantum* using feline serum or plasma, by enzyme immunoassay, in microplate. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzymatic. This product is an immunoenzymatic assay in solid phase, based on the principle of immunocapture for the qualitative detection of IgG antibodies against *Leishmania infantum* in feline serum or plasma. Antibodies against *L. infantum* present in the sample, bind to the recombinant promastigote antigen of *L. infantum* coated on the microplate, forming immobilized antigen-antibody complexes. After the initial incubation, the microplate is washed to eliminate unbound materials. The Conjugate, formed by the anti-feline IgG antibody conjugated to peroxidase, binds the immobilized antigen-antibody complex on the plate. Then, the microplate is washed and incubated with Substrate. The intensity of the blue color produced by the addition of the Substrate is proportional to the number of antibodies against *Leishmania* present in the sample. The Stop Solution is added to stop the reaction, promoting a color change to yellow, measured in a microplate reader.

REAGENTS

- 1- Sensitized Plate** - Store between 2 and 8 °C. Contains: Polystyrene plate, divided into 12 strips of 8 cavities each, impregnated with recombinant *L. infantum* antigen.
- 2- Concentrated Conjugate (100x)** - Store between 2 and 8 °C. Contains: Peroxidase-linked anti-feline IgG antibody solution.
- 3- Concentrated Wash Solution (20x)** - Store between 2 and 8 °C. Contains: Buffer Solution, surfactant and preservative.
- 4- Diluent - Store between 2 and 8 °C. Contains:** Buffer solution, surfactant, stabilizer and preservative.
- 5- Substrate TMB** - Store between 2 and 8 °C. Contains: Solution containing Tetramethylbenzidine (TMB <1.0mg/mL), Citric Acid Solution <5% and Urea Peroxide <1%.
- 6- Stop Solution** - Store between 2 and 8°C. Contains: 1M Sulfuric Acid.
- 7- Negative Control** - Store between 2 and 8 °C. **Potentially infectious.**
- 8- Positive Control** - Store between 2 and 8 °C. **Potentially infectious.**
- 9- Plate Sealers**

PRESENTATIONS

Componentes	Apresentação			
	1	2	3	4
	96 wells	192 wells	480 wells	960 wells
1- Sensitized Plate	1 unidade (96 wells)	2 unidades (192 wells)	5 unidades (480 wells)	10 unidades (960 wells)
2- Concentrated Conjugate (100X)	1 x 300 µL	1 x 350 µL	1 x 650 µL	1 x 1300 µL
3- Concentrated Wash Solution (20X)	1 x 50 mL	2 x 50 mL	4 x 50 mL	6 x 50 mL
4- Diluent	1 x 60 mL	2 x 60 mL	5 x 60 mL	10 x 60 mL
5- Substrate TMB	1 x 15 mL	2 x 15 mL	5 x 15 mL	10 x 15 mL
6- Stop Solution	1 x 15 mL	2 x 15 mL	5 x 15 mL	10 x 15 mL
7- Negative Control	1 x 1 mL	1 x 2 mL	1 x 4 mL	1 x 6 mL
8- Positive Control	1 x 1 mL	1 x 2 mL	1 x 4 mL	1 x 6 mL
9- Plate Sealers	3 units	6 units	15 units	30 units

EQUIPMENT AND OPERATIONAL INPUTS

Materials contained in the kit:

- Reagents described in the previous item.
- Usage Instructions (manual).

Necessary materials not contained in the kit:

- 1-** Pipettes capable of dispensing volumes of 5 to 500 µL with a variation coefficient of less than 1.5%.
- 2-** Repipettor for repetitive pipetting of volumes of 500 µL with coefficient of variation less than 1.5% or multichannel pipette (optional).
- 3-** Microplate washer (optional).
- 4-** ELISA reader with absorbance capacity at 450 and 630 nm wavelength.
- 5-** Absorbent paper to dry the wells.
- 6-** Stopwatch or watch.
- 7-** Bottle to store the Washing Solution after dilution.
- 8-** Distilled or deionized water.
- 9-** Quality Control Tools.
- 10-** Incubator at 37 °C ± 2 °C.

TRANSPORTATION AND STORAGE CONDITIONS

The product storage temperature should be 2 to 8 °C. Transport can be carried out at room temperature 30 °C for up to 5 days. Keep away from light and avoid humidity.**Do not freeze.**

SPECIAL CARE

- 1-** For professional *in vitro* diagnostic use only.
- 2-** Strictly follow the proposed methodology to obtain exact results.
- 3-** The sachet containing the microplate must be opened only after reaching room temperature. Replace the unused microwell strips in the sachet, seal and store at 2 to 8 °C.
- 4-** The water used to clean the material must be fresh and free from contaminants.
- 5-** Saturated deionizing columns release alkaline water, various ions, and oxidizing and reducing agents that can significantly alter the results.
- 6-** The Stop Solution contains Hydrochloric Acid which is a strong acid. Therefore, handle it with due care.
- 7-** The handling of any product that contains serum is potentially capable of transmitting disease. Therefore, it is necessary to take due care of biosafety in the handling of these products.
- 8-** Pipette the reagents always in the same order to minimize the difference in reaction time between the wells.
- 9-** As a protection measure, the plate must be covered during the reaction.
- 10-** Make sure that the bottom of the cavity is clean and dry and that there are no bubbles on the liquid surface before reading the plate. Do not allow the wells to dry out during the test.
- 11-** Do not expose the reagents, especially the Substrate, to strong light or Hypochlorite vapors during storage or incubation steps.
- 12-** We recommend applying local, state and federal environmental protection standards so that the disposal of reagents and biological material is carried out in accordance with current legislation.
- 13-** To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Safety Information Sheet for Chemical Products) available on the website www.bioclin.com.br or upon request by the SAC (Service of Chemicals) Customer Service) of Quibasa.
- 14-** Do not use the product in case of damage to the packaging.
- 15-** It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or Plasma (EDTA or Heparin)

Hemolyzed or highly lipemic samples should not be used. The samples can be kept refrigerated, between 2 and 8 °C, for a maximum period of 5 days. If the samples cannot be analyzed within 5 days, they can be stored for up to 30 days at a temperature of -20 °C.

PROCESS DESCRIPTION

Stability After Open

The results of the stability test show that the VETLISA Feline Leishmaniasis IgG kit is stable after being opened for at least 30 days. This stability may vary according to the test conditions and the environment. Therefore, it is suggested to monitor the performance of the product using internal controls of the kit and the criteria for validating the technique.

ATTENTION: The Positive and Negative Controls are ready to use.

Preparation of Work Reagents

WASHING SOLUTION

Dilute the contents of vial No. 3 (Concentrated Washing Solution) in 1000 mL of distilled or deionized water. After preparation, the solution can be stored at 2 to 30 °C for at least 30 days. It can be stored at room temperature. If crystallization occurs, heat to 37 °C until dissolved.

SUBSTRATE

The Substrate is ready for use.

CONJUGATE SOLUTION

Dilute the Concentrated Conjugate (Reagent No. 2) in the ratio of 1:101 in Diluent (Reagent No. 4). Prepare the solution at the moment of performing the test. To perform an assay using all cavities in the kit, mix 110 µL of Concentrated Conjugate in 11 mL of Diluent.To perform an assay using 8 cavities (1 strip), mix 10 µL of the Concentrated Conjugate in 1 mL of Diluent. **IMPORTANT:** The diluted conjugate solution cannot be stored. Therefore, prepare only the amount necessary to perform the test.

SAMPLE DILUTION

In a test tube, dilute 15 µL of the sample in 300 µL of Diluent (Reagent No. 4), if the test is going to be performed in duplicate. Cap the tube and gently vortex or manually mix by inversion. Dilutions cannot be stored.

Technique

For use in automatic equipment, consult the Customer Service Department (SAC)

- Before starting the assay, place all reagents, Controls and Samples to stabilize at room temperature (15 - 30 °C) for at least 30 minutes. Return the unused strips to the original sealed packaging.
- 1-** Select the cavities to be used considering: Controls and Samples (it is recommended to test in duplicate).
- 2-** Select the first cavity for the Blank (OPTIONAL).
- 3-** Pipette 100 µL of the Negative Control and the Positive Control into the previously determined cavities. **Note: The controls are ready to use, there is no need to dilute them.**
- 4-** Pipette 100 µL of the previously diluted samples into the previously determined cavities. In the Blank cavity (OPTIONAL), if you have chosen, pipette only 100 µL of the diluent.
- 5-** Homogenize gently for ± 10 seconds. Cover the cavities with plate sealer.
- 6-** Incubate for 30 minu tes at 37 °C ± 2 °C.
- 7-** Discard the content of the cavities by aspiration (washing machine) or by decantation (manual). Use approximately 300 µL/ cavity of previously diluted Wash Solution, for a total of five (5) wash cycles. Shake for three seconds after each wash cycle. To guarantee the drying of the plate, at the end of the wash, tap the plate for a few seconds on absorbent paper.
- IMPORTANT:** Poor washing/drying can cause inadequate results.
- 8-** Pipette 100 µL of the Conjugate previously diluted into each cavity, even in the Blank cavity.
- 9-** Homogenize gently for ± 10 seconds. Cover the cavities with plate sealer.
- 10-** Incubate for 30 minutes at a 37 °C ± 2 °C.
- 11-** Repeat item 7.
- 12-** Add 100 µL of Substrate TMB into each cavity.
- 13-** Homogenize gently for ± 3 seconds. Cover the cavities with plate sealer.
- 14-** Incubate for 10 minutes at a 37 °C ± 2 °C, protected from light.
- 15-** Remove the plate sealer from the cavities.
- 16-** Pipette 100 µL of Stop Solution into all cavities.
- 17-** Homogenize gently for ± 3 seconds
- 18-** Read the absorbance in an ELISA reader on a dual filter 450 nm (primary filter) / 630nm (secondary filter) within a maximum of 10 minutes after adding the Stop Solution.

VETERINARY USE

Technique Verification

Check if the results obtained for reading the Blank and the Controls are compatible with the values presented below:

ITEM	ABSORBANCE (DUAL FILTER)
Blank	< 0,100
Negative Control	0,050 to 0,250
Positive Control	> 1,000

If the values are outside the expected values, the technique should be repeated.

Calculations

Calculate Cut Off according to the following formula:

Cut Off = Average Negative Control Absorbance + 0.350

Example:

ITEM	ABSORBÂNCIA
Negative Control	A1 = 0,107
	A2 = 0,112
Cut Off = Average Negative Control Absorbance 0.350	(0,107 + 0,112) / 2 + 0,350 = 0,459

Calculate the Index by dividing the Sample's absorbance by the Cut Off value.

Example:

ITEM	ABSORBÂNCIA
Sample	1,093
Cut Off Value	0,459
Index = Samples / Cut Off Value	1,093 / 0,459 = 2,381

Note: The data presented in the examples are for illustration only and cannot be used to calculate the results.

INTERPRETATION OF THE RESULTS

RESULTS	INDEX
Negative	< 0.9
Positive	> 1.1
Undetermined	Between 0,9 - 1,1

Note: In the case of an undetermined result, the sample must be re-analyzed. Samples that obtain repeatedly indeterminate results should be retested using an alternative method. If the results remain indeterminate, a new sample should be collected in two weeks. The result of the last collected sample must prevail. The results provided by this kit must be interpreted by the responsible veterinarian medical professional, and it is not the only criterion for determining the diagnosis and / or treatment of the patient.

PROCEDURE LIMITATIONS

The interpretation of a diagnostic test should not be established on the basis of a single test. Other confirmatory tests must be included before a sample is considered positive. A negative result does not exclude the possibility of exposure. All results must be interpreted in conjunction with other available clinical information, before the descriptive diagnosis of the disease.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to note that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. For that, it is recommended to use controls, which allow to evaluate the precision and accuracy of the dosages.

PRODUCT PERFORMANCE

Accuracy
REPEATABILITY
Repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

REPEATABILITY	Sample		
	1	2	3
Mean	0.311	1.028	1.713
Standard Deviation	0.034	0.103	0.172
Coefficient of Variation (%)	10.95%	10.01%	10.02%

REPRODUTIBILITY
Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

REPRODUCTIBILITY	Sample		
	1	2	3
Mean	0.303	1.057	1.681
Standard Deviation	0.031	0.082	0.146
Coefficient of Variation (%)	10.16	7.76	8.66

Clinical Sensitivity and Specificity
The VETLISA FELINE LEISHMANIASIS IgG was used to analyze clinical samples in comparison with another method of indirect immunofluorescence. The results show that the clinical sensitivity of the kit is 97.6% and the clinical specificity is >99%.

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Clinical Sensitivity: 97.6% (40/41)
Clinical Specificity:> 99% (132/132)
Accuracy: 99.4% [(132 + 40) / (132 + 41)]

CLINICAL SIGNIFICANCE
Leishmaniasis is an infectious and zoonotic disease caused by the protozoan Leishmania infantum, which parasites cells of the phagocytic system. Transmission occurs through the bite of the vector Lutzomyia longipalpis, also popularly known as the straw mosquito. The increasing number of cases of feline leishmaniasis has drawn attention to the role of felines in the epidemiology of the disease. The incubation period for leishmaniasis is long, with an average of 3 to 7 months, but it can vary from 3 months to several years. The main clinical manifestations are generalized lymphadenopathy, exfoliative dermatitis, progressive weight loss and keratoconjunctivitis. In cats, animals are commonly asymptomatic and bone marrow is an important reservoir of the parasite. Often, clinical signs are mistaken for sporotrichosis and the differential diagnosis must be made. The serological diagnosis by the detection of antibodies is an effective and safe tool. In addition, the diagnosis of leishmaniasis is essential to control the reservoirs and prevent transmission to humans.

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QUALITY ASSURANCE
Before being released for consumption, all Bioclin reagents are tested by the Quality Control Department. The quality of the reagents is ensured until the expiration date mentioned on the presentation packaging, as long as they are stored and transported under the appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda**
Rua Teles de Menezes. 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Tel.: (31) 3439.5454
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Industria Brasileira

CONSUMER SERVICE
Customer Support Service
Tel.: 0800 0315454
E-mail: sac@bioclin.com.br

Product licensed in the Ministry of Agriculture, Livestock and Supply since 03/09/2020 with the number 10.353/2020.

Technical manager: Dra. Camila Eckstein (CRMV/MG 20611)

Review: September/2021

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR «>» TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER