FINALIDADE

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PURPOSE

Rapid immunochromatographic test for qualitative determination of canine parvovirus antigen in samples of feces.

PRINCIPLE OF ACTION

Methodology: Immunochromatography.
Parvovirus Ag VET FAST is an *in vitro* diagnostic test for rapid and qualitative detection of parvovirus antigens in dog feces. The cassette consists of a nitrocellulose membrane, impregnated in the test line with antibody anti-canine parvovirus. The conjugate consists of antibody anti-canine parvovirus, conjugated with colloidal gold particles. When the sample is added to the cassette, it migrates by capillary action until it reaches the conjugate. In this region, the antigens present in the sample interact with the conjugate, forming the complex: antigen – conjugate.
The sample continues to migrate through the immunochromatographic strip until reaching the test line, where the parvovirus antigens also interact with the antibody immobilized in the test line, forming the complex: immobilized antibody - antigen - conjugate. This shows the formation of a pink-red line on the test line; indicating positive result. Sample excess continues to migrate through the test strip until it reaches the control line, which is composed of impregnated anti-conjugate antibody. Thus, the formation of a second line (control line) is observed; indicating that the test yielded a valid result.

REAGENTS

Reagent Number 1 (R1): Cassette - Store between 15 and 30°C. **Do not freeze.** It contains:

- A)** 01 plastic cassette with 01 test strip containing:
- Control Region (C): Antibody anti-mouse antibody immobilized.
- Test Region (T): Immobilized monoclonal anti-canine parvovirus antibody.
- Conjugate: Monoclonal mouse antibody anti-canine parvovirus, colloidal gold conjugate.

B) 01 sachet of silica

Reagent Number 2 (R2): Diluent - Store between 15 and 30°C. **Do not freeze.** Contains: Buffer solution and preservative.

Reagent 3 (R3): Swab - Sterile flexible cotton swab.

Collection Accessory:

- Disposable plastic pipette.

PRESENTATIONS

Reagent	Presentation	Number of Cassettes for Package	Diluent
Cassette + Diluent	VET007-1	10	10 x 1,0 mL
	VET007-2	25	25 x 1,0 mL
	VET007-3*	10	10 x 1,0 mL
	VET007-4*	25	25 x 1,0 mL
	VET007-5	5	5 x 1,0mL
	VET007-6	1	1 x 1,0mL

* Individualized presentation.

EQUIPMENTS AND OPERATIONAL INPUTS

Clock or stopwatch. They are in the specialized market of articles for Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperature of the product must be between 15 and 30°C. Keep away from direct sunlight. **Do not freeze.** The cassette is sensitive to air humidity. Therefore, it should always be kept in its inviolate primary packaging (aluminized sachet) until the time of use.

SPECIAL CARE


- 1- Product for *in vitro* diagnosis, for veterinary use only.**
- 2- Strictly follow the methodology proposed for obtaining results.
- 3- Cassettes (test devices) do not contain contaminants; but since the samples are biological, the professional must strictly follow the safety rules and routines. The use of disposable gloves and other personal protective equipment is essential.
- 4- Use clean, dry containers to package biological samples.
- 5- Open the aluminum envelop to remove the cassette only at the time of use.
- 6- The kit cannot be frozen or exposed to high temperatures, since this will cause deterioration of its components.
- 7- The samples should reach room temperature before starting the test.
- 8- Do not use reagents from different lots or kits.
- 9- Do not use kit beyond its expiration date. The expiry date is stated on the product label.
- 10- Kit components are disposable and must not be reused.
- 11- The kit materials that had contact with the sample of the patient must have a special destination. The disposal of sterile flexible cotton swab, the plastic pipettes, the cassettes (test devices) and diluent vials should be disposed of in organic waste.
- 12- Secondary packaging (box). Instruction for use, envelope aluminum and silica sachet can be disposed of in trash common.
- 13- We recommend to apply the local, state and federal environmental protection so that the disposal of kit components and biological material is done in accordance with the current legislation.
- 14- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or through request by the Customer Service Department (SAC) of the Quibasa.

SAMPLES

Perform the test with a sample of dog feces.
Test the sample immediately after collection.
If necessary, samples can be stored between 2 to 8°C for up to 2 days or at -20°C for longer periods.

DESCRIPTION OF PROCESS

- 1- The sample should be at room temperature before starting the test.
- 2- Using the swab, collect a sample from 4 different points of feces.
Note: In samples with liquid or pasty consistency, ensure that excess feces are not absorbed by the swab. Sampling directly from the animal's rectum must be carried out only in the absence of a sample in a collection bottle and attention must be paid to the amount of feces sampled, which when insufficient can impair the accuracy of the results.
- 3- Insert the swab into the diluent tube, dissolving the sample.
Note: Use rotating swab movements inside the tube and rubbing it against the wall for 10 seconds to maximize sample dissolution in the diluent.
- 4- Wait for the larger particles to settle at the bottom of the tube (we recommended 5 minutes of decanting).
Note: Excess particulate material in the sample diluent supernatant can compromise the run.
- 5- Use the pipette that comes with the cassette to collect the dissolved sample.
- 6- Add 4 drops of sample (80 µL), slowly, into the sample application well of the cassette.
- 7- Wait for the formation of the lines. Interpret the results between 5 and 10 minutes. Do not interpret after 10 minutes.

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Reactive Test: Formation of two lines: one in the control region (C) and one in the test region (T). Do not interpret after 10 minutes.	
Not Reactive Test: Line formation only in control region (C). Do not interpret after 10 minutes.	
Invalid Test: The absence of line formation in the control region (C), with or without line formation in the test region (T), indicates procedure error or deterioration of the cassette. In this case, repeat the test using a new cassette.	

TEST LIMITATIONS

This product is intended to detect parvovirus infections in dog. A negative result can be observed if the amount of analyte present in the sample is below the detection limit of the test. This can happen, depending on disease period in which the sample was collected, the type of sample used in the test and its storage conditions. A negative result may also occur in case of inefficiency in the sample extraction step in the diluent.
Therefore, for the diagnosis of the disease, the professional must take into account not only the result of this test; but also, clinical signs, animal health conditions, and other diagnostic results that may be available.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal program of quality control, where procedures, standards, limits and tolerance for variations are clearly established. It is important to note that all measuring systems characteristic analytical variability, which should be monitored by the laboratories themselves. Therefore, it is the use of controls, which allow the assessment of the precision and accuracy of dosages.

PRODUCT PERFORMANCE
QUALITY CONTROL

CLINICAL SENSITIVITY AND SPECIFICITY

This product has been tested comparing with other methods. A total of 106 samples characterized by parvovirus were analyzed; 47 samples were diagnosed as positive and 59 were diagnosed as negative. The results show that the clinical sensitivity of the product is >99% and the clinical specificity is 94,91%.

METHOD		REFERENCE		Total
Cassette Parvovirus Ag	Result	Reactive	Not reactive	
	Reactive	47	3	50
	Not reactive	0	56	56
Total result		47	59	106

Clinical sensitivity: (47/47): >99%
Clinical Specificity: (56/59): 94,91%
Precision: (103/106): 97,17%

ANALYTICAL SENSITIVITY

The product was tested with samples of feces, which contained a determined canine parvovirus load. The product was able to identify with positive result all the concentrations tested with a value equal to or greater than 10 HAU/mL. Therefore, the analytical sensitivity (or detection limit) of 10 HAU/mL was established.

REPEATABILITY

Repeatability was calculated from 10 successive determinations, using 2 different samples, obtaining the following results:

Samples	N° of Repetitions	Expected Result	Found Result
Samples 1	10	Reactive	Reactive
Samples 2	10	Reactive	Reactive

REPRODUCIBILITY

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 2 different samples, obtaining the following results:

Samples	N° de Days	N° of Repetitions	Expected Result	Found Result
Samples 1	03	10	Reactive	Reactive
Samples 2	03	10	Reactive	Reactive

CLINICAL SIGNIFICANCE

Canine parvovirus is one of the main causes of dog mortality. The disease is caused by canine parvovirus, which has high rates of morbidity and mortality. Puppies are the most affected, due to the immaturity of the immune system and flaws in vaccine protocols. The incubation period varies between 3 to 7 days and the main clinical changes observed are intense vomiting and pasty and fetid diarrhea, containing variable amounts of blood, fever, lethargy, tachycardia, dehydration, weight loss and pale mucous membranes. Transmission can occur through direct contact with the sick dog or contaminated secretions and objects. Due to the rapid course of the disease, detection of viral antigen favors early diagnosis even before the appearance of clinical signs.

BIBLIOGRAPHIC REFERENCES

- 1- DECARO N, BUONAVOGLIA C. Canine parvovirus--a review of epidemiological and diagnostic aspects, with emphasis on type 2c. Vet Microbiol. 2012;155:1-12.
2- PRITTIE J. Canine parvoviral enteritis: a review of diagnosis, management, and prevention: J Vet Emerg Crit Care 2004;4:167-176.
3- MIRANDA C, THOMPSON G. Canine parvovirus: the worldwide occurrence of antigenic variants. Journal of General Virology, v.97, p.2043-2057, 2016.
4- QUIBASA: Data from the Research and Development Department.

QUALITY GUARANTEE

Before being released for consumption, all Quibasa reagents are tested by the Quality Control Department. The quality of the reagents is assured until the validity date mentioned on the packaged, as long as they are stored and transported under the appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Tel.: (31) 3439.5454
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Industria Brasileira

CONSUMER SERVICE

Customer Support Service
Tel.: 0800 0315454
E-mail: sac@bioclin.com.br

Product licensed in the Ministry of Agriculture, Livestock and Supply since 08/11/2019 with the number 10.279/2019.

Technical manager: Dra. Camila Eckstein (CRMV/MG 20611)

Review: September/2021

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MANUFACTURED BY

BATCH CODE

CONTROL

DATE OF MANUFACTURE

POSITIVE CONTROL

USED BY
(last day of month)

NEGATIVE CONTROL

TEMPERATURE LIMITATION
(store at)

BIOLOGICAL RISK

CONTAINS SUFFICIENT
FOR <N> TESTS

INFLAMMABLE

CONSULT INSTRUCTIONS
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC DEVICE

POISON

EUROPEAN AUTHORIZED
REPRESENTATIVE

CE MARK

KEEP AWAY
FROM SUNLIGHT

DO NOT USE IF
PACKAGE IS
DAMAGED