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
INTERPRETATION OF THE RESULT


IgM Reagent Test: Formation of a red line on the test region M (1) and another line in the control region (C) in first 10 to 15 minutes. Do not interpret after 15 minutes.

IgG Reagent Test: Formation of a red line on the test region G (2) and another line in the control region (C) in the first 10 to 15 minutes. Do not interpret after 15 minutes.


IgM and IgG Reagent Test: Formation of two red lines in the test region (M = 1 and G = 2) and a red line in the control (C) in the first 10 to 15 minutes. Do not interpret after 15 minutes.

Non-Reagent Test: Formation of a red line in the region control (C) and complete absence of a red line in the regions tests (M = 1 and G = 2). Do not interpret after 15 minutes.

Invalid Result: Complete absence of a line in the control region (C) with or without a red line in the test regions (M = 1 and G = 2). Test the sample again.

Notes:

- 1- The line in the control region can appear before 10 minutes incubation of the test, and this does not mean that results do not reagents can be interpreted ahead of time.
- 2- The results should not be interpreted after 15 minutes to avoid incorrect results.
- 3- The results provided by this kit must be interpreted responsible medical professional, and it is not the only criterion for determining the diagnosis and/or treatment of the patient.

CROSS REACTIVITY

1349 samples negative for COVID-19 and positive for influenza, rinovirus and/or respiratory syncytial virus, were subjected to the cross-reactivity test. 2 false positive results were observed for influenza, 2 for rinovirus and 2 for respiratory syncytial virus. Despite the results found, the possibility of other cross-reactivities cannot be completely ruled out. The final diagnosis should consider the patient's clinical data along with other laboratory data.

REPRODUCIBILITY

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 5 different samples, obtaining the following results:

Samples	Nº of Days	Repetition Nº	Expected Result	Found Result
Human Blood 01	03	10	100% Non-Reagent	100% Non-Reagent
Human Blood 02	03	10	100% Reagent IgM	100% Reagent IgM
Human Serum 01	03	10	100% Reagent IgM/IgG	100% Reagent IgM/IgG
Human Serum 02	03	10	100% Non-Reagent	100% Non-Reagent
Human Plasma 01	03	10	100% Reagent IgM/IgG	100% Reagent IgM/IgG

DIAGNOSTIC SIGNIFICANCE

Coronaviruses (CoV) belong to the Coronaviridae family, and are widely distributed infecting humans and other mammals. In humans, the most common symptoms presented are fever, cough, dyspnoea, myalgia or fatigue. It can cause diseases ranging from the common cold to more serious illnesses, such as the Middle East Respiratory Syndrome (MERSCoV) and the Severe Acute Respiratory Syndrome (SARS-CoV). In December 2019, a series of pneumonia cases of unknown cause appeared in Wuhan in China, with clinical symptoms very similar to viral pneumonia, after complete sequencing of the respiratory samples it was shown to be Coronavirus infection, which was called the new coronavirus (2019-nCoV). In approximately 150 countries on all continents, cases of infection have already been reported.

BIBLIOGRAPHIC REFERENCES

- 1- CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
- 2- Corman VM, Landt O, *et al.* Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
- 3- Huang C, Wang Y, *et al.* Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb 15;395(10223):497-506.
- 4- ANVISA. RDC N°107, DE 5 DE SETEMBRO DE 2016. DOU n°172, 06 de setembro de 2016. Ministério da Saúde.
- 5- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.