Bioclin

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[REF] K136

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FINALIDADE

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- Röhle G, Siekmann L. Quality assurance of quantitative determination. In: Thomas L, editor. Clinical laboratory diagnostics. 1st ed. Frankfurt:TH-Books Verlagsgesellschaft; 1998. p. 1393-1401.
- Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).
- National Committee For Clinical Laboratory Standards: Colecion, Transport And Processing Of Blood Specimens For Coagulation Test And Performance Of Coagulation Assays. Approved Guideline Fourth Edition. NCCLS document H12-A4, 2003.
- International Committee For Standardization In Haematology And International Committee On Trombosis And Haemostasis . Amer. J. Clin. Path . 88:779 1985.
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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- Röhle G, Siekmann L. Quality assurance of quantitative determination. In: Thomas L, editor. Clinical laboratory diagnostics. 1st ed. Frankfurt:TH-Books Verlagsgesellschaft; 1998. p. 1393-1401.
- Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).
- National Committee For Clinical Laboratory Standards: Coleccion, Transport And Processing Of Blood Specimens For Coagulation Test And Performance Of Coagulation Assays. Approved Guideline Fourth Edition. NCCLS document H12-A4, 2003.
- International Committee For Standardization In Haematology And International Committee On Trombosis And Haemostasis . Amer. J. Clin. Path . 88:779 1985.
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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USAGE INSTRUCTIONS

FUNCTION

Method for quantitative determination of Fibrinogen in biological samples of citrated plasma. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Coagulometry

Consists on the quantitative determination of fibrinogen in citrated plasma, after the addition of a standard quantity of thrombin in a diluted citrated plasma sample and measurement of the time of coagulation. When the concentration of the thrombin is high, the coagulation time of the diluted citrated plasma is inversely proportional to the fibrinogen concentration. To achieve the fibrinogen concentration, the coagulation time of analyzed plasma is compared with the time of coagulation of a plasma with known concentration of fibrinogen.

REAGENT

Number 1 - Buffer - Store between 2 and 8°C. Contains: Imidazole 50 mmol/L, Sodium Azide 0.05% (Slightly turbid aspect).

Number 2 - Thrombin - Store between 2 and 8°C. Contains: lyophilized Thrombin 100 units NIH/mL, Albumin 0.9%, TRIS Buffer 123 mmol/L.

Number 3 - Reference Plasma - Store between 2 and 8°C. Contains: Buffer HEPES 25 mmol/L, Plasma lyophilized. **Potentially infectious.**

ATT: refer to the concentration on the flask label.

PRESENTATION

Presentation	Reagent 1	Reagent 2	Reagent 3
1	1 x 60 mL	5 x 1 mL	1 mL
2	2 x 60 mL	10 x 1 mL	1 mL
3	4 x 60 mL	20 x 1 mL	1 mL
4	60 mL	5 mL	1 mL
5	2 x 60 mL	1 x 10 mL	1 mL
6	2 x 60 mL	2 x 5 mL	1 mL
7	4 x 60 mL	2 x 10 mL	1 mL
8	4 x 60 mL	4 x 5 mL	1 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Coagulometer, water bath at 37°C, watch or stopwatch, pipettes, test tubes, Biocontrol N Coagulation and Biocontrol P Coagulation. These items are found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be 2 to 8°C. Transportation, at temperatures up to 30°C, should not exceed 5 days. Protect from light and avoid moisture.

SPECIAL CARE

1- For *in vitro* diagnostic use only.

2- Follow strictly the proposed methodology to obtain results accurate.

3- The water used in material cleansing must be recent and free of contaminants.

4- The reagent must be handle cautiously because it is **potentially infectious.**

5- The water level on water bath must be superior to the level of the reagents in the test tubes.

6- The reagent N° 1 contain Sodium Azide, irritating to mucous and skin. Handle with care.

7- Avoid the heating of the **reagent N° 2.**

8- The reagent was tested for anti-HIV, anti-HCV antibodies and antigen HBs using methods of recent latest generation and present negative results. The risk of infection can not be excluded and the reagent must be handle with the same caution taken for the patient sera. **Potentially infectious.**

9- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

10- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Advisory Service Customer) of Quibasa.

11- Do not use the product in case of damaged packaging.

SAMPLES

Material: Plasma

Anticoagulant: Use as anticoagulant citrate Bioclin, as written in instructions manual, or 3.8% Sodium Citrate (0.130 M) at a ratio of 9 parts blood to 1 part anticoagulant (Example: 4.5 mL blood + 0.5 mL of anticoagulant). **Do not EDTA or Heparin.**

Collection: Blood should be obtained by venipuncture, avoiding hemolysis, prolonged tourniquet, blistering and aspiration tissue fluid.

Preparation: Blood should be mixed with the anticoagulant soon after collection.

Centrifuge immediately at 3000 rpm for 15 minutes.

Remove the plasma without cellular pipette red or yellow layer.

The samples should be tested in less than 4 hours. If the test can not be done in this period the plasma should be frozen for at maximum 2 week at -20°C.

Do not use clotted samples.

PROCESS DESCRIPTION
TECHNIQUE

Bioclin recommends, for use of the kit, to use the Biocontrol N Coagulation and Biocontrol P Coagulation kits as the control serum.

Carry out tests and controls in duplicate.

Reagent 1, Buffer, was developed to facilitate the detection of clots by coagulometers with an optical detection system.

Homogenize before use.

Preparation of the Reference Plasma:

Carefully open the flask and add 1 mL of deionized water. Homogenized. Let the flask resting for 30 minutes. Each 10 minutes, agitate carefully the flask with circular moves, to avoid the adherence of the material to the flask's walls. Avoid foam formation. The reconstituted Control is stable for 8 hours at 2 to 8°C. After reconstitution, freezing is not recommended.

Dilution:

1- Prepare the dilutions to assemble the **calibration curve** according to the table 1 and carry out the technique from item 3 on, considering each dilution as sample.

2- Dilute **sample and control** 1:10 (0.1 mL of the sample + 0.9 mL of R1 Buffer). Carry out the technique from item 3 on.

Procedure:

3- Incubate 0.2 mL of the diluted sample in water bath at 37°C during 4 to 6 minutes

4- Add 0.1 mL of the R2 (Thrombin) **not heated**, simultaneously start the chronometer and register the time of clot formation.

CALIBRATION CURVE

Measure in duplicate the time of coagulation of each dilution of the R3 (Reference Plasma) and obtain the average of the time.

Plot a graph in a log-log applying the concentration values on the axis x (abscissa) and the coagulation time on the axis Y (ordinate).

After knowing the patients coagulation time and control, place the result in the calibration curve and read the value of the fibrinogen for the time found.

Table 1 - Calibration Curve

R3 - Reference Plasma		R1 Buffer (mL)	Dilution Factor	Fibrinogen = R3 Concentration (mg/dL) x Dilution Factor	Time (sec.)
Dilution	Volume (mL)				
1:5	0,1	0,4	10/5 = 2.0	-- X 2.0	--
1:10	0,1	0,9	10/10 = 1.0	-- X 1.0	--
1:15	0,1	1,4	10/15 = 0.67	-- X 0.67	--
1:20	0,1	1,9	10/20 = 0.5	-- X 0.5	--
1:30	0,1	2,9	10/30 = 0.33	-- X 0.33	--

The R3 - Reference Plasma diluted 1:10 represents 100% of the value indicated on the flask's label, which means that the value indicated on the related label to the plasma concentration already diluted 1:10.

RESULTS

If the coagulation time of the sample is lower than the coagulation time of the dilution 1:5 from the curve, the test shall be repeated diluting the sample 1:30, read the result at the calibration curve and multiply it by 3.

If the coagulation time of the sample is higher than the coagulation time of the dilution 1:30 from the curve, the test shall be repeated diluting the sample 1:3, read the result at the calibration curve and multiply it by 0.3.

If the coagulation time of the dilution 1:3 is higher than the higher time of coagulation from the curve, the result shall be reported as lower than the lower concentration from the calibration curve.

INTERFERENCES

Samples with icterus, hemolyzed and lipemia.

Use of plasma with Heparin and EDTA or Sodium Oxalate.

High levels of products with Fibrinogen degradation or fibrin can extend the coagulation time specially when the fibrinogen levels are lower than 150 mg/dL.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The Plasma Reference of Fibrinogen kit is traceable to the reference material NIBSC 98/612 (International Standard for Fibrinogen Plasma - WHO International Standard).

REFERENCE VALUES

The interval shall be used only as a guideline. We recommend that each laboratory establishes, in the watched population, their own reference intervals, using a normal plasmas pool or a reference plasma.

Fibrinogen: 150 to 350 mg/dL .

Conversion: conventional units (mg/dL) x 0.01 = Units SI (g/L)

PRODUCT PERFORMANCE
ACCURACY

Comparison of Methods

The Fibrinogen Kit was compared with another method for evaluating Fibrinogen commercially available. Were performed 07 analysis and the results were evaluated. The linear equation obtained was $Y = 1.001X - 0.408$ and the coefficient of correlation 0.999. With these results, we can conclude that the kit has good methodological specificity.

PRECISION

Repeatability

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	170.40	257.55	332.40
Standard Deviation (mg/dL)	0.60	0.69	0.68
Coefficient of Variation (%)	0.35	0.27	0.20

Reproducibility

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	170.40	257.55	332.35
Standard Deviation (mg/dL)	0.05	0.10	0.05
Coefficient of Variation (%)	0.03	0.04	0.02

SENSITIVITY

The sensitivity of the kit was determined by the last calibration point (1:30).

DIAGNOSTIC SIGNIFICANCE

Fibrinogen is a glycoprotein present in plasma and platelet granules, which is the coagulation factor (Factor I) found in the highest plasma concentration (150-350 mg/dL).

Fibrinogen originates in the liver and is converted to fibrin during coagulation. Because fibrin is required for clot formation, fibrinogen deficiency can produce mild to severe bleeding disorders.

Also known as acute phase protein, fibrinogen can occur at high levels in patients with acute infection, collagen disorders, neoplasms, nephrosis or hepatitis. High levels may also occur during the third trimester of pregnancy and in women taking oral contraceptives.

Depressed fibrinogen levels may indicate congenital afibrinogenemia, hypofibrinogenemia, DIC, fibrinolysis, severe liver disease, metastatic cancer, or bone marrow lesions. Trauma or obstetric complications can also cause low levels.

Fibrinogen levels below 100 mg/dL prevent an interpretation of coagulation tests that have a fibrin clot as an end point.

BIBLIOGRAPHIC REFERENCES

- Röhle G, Siekmann L. Quality assurance of quantitative determination. In: Thomas L, editor. Clinical laboratory diagnostics. 1st ed. Frankfurt:TH-Books Verlagsgesellschaft; 1998. p. 1393-1401.
- Biosafety in Microbiological and Biomedical. Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).
- National Committee For Clinical Laboratory Standards: Colecion, Transport And Processing Of Blood Specimens For Coagulation Test And Performance Of Coagulation Assays. Approved Guideline Fourth Edition. NCCLS document H12-A4, 2003.
- International Committee For Standardization In Haematology And International Committee On Trombosis And Haemostasis . Amer. J. Clin. Path . 88:779 1985.
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

 QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA resgistration for the Fibrinogen Kit: 10269360218

Review: May/2022

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER