TROPONINA I

REF K085

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FINALIDADE

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- Mair J., Genser N., Morandell D., et al. Cardiac troponin-I in the diagnosis of myocardial injury and infarction. Clin. Chem. Acta 1996;245:19-38;
- Vallins W.J., Brand N.J., Dabhade N., et al. Molecular cloning of human cardiac troponin-I using polymerase chain reaction. FEBS Lett 1990;270:57-61;
- Bodor G.S., Porter S., Landt Y., et al. Development of monoclonal antibodies for an assay of cardiac troponin-I and preliminary results in suspected cases of myocardial infarction. Clin. Chem. 1992;38:2203-14;
- Zaninotto M., Altinier S., Lachin M., et al. Fluoroenzymometric method to measure cardiac troponin-I in the sera of patients with myocardial infarction. Clin. Chem. 1996;42:1460-6;
- OH S., Foster K., Datta P., et al. Use of a dual monoclonal solid phase and polyclonal detector to create an immunoassay for the detection of human cardiac troponin-I. Clinical Biochemistry 2000;33(4);255-262.
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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1. Mair J., Genser N., Morandell D., et al. Cardiac troponin-I in the diagnosis of myocardial injury and infarction. Clin. Chem. Acta 1996;245:19-38;
2. Vallins W.J., Brand N.J., Dabhade N., et al. Molecular cloning of human cardiac troponin-I using polymerase chain reaction. FEBS Lett 1990;270:57-61;
3. Bodor G.S., Porter S., Landt Y., et al. Development of monoclonal antibodies for an assay of cardiac troponin-I and preliminary results in suspected cases of myocardial infarction. Clin. Chem. 1992;38:2203-14;
4. Zaninotto M., Altinier S., Lachin M., et al. Fluoroenzymometric method to measure cardiac troponin-I in the sera of patients with myocardial infarction. Clin. Chem. 1996;42:1460-6;
5. OH S., Foster K., Datta P., et al. Use of a dual monoclonal solid phase and polyclonal detector to create an immunoassay for the detection of human cardiac troponin-I. Clinical Biochemistry 2000;33(4);255-262.
6. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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USAGE INSTRUCTIONS

FUNCTION

Immunochromatographic method for rapid determination and qualitative cTnI (Human Cardiac Troponin I), in biological samples of total blood, serum or plasma (Heparim or EDTA). For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Immunochromatography

The method uses Monoclonal Antibody mouse anti-cTnI, which reacts with Troponin found in samples of serum, plasma or whole blood. The samples migrate by capillarity through the chromatographic membrane. Positive samples cTnI will form a red line in the region where the Monoclonal Antibodies anti-cTnI are immobilized. The samples remain absorbed by the membrane to the region of the Control Antibody, as formation of another line, confirming the correct processing of test.

REAGENTS:

Number 1 (R1): Cassette - Store between 2 and 30°C. **Do not freeze.** Contains:

- A) 01 plastic (device) cassette with 01 test strip containing:
- Control Region: Polyclonal Antibody goat anti-IgG from mouse.
 - Region Test: Monoclonal antibodies mouse anti-cTnI.
 - Conjugate: Region B Monoclonal antibodies mouse anti-cTnI conjugated with Colloidal Gold.

- B) 01 sachet of silica;

PRESENTATION

The kit includes the following presentations:

Cassette per Package
1
5
10
20
25
30
40
50
100

EQUIPMENTS AND OPERATIONAL INPUTS

Watch or stopwatch. They can be found at markets specialized in laboratories of clinical analysis.

TRANSPORTATION AND STORAGE CONDITIONS

Transport can be carried out at temperatures up to 45 °C for a maximum of 3 weeks. The storage temperature should be between 2 and 30 °C. Avoid humidity. **Do not freezer.**

SPECIAL CARE

- 1- **For *in vitro* diagnostic use only.**
- 2- Strictly follow the methodology proposed to obtain accurate results.
- 3- The cassettes (device) do not show contaminants, but how the samples are biological, safety routines should be followed accurately and the use of disposable gloves is essential.
- 4- Use clean and dry containers for storage of samples.
- 5- Open the envelope of aluminum, to remove the cassette (device), only in time of use.
- 6- The kit should not be frozen or exposed to high temperatures, as this will cause deterioration.
- 7- The reagents and the sample must be equilibrated environment before starting the test.
- 8- It is interesting to conduct the test immediately after sample collection.
- 9- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 10- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 11- Do not use the product in case of damaged packaging.
- 12- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Total Blood: Collect blood by venopuncture or with the aid of a lancet, using the appropriate anti-coagulant (Heparin or EDTA). Total blood can be used to immediate testing or may be stored between 2 and 8°C for up to two days. **Do not freeze total blood.**

Serum/Plasma: Samples of serum and plasma are stable for 2 days when kept between 2 and 8°C and for 03 months at -20°C and should be thawed for use only one time. Hemolyzed and lipemic samples should be avoided because they may cause a false positive result. Sera and plasmas containing precipitate may give inconsistent results.

Prolonged storage or high temperature can deteriorate the cTnI protein. To obtain the plasma, collect Heparin or EDTA as anticoagulant.


PROCESS DESCRIPTION

The stability of the Troponin I cassette (after opening the aluminized sachet) is 2 hours. This stability may vary according to the environmental conditions after opening the aluminized sachet.

TECHNIQUE

- 1- The sample should be at room ambient temperature before starting the test.
- 2- Remove the cassette from its protective packaging and appropriately identify it so.


- 3- Add 100 µL of serum, plasma or whole blood on the center of the application window of the cassette.


- 4- Wait for the formation of lines. Interpret results between 15 and 30 minutes. Do not take a reading test after 30 minutes of reaction.

RESULTS INTERPRETATION



Reagent Test: Formation of two lines, one in the control region (C) and one in the test region (T) between 15 and 30 minutes. Do not interpret after 30 minutes. The appearance of two lines indicates the presence of cTnI in the sample. The color of the line increases with increasing test concentration of cTnI.

Non-Reagent Test: Formation of a red line in the control region (C) and complete absence of line in the test region (T) between 15 and 30 minutes. Do not interpret after 30 minutes.

Invalid Testing: The absence of line formation in the control region (C) indicates an error in the procedure or deterioration of the cassette. In this case, repeat the test using a new cassette.

PROCEDURE LIMITATIONS

- 1- The line in the control region appears before the 15 minute test incubation, and this does not mean that non-reagent results can be interpreted ahead of time.
- 2- The reagent result will not change after 30 minutes, but the result must not be interpreted after this time to avoid incorrect results.
- 3- The results after 30 minutes of reaction must be considered invalid.
- 4- It is important to use the correct volume of sample, since lower or higher volumes may determine erroneous results.
- 5- The reading time of the reaction must be followed according to established technique in order to avoid false interpretation of the results.

6- The Troponin I kit, for cTnI, is a specific test to aid in the confirmation of a myocardial infarction. Reagent results should be confirmed with testing and other clinical information available. Results isolates the test could not be used for diagnosis definitive.

7- The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

INTERFERENCES

The Troponin I kit was evaluated for possible interference for samples with visible hemolysis, lipemia and icterus. Human hemoglobin, human Albumin, Bilirubin or Triglycerides were added to samples with different concentrations of Troponin, using samples are not add as a standard. No significant interferences were observed in 20 samples with positive or negative Troponin.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

PRODUCT PERFORMANCE

ACCURACY

Comparison of Methods, Methodological Specificity and Sensitivity

To establish the sensitivity and specificity of the Troponin I kit, another kit qualitative commercially available was used for comparison. It was showing an agreement of 99.31%. The observed sensitivity of the kit was 99.05% and specificity of 99.46%.

Tests of cross-reactivity of the kit was evaluated cTnI using homologous proteins Troponin. Were employed Human Cardiac Troponin T (cTnT), Troponin Human Skeletal (STnI) and Human Cardiac Troponin C (cTnC). No reactivity was observed.

PRECISION

Repeatability

The repeatability was calculated from 20 successive determinations, using 5 different sample, obtaining the following results:

Samples	N° of Repetitions	Expected Results	Found Results
Human Blood 01	20	100% Non-Reagent	100% Non-Reagent
Human Blood 02	20	100% Reagent	100% Reagent
Human Serum 01	20	100% Reagent	100% Reagent
Human Serum 02	20	100% Non-Reagent	100% Non-Reagent
Human Plasma 01	20	100% Non-Reagent	100% Non-Reagent

Reproducibility

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 5 different samples, obtaining the following results:

Samples	N° of Days	N° of Repetitions	Expected Results	Found Results
Human Blood 01	03	20	100% Non-Reagent	100% Non-Reagent
Human Blood 02	03	20	100% Reagent	100% Reagent
Human Serum 01	03	20	100% Non-Reagent	100% Non-Reagent
Human Serum 02	03	20	100% Reagent	100% Reagent
Human Blood 01	03	20	100% Reagent	100% Reagent

DETECTION LIMIT

The test can detect cTnI concentrations greater than or equal to 0.5 ng/mL.

DIAGNOSTIC SIGNIFICANCE

The serum levels of cardiac enzymes and iso-enzymes are essential for the diagnosis or exclusion of myocardial damage. Cardiac Troponin I is specific for cardiac tissue and is detected in serum only in cases of damage to the myocardium. She allows early identification of patients with symptoms suggestive of ischemia, allowing the identification of patients from 4 to 6 hours until 6 to 10 days after infarction.

You can also identify false elevations of CK-MB as in rhabdomyolysis. Troponin I increases its value until forty times in the first peak, whereas CK-MB only six to nine times.

BIBLIOGRAPHIC REFERENCES

1. J. Mair, N. Genser, Morandell D., et al. Cardiac troponin-I in the diagnosis of myocardial injury and infarction. Clin. Chem. Acta 1996, 245:19-38;
2. Vallins WJ, Brand NJ, Dabhade N, et al. Molecular cloning of human cardiac troponin-I using polymerase chain reaction. FEBS Lett 1990, 270:57-61;
3. Bodor GS, Porter S., Landt Y., et al. Development of monoclonal antibodies for an assay of cardiac troponin-I and preliminary results in Suspected cases of myocardial infarction. Clin. Chem. 1992, 38:2203-14;
4. M. Zaninotto, S. Altinier, M. Lachin, et al. Fluoroenzymometric method to measure cardiac troponin-I in the sera of Patients with myocardial infarction. Clin. Chem. 1996, 42:1460-6;
5. S. OH, K. Foster, P. Datta, et al. Use of a dual monoclonal solid phase and polyclonal detector to create an immunoassay for the detection of human cardiac troponin-I. Clinical Biochemistry 2000, 33 (4), 255-262.
6. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE
Customer Advisory Service
Phone: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for Troponin I kit: 10269360147

Review: May/2022

UNIVERSAL SYMBOLOGY

 CATALOG NUMBER

 LOT NUMBER

 MANUFACTURING DATE

 VALIDITY DATE
(last day of the month)

 TEMPERATURE LIMIT
(store)

 CONTENT IS SUFFICIENT
FOR <N> TEST

 SEE INSTRUCTIONS
FOR USE

 IN VITRO DIAGNOSTIC
PRODUCT

 KEEP AWAY
FROM SUNLIGHT

 DO NOT REUSE

 CAUTION

 MADE BY

 CONTROL

 POSITIVE CONTROL

 NEGATIVE CONTROL

 BIOLOGICAL RISK

 FLAMMABLE

 CORROSIVE

 TOXIC

 DO NOT USE IF
PACKAGE IS
DAMAGED

 PRODUCT
STERILIZED

 DANGER