# Bioclin

## LIPASE

REF K025

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### FINALIDADE

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### USAGE INSTRUCTIONS

#### FUNCTION

Method for determination of lipase activity in biological samples of serum or plasma (Heparin). Colorimetric test, only for *in vitro* diagnostic use.

#### PRINCIPLE OF ACTION

**Methodology:** Colorimetric

The Lipase acts hydrolyzing the Grease Acids Glycerol esters of long chains (Triacylglycerol) in Diglycerides, Monoglycerides and free Grease Acids. During the reaction the substrate, in a covered and stable environment, acquire the emulsification form (micelles) forming interfaces (water-lipids) necessary to the action of Lipase who in the presence of Dithionitrobenzoic (DTNB) form a chromogen of yellow color proportional to the serum Lipase.

#### REAGENTS

**Number 1 - Buffer** - Store between 2 and 8°C. Contains: Tris (Hydroxymethyl) Aminemethane < 200 mmol/L and preservative.

**Number 2 - Enzymatic Inhibitor** - Store between 2 and 8°C. Contains: PMSF < 15 mmol/L and solubilizer.

**Number 3 - Color Reagent** - Store between 2 and 8°C. Contains: DTNB < 6 mmol/L, Sodium Acetate < 200 mmol/L and preservative.

**Number 4 - Substrate** - Store between 2 and 8°C. Contains: 2-3-Dimercapto-1-Propanol Tributirato < 50 mmol/L, solubilizer and surfactant.

**Number 5 - Acetone** - Store between 2 and 8°C. Contains: Acetone < 20 mol/L.

#### PRESENTATION

| Reagents | Volume |
|----------|--------|
| Nº 1     | 40 mL  |
| Nº 2     | 1 mL   |
| Nº 3     | 4 mL   |
| Nº 4     | 2 mL   |
| Nº 5     | 80 mL  |

#### EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, water-bath at 37°C, watches and stopwatches, centrifuge, pipettes test tubes, Biocontrol N and Biocontrol P Bioclin. They can be found at market specialized on Laboratories of Clinical Analysis.

#### TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature is between 2 and 8°C. The transport at temperatures between up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

#### SPECIAL CARE

**1- For *in vitro* diagnostic use only.**

**2-** Strictly follow the methodology proposed to obtain exact results.

**3-** Water used in material cleaning must be recent and free of contaminants.

**4-** The water level in water-bath must be superior to the reagents level in test tubes.

**5-** The colorimetric cuvettes or the spectrophotometer must be strictly dry.

**6-** It is of great importance, for the good test development, a strict control of time and temperature.

**7-** The reagent Nº 1 and 3 are toxic and irritating to skin and mucous membrane. Handle the reagents carefully.

**8-** We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

**9-** To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.

**10-** Do not use the product in case of damaged packaging.

**11-** It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

#### SAMPLE

Serum obtained free of hemolysis, plasma collected with Heparin. Other anticoagulants (EDTA, Citrate, Oxalate) inhibit the lipase activity. Serum samples are stable for 7 days between 2 and 30°C and 12 months when frozen at -20°C. Plasma samples are stable for 7 days between 2 and 30°C and 2 months when frozen at -20°C.

#### PROCESS DESCRIPTION

##### TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Mark 2 test tubes with letters C (Control), A (Sample) and proceed as follows:

|              | Control | Sample |
|--------------|---------|--------|
| Reagent Nº 1 | 1.0 mL  | 1.0 mL |
| Sample       | 50 µL   | 50 µL  |
| Reagent Nº 2 | ---     | 50 µL  |
| Reagent Nº 3 | 100 µL  | 100 µL |

Place the tubes in water-bath at 37°C for 2 minutes to equilibrate the temperature. Then add:

|              |     |        |
|--------------|-----|--------|
| Reagent Nº 4 | --- | 0.1 mL |
|--------------|-----|--------|

Homogenize and incubate at 37°C for approximately 30 minutes. Then add:

|              |        |        |
|--------------|--------|--------|
| Reagent Nº 5 | 2.0 mL | 2.0 mL |
|--------------|--------|--------|

Homogenize well and let it rest for 3 minutes at room temperature.

Transfer the contents of the Control tubes and Sample for 2 centrifuge tubes. Centrifuge at 3.500 rpm for 5 minutes. Transfer the translucent supernatants to dry cuvettes of the colorimeter or spectrophotometer and determine the absorbencies of Control and Sample tubes at 410 nm (400 - 415 nm), matching the zero with the distilled or deionized water. The formed color is stable for 30 minutes.

#### CALCULATIONS

$$\text{Lipase (U/L)} = \frac{\text{Sample Abs.} - \text{Control Abs.} \times 1.000}{7}$$

**Example:**

Sample Absorbance: 0.131

Control Absorbance: 0.102

$$\text{Lipase (U/L)} = \frac{0.131 - 0.102 \times 1.000}{7} = 4.14 \text{ U/L}$$

#### PROCEDURE LIMITATIONS

The method is indicated only for manual determination of serum Lipase.

#### INTERFERENCES

Some drugs such as: cholinergic, codeine, indomethacin, meperidine and morphine can cause an increase in the serum Lipase. However the calcium ions can interfere reducing the values of Lipase.

#### INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

#### REFERENCE VALUES

The reference values in U/L for this method were obtained through the determination of lipase in healthy populations of male and female.

Serum: 2 to 15 U/L

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

## PRODUCT PERFORMANCE

### ACCURACY

#### Recovery

The recovery analysis was made with 05 sample determinations. The accuracy was calculated and was found in a good concordance with the reference values, obtaining a recovery between 95 and 107%.

#### Comparison of Methods

The Lipase kit was compared with other methods for determination of Lipase dosage commercially available. 07 analyzes were performed and the results were evaluated. The linear equation obtained was  $Y = 1.030X - 0.222$  and correlation coefficient 0.999. With these results, we can conclude that the kit shows good methodological specificity.

### PRECISION

#### Repeatability

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

| Repeatability                | Sample 1 | Sample 2 | Sample 3 |
|------------------------------|----------|----------|----------|
| Average Concentration (U/L)  | 4.22     | 13.34    | 17.91    |
| Standard Deviation (U/L)     | 0.10     | 0.11     | 0.08     |
| Coefficient of Variation (%) | 2.25     | 0.82     | 0.46     |

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The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

| Reproducibilidad             | Sample 1 | Sample 2 | Sample 3 |
|------------------------------|----------|----------|----------|
| Average Concentration (U/L)  | 4.25     | 13.29    | 17.90    |
| Standard Deviation (U/L)     | 0.03     | 0.04     | 0.01     |
| Coefficient of Variation (%) | 0.72     | 0.33     | 0.06     |

### SENSITIVITY

The sensitivity was calculated from 20 determinations of a sample free of lipase. The average found was 1.08 U/L with standard deviation of 0.01 U/L. The sensitivity, which indicates the methods detection limit, corresponds the average plus 3 times the standard deviation and is equal to 1.11 U/L.

### LINEARITY

The reaction is linear up to 300 U/L. For higher values, carry out a new determination, reducing the incubation time of the second step, after the addition of Reagent N° 4. Multiply the results by 30 and divide by the used incubation time, in minutes.

### DIAGNOSTIC SIGNIFICANCE

Lipase is an enzyme produced almost exclusively by pancreas.

In the presence of pancreatitis, the lipase is released at high concentrations in the blood, persisting for much longer than the amylase, allowing confirmation of late diagnosis disease. The values of lipase in pancreatitis can reach 5 to 10 times the upper limit of reference, while the other diseases lipase levels are always below 3 times the values of reference.

Lipase increases in acute pancreatitis, chronic pancreatitis (cancer of pancreas), obstruction of the ducts by calculation or neoplasms. Since the excretion of lipase is via the kidneys, kidney failure in general, lipase levels are high. In intestinal infarction, in acute cholangitis and in small bowel obstruction increased values of lipase can also be observed.

The lipase may be reduced physiologically during pregnancy, in the course of tuberculosis and other infectious diseases.

### BIBLIOGRAPHIC REFERENCES

- 1 - CHERRY, I.S., CRANDAL, L.A: Am J. Physiol.: 100:266, 1932.
- 2 - RICK; W.: Klin. Chem. Klin. Biochem.: 7: 530, 1969.
- 3 - FRIED. R., HOEFMAYER, J.: Z. Klin Chem. Klin Biochem.: 11:89, 1973.
- 4 - ROE. J.H., BYLER, R.E.: Anal. Biochem. 6:451, 1963.
- 5 - TIETZ, N.W., FIERECK, E.ª: Clin. Chim. Acta. 13:352, 1966.
- 6 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

### QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

### QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone.: +55 31 3439.5454 | E-mail: [bioclin@bioclin.com.br](mailto:bioclin@bioclin.com.br)  
CNPJ: 19.400.787/0001-07 - Made in Brazil

### CUSTOMER SERVICE

Customer Advisory Service  
Phone.: 0800 0315454 | E-mail: [sac@bioclin.com.br](mailto:sac@bioclin.com.br)

ANVISA registration for Lipase kit: 10269360102

Review: May/2022

## UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE  
(last day of the month)

TEMPERATURE LIMIT  
(store)

CONTENT IS SUFFICIENT  
FOR <N> TEST

SEE INSTRUCTIONS  
FOR USE

IN VITRO DIAGNOSTIC  
PRODUCT

KEEP AWAY  
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF  
PACKAGE IS  
DAMAGED

PRODUCT  
STERILIZED

DANGER