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REF K222

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FINALIDADE

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PELIGRO

PROCEDURE LIMITATIONS

This process can only be accomplished with a spectrophotometer cuvette at 37°C, due to the small time interval reaction. **The methodology described should not be used for extremely lipemic samples.**

INTERFERENT

Triglyceride concentrations up to 900 mg/dL, Bilirubin up to 5 mg/dL and Hemoglobin up to 180 mg/dL do not interfere significantly with the reaction. Samples with Bilirubin, Hemoglobin and Triglycerides in concentrations greater than those above mentioned can be diluted in NaCl 150 mmol/L (0.85%) before performing the tests to reduce interferences. The result obtained should be multiplied by the dilution factor.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard and the BIOCAL calibrator are traceable to the reference material NIST (National Institute of Standards and Technology) SRM 914, making it equivalent to the IDMS method (Isotope Dilution Mass Spectrometry), as recommended by the NKDEP (National Kidney Disease Education Program).

REFERENCE VALUES

The reference values in mg/dL, for this method were obtained through the determination of uncompensated Creatinine in healthy populations of male and female.

Serum/Plasma ⁵	Adult (Male)	0.70 - 1.20 mg/dL
	Adult (Female)	0.50 - 1.00 mg/dL
	Newborn	0.31 - 0.92 mg/dL
	2 weeks - 1 year	0.16 - 0.39 mg/dL
	1 - 3 years	0.17 - 0.35 mg/dL
	3 - 5 years	0.26 - 0.42 mg/dL
	5 - 7 years	0.29 - 0.48 mg/dL
	7 - 9 years	0.34 - 0.55 mg/dL
	9 – 11 years	0.32 - 0.64 mg/dL
Urine	11 - 13 years	0.42 - 0.71 mg/dL
	13 - 15 years	0.46 - 0.81 mg/dL
	Adult (Male)	1500 - 2500 mg/24h or 21 - 26 mg/kg/24h
	Adult (Female)	800 - 1500 mg/24h or 16 - 22 mg/kg/24h
Creatinine clearance	Children 2 - 3 years	6 - 22 mg/kg/24h
	Children > 3 years	12 - 30 mg/kg/24h
	Adult (Male)	97 - 137 mL/min/1.73m²
	Adult (Female)	88 - 128 mL/min/1.73m²
Creatinine clearance	Children	70 - 140 mL/min/1.73m²

mg/kg weight = mg/24h divided by body weight

To convert the values to mg/dL to mmol/L (SI), multiply by 0.0884. These values should be used as guidelines, each laboratory should create their own range of reference values, according to the population served. The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

ACCURACY

Comparison of Methods

The Creatinine Automated kit was compared with other methods commercially available for measurement of Creatinine. 40 analyses were performed and the results were evaluated. The linear equation was obtained $y = 0.9913x - 0.0248$ and correlation coefficient 0.9991. With these results, we can conclude that the kit shows good methodology specificity.

PRECISION

Repeatability

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Samples		
	1	2	3
Average Concentration	1.01	1.94	5.55
Standard Deviation	0.01	0.01	0.04
Variation Coefficient (%)	1.20	0.74	0.72

Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Samples		
	1	2	3
Average Concentration	1.02	1.95	5.60
Standard Deviation	0.01	0.02	0.05
Variation Coefficient (%)	1.34	0.97	0.87

SENSITIVITY

Sensitivity was calculated from 40 determinations of a sampel free of Creatinine. The average found was -0.009 mg/dL, with a standard deviation 0.011 mg/dL. The sensitivity, which indicates the method detection limit, corresponds the average plus 3 times the standard deviation, and is equal to 0.034 mg/dL.

LINEARITY

The reaction color is linear to the concentration of 12 mg/dL. For higher values dilute the sample with Sodium Chloride 0.85% and repeat the determination. Multiply the value obtained by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

Creatinine is synthesized primarily in the liver and kidneys, is produced and excreted at a constant rate independent of diet, degree of hydration and protein metabolism. With renal failure, the laboratory data provide insurance rates as the ability kidneys to excrete, secrete and reabsorb creatinine.

Creatinine is elevated in acute and chronic renal failure, in urinary tract obstruction, congestive heart failure, in dehydration, shock and in diabetes *mellitus*.

Decreased levels may be observed in muscular dystrophies, malnutrition, reduced muscle mass and severe renal disease.

BIBLIOGRAPHIC REFERENCES

- 1 - SLOT, C.; Sand. J. Clin. Lab. Invest. 17:381, 1965.
- 2 - HENRY, R. J.: Clinical Chemistry - Principles and Technics, 4 Ed., Harper and Row, New York, 1964.
- 3 - TODD,SANFORD,DAVIDSOHN;Diagnósticos Clínicos, 16 Ed., Editora Manole, 1983.
- 4 - CARL, A. B. and EDWARD, R. A.: Tietz Textbook of Clinical Chem 2nd ed., 1994, 1531-1539.
- 5 - CERIOTTI, F. et al. Reference intervals for serum creatinine concentrations: Assessment of available data for global application. Clinical Chemistry, 2008.
- 6 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda


Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for Creatinine Automated kit: 10269360328

Review: May/2022

REF

CATALOG NUMBER

LOT

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER