DENGUE BIO

REF K105

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FINALIDADE

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- Halstead SB, Selective primary health care: strategies for control of disease in the developing world: XI, Dengue. Rev. Infect. Dis. 1984; 6:251-264.
- Halstead SB, Pathogenesis of dengue: challenges to molecular biology. Science 1988; 239:476-481
- Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). Japanese J. Trop. Med. Hygiene 1994; 22: 9-12
- Lam SK. Dengue haemorrhagic fever. Rev. Med. Micro. 1995; 6:39-48
- Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.
- Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. Clin Diagn Lab Immunol. 2003 Jul; 10(4): 725-8.
- Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. Wien Med Wochenschr (in German). 1997; 147(19-20): 463-4
- Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus infections. Microbiol Immunol. 1994; 38(12): 951-5.]
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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1. Halstead SB, Selective primary health care: strategies for control of disease in the developing world: XI, Dengue. Rev. Infect. Dis. 1984; 6:251-264.
2. Halstead SB, Pathogenesis of dengue: challenges to molecular biology. Science 1988; 239:476-481.
3. Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). Japanese J. Trop. Med. Hygiene 1994; 22: 9-12.
4. Lam SK. Dengue haemorrhagic fever. Rev. Med. Micro. 1995; 6:39-48 .
5. Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.
6. Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. Clin Diagn Lab Immunol. 2003 Jul; 10(4): 725-8.
7. Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. Wien Med Wochenschr (in German). 1997; 147(19-20): 463-4.
8. Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus infections. Microbiol Immunol. 1994; 38(12): 951-5.].
9. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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FUNCTION

Quick Immunochromatographic. Test for determination of antibodies IgM and IgG for the Dengue virus in biological samples whole blood, serum or plasma. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Immunochromatography

The Dengue Bio kit is an immunochromatographic test who detects separately classes antibodies IgM and IgG for the Dengue for virus. The kit contains antibodies of goat anti-human IgG in the test line for IgG, and goat anti-human IgM. In test line for IgM. After the sample adding, antigens of the Dengue virus conjugated to gold colloidal particles, present in the kit, bind to the sample antibodies and migrate by capillarity through the membrane. In case of reagent samples for IgG, it will form a red line in the IgG test area. Reagent samples for IgM will form a red line in the IgM test area. The mixture of the antigen-sample continues to run through the membrane until the control line, where occurs a formation of a red line, for the reagent and non-reagent samples. The red line formation in the control area demonstrates that the kit is adequate for use.

REAGENTS

Number 1: Cassette: Store between 15 and 30°C. **Do not freeze.** Contains:

- A) 01 plastic cassette (device) with 01 test strip containing:
- C control region: rabbit IgG anti goat antibody.
 - G test region (T1): goat anti-human IgG antibody.
 - M test region (T2): goat antihuman IgM antibody.
 - Conjugate: Dengue virus antigen (I, II, III and IV) conjugated with the gold colloidal particles.

B) 01 sachet of silica

Number 2: Diluent - Store between 15 and 30°C. **Do not freeze.** Contains: Tris buffer 100 mmol/L pH 7.4 and preservative.

PRESENTATION

Reagent	Number of Cassettes for Package	Diluent
Cassette + Diluent	10	1.00 mL
	25	2.50 mL
	50	2 x 2.50 mL
	100	4 x 2.50 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Watches and stopwatches. These items are found in the specialized market of articles for Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperature should be between 15 and 30°C. Protect from light and avoid moisture.**Do not freeze.**

SPECIAL CARE

- 1- **For *in vitro* diagnostic use only.**
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- The cassettes do not present contaminant substances, but once the samples are biological, the safety routine must be followed strictly and the use of disposable gloves is indispensable.
- 4- Use clean and dry recipients to settle the samples.
- 5- The reagent N° 2 contain sodium azide, irritating to skin and mucous membranes. Handle the reagent with care.

6- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

7- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

8- Do not use the product in case of damaged packaging.

9- It is essential that the instruments used are properly calibrated and subjected to periodic maintenance.

SAMPLE

Whole Blood: shall be stored between 2 and 8°C and tested In 24 hours after collection.

Plasma or serum: shall be separate from blood interface immediately after centrifugation, and must occur instantly after the collection. The plasma or serum can be stored for 3 days between 2 and 8°C and 6 months at -20°C.

PROCESS DESCRIPTION

The stability of the Dengue BIO cassette (after opening the aluminized sachet) is 2 hours. This stability may vary according to the environmental conditions after opening the aluminized sachet.

- 1- The sample must be in a temperature between 15 and 30°C before begin the test;
- 2- Remove the cassette from is protective packaging and appropriately identify it so.
- 3- Transfer 10 µL of the total blood, plasma or serum inside the sample window;
- 4- Hold the flask of Diluent (Reagent N° 2) vertically and apply 4 drops (100 µL) of diluent in the sample well;

5- Wait the formation of the lines after 15 to 20 minutes. **Do not interpret it after 20 minutes.**

RESULTS INTERPRETATION

IgG Reagent Test: a red line forms in the test area G (T1) and another line in the control area (C) in the first 15 to 20 minutes. Do not interpret after 20 minutes.


IgM Reagent Test: a red line forms in the test area M (T2) and another line in the control area (C) in the first 15 to 20 minutes. Do not interpret after 20 minutes.


IgM and IgG Reagent Test: two red line forms in the test area (G = T1 and M = T2) and a red line in the control area (C) in the first 15 to 20 minutes. Do not interpret after 20 minutes.

Non-Reagent Test: a red line forms in the control area (C) and complete absence of red lines in the test areas (G = T1 and M = T2). Do not interpret after 20 minutes.

Invalid Test: complete absence of red lines in the control area (C) and test regions (G = T1 and M = T2).

NOTES:

- 1- The non-reagent results must only be read in 20 minutes, and not before or after this time.
- 2- The line in the control area must show before 10 minutes of incubation test, and that does not mean the non-reagent results can be interpreted before the time.
- 3- The results must not be interpreted after 20 minutes in order to avoid incorrect results.
- 4- The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

INTERFERENCES

Cross reactivity with other virus of the flavivirus can be common, according to bibliographic references, however there was no evidence in the product developing tests.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

PROCEDURE LIMITATIONS

- 1- The Bioclin's Bio Dengue kit can only be used for the *in vitro* diagnosis.
- 2- The Bioclin's Bio Dengue Kit determines only the the presence of antibodies against the dengue virus (I, II, III and IV) and shall not be used isolated as a diagnosis criteria. Reagent results should be evaluated with other clinical information available. The isolated results of the test must not be used for a definitive diagnosis.
- 3- A non-reagent result does not extinct the possibility of infection by the Dengue virus (I, II, III and IV). Scientific studies reveal that infected people can present an immunological window of 5 to 10 days for IgM and IgG detection, respectively. The patient shall be evaluated again during this period.
- 4- Immunosupressed patients' results should be carried out with caution.
- 5- In every diagnosis test, the results must be interpreted within other patients clinical data.
- 6- If the non-reagent result continues and the clinical symptoms persists, an additional test can be carried out by another method. A non-reagent result does not extinct, in any way, the possibility of infection by the dengue virus.

PRODUCT PERFORMANCE

ACCURAY

Comparison of Methods, Methodology and Specificity

To establish the sensitivity and specificity of the Dengue Bio kit in relation to other Dengue test, 221 known clinical samples were studied. Other qualitative kit available in the market was used to compare with the kit and for determination of relative sensitivity and specificity. The found clinical sensitivity was 96% and the specificity found was of 95.8%.

PRECISION

Repeatability

The repeatability was calculated from 20 successive determinations, using 5 samples with different concentrations, obtaining the following results:

Samples	Repetition Nº	Expected Result	Found result
Blood 01	20	100% Reagent IgG	100% Reagent IgG
Blood 02	20	100% Non-Reagent	100% Non-Reagent
Serum 01	20	100% Non-Reagent	100% Non-Reagent
Serum 02	20	100% Reagent IgM	100% Reagent IgM
Plasma 01	20	100% Reagent IgG/IgM	100% Reagent IgG/IgM

Reproducibility

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Samples	Nº of Days	Repetition Nº	Expected Result	Found Result
Blood 01	3	20	100% Non-Reagent	100% Non-Reagent
Blood 02	3	20	100% Reagent IgM	100% Reagent IgM
Serum 01	3	20	100% Reagent IgG/IgM	100% Reagent IgG/IgM
Serum 02	3	20	100% Non-Reagent	100% Non-Reagent
Plasma 01	3	20	100% Reagent IgG/IgM	100% Reagent IgG/IgM

DIAGNOSTIC SIGNIFICANCE

The Dengue virus is a flavivirus transmitted by the *Aedes aegypti* mosquito. It is distributed in the tropical and subtropical areas worldwide, causing over 100 million infections yearly. The classic Dengue is characterized by high fever (normally between 38 and 40°C of suddenly beginning, low appetite, minor illnesses, caphaleia, muscle and eyes sore. Hemorrhagic Dengue can cause gingivorrhagy and epistaxis, internal bleeding and disseminated intravascular coagulation, with many vital organs damage, occasioning in death. The primary infection of dengue causes a increase in the IgM antibodies after 3 to 5 days in the beginning of the fever. IgM antibodies normally remain in the circulation for 30 to 90 days. Patients of endemic regions can present the secondary infections, that result in high levels of IgG antibody, isolated or simultaneously as a response of IgM. Therefore, an specific detection of antibodies igM or IgG can help to difference the primary or secondary.

BIBLIOGRAPHIC REFERENCES

- 1.Halstead SB, Selective primary health care: strategies for control of disease in the developing world: XI, Dengue. Rev. Infect. Dis. 1984; 6:251-264
- 2.Halstead SB, Pathogenesis of dengue: challenges to molecular biology. Science 1988; 239:476-481
- 3.Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). Japanese J. Trop. Med. Hygiene 1994; 22: 9-12.
- 4.Lam SK. Dengue haemorrhagic fever. Rev. Med. Micro. 1995; 6:39-48.
- 5.Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.
- 6.Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. Clin Diagn Lab Immunol. 2003 Jul; 10(4): 725-8.
7. Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. Wien Med Wochenschr (in German). 1997; 147(19-20): 463-4.
- 8.Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus infections. Microbiol Immunol. 1994; 38(12): 951-5.].
- 9.QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca

CEP 31565-130 - Belo Horizonte - MG - Brasil

Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service

Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for Dengue Bio kit: 10269360174

Review: May/2022

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

CONTROL +
POSITIVE CONTROL

CONTROL -
NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

STERILE R
PRODUCT
STERILIZED

DANGER