TRANSAMINASE AST (TGO) CINÉTICA

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REF K048

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FINALIDADE

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- 1 - The Committee on Enzymes of the Scandinavian Society for Clinical Chemistry and Clinical Physiology Scand J. Clin Lab Invest., 1974, 33, 291 - 306.
- 2 - BERGMEYER, Bowers and cols., Clin. Chim., Acta., 1976, 70, 19 - 42, 1977, 21 - 22.
- 3 - BERGMEYER, HV.; SCHEIBE, P.; WAHLEFELD, A.W., Clin., Chem., 1978, 24, 58 - 73.
- 4 - Expert Panel on Enzymes of the International Federation of Clinical Chemistry: Part 3. Revised IFCC Method for Aspartate Aminotransferase, Clin., Chem., 1978, 24, 720.
- 5 - Scandinavian Committee on Enzymes Scand J. Clin Lab Invest, 1981, 41, 107 - 116.
- 6 - BURTIS; CARL, A.; ASHWOOD; EDWARD, R., Clin. Chem., Tietz TEXT BOOK of; 2a ed., 1.986, 788 - 797.
- 7 - PESCE, A., J.; KAPLAN, L., A., Methods in Clin. Chem., 1.987.
- 8 - SOLDIN, S.J., BRUGNARA, C., WONG, E.C.: Pediatric Reference Intervals, 5.ed. Washington: AACC Press, 2005.p.3-4.
- 9 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31
- 10 - ELLIS S. Benson M.D. & D. B. Tonks PH.D. (1963) Laboratory medicine, Postgraduate Medicine, 34:4, A-58-A-70, DOI: 10.1080/00325481.1963.11694890
- 11 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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- 1 - The Committee on Enzymes of the Scandinavian Society for Clinical Chemistry and Clinical Physiology Scand J. Clin Lab Invest., 1974, 33, 291 - 306.
- 2 - BERGMAYER, Bowers and cols., Clin. Chim., Acta., 1.976, 70, 19 - 42, 1977, 21 - 22.
- 3 - BERGMAYER, HV.; SCHEIBE, P.; WAHLEFELD, A.W., Clin., Chem., 1.978, 24, 58 - 73.
- 4 - Expert Panel on Enzymes of the International Federation of Clinical Chemistry: Part 3. Revised IFCC Method for Aspartate Aminotransferase, Clin.,Chem., 1.978, 24, 720.
- 5 - Scandinavian Committee on Enzymes Scand J. Clin Lab Invest. 1.981, 41, 107 - 116.
- 6 - BURTIS; CARL, A.; ASHWOOD; EDWARD, R., Clin. Chem., Tietz TEXT BOOK of; 2a ed., 1.986, 788 - 797.
- 7 - PESCE, A., J.; KAPLAN, L., A., Methods in Clin. Chem., 1.987.
- 8 - SOLDIN, S.J., BRUGNARA, C., WONG, E.C.: Pediatric Reference Intervals, 5.ed. Washington: AACC Press, 2005.p.3-4.
- 9 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31
- 10 - ELLIS S. Benson M.D. & D. B. Tonks PH.D. (1963) Laboratory medicine, Postgraduate Medicine, 34:4, A-58-A-70, DOI: 10.1080/00325481.1963.11694890
- 11 - QUIBASA: Datos del Departamento de Investigación y Desarrollo.

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REF K048

USAGE INSTRUCTIONS

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:
SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br



FUNCTION

Test for determination of Aspartate Amino Transferase (AST or TGO) in biological sample (serum or plasma), through kinetic test. Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Kinetic (UV)

AST catalyzes the transfer of Amine groups from Aspartate to Alfa - Ketoglutarate, leading to the formation of Glutamate and Oxalacetate. Oxalacetate in the presence of MDH reacts with NADH, reducing to Malate and NADH oxidizes to NAD +. The oxidation rate is proportional to the activity of AST in the sample.

REAGENTS

Number 1 - Substrate - Store between 2 e 8 °C. Contain Buffer < 200 mmol/L, Sodium Aspartate < 450 mmol/L, D-Lactate Dehydrogenase < 5 KU, Malate Dehydrogenase < 5 KU, chelator, surfactant and preservative.

Number 2 - Coenzyme - Store between 2 e 8 °C. Contém: Tris Buffer < 300 mmol/L, Alpha Ketoglutaric Acid < 180 mmol/L, NADH < 5 mmol/L, surfactant and preservative.

PRESENTATION

Presentation	Reagent Nº 1	Reagent Nº 2
K048-6-VET	2x40 mL	2x10 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Thermostatic spectrophotometer, pipettes and tips, test tubes, clock or stopwatch, Biocal, Biocontrol N e Biocontrol P Bioclin. These items are found in the specialized market for articles for Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8 °C. The transport, at temperatures up to 30 °C, should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must be recent and free of contaminants.

4- Saturated deionizing columns release alkaline water, various ions and oxidizing and reducing agents, which can significantly alter the results.

5-Visible hemolysis can cause slight variations in results.

6- We recommend applying local, state and federal environmental protection standards so that the disposal of reagents and biological material is carried out in accordance with current legislation.

7- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or

upon request by the SAC (Customer Advisory Service) of Quibasa.

8-Do not use the product in case of damage to the packaging.

9- It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

10-Reagents No. 1 and 2 contain sodium azide and should be handled with care.

11- It is important, for the good performance of the test, a strict control of time, temperature and pH.

12- Lipemic and icteric samples increase the absorbance by 340 nm. In this case, dilute the sample 1: 2 with physiological saline multiply the result obtained by 2.

13- Samples stored without centrifugation for long periods may have their pyruvate concentration increased, which may alter the results. See more information in Interferents.

14- Absorbance readings below 0.800 of the Working Reagent indicate loss of the same. In this case, do not use the reagent.

SAMPLES

Serum or Plasma: collected with EDTA or Heparin, obtained free from hemolysis. The enzyme is stable for 3 days at 2 to 8 °C and 3 months at -20 °C.

PROCESS DESCRIPTION

PREPARATION OF WORKING REAGENT

Mix 4 parts of Reagent Nº 1 with 1 part of Reagent Nº 2. Working Reagent is stable during 72 hours between 15 and 30 °C and 14 days between 2 and 8 °C.

The stability of the Transaminasa AST (TGO) product Kinetics after opening (In use / On board - Product installed in equipment) is 30 days. This stability can vary depending on the condition of the problem. Therefore, it is suggested to monitor the performance of the product using control numbers.

REACTION CONDITIONS

It is indispensable to use a thermostated cuvette at 37 °C, 1 cm optical path and reading at 340 nm (334 - 365 nm).

TECHNIQUE

Bioclin recommends using the Biocal kit as a calibrator and the Bioclin Biocontrol N and P kits as a control serum.

Add 100 µL of Sample to 1.0 mL of Working Reagent, mix and transfer to a thermostated cuvette at 37 °C and wait 1 minute. For equine samples, it is suggested to previously dilute the sample in a 1:2 ratio in saline (50 µL of sample and 50 µL of saline solution sample - 0.85-0.9% sodium chloride solution).

Take the initial reading, simultaneously starting the stopwatch. Repeat the readings after 1, 2 and 3 minutes. Calculate the average of the absorbance differences per minute (ΔAbs./min.) and use it to calculate the result. For samples that have been diluted, the result obtained should be multiplied by 2.

DESCRIPTION OF THE CALCULATIONS

AST - Results expressed in U/L

340 nm = ΔAbs./min. x 1746

334 nm = ΔAbs /min. x 1780

365 nm = ΔAbs/min. x 3235

PROCESS LIMITATIONS

The specifications below refer to semi-automatic equipment:

The kinetic method is based on molar absorptivity, for this reason, the readings must be performed on a spectrophotometer that fulfills the following conditions:

- Wavelength 340 nm
- Semi-trajectory of the 10 nm passband
- Spurious light less than 0.5%
- Thermostatized 1cm bucket

INTERFERENT

No interference was observed by Bilirubin until 19 mg/dL, Hemoglobin until 45 mg/dL, Triglycerides until 650 mg/dL and Pyruvate until 0.2 mmol/L.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and ccuracy of the dosages.

TRACEABILITY

The kit calibration can be made using the theoretical calibration factor based on the molar absorptivity of the NADH, or through the BIOCAL calibrator. Bioclin recommends the usage of the BIOCAL calibrator which is traceable to the reference material ERM-AD457.

REFERENCE VALUES

Reference values were determined in serum samples of dogs, cats, healthy cattle and horses, males and females, adults without predilection of breeds.

Specie	Reference Value
Cat	8 - 59 U/L
Dog	16 - 71 U/L
Cattle	49 - 109 U/L
Horses	212 - 479 U/L

These values should be used as a guide, and each laboratory should create its range of reference values, according to the population served. The results provided by this kit must be interpreted by the responsible medical professional, not being the only criterion for determining diagnosis and / or treatment of the patient.

PRODUCT PERFORMANCE

Accuracy

Comparison of Methods

The Transaminase AST (TGO) Kinetic kit was compared with another commercially available method for measuring aspartate aminotransferase. 40 analyzes were carried out and the results were evaluated. The linear equation obtained was $Y = 0.9765X + 1.0283$, with a linear correlation coefficient equal to 0.9999. With these results, it can be concluded that the kit has good methodological specificity.

Precision

Repeatability

The repeatability was calculated from 10 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average Concentration(U/L)	35.28	41.92	58.86
Standard Deviation (U/L)	0.52	0.06	0.11
Coefficient of Variation (%)	1.46	0.13	0.19

Reproducibility

Reproducibility or was calculated from the analysis of 3 samples of different levels during 15 days, not necessarily consecutive, totaling 60 determinations for each sample, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Average Concentration(U/L)	40.34	45.71	58.53
Standard Deviation (U/L)	0.4	0.5	0.6
Coefficient of Variation (%)	1.01	1.01	1.02

SENSITIVITY

Sensitivity was calculated from 40 determinations of a sample free from Aspartate aminotransferase. The average found was 2.10 U/L with a standard deviation of 0.585 U/L. The sensitivity, which indicates the detection limit of the method, corresponds to 3 times the standard deviation, and is equal to 1.756.

LINEARITY

The reaction is linear up to a concentration of 400 U/L. For samples with values above 400 U/L, it is recommended to dilute the sample with 0.85% Sodium Chloride, repeat the dosage and multiply the result obtained by the dilution factor. For example, if the sample is diluted 1:10 and the result obtained is 80 U/L, multiply by 10 (80 x 10) U/L and consider the activity value of 800 U/L as the final result.

DIAGNOSTIC SIGNIFICANCE

The enzyme Aspartate Amino Transferase (AST) is located in the cytosol and mitochondria of skeletal muscle

cells, liver and cardiac muscle. An increase in serum AST indicates lesions in the liver and/or muscle tissue.

In cattle and horses, an increase in AST is common due to hepatic lipidosis, liver abscesses, bacterial (mainly leptospirosis in cattle) and parasitic (*Fasciola hepatica*) infections. An increase in AST may also occur due to musculoskeletal injuries.

In cats, an increase in AST may occur in cases of pyogranulomatous hepatitis secondary to Feline Infectious Peritonitis (FIP) and hyperthyroidism.

In cases of *in vivo* hemolysis, a false increase in serum AST may occur due to the release of intraerythrocytic AST. In horses, physical activity a few hours before collection may result in the release of intramuscular AST, resulting in a physiological increase in the enzyme without the occurrence of injury. The use of corticosteroids and anticonvulsants may cause liver injury and consequently an increase in AST.

BIBLIOGRAPHIC REFERENCES

- 1 - The Committee on Enzymes of the Scandinavian Society for Clinical Chemistry and Clinical Physiology Scand J. Clin Lab Invest., 1974, 33, 291 - 306.
- 2 - BERGMAYER, Bowers and cols., Clin. Chim., Acta., 1.976, 70, 19 - 42, 1977, 21 - 22.
- 3 - BERGMAYER, HV.; SCHEIBE, P.; WAHLEFELD, A.W., Clin., Chem., 1.978, 24, 58 - 73.
- 4 - Expert Panel on Enzymes of the International Federation of Clinical Chemistry: Part 3. Revised IFCC Method for Aspartate Aminotransferase, Clin.,Chem., 1.978, 24, 720.
- 5 - Scandinavian Committee on Enzymes Scand J. Clin Lab Invest, 1.981, 41, 107 - 116.
- 6 - BURTIS; CARL, A.; ASHWOOD; EDWARD, R., Clin. Chem., Tietz TEXT BOOK of: 2a ed., 1.986, 788 - 797.
- 7 - PESCE, A., J.; KAPLAN, L., A., Methods in Clin. Chem., 1.987.
- 8 - SOLDIN, S.J., BRUGNARA, C., WONG, E.C.: Pediatric Reference Intervals, 5.ed. Washington: AACC Press, 2005.p.3-4.
- 9 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31
- 10 - ELLIS S. Benson M.D. & D. B. Tonks PH.D. (1963) Laboratory medicine, Postgraduate Medicine, 34:4, A-58-A-70, DOI: 10.1080/00325481.1963.11694890
- 11 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone: 0800 0315454 | E-mail: customerservice@bioclin.com.br

Registration Number: Product exempt from registration with the Ministry of Agriculture, Livestock and Supply.

Technical Responsible: Dr. Camila Eckstein (CRMV/ MG 20.611)

Review: June/2025

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MADE BY

LOT NUMBER

CONTROL

MANUFACTURING DATE

POSITIVE CONTROL

VALIDITY DATE
(last day of the month)

NEGATIVE CONTROL

TEMPERATURE LIMIT
(store)

BIOLOGICAL RISK

CONTENT IS SUFFICIENT
FOR <N> TEST

FLAMMABLE

SEE INSTRUCTIONS
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC
PRODUCT

TOXIC

KEEP AWAY
FROM SUNLIGHT

DO NOT USE IF
PACKAGE IS
DAMAGED

DO NOT REUSE

PRODUCT
STERILIZED

CAUTION

DANGER