TRANSAMINASE ALT (TGP) CINÉTICA

REF K049

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FINALIDADE

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REF **K049**

INSTRUCTIONS FOR USE

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:
SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

FUNCTION

Method for determination of Alanine Amino Transferase (ALT) in biological samples of serum or plasma (EDTA or Heparin). Kinetic test, for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Kinetic (UV)

Kinetic determination (UV) of ALT, according to the reaction:

ALT catalyzes the transfer of Amino groups of Alanine for α - Ketoglutarate, leading to formation of Pyruvate and Glutamate. Pyruvate in the presence of LDH reacts with NADH, reducing the Lactate and oxidizes NADH to NAD⁺. The rate of oxidation is proportional to ALT activity in the sample.

REAGENTS

Number 1 - Substrate - Store between 2 and 8 °C. Contains: Tris Buffer < 200 mmol/L (pH 7.8), LDH 2400 U/L, L-Alaninane < 500 mmol/L, stabilizer and preservative.

Number 2 - Coenzyme - Store between 2 and 8 °C. Contains: Tris Buffer < 200 mmol/L (pH 10), Alpha-Ketoglutarate < 100 mmol/L, NADH < 5 mmol/L and preservative.

PRESENTATION

Presentation	Reagent N° 1	Reagent N° 2
K049-6-VET	2 x 40 mL	2 x 10 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Thermostated spectrophotometer, pipettes, watch or stopwatches, test tubes, Biocontrol N, Biocontrol P Bioclin and Biocal Bioclin. These items are found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures between 15 and 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning and preparation of the working reagent must be recent and free of contaminants.

4- It is important, for the good development of the test, a rigorous control of time, temperature and pH.

5- Reagent N° 1 contains Sodium Azide, and should be handled with care.

6- Lipemic and icteric samples increase the absorbance at 340 nm. In this case, you should dilute the sample 1:2 with Saline Solution. Multiply the results by 2.

7- Readings of absorbance less than 0.800 of Reagent Working indicate loss of it. In this case, do not use the reagent.

8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and

biological material can be made in accordance with current legislation.

9- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

10- Do not use the product in case of damaged packaging.

11- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or plasma with EDTA or heparin, obtained free of hemolysis. The sample is stable for 3 days between 2 and 8°C, and 3 months at -20°C.

PROCESS DESCRIPTION

The calibration stability of the Transaminase ALT Kinetic kit installed on refrigerated equipment is at least 12 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

PREPARATION OF WORKING REAGENT

Mix 4 parts of Reagent N° 1 with 1 part of Reagent N° 2. Working Reagent is stable during 72 hours between 15 and 30 °C and 14 days between 2 and 8 °C.

REACTION CONDITIONS

Is indispensable condition the use of thermostated cuvette at 37°C, 1 cm optical path and reading at 340 nm (334 - 365 nm).

TECHNIQUE

Bioclin recommends using the kit Biocal Bioclin as calibrator and as control serum, Biocontrol N and P Bioclin Kits.

Add 100 μ L of Sample to 1.0 mL of Working Reagent, mix and transfer to a thermostated cuvette at 37 °C and wait for 1 minute. Make the initial reading, simultaneously starting the timer. Repeat readings after 1, 2 and 3 minutes. Calculate the average differences in absorbance per minute ($\Delta A/\text{min.}$) and use it to calculate the result.

CALCULATIONS

ALP (U/L) 340 nm = $\Delta A/\text{min.} \times 1746$

334 nm = $\Delta A/\text{min.} \times 1780$

365 nm = $\Delta A/\text{min.} \times 3235$

Results are expressed in U/L.

For an average variation in absorbance ≥ 0.15 in 340 nm and 334 nm or ≥ 0.080 in 365 nm, repeat the determination diluting the sample with NaCl 0.85%. Multiply the results obtained by the dilution factor.

PROCESS LIMITATIONS

The kinetic method is based on the molar absorptive. By this reason, the readings must be conducted in a spectrophotometer that satisfies the following conditions:
Wavelength 340 nm
Semi trajectory of the pass band 10 nm
Stray light less than 0.5%
1cm thermostated cuvette

INTERFERENT

No interference was observed by Bilirubin until 19 mg/dL, Hemoglobin until 180 mg/dL and Triglycerides until 650 mg/dL. The use of anticonvulsant drugs such as primidone and phenobarbital can elevate ALT activity.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit calibration can be made using the theoretical calibration factor, based on the molar absorptivity of the NADH, or through the BIOCAL calibrator. Bioclin recommends the usage of the BIOCAL calibrator which is traceable to the reference material ERM-AD454.

REFERENCE VALUES

The reference values were determined in serum samples of healthy dogs, cats, cattle and horses, males and females, adults and without breed predilection.

Species	Value
Dogs	6 - 80 U/L
Cats	21 - 74 U/L
Cattle	13 - 39 U/L
Horses	4 - 17 U/L

These values should be used as a guideline and each laboratory should establish its own reference value, according to the population served. The results provided by this kit must be interpreted by the veterinary medical professional, and are not the only criterion for determining the diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

ACCURACY

Recovery

The recovery analysis was performed with 05 determinations of sample. Accuracies were calculated and were found in good agreement with the reference values, obtaining a recovery between 95% and 102%.

Comparison of Methods

The Transaminase ALT Kinetic kit was compared with another method commercially available to measure Alanine Amino Transferase. 42 analyzes were conducted and the results were evaluated. The linear equation obtained was $Y = 0.9693X + 1.33$, with correlation coefficient 0.9995. With these results we can conclude the kit shows good methodological specificity.

PRECISION

Repeatability

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average Concentration (U/L)	34.05	146.95	90.55
Standard Deviation (U/L)	0.22	0.64	0.50
Coefficient of Variation (%)	0.65	0.43	0.56

Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Average Concentration (U/L)	34.11	144.88	89.83
Standard Deviation (U/L)	0.47	1.97	1.11
Coefficient of Variation (%)	1.36	1.36	1.24

SENSITIVITY

Sensitivity of Transaminase ALT Kinetic kit was calculated from 40 determinations of a sample free of Alanine Amino Transferase. The average found, which indicates the method detection limit, was of 0.050 U/L with a standard deviation of 0.316 U/L. The sensitivity, which indicates the method quantification limit, corresponds the average plus 3 times the standard deviation and is equal to 0.998 U/L.

LINEARITY

This reaction is linear up to the concentration of 400 U/L. For samples with values higher than 400 U/L it is recommended to dilute the sample with Sodium Chloride 0.85%, repeat the dosage and multiply the results obtained by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

Alanine Amino Transferase (ALT) is an enzyme present in many tissues. Damage to the cell membrane leads to the release of this enzyme, which reaches the blood stream and raises the serum level of ALT. The elevation of serum ALT has an acute course (generally proportional to the intensity of the lesion found), and its baseline return occurs in approximately 14 days after removal of the cause. In chronic diseases, the increase in serum ALT is discreet, and may not be evident in the biochemical evaluation of this enzyme. The increase in the ALT enzyme occurs mainly with the occurrence of liver disease, muscle (mainly in large animals) or hyperthyroidism (mainly in cats).

Intravascular hemolysis can cause a false increase in ALT, due to the release of the enzyme present inside the red blood cells. The use of some anticonvulsant drugs, corticosteroids, antibiotics (tetracycline in cats), acetaminophen and caparsolate (in dogs)

can promote changes in the ALT result.

When used concomitantly with the evaluation, their effect on this alteration must be considered.

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QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

Registration Number: Product exempt from registration with the Ministry of Agriculture, Livestock and Supply.

Technical Responsible: Dr. Camila Eckstein (CRMV/MG 20.611)

Review: June/2025

UNIVERSAL SYMBOLOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER