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REF **K093**

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FINALIDADE

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REF K093

USAGE INSTRUCTIONS

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:
SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

FUNCTION

Method for quantitative determination of Lipase in biological samples of serum or plasma heparinized. Colorimetric test, for *in vitro* diagnostics use only.

PRINCIPLE OF ACTION

Methodology: Enzymatic Colorimetric

The Lipase chromogenic substrate 1.2 - O - Dilauril - Rac - Glicero - 3 - Acid - Glutaric - (6 - Methylresorufin) - Ester is cleaved by the catalytic action of alkaline lipase forming 1.2 - O - Dilauril - Rac - Glycerol and a intermediate unstable, Glutaric Acid - (6 - Methylresorufin) Ester. This spontaneously decompose to form Glutaric Acid in alkaline solution and Methylresorufin. The intensity of red color is formed directly proportional to the Lipase activity and can be determined photometrically.

1.2 - O - Dilauril - Rac - Glicero - 3 - Glutaric Acid - (6 - Methylresorufin) Ester $\xrightarrow{\text{Lipase}}$
1.2 - O - Dilauril - Rac - Glycerol + Glutaric Acid - (6 - Methylresorufin) Ester

Glutaric Acid - (6 - Methylresorufin) Ester $\xrightarrow{\text{Spontaneous decomposition}}$
Glutaric Acid + Methylresorufin

REAGENTS

Number 1 – Buffer Reagent – Store between 2 and 8°C. Contains: Buffer BICIN (N,N-Bis(2-Hydroxyethyl)-Glycine) 50 mmol/L (pH 8.0), Colipase > 0.9 mg/L, Na-Deoxycolate 1.6 mmol/L, Calcium Chloride 10 mmol/L, detergent and preservative.

Number 2 – Substrato – Store between 2 and 8°C. Contains: Buffer Tartrate 10 mmol/L (pH 4.0), 1.2-O-Dilauril-Rac-Glicero-3-Ácido Glutaric-(6-Methylresorufin) Ester 0.27 mmol/L, Taurodeoxicolato 8.8 mmol/L, detergent and preservative.

PRESENTATION

Presentation	Reagent N° 1	Reagent N° 2
K093-12-VET	2 x 15 mL	2 x 10 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Thermostated spectrophotometer, pipettes and tips, test tubes, watch or stopwatch, Biocal Bioclin, Biocontrol N and Biocontrol P Bioclin. These items are found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transportation temperature should be between 2 and 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid humidity. **Do not freeze.**

SPECIAL CARE

- 1- For *in vitro* diagnostic use only.**
- Strictly follow the methodology proposed to obtain exact results.
- Water used in the cleansing of the material should be recent and free of contaminants.
- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- Reagent N° 2 must be protected from light.
- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- Do not use the product in case of damaged packaging.
- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or heparinized plasma can be uses. Serum samples are stable for 7 days between 2 and 30°C and 12 months when frozen at -20°C. Plasma samples are stable for 7 days between 2 and 30°C and 12 months when frozen at -20°C.

PROCESS DESCRIPTION

The calibration stability of the Lipase Automated kit installed on refrigerated equipment is at least 28 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

TECHNIQUE

For the use of Lipase Automated, use as calibrator the Biocal Bioclin kit and as control serum the Biocontrol N and Biocontrol P Bioclin kits.

This kit is for use in automated biochemical analyzers. **Check the application sheet for this equipment on www.bioclin.com.br or through SAC.**

INTERFERENT

Studies have been performed to determine the level of interference of different compounds: Conjugate Bilirubin and non conjugate do not produce interference up to 60 mg/dL; Hemoglobin does not produce significant interference up to 500 mg/dL; Lipids do not produce significant interference up to 1000 mg/dL; Rare cases of gammopathy, particularly the IgM type (Waldenström macroglobulinemia) may cause unreliable results.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present an analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The calibration kit must be carried on using the BIOCAL calibrator, which is traceable to the reference material BCR-694 (Community Bureau of Reference).

REFERENCE VALUES

The reference values were determined in serum of healthy dogs, cats, cattle and horses, males and females, adults and mixed breed.

Specie	Reference Value
Dog	19.33 - 106.33 U/L
Cat	2.22 - 21.48 U/L
Cattle	3 - 10 U/L
Horses	3.43 - 11.24 U/L

These values should be used as guidelines and each laboratory should establish its own reference value, according to the population served. The results provided by this kit must be interpreted by the veterinary medical professional, and are not the only criterion for determining the diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

ACCURACY

Comparison of Methods

Lipase Automated kit was compared with other method commercially available for dosage of Lipase. 42 analyzer were performed and the results were evaluated. The linear equation obtained was $Y = 0.8219X + 3.2388$ and correlation coefficient 0.9978. With these results, we can conclude that the kit shows good methodological specificity.

PRECISION

Repeatability

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average Concentration (U/L)	50.30	139.20	11.45
Standard Deviation (U/L)	0.85	1.16	0.81
Coefficient of Variation (%)	1.70	0.83	7.12

Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Average Concentration (U/L)	51.28	144.25	12.78
Standard Deviation (U/L)	0.91	6.93	1.17
Coefficient of Variation (%)	1.77	4.80	9.17

SENSITIVITY

Sensitivity of Lipase Automated kit was calculated from 40 determinations of a sample free of Lipase. The average found, which indicates the method detection limit, was of 2.063 U/L with a standard deviation of 0.046 U/L. The sensitivity, which indicates the method quantification limit, corresponds the average plus 3 times the standard deviation and is equal to 2.201 U/L.

LINEARITY

Test is linear up to 300 U/L.

DIAGNOSTIC SIGNIFICANCE

Lipase is an enzyme produced by the pancreas and whose function is to hydrolyze triglycerides. Damage to pancreatic cells in cases of pancreatitis causes the release of enzymes present in these cells, such as lipase, which reaches the bloodstream. Thus, the increase in serum lipase is a strong indication of pancreatic injury but can also indicate gastrointestinal, hepatic disease

or azotemia, but the mechanism of lipemia elevation in these changes is still unknown. The elevation of lipase in cases of pancreatic alteration begins 24 hours after the injury and reaches its peak between 2 to 5 days.

In calves, pancreatic lipase is little active until the second week of life, however the activity of this enzyme triples its activity in these animals at eight weeks of age. In ruminants, salivary lipase has marked activity in calves and decreases its activity significantly after weaning.

BIBLIOGRAPHIC REFERENCES

- 1 - Steinberg WM, Goldstein SS, Davies ND et al. Diagnostics assays in acute pancreatitis [Review]. Ann Intern Med 1985; 102:576-580.
- 2 - Panteghini M et al. Diagnostic value of four assays for lipase determination in serum: A comparative reevaluation. Clin Biochem 1991; 24:497-503.
- 3 - TIETZ, N. W., Chlinical Guide to Laboratory Tests., 3ª edição. Philadelphia, PA: WB Saunders, 1995: 865.
- 4 - Neumann U et al. New substrates for the optical determination of lipase. EP 207252 (1987).
- 5 - Gargouri Y, Julien R, Bois A et al. Studies on the detergent inhibition of pancreatic lipase activity. J of Lipase Research 1983; 24:1336-1342.
- 6 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone: 0800 0315454 | E-mail: sac@bioclin.com.br

Registration Number: Product exempt from registration with the Ministry of Agriculture, Livestock and Supply.

Technical Responsible: Dr. Camila Eckstein (CRMV/MG 20.611)

Review: June/2025

UNIVERSAL SYMBOLOGY