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REF **K106**

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- 1 - Thomas L ed. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft, 1998:192-202.
- 2 - Tolman KG, Rej R. Liver function. In: Burtis CA, Ashwood ER, editors, Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 1125-77.
- 3 - Rand RN, di Pasqua A. A new diazo method for the determination of bilirubin. Clin Chem 1962; 6:570-8.
4. Kaneko JJ, Harvey JW, Bruss ML. Clinical Biochemistry of Domestic Animals. San Diego: Academic Express, 6. ed., 1998. 936 p.
5. EClínPahth. Cornell University College of veterinary Medicine. 2013.
6. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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1. Thomas L ed. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft, 1998:192-202.
2. Tolman KG, Rej R. Liver function. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 1125-77.
3. Rand RN, di Pasqua A. A new diazo method for the determination of bilirubin. Clin Chem 1962; 6:570-8.
4. Kaneko JJ, Harvey JW, Bruss ML. Clinical Biochemistry of Domestic Animals. San Diego: Academic Express, 6. ed., 1998. 936 p.
5. EclinPanth. Cornell University College of veterinary Medicine. 2013.
6. QUIBASA: Informaciones del Departamento de Investigación y Desarrollo.

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**BILIRUBIN TOTAL AUTOMATED****VETERINARY USE**REF **K106**TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST,  
CONTACT CUSTOMER ADVISORY SERVICE:**INSTRUCTIONS FOR USE**

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

**FUNCTION**

Method for determination of Total Bilirubin in biological samples serum or plasma (EDTA or Heparin). Colorimetric test, for *in vitro* diagnostic use only.

**PRINCIPLE OF ACTION**

**Methodology:** Colorimetric Test - DPD

Total Bilirubin, through the coupling reaction with Dichlorophenyl Diazonium Tetrafluoroborate, forms an azocompound, a red colored complex with maximum absorption at 546 nm. The intensity of the color formed is directly proportional to the concentration of Total Bilirubin in the sample.

**PRESENTATION**

Presentation	Reagent N° 1	Reagent N° 2
K106-1-VET	1 x 40 mL	1 x 10 mL

**REAGENTS**

**Reagent Number 1 - Buffer** - Store between 2 and 8 °C. Contains: Sulfamic Acid < 100 mmol/L, stabilizer, surfactants and preservative.

**Reagent Number 2 - Color Reagent** - Store between 2 and 8 °C. Contains: Dichlorophenyl Diazonium Tetrafluoroborate < 10 mmol/L, stabilizer and surfactant.

**EQUIPMENTS AND OPERATIONAL INPUTS**

To perform this technique it is required the Biochemistry Automated Equipment, Biocal kit and Biocontrol N and P Bioclin.

**TRANSPORTATION AND STORAGE CONDITIONS**

The storage temperature should be between 2 to 8 °C. The transport at temperatures up to 30 °C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

**SPECIAL CARE**

**1- For *in vitro* diagnostic use only.**

**2- Reagent N° 2** must be kept protected from light.

**3- Reagent N° 1** may show turbidity when subjected to low temperatures. If this occurs, allow the kit to stabilize at room temperature (15 to 30 °C).

**4- Water** used in cleaning must be recent and free of contaminants.

**5- Saturated deionizer columns** release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

**6- Hemolyses**, even if discreet, may interfere in dosage.

**7- Reagent N° 1** contains Sodium Azide, irritant to skin and mucous. Handle reagent with care.

**8- We recommend** applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

**9-** To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.

**10-** Do not use the product in case of damaged packaging.

**11-** It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

**SAMPLES**

Serum or plasma (EDTA or Heparin) free of hemolysis.

**Sample must be protected from light.**

The samples is stable for 3 days between 2 and 8°C and 3 months at -20°C.

Sample for therapeutic control should always be collected at the same hour.

**PROCESS DESCRIPTION**

The calibration stability of the Bilirubin Total Automated kit installed on refrigerated equipment is at least 30 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

**PREPARATION OF REAGENTS**

Reagents are ready for use.

**TECHNIQUE**

For the use of the kit, use as calibrator the Biocal Bioclin kit and as control serum the Biocontrol N and P Kit Bioclin.

This kit is for use in automated biochemical analyzers.

**Check the application sheet for this equipment on [www.bioclin.com.br](http://www.bioclin.com.br) or through SAC.**

**INTERFERING**

Hemolysis, even if slight, interferes with the bilirubin dosage. No interference was observed by triglycerides up to 1500 mg/L. Exposing the sample to light may lead to falsely diminished results.

**INTERNAL QUALITY CONTROL**

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present an analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

**TRACEABILITY**

Bioclin recommends the usage of the BIOCAL calibrator which is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 916A.

**REFERENCE VALUES**

The reference values were determined in serum samples of healthy populations of dogs, cats, cattle and horses, males and females, adults and no breed preference.

Specie	Reference Value
Dogs	0.046 - 0.257 mg/dL
Cats	0.05 - 0.29 mg/dL
Cattle	0.07 - 0.7 mg/dL
Horses	0.067 - 1.54 mg/dL

To convert the values from mg/dL into mmol/L (SI), multiply the results obtained for 0.0171.

These values must be used as a guide and each laboratory must establish its own reference value, according to the population served. The results provided by this kit must be interpreted by a veterinary professional, and is not the only criterion to determine the diagnosis and/or treatment of the patient.

**PRODUCT PERFORMANCE ACCURAY****Comparison of Methods**

Bilirubin Total Automated kit was compared with other commercially available methods for measurement of Total Bilirubin. 42 analyses were performed and the results were evaluated. The linear equation obtained was  $Y = 0.9983X - 0.0654$  and the correlation coefficient 0.9955. With these results, we can conclude that the kit shows good methodological specificity.

## PRECISION

### Repeatability

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
<b>Average Concentration (mg/dL)</b>	1.24	2.06	0.78
<b>Standard Deviation (mg/dL)</b>	0.02	0.02	0.01
<b>Coefficient of Variation (%)</b>	1.45	0.97	1.69

### Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
<b>Average Concentration (mg/dL)</b>	1.27	2.05	0.80
<b>Standard Deviation (mg/dL)</b>	0.04	0.05	0.03
<b>Coefficient of Variation (%)</b>	3.16	2.50	3.51

## SENSITIVITY

The sensitivity was calculated from 40 determinations of a sample free of Bilirubin Total. The average found was 0.044 mg/dL with a standard deviation of 0.005 mg/dL. The sensitivity, which indicates the method detection limit, correspond the average plus 3 times the standard deviation, is equal to 0.058 mg/dL.

## LINEARITY

The reaction is linear up to the concentration of 30 mg/dL. For higher values dilute the sample with NaCl 0.85% and repeat the dose. Multiply the resulting number by the dilution factor used.

## DIAGNOSTIC SIGNIFICANCE

Bilirubin is the result of the degradation of hemoglobin present in red blood cells by the reticuloendothelial system. After the oxidation and reduction process, hemoglobin is transformed into bilirubin, which is conjugated in hepatocytes and eliminated in the intestinal lumen. Total bilirubin is made up of direct bilirubin (which has been conjugated) and indirect bilirubin (not conjugated). Changes in normal bilirubin concentration are strong indicators of hepatobiliary changes and other causes of jaundice.

The hepatobiliary changes can be of origin:

- Prehepatic: resulting from hemolytic processes in which the conjugated fraction of bilirubin exceeds the conjugated fraction.
- Hepatic: occurs due to failure in the process of conjugation of hemoglobin in the hepatocytes in cases of failure in the functioning of these cells.
- Post-hepatic: It occurs due to defects in the elimination of conjugated bilirubin. This is usually due to partial or

total obstruction in the biliary tract.

In cattle, hemolytic anemia caused by anaplasmosis is an important cause of hyperbilirubinemia due to hemolytic anemia caused by this pathogen. Another important cause of the increase in total circulating bilirubin is hepatic lipidosis and copper intoxication.

In the caballos, the elevation of bilirubin can occur for prolonged periods of time, and therefore it can occur secondary to other pathologies when the animal stops feeding. Neonatal isoerythrolysis and pathologies that cause acute hepatitis are common causes of bilirubin elevation in this species.

In dogs and cats, the elevation of bilirubin is associated with hemolytic, inflammatory and infectious liver diseases and biliary obstructions of various causes. In dogs, it is common that the elevation of bilirubin in the blood occurs in an acute way, and that after 2 to 3 days there is a renal compensation between the rhythm of production and excretion of bilirubin, which tends to return the serum concentration normal after this period. Unlike dogs, this compensatory mechanism does not occur in cats and posthepatic hyperbilirubinemia may persist until the cause is eliminated.

## BIBLIOGRAPHIC REFERENCES

1. Thomas L ed. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft, 1998:192-202.
2. Tolman KG, Rej R. Liver function. In: Burtis CA, Ashwood ER, editors, Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 1125-77.
3. Rand RN, di Pasqua A. A new diazo method for the determination of bilirubin. Clin Chem 1962; 6:570-8.
4. Kaneko JJ, Harvey JW, Bruss ML. Clinical Biochemistry of Domestic Animals. San Diego: Academic Express, 6. ed., 1998. 936 p.
5. EClínPahth. Cornell University College of veterinary Medicine. 2013.
6. QUIBASA: Data from Research and Development.

## QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

### QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

## CUSTOMER SERVICE

Customer Advisory Service  
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

**Registration Number:** Product exempt from registration with the Ministry of Agriculture, Livestock and Supply.

**Technical Responsible:** Camila Eckstein (CRMV/MG 20.611).

**Review:** June/2025

## UNIVERSAL SYMBOLOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER