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1. Thomas L ed. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft, 1998:192-202.
2. Tolman KG, Rej R. Liver function. In: Burtis CA, Ashwood ER, editors, Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 1125-77.
3. Rand RN, di Pasqua A. A new diazo method for the determination of bilirubin. Clin Chem 1962; 6:570-8.
4. Kaneko JJ, Harvey JW, Bruss ML. Clinical Biochemistry of Domestic Animals. San Diego: Academic Express, 6. ed., 1998. 936 p.
5. ECLinPanth. Cornell University College of veterinary Medicine. 2013.
6. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento

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1. Thomas L ed. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft, 1998:192-202.
2. Tolman KG, Rej R. Liver function. In: Burtis CA, Ashwood ER, editors, Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 1125-77.
3. Rand RN, di Pasqua A. A new diazo method for the determination of bilirubin. Clin Chem 1962; 6:570-8.
4. Kaneko JJ, Harvey JW, Bruss ML. Clinical Biochemistry of Domestic Animals. San Diego: Academic Express, 6. ed., 1998. 936 p.
5. EClinPath. Cornell University College of veterinary Medicine. 2013.
6. QUIBASA: Informaciones del Departamento de Investigación y Desarrollo.

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**BILIRUBIN DIRECT AUTOMATED**REF **K107****INSTRUCTIONS FOR USE****FUNCTION**

Method for determination of Direct Bilirubin in biological samples serum or plasma (EDTA or Heparin). Colorimetric test, for in vitro diagnostic use only.

**PRINCIPLE OF ACTION**

**Methodology:** Colorimetric Test - DPD

Direct Bilirubin, through the coupling reaction with Dichlorophenyl Diazonium Tetrafluoroborate forms an azo compound, a red colored complex with maximum absorption at 546 nm. The intensity of color formed is directly proportional to the concentration of Direct Bilirubin in the sample.

**REAGENTS**

**Number 1 - Buffer** - Store between 2 and 8 °C. Contains: Sulfamic Acid < 200 mmol/L and stabilizers.

**Number 2 - Color Reagent** - Store between 2 and 8 °C. Contains: Dichlorophenyl Diazonium Tetrafluoroborate < 5 mmol/L and stabilizer.

**PRESENTATION**

Presentation	Reagent N° 1	Reagent N° 2
K107-1-VET	1 x 40 mL	1 x 10 mL

**EQUIPMENTS AND OPERATIONAL INPUTS**

To perform this technique it is required the biochemistry automated equipment, Biocal kits and Biocontrol N and P Bioclin.

**TRANSPORTATION AND STORAGE CONDITIONS**

The storage temperature should be between 2 to 8 °C. The transport at temperatures up to 30 °C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

**SPECIAL CARE**

- 1- For *in vitro* diagnostic use only.**
- 2- The reagent N° 2 must be kept protected from light.**
- 3- Water used in material cleaning must be recent and free of contaminants.**
- 4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.**
- 5- Calibration must be repeated periodically to verify any alterations in colorimetric responses or spectrophotometer.**
- 6- Hemolysis, even if discreet, interfere in dosage.**
- 7- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.**
- 8- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.**
- 9- Do not use the product in case of damaged packaging.**
- 10- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.**

**SAMPLES**

Serum or plasma (EDTA or Heparin) free of hemolysis. **Sample must be protected from light.** The samples is stable for 3 days between 2 and 8 °C and 3 months at -20 °C. Sample for therapeutic control should always be collected at the same hour.

**PROCESS DESCRIPTION**

The calibration stability of the Bilurubin Direct Automated kit installed on refrigerated equipment is at least 30 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

**PREPARATION OF REAGENTS**

Reagents are ready for use.

**TECHNIQUE**

For the use of the kit, use as calibrator the Biocal Bioclin kit and as control serum the Biocontrol N and P Bioclin kit.

This kit is for use in automated biochemical analyzers only. **Check the application sheet for this equipment on the website [www.bioclin.com.br](http://www.bioclin.com.br) or through SAC 0800 031 5454.**

**INTERFERENCES**

The hemolysis, even if it is mild, interferes with the dosage of Bilirubin. No interference will be observed up to 1000 mg/L of Triglycerides. Exposing the sample to light may lead to falsely reduced bilirubin results.

**INTERNAL QUALITY CONTROL**

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present an analytical variety, and it must be monitored by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

**TRACEABILITY**

Bioclin recommends the usage of the BIOCAL calibrator which is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 916A.

**REFERENCE VALUES**

Reference values were established in serum samples from healthy populations of dogs, cats, cattle and horses, males and females, adults and without racial predilection.

Specie	Value
Dog	0 - 0.06 mg/dL
Cat	0 - 0.04 mg/dL
Cattle	0.02 - 0.7 mg/dL
Horse	0.13 - 0.547 mg/dL

**VETERINARY USE**

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / [sac@bioclin.com.br](mailto:sac@bioclin.com.br)

To convert the values to mg/dL to mmol/L (SI), multiply the results obtained by 0.0171.

These values should be used only as a guide. Each laboratory should establish its own reference range, according to the population served. The results provided by this kit must be interpreted by the veterinary responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

**PRODUCT PERFORMANCE ACCURACY****Comparison of Methods**

Bilirubin Direct Automated kit was compared with other methods for determination of Direct Bilirubin dosage commercially available measurement. 42 analyzes were performed and the results were evaluated. The equation linear obtained was  $Y = 1.0228X + 0.0034$  and the correlation coefficient 0.9916. With these results, we can conclude that the kit shows good methodological specificity.

**PRECISION****Repeatability**

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	0.42	0.69	0.58
Standard Deviation (mg/dL)	0.01	0.01	0.01
Coefficient of Variation (%)	1.50	1.13	2.14

## Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
<b>Average Concentration (mg/dL)</b>	0.43	0.70	0.60
<b>Standard Deviation (mg/dL)</b>	0.01	0.03	0.02
<b>Coefficient of Variation (%)</b>	1.51	4.75	3.43

## SENSITIVITY

The sensitivity was calculated from 40 determinations of a sample free of Direct Bilirubin. The average found was 0.055 mg/dL, with deviation standard of 0.005 mg/dL. The sensitivity, which indicates the method limit detection, corresponds the average plus 3 times the standard deviation, and is equal to 0.070 mg/dL.

## LINEARITY

The reaction is linear to the concentration of 12 mg/dL. For higher values dilute the sample with NaCl 0.85% and repeat the dose. Multiply the resulting number by the dilution factor used.

## DIAGNOSTIC SIGNIFICANCE

The elimination of red blood cells in the mononuclear phagocytic system results in the release of bilirubin, which is transported to the liver where it undergoes the conjugation process that is soluble in water. Conjugated or direct bilirubin is the metabolized fraction that is excreted in the bile, what makes its serum concentration and normal urine low in most species, except in dogs where it is common to detect it in the urine (in low concentrations). In horses and cattle, prolonged stay can promote a physiological increase in bilirubin with clinical presentation of jaundice without hepatic changes.

Hemolysis is considered to be the main cause of the increase in bilirubin. In the case of severe hemolysis in a short period of time, it is common that the predominant bilirubin is indirect. In small animals and in foals with neonatal isoerythrolysis (prolonged hemolysis) the appearance of hepatobiliary alterations with direct and indirect bilirubin accumulation is frequent. High concentrations of bilirubin for a prolonged period can cause kernicterus (an action resulting from the deposition of bilirubin in the central nervous system and which can cause death in an animal). Hepatic changes can result in an increase in the concentration of direct and indirect bilirubin, and generally occur simultaneously with an increase in other enzymes such as alkaline phosphatase and GGT.

Cholestasis is an important cause of increased bilirubin and can be characterized as intrahepatic (direct change in hepatocytes) or extrahepatic (eg, obstruction by neoplasms, gallstones or parasites such as Fasciola hepatica in cattle and horses). Anatomical or hereditary defects can also occur in the conjugation process and lead to an increase in the concentration of serum bilirubin.

## BIBLIOGRAPHIC REFERENCES

1. Thomas L ed. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft, 1998:192-202.

2. Tolman KG, Rej R. Liver function. In: Burtis CA, Ashwood ER, editors, Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 1125-77.
3. Rand RN, di Pasqua A. A new diazo method for the determination of bilirubin. Clin Chem 1962; 6:570-8.
4. Kaneko JJ, Harvey JW, Bruss ML. Clinical Biochemistry of Domestic Animals. San Diego: Academic Express, 6. ed., 1998. 936 p.
5. EclinPalth. Cornell University College of veterinary Medicine. 2013.
6. QUIBASA: Data from Reseach and Development.

## QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

### QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

## CUSTOMER SERVICE

Customer Advisory Service  
Phone: 0800 0315454 | E-mail: sac@bioclin.com.br

**Registration Number:** Product exempt from registration with the Ministry of Agriculture, Livestock and Supply.

**Technical Responsible:** Camila Eckstein (CRMV/MG 20.611).

**Review:** June/2025

## UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MADE BY

LOT NUMBER

CONTROL

MANUFACTURING DATE

POSITIVE CONTROL

VALIDITY DATE  
(last day of the month)

NEGATIVE CONTROL

TEMPERATURE LIMIT  
(store)

BIOLOGICAL RISK

CONTENT IS SUFFICIENT  
FOR <N> TEST

FLAMMABLE

SEE INSTRUCTIONS  
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC  
PRODUCT

TOXIC

KEEP AWAY  
FROM SUNLIGHT

DO NOT USE IF  
PACKAGE IS  
DAMAGED

DO NOT REUSE

PRODUCT  
STERILIZED

CAUTION

DANGER