FOSFATASE ALCALINA IFCC

REF **K224**

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FINALIDADE

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INSTRUCTIONS FOR USE

FUNCTION

Method for determination of Alkaline Phosphatase. In serum or plasma (Heparin) samples. Kinetic test, only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Kinetic IFCC.

The Alkaline Phosphatase catalyzes the transfer of Phosphate group from the substrate p-Nitrophenylphosphate (pNPP) to 2-Amino-2-Methyl-1-Propanol (AMP), releasing p-Nitrophenol according to the equation below:

REAGENTS

Number 1 - Buffer - Store between 2 and 8°C. Contains: 2-Amino 2-Methyl Propanol <1.0 mol/L, activators, stabilizer and preservative.

Number 2 - Substrate - Store between 2 and 8 °C. Contains: p-NFF Substrate p-Nitrophenylphosphate.

PRESENTATION

Presentation	Reagent N° 1	Reagent N° 2
K224-2-VET	2 x 40 mL	2 x 10 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Thermostated spectrophotometer, pipettes, atch or stopwatch, test tubes, Biocontrol N, Biocontrol P Bioclin and Biocal Bioclin They can be found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures between 15 and 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must to be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers, that may alter the results significantly.

5- It is important, for the good development of the test, a rigorous control of time and temperature.

6- We recommend applying the local, state and federal rules

for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

7- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

8- Do not use the product in case of damaged packaging.

9- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or plasma (Heparin) free of hemolysis. The analyte is stable for 7 days if kept in temperature between 2 and 8°C and 30 days if kept at 10°C negative.

PROCESS DESCRIPTION

The calibration stability of the Alkaline Phosphatase IFCC kit installed on refrigerated equipment is up to 3 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

PREPARATION OF WORKING REAGENT

Mix four (4) parts of Reagent N° 1 with one (1) part of Reagent N° 2. Working Reagent is stable during 14 days, if kept between 2 and 8 °C. Store away from light.

REACTION CONDITIONS

It is indispensable conditions to use a thermostated cuvette at 37°C, 1 cm optical path and reading at 405 nm.

TECHNIQUE

Bioclin recommends using the kit Biocal Bioclin as calibrator and as control serum, Biocontrol N and P Bioclin Kits. Add 20 µL of Sample to 1.0 mL of Working Reagent, mix and transfer to a thermostated cuvette at 37 °C and wait for 1 minute. Make the initial Reading simultaneously starting the stopwatch. Repeat readings after 1, 2 and 3 minutes. Calculate the mean differences in absorbance per minute (Abs/min.) and use it to calculate the result.

CALCULATIONS

ALP (U/L) = Abs/min. x 2757

Results are expressed as U/L.

LIMITATIONS OF THE PROCEDURE

The kinetic method is based on the molar absorptivity. By this reason, the readings must be conducted in a spectrophotometer that satisfies the following conditions:

Wavelength 405 nm

Semi trajectory of the pass band 10 nm

Stray light less than 0.5%

1 cm thermostated cuvette.

INTERFERENCES

Do not use hemolyzed or lipemic samples, as the results obtained may be falsely high.

EDTA, citrate and oxalate anticoagulants should be avoided.

The use of anticonvulsant drugs (phenobarbital and primidone) and corticosteroids can affect AF activity.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit calibration can be made using the theoretical calibration factor based on the molar absorptivity of the p-Nitrophenol, or through the Biocal calibrator. Bioclin recommends the usage of the Biocal calibrator which is also traceable to the molar absorptivity of the p-Nitrophenol.

REFERENCE VALUES

The reference values were determined in serum samples of healthy dogs, cats, cattle and horses, males and females, adults without predilaction of breed.

Specie	Reference Value
Dog	15 - 124.6 U/L
Cat	14.5 - 96 U/L
Cattle	28.3 - 114.6 U/L
Horses	74.6 - 406 U/L

These values should be used as a guideline and each laboratory should establish its own reference value, according to the population served. The results provided by this kit must be interpreted by the veterinary medical professional, and are not the only criterion for determining the diagnosis and/or treatment of the patient.

VETERINARY USE

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

PRODUCT PERFORMANCE QUALITY CONTROL ACCURACY

Comparison Of Methods And Specificity Methodology

The Alkaline Phosphatase IFCC was compared with other commercially available methods for measurement of Phosphatase. 42 analyzes were performed and the results were evaluated. The linear equation obtained was $y = 1.1027x - 0.2651$ and the correlation coefficient 0.9999. With these results, we can conclude that the kit shows good methodological specificity.

PRECISION

Repeatability

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average Concentration (U/L)	153.1	145.9	327.9
Standard Deviation (U/L)	1.7	1.8	2.6
Coefficient of Variation (%)	1.14	1.22	0.78

Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Average Concentration (U/L)	154.5	147.4	328.9
Standard Deviation (U/L)	1.7	1.4	0.9
Coefficient of Variation (%)	1.13	0.92	0.26

SENSITIVITY

The sensitivity was calculated from 40 determinations of a sample free the Phosphatase. The standard deviation was 0.83 U/L. The sensitivity, that indicates the method detection limit, corresponds to 3 times the standard deviation is equal to 2.48 U/L.

LINEARITY

The reaction is linear up to the concentration of 1500 U/L. For sample with value above 1500 U/L, we recommend diluting the sample with NaCl 0.85%. Multiply the results obtained by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

Alkaline Phosphatase (AP) is present in high concentrations in bones, liver, intestine, kidneys and placenta. The increase in the concentration of this enzyme is associated with cholestasis, increase in non-osteolytic and steroid osteoblastic activity (in the canine species). In felines, the increase in FA occurs in situations similar

to what occurs in dogs (liver diseases, cholangitis and cholangiohepatitis, bile duct obstructions and tumors), however, in cats the elevation of this enzyme occurs less markedly and for shorter periods, due to the lower amount of FA in the hepatocytes of this species. In cattle and horses, FA has lower specificity and generally other enzymes such as bilirubin are elevated earlier in the case of cholestasis.

Physiological conditions in which a non-pathological increase in alkaline phosphatase can occur include: Age (young animals can have AF up to 10 times that of an adult animal), Breeds (Siberian Huskie and Scottish Terrier dogs are predisposed to Beningian hyperphosphatemia, reaching up to 700 U/L) and Stress.

BIBLIOGRAPHIC REFERENCES

- 1 - International Federation of Clinical Chemistry (1978). Section N° 2. Provisional recommendations on the theory of reference values. Part. I. The concept of reference values. Expert panel on the Theory of Reference Values – Clin. Chim. Acta 87/3: 459F (1978).
- 2 - Spandrio, L. et al – Considerazioni sulla determinazione della fosfatasi alcalina del siero con método cinético – LAB2/2:123 (1975).
- 3 - Gardner, M. D. et al – Frequency distribution and 'reference values' of plasma alkaline phosphatase (EC3.1.3.1) activity in the adult population of Scottish new town – J. Clin. Pathol. 31/12:1202 (1978).
- 4 - Demaría, I; Setta, F.; Lorenzo, L. – VIII Congreso Argentino de Bioquímica, Octubre 1990 – Revista de la Asociación Bioquímica Argentina 54/3 (1990).
- 5 - Kaneko JJ, Harvey JW, Bruss ML. Clinical Biochemistry of Domestic Animals. San Diego: Academic Express, 6. ed., 1998. 936 p.
- 6 - EClinPahth. Cornell University College of veterinary Medicine. 2013.
- 7 - Burtis, C. A., Ashwood, E. R., & Bruns, D. E. (2008). Tietz fundamentos de química clínica. B. G. Sawyer (Ed.). Rio de Janeiro: Elsevier.
- 8 - Nestor DD, Holan KM, Johnson CA, Schall W, Kaneene JB. (2006). Serum alkaline phosphatase activity in Scottish Terriers versus dogs of other breeds. Journal of the American Veterinary Medical Association, 228(2), 222–224.
- 9 - Lawler DF, Keltner DG, Hoffman WE, et al. Benign familial hyperphosphatemia in Siberian Huskies. Am J Vet Res 1996;57:612–617.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP: 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

Registration Number: Product exempt from registration with the Ministry of Agriculture, Livestock and Supply.

Technical Responsible: Dr. Camila Eckstein (CRMV/ MG 20.611)

Review: June/2025

UNIVERSAL SYMBOLOLOGY