AMILASE G7 IFCC

REF K242

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FINALIDADE

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1-Schumann G, Klauke R, Canalias F. , et al. "IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 °C: International Federation of Clinical Chemistry and Laboratory Medicine (IFCC): Scientific Division, Committee on Reference Systems for Enzymes (C-RSE): Part 8. Reference procedure for the measurement of catalytic concentration of α -amylase: [α -Amylase: 1,4- α -Dglucan 4-glucanohydrolase (AMY), EC 3.2.1.1]". Clin Chem Lab Med, 44(9), 2006.

2-WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.

3-Benson E, Tonks DB. Laboratory medicine: Quality Control System in Clinical Chemistry Laboratories. Postgraduate Medicine. 34:4. 1963.

4-EClinPalth. Cornell University College of veterinary Medicine. 2013.

5- Kaneko JJ, Harvey JW, Bruss ML. Clinical Biochemistry of Domestic Animals. San Diego: Academic Express, 6. ed., 1998. 936 p.

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INSTRUCTIONS FOR USE



FUNCTION

Test to determine Alpha - Amylase activity in biological samples (serum, plasma, urine), through kinetic test. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: IFCC Kinetics

Alpha-Amylase cleaves substrate 4 - Ethidylene - NP-G7 (4NP-G7) in free oligosaccharides (G5, G4 and G3); and 4NP-G5, 4NP-G4 and 4NP-G3. The NP and G5, G4 and G3 fragments are hydrolyzed by the combined action of amylase and Alpha-Glycosidase, producing Glucose and free p-Nitrophenol (NP), which leads to an increase in absorbance. The rate of NP formation is proportional to the total alphaamylase activity present in the sample.

4NP - Nitrofenol

G - Glicose

REAGENTS

Number 1 - Enzyme Buffer - Store at 2 to 8 °C. Contains: Buffer > 10 mmol/L, Sodium Chloride > 70 mmol/L, Magnesium Chloride > 1 mmol/L, Alpha Glycosidase < 50 KU/L, surfactant and stabilizer.

Number 2 - Substrate - Store at 2 to 8 °C. Contains: Buffer > 10 mmol/L, 4 - Ethidylene-NP-G7 < 10 mmol/L, surfactant and preservative.

PRESENTATION

Presentation	Reagent 1	Reagent 2
K242-2-VET	1 x 40 mL	1 x 10 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Thermostatic spectrophotometer, pipettes, clock or stopwatch, test tubes, Biocontrol N, Biocontrol P Bioclin and Biocal Bioclin. These items are found in the specialized market of articles for Clinical Analysis Laboratories.

TRANSPORATION AND STORAGE CONDITIONS

The storage temperature should be 2 to 8 °C. Transport at temperatures up to 30 °C should not exceed 5 days. Keep away from light and avoid humidity. **Do not freeze**

SPECIAL CARES

1-For *in vitro* diagnostic use only.

2-Strictly follow the proposed methodology to obtain exact results.

3-The water used to clean the material must be fresh and free from contaminants.

4-Saturated deionizing columns release alkaline water, various ions and oxidizing and reducing agents, which can significantly alter the results.

5-Visible hemolysis can cause slight variations in results.

6-Reagent N° 1 and 2 are easily contaminated by saliva. Therefore, we recommend the use of masks, pears and automatic pipettes when handling this product.

7-We recommend applying local, state and federal environmental protection standards so that the disposal of reagents and biological material is carried out in accordance with current legislation.

8-To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

9-Do not use the product in case of damage to the packaging.

10-It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum obtained free of hemolysis, plasma collected with heparin. Other anticoagulants (EDTA, Citrate, Oxalate) inhibit the activity of alpha amylase. Stable for 7 days between 15 and 30 °C, 2 months between 2 and 8 °C and 1 year at -20 °C.

Urine. Stable for 2 days at 15 to 30 °C, 10 days at 2 to 8 °C and 3 weeks at -20 °C.

PROCESS DESCRIPTION

The calibration stability of the Amylase G7 IFCC kit installed in refrigerated equipment is up to 10 days. This stability may vary according to the conditions of the reagent, the equipment and the environment. Therefore, it is suggested to monitor the product's performance using control.

WORKING REAGENT PREPARATION

Mix four (4) parts of Reagent N° 1 with one (1) part of Reagent N° 2. The working reagent is stable for 14 days at 2 to 8 °C. Store away from light.

REACTION CONDITIONS

It is essential to use a thermostated cell at 37 °C, an optical path of 1 cm and reading at 405 nm.

TECHNIQUE

Bioclin recommends, for product use, to use the Biocontrol N and P Bioclin kits as control serum.

Add 20 µL of Sample to 1.0 mL of the Working Reagent, mix and transfer to a thermostated cell at 37 °C and wait 2 minutes. Take the initial reading, simultaneously triggering the stopwatch. Repeat the readings after 1, 2 and 3 minutes. Use the same technique to determine the reaction blank, using water as a sample.

Calculate the average of the absorbance differences per minute (Abs/min.) and use it to calculate the result.

CALCULATIONS

Serum or Plasma

Using Factor

Abs/min. = Average (Abs/min Sample) - Average (Abs/min Water)

Alpha Amylase Activity (U/L) = Abs/min x Factor

Factor 5,040

Using Calibrator

Alpha Amylase Activity (U/L) =

$\frac{\text{Average (Abs/min. Sample)} - \text{Average (Abs/min Water)}}{\text{Average (Abs/min. Calibrator)} - \text{Average (Abs/min Water)}} \times \text{Cal Conc.}$

Average (Abs/min. Calibrator) - Average (Abs/min Water)

Cal Conc. = Calibrator Concentration in U/L

Urine

Using Factor

Abs/min. = Average (Abs/min Sample) - Average (Abs/min Water)

Alpha Amylase Activity (U/L) = Abs/min x Factor

Factor 10,080

Using Calibrator

Alpha Amylase Activity (U/L) =

$\frac{\text{Average (Abs/min. Sample)} - \text{Average (Abs/min Water)}}{\text{Average (Abs/min. Calibrator)} - \text{Average (Abs/min Water)}} \times \text{Cal Conc.}$

Average (Abs/min. Calibrator) - Average (Abs/min Water)

Cal Conc. = Calibrator Concentration in U/L

Urinary Amylase (U/h) = $\frac{\text{U/L} \times \text{V Amylase}}{\text{H} \times 1000}$

V = Urinary Volume in mL

H = Number of hours of urinary collection

Example

Urinary Amylase (value found in the dosage) = 90 U/L

V = 200 mL

H = 2

Urinary Amylase (U/h) = $(90 \times 200) / (2 \times 1000) = 9 \text{ U/h}$

LIMITATIONS OF THE PROCESS

The specifications below refer to semi-automatic equipment:

The kinetic method is based on molar absorptivity, for this reason, readings must be performed on a spectrophotometer that meets the following conditions:

Wavelength 405 nm

Semi-path of the 10 nm passband

Spurious light less than 0.5%

Thermostatic 1cm bucket

INTERFERENT

Saliva contains Amylase. Therefore, you should not pipette by mouth and avoid contact of the reagent with the skin. Lipemia (Triglycerides up to 1000 mg/dL) and Bilirubin (up to 50 mg/dL) do not interfere with the methodology. Hemoglobin up to 600 mg/dL and Ascorbic Acid up to 45 mg/dL do not interfere with the results. Some drugs can interfere with laboratory results by increasing serum Amylase levels (morphine, meperidine, codeine, thiazide diuretics) or decreasing, as in cases of barbiturate poisoning.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to note that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. For that, it is recommended to use controls, which allow to evaluate the precision and accuracy of the dosages.

TRACEABILITY

The calibration of the kit can be done using the theoretical calibration factor, based on the molar absorptivity of 4 - Ethidylene-NP-G7, or through the BIOCAL calibrator. Bioclin recommends using the BIOCAL calibrator, which is traceable to the Reference method described in the IFCC: SCHUMANN, G., et al. "IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 °C:

International Federation of Clinical Chemistry and Laboratory Medicine (IFCC): Scientific Division, Committee on Reference Systems for Enzymes (C-RSE): Part 8. Reference procedure for the measurement of catalytic concentration of α -amylase: [α -Amylase: 1,4- α -D-glucan 4-glucanohydrolase (AMY), EC 3.2.1.1]". Clinical Chemistry and Laboratory Medicine (CCLM), vol. 44, no 9, janeiro de 2006.

REFERENCE VALUES

The reference values were determined in healthy dogs, cats, cattle and horses, males and females, adults and with no predilection for breed.

Specie	Reference Value
Dogs	562.5 - 2281 U/L
Cats	787 - 2609 U/L
Cattle	9.0 - 87 U/L
Horses	3 - 8 U/L

These values should be used as a guideline and each laboratory should establish its own reference value, according to the population served. The results provided by this kit must be interpreted by the veterinary medical professional, and are not the only criterion for determining the diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE ACCURACY

Comparison of Methods

The Amylase G7 IFCC kit was compared with another commercially available method for measuring Amylase. 40 analyzes were performed and the results were evaluated. The linear equation obtained was $Y = 0.9138X - 3.3580$, with a linear correlation coefficient equal to 1.000. With these results, it can be concluded that the kit has good methodological specificity.

PRECISION

Repeatability

Repeatability was calculated from the analysis of 10 repetitions of 3 samples of different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average Concentration (U/L)	312	751	1467
Standard Deviation (U/L)	2.3	3.9	12.8
Coefficient of Variation (%)	0.74	0.52	0.88

Reproducibility

Reproducibility or was calculated from the analysis of 3 samples of different levels during 15 days, not necessarily consecutive, totaling 60 determinations for each sample, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Average Concentration (U/L)	156	378	753
Standard Deviation (U/L)	2.7	6.0	9.9
Coefficient of Variation (%)	1.72	1.59	1.32

SENSITIVITY

Sensitivity was calculated from 40 determinations of a sample free from the presence of amylase. The standard deviation was 1.73 U/L. The sensitivity, which indicates the detection limit of the method, corresponds to 3 times the standard deviation 5.2 U/L.

LINEARITY

The reaction is linear up to 2000 U/L. For values above this value, dilute the sample 1:10 in physiological saline and multiply the result found by the dilution factor (10). For example, if the sample is diluted 1:10 and the result obtained is 80 U/L, multiply by 10 (80 x 10) U/L and consider the final value of the activity value of 800 U/L.

DIAGNOSTIC MEANING

Amylase is an enzyme found in high concentrations in the pancreas, but can also be found in other organs such as the kidneys, intestines and liver. Amylase is found in high concentrations in human and swine saliva, but it is not relevant in other animal species.

Amylase is produced and stored inside the vesicles of pancreatic cells, which prevents its direct contact with the cell cytoplasm and adjacent cells. Acute injuries to the pancreatic cell promote the release of the contents of these vesicles, leading to contact of amylase with the bloodstream, leading to marked increase in this enzyme in the serum within 48 hours after the injury. In dogs, increased serum amylase levels can also occur in cases of kidney disease, diabetes *mellitus* and tumors. Chronic renal failure and reduced glomerular filtration rate can lead to increased amylase concentration, as well as in cases of inflammatory and obstructive intestinal changes. The use of corticosteroids such as prednisone and dexamethasone can result in a false reduction in amylase activity in this species. In felines in cases of acute pancreatitis are rare, but can occur in cases of

abdominal trauma, hypercalcemia (increased calcium in the blood), poisoning, ischemia, cholangiohepatitis (inflammation of the liver and bile ducts), inflammatory bowel disease, parasitic, bacterial, fungal and viral infections that directly or indirectly affect pancreatic tissue, and conditions in other abdominal organs can lead to the occurrence of pancreatitis due to proximity. In cattle, the main cause of pancreatitis is the presence of parasites in the pancreatic duct, such as *Eurytrema* sp., and lesions in adjacent tissues. In horses, pancreatitis generally occurs secondary to gastrointestinal, hepatic or endocrine pathological changes.

BIBLIOGRAPHIC REFERENCES

1-Schumann G, Klauke R, Canalias F. , et al. "IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 °C: Part 8. Reference procedure for the measurement of catalytic concentration of α -amylase: [α -Amylase: 1,4- α -Dglucan 4-glucanohydrolase (AMY), EC 3.2.1.1]". Clin Chem Lab Med, 44(9), 2006.

2-WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.

3-Benson E, Tonks DB. Laboratory medicine: Quality Control System in Clinical Chemistry Laboratories. Postgraduate Medicine. 34:4. 1963.

4-EClinPanth. Cornell University College of veterinary Medicine. 2013.

5- Kaneko JJ, Harvey JW, Bruss ML. Clinical Biochemistry of Domestic Animals. San Diego: Academic Express, 6. ed., 1998. 936 p.

6 - QUIBASA: Data from the Research and Development Department

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Quality Control Department. The quality of the reagents is ensured until the expiration date mentioned on the presentation packaging, provided they are stored and transported under the appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.
br. CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

Registration Number: Product exempt from registration with the Ministry of Agriculture, Livestock and Supply.

Technical Responsible: Dra. Camila Eckstein (CRMV/ MG 20.611)

Review: October//2025

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N>- TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER