COLESTEROL MONORREAGENTE

REF **K306**

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FINALIDADE

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INSTRUCTIONS FOR USE

FUNCTION

Method for determination of Cholesterol in biological samples of serum or plasma (Heparin). Enzymatic colorimetric test, for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Enzymatic Colorimetric - COD - PAP (Trinder)

The enzymatic determination of Cholesterol is made accordingly with the following reactions:

The intensity of formed cherry color is directly proportional to the concentration of Cholesterol in the sample.

REAGENTS

Number 1 - Enzymatic Reagent - Store between 2 and 8°C. Contains: Buffer < 150 mmol/L, Phenol < 9 mmol/L, 4-Aminoantipyrine < 3 mmol/L, Cholesterol Oxidase < 5 KU, Lipoprotein Lipase < 5 KU, Peroxidase < 5 KU, chelator, surfactants, stabilizers and preservative.

Number 2 - Standard - Store between 2 and 8°C. Contains: Cholesterol 200 mg/dL, solubilizers, surfactant and stabilizer.

PRESENTATION

PRESENTATION	REAGENT Nº 1	REAGENT Nº 2
K306-2-VET	2 x 100 mL	3 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, water bath (37°C), watch or stopwatch, pipettes, test tubes, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For *in vitro* diagnostic use only.

- Strictly follow the methodology proposed to obtain exact results.
- Water used in material cleaning must be recent and free of contaminants.
- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- The Sodium Hypochlorite is a contaminating agent that can significantly alter the results, so the materials used to perform the tests must be properly washed and free from this type of residue.
- Water level in water bath should be above the level of reagents in test tubes.
- Plasma collected with anticoagulants such as oxalate, EDTA or citrate, generate wrongly lowered results.

8- Handle reagents with care. Reagent Nº 1 contains Sodium Azide, irritant to skin and mucous.

9- Slightly pink color development in the Reagent Nº 1 does not damage their quality, and stability as long as the corresponding Blank and Standard dosages are used.

10- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

11- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

12- Do not use the product in case of damaged packaging.

13- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum obtained free of hemolyses or plasma collected with Heparin. Cholesterol is stable in plasma or serum for up to 7 days between 2 and 8 °C and 3 months at -20 °C. ⁷

PROCESS DESCRIPTION

The calibration stability of the Cholesterol Monoreagent kit installed on refrigerated equipment is at least 30 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Mark three test tubes: B (Blank), A (Sample), P (Standard), and proceed as follows:

	Blank	Standard	Sample
Sample	--	--	10 µL
Reagent Nº2	--	10 µL	--
Reagent Nº 1	1.0 mL	1.0 mL	1.0 mL

Homogenize well and put it in water bath at 37°C for 10 minutes. Read the absorbance of Sample and Standard at 500 nm (490 - 540 nm), matching the zero with the Blank. The color is stable for 30 minutes.

CALCULATIONS

$$\text{Cholesterol (mg/dL)} = \frac{\text{Sample Absorbance}}{\text{Standard Absorbance}} \times 200$$

Since reaction follows the Lambert-Beer Law, the Calibration Factor can be used.

$$\text{Calibration Factor} = \frac{\text{Standard Concentration (200 mg/dL)}}{\text{Standard Absorbance}}$$

$$\text{Cholesterol (mg/dL)} = \text{Sample Absorbance} \times \text{Calibration Factor}$$

The results shall be express as mg/dL.

VETERINARY USE

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

INTERFERENCES

The equilibrium of the assay is affected by some interference substances such as the Ascorbic Acid (even in small concentrations), Hemoglobin higher than 150 mg/dL and Bilirubin higher than 20 mg/dL. Failure to fast (approximately 12 hours for dogs and cats) and the use of corticosteroids can lead to falsely increased results.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 911C.

REFERENCE VALUES

The reference values are determined in healthy populations of dogs, cats, cattle and horses, males and females, adults and without breed predilection.

Specie	Reference Value
Dogs	114 - 366 mg/dL
Cats	65.5 - 220.5 mg/dL
Cattle	46.54 - 236.5 mg/dL
Horses	55.33 - 102.7 mg/dL

These values must be used as a guide, each laboratory must establish its own reference value, according to the population served. The results provided by this kit are to be interpreted by the professional veterinarian, but no single criterion is used to determine the patient's diagnosis and/or treatment. To convert the values from mg/dL into mmol/L (SI) multiply by 0.026.

PRODUCT PERFORMANCE**ACCURACY****Recovery**

The recovery analysis was made through 5 determinations of samples. Accuracies were calculated and were found in good concordance to the reference values, obtaining a recuperation between 95% and 105%.

Comparison of Methods

The Cholesterol Monoreagent kit was compared with another method commercially available to measure the dosage of Cholesterol. 42 analyzes were performed and results were evaluated. The linear equation obtained was $Y = 1.0118X - 0.5703$ with correlation coefficient 0.9932. With these results we can conclude that the kit shows good methodological specificity.

PRECISION**Repeatability**

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	102.75	207.40	386.73
Standard Deviation (mg/dL)	2.82	1.08	3.10
Correlation Coefficient (%)	2.74	0.52	0.80

Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	100.48	203.86	385.53
Standard Deviation (mg/dL)	2.23	5.60	2.26
Correlation Coefficient (%)	2.22	2.75	0.59

SENSITIVITY

Sensitivity was calculated through 40 determinations of a sample free of cholesterol. The average found, which indicates the limit of detection, was -1.99 mg/dL with a standard deviation of 0.89 mg/dL, the sensitivity, which indicates the limit of quantification from the method, corresponds to the average plus 3 times the standard deviation, and is equal to 0.67 mg/dL.

LINEARITY

Reaction is linear up to 500 mg/dL. For samples with higher values than 500 mg/dL, dilute the sample with Sodium Chloride 0.85%, repeat the dosage and multiply the results by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

Cholesterol is a precursor to cholesteryl esters, bile acids and steroid hormones. In cases of nephrotic syndrome, hypothyroidism (in dogs), cholestasis, diabetes mellitus, hyperadrenocorticism, pancreatitis and negative energy balance, cholesterol elevation is common. Severe hemolysis, absence of fasting and use of corticosteroids can lead to falsely increased cholesterol levels. Conditions of hypocholesterolemia can occur due to: genetic defects in the production of apolipoprotein, digestive problems that result in malabsorption of precursors, reduction in production due to liver damage, metabolic alterations that affect cholesterol

metabolism, and an increase in the ingestion of lipoproteins. Severe jaundice can result in a false reduction in cholesterol concentration.

BIBLIOGRAPHIC REFERENCES

- 1 - Allain CC et al. Clin. Chem., 1974. 20-470.
- 2 - Tonks DB. Qualyt Control in Clinical laboratories, 1983.
- 3 - Trinder P. Ann. Clin. Biochem., 1969, 6-24.
- 4 - Henry JB. Diagnosis and Management by Laboratory Methods, 19^a ed., 1996.
- 5 - Carl AB, Edwards A. Tiedz Textbook of Clinical Chem. 2nd ed., 1994, 1002-1081.
- 6 - Reunião Conjunta – Laudos Laboratoriais. SBC/DA, SBAC, SBPC/ML, SBBM; 15 de Outubro de 2013.
- 7 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31
- 8 - Faludi AA, Izar MCO, Saraiva JFK, Chacra APM, Bianco HT, Afune Neto A et al. Atualização da Diretriz Brasileira de Dislipidemias e Prevenção da Aterosclerose – 2017. Arq Bras Cardiol 2017; 109(2Supl.1):1-76
- 9 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

**QUIBASA QUÍMICA BÁSICA Ltda**

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone: 0800 0315454 | E-mail: sac@bioclin.com.br

Registration Number: Product exempt from registration with the Ministry of Agriculture, Livestock and Supply.

Technical Responsible: Dr. Camila Eckstein (CRMV/MG 20.611)

Review: June/2025

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MADE BY

LOT NUMBER

CONTROL

MANUFACTURING DATE

POSITIVE CONTROL

VALIDITY DATE
(last day of the month)

NEGATIVE CONTROL

TEMPERATURE LIMIT
(store)

BIOLOGICAL RISK

CONTENT IS SUFFICIENT
FOR <N> TEST

FLAMMABLE

SEE INSTRUCTIONS
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC
PRODUCT

TOXIC

KEEP AWAY
FROM SUNLIGHTDO NOT USE IF
PACKAGE IS
DAMAGED

DO NOT REUSE

PRODUCT
STERILIZED

CAUTION

DANGER