## GLICOSE MONORREAGENTE

REF K310

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
## INSTRUCTIONS FOR USE

### FUNCTION

Method for determination of glucose in biological samples of serum or plasma (EDTA, Heparin, Fluoride or Citrate), urine or cavity fluids. Colorimetric enzymatic test. For *in vitro* diagnostic only.

### PRINCIPLE OF ACTION

**Methodology:** Colorimetric Enzymatic - GOD - PAP (Trinder)  
The Glucose is enzymatically oxidized by Oxidase Glucose (GOD) according to the following reaction:

The Peroxide of Hydrogen in the presence of Peroxidase (POD) reacts with 4 - Aminoantipyrine and Phenol, forming a cherry red chromogen which color intensity is proportional to the glucose concentration.

### REAGENTS

**Number 1 - Enzymatic Reagent** - Store between 2 and 8°C. Contains: Tampon < 36 mmol/L, Phenol < 20 mmol/L, 4 - Aminoantipyrin < 5 mmol/L, Oxidase Glucose > 10.000 U/L, Peroxidase > 700 U/L, stabilizers, surfactant and preservative.  
**Number 2 - Standard** - Store between 2 and 8 °C, well sealed. Contains: Glucose 100,0 mg/dL and preservative.

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### EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, water-bath at 37°C, watches and stopwatches, pipettes, test tubes, Biocontrol N and Biocontrol P Bioclin. These items are found at markets in the specialized market of clinical analysis laboratories.

### TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 and 8°C. The transport temperature up to 30°C, must not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

### SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- Water used in material cleaning must be recent and free of contaminants.
- 4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- 5- The Sodium Hypochlorite is a contaminating agent that can significantly alter the results, so the materials used to perform the tests must be properly washed and free from this type of residue.
- 6- Water level in water bath should be above the level of reagents in test tubes.

7- Reagents N° 1 and N°2 contains Sodium Azide which can be irritating to skin and mucous membranes. Handle with care.

8- The development of a pinkish color in reagent N° 1 do not interfere in reagents quality and stability, once the correspondent Blank be used and dosages in regular bases.

9- Determine the factor periodically and in every batch.

10- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

11- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

12- Do not use the product in case of damaged packaging.

13- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

### SAMPLE

Serum, plasma (EDTA, Heparin, Fluoride or Citrate), urine and cavity fluids.

For serum collection, blood without the addition of anticoagulant should be centrifuged immediately after collection.

The use of the anticoagulant Fluoride Bioclin is recommended because it inhibits glycolysis. The use of fluoride in inadequate concentrations can cause the lysis of red blood cells and change the concentration of glucose in the sample.

Glucose is stable for up to 8 hours in whole blood samples treated with an antiglycolytic and kept at 4°C. In serum and plasma samples kept at 2 - 8 °C glucose is stable for up to 3 days.

Urine should be stored at 2 to 8° for a maximum of 24 hours. Keeping the urine at room temperature can result in the loss of up to 40% of the glucose in the sample within 24 hours. Cavity fluids must be analyzed immediately after collection as they are influenced by constituents that may be present in the sample.

### PROCESS DESCRIPTION

The calibration stability of the Glucose Monoreagent kit installed on refrigerated equipment is at least 30 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

### TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Mark 3 tests tubes with the letters B (Blank), A (Sample), P (Standard) and proceed as follows:

	Blank	Standard	Sample
Sample	--	--	10 mL
Reagent N° 2	--	10 mL	--
Reagent N° 1	1.0 mL	1.0 mL	1.0 mL

Homogenize well and place in water-bath at 37°C for 10 minutes.

Read the sample and standard absorbance at 505 nm (490 - 550 nm), matching the zero with the blank. The color is stable for 30 minutes.

### CALCULATIONS

$$\text{Glucose (mg/dL)} = \frac{\text{Sample Absorbance} \times 100}{\text{Standard Absorbance}}$$

#### Example:

$$\begin{aligned} \text{Sample Absorbance} &= 0.347 \\ \text{Standard Absorbance} &= 0.350 \end{aligned}$$

$$\text{Glucose (mg/dL)} = \frac{0.347}{0.350} \times 100 = 99$$

The reaction follows the Lambert-Beer Law, the calibration factor can be used.

$$\text{Calibration Factor} = \frac{\text{Standard Concentration (100 mg/dL)}}{\text{Standard Absorbance}}$$

$$\text{Glucose (mg/dL)} = \text{Absorbance Sample} \times \text{Calibration Factor}$$

#### Example:

$$\text{Calibration Factor} = \frac{100}{0.350} = 286$$

$$\text{Glucose (mg/dL)} = 0.347 \times 286 = 99$$

The results are express in mg/dL.

### PROCEDURE LIMITATIONS

The proposed method is not indicated for glucose dosage in urine.

### INTERFERENCES

Samples with concentrations up to 20 mg/dL of Bilirubin, 750 mg/dL of Triglycerides and 160 mg/dL of Hemoglobin do not produce significant interference.

In cases where the action of the sample interferers is verified, it is recommended to perform the Sample Blank to minimize the impact of the interferent on the measurement. Proceed as follows:

Mark 1 tube as a Sample Blank and place 1.0 mL of 0.85% Sodium Chloride with 10 µL of the Sample. Determine

## VETERINARY USE

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

its absorbance in 490 - 550 nm, zeroing with distilled or deionized water. Subtract the absorbance thus obtained from the absorbance of the Sample tube. Calculate the concentration by multiplying the result by the Calibration Factor.

Sodium fluoride, despite inhibiting glycolysis and stabilizing glucose levels in the sample, has a hypotonic action that can promote red blood cell lysis if used in a proportion different from that indicated for its application.

The use of drugs such as xylasin, detomidine, propranolol, ketamine and megestrol acetate can cause an increase in blood glucose due to their pharmacological action.

### INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

### TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 917.

### REFERENCE VALUES

The reference values were determined in serum samples of healthy dogs, cats, cattle and horses, males and females, adults and with no predilection for breed.

Species	Value
Dogs	68.9 - 124.5 mg/dL
Cats	69.0 - 198.5 mg/dL
Cattle	47.6 - 83.7 mg/dL
Horses	57.3 - 100.3 mg/dL

These values should be used as a guideline and each laboratory should establish its own reference value, according to the population served. The results provided by this kit must be interpreted by the veterinary medical professional, and are not the only criterion for determining the diagnosis and/or treatment of the patient.

**PRODUCT PERFORMANCE****ACCURACY****Recovery**

The recovery analysis was made with 5 sample determinations. The accurate was calculated and was found in a good parity with the reference values, obtaining a recovery between 94% and 104%.

**Comparison of Methods**

Glucose Monoreagent kit was compared with other commercially available methods for measurement of glucose. For this test, 42 analyses are performed and the results were evaluated. The linear equation obtained was  $Y = 1.0053X + 0.4653$  and the correlation coefficient 0.9982. With these results, we can conclude that the kit shows good methodological specificity.

**PRECISION****Repeatability**

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	240.75	115.01	73.59
Standard Deviation (mg/dL)	3.72	4.15	1.92
Coefficient of variation (%)	1.55	3.61	2.61

**Reproducibility**

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	242.29	114.58	74.06
Standard Deviation (mg/dL)	2.56	0.60	0.56
Coefficient of variation (%)	1.06	0.52	0.76

**SENSITIVITY**

The sensitivity was calculated from 40 determinations of a sample free of glucose. An average value, that indicates the method detection limit, of 0.55 mg/dL, with standard deviation of 0.25 mg/dL. The sensitivity that indicates the method quantification limit corresponds the average plus three times the standard deviation and is equal to 1.31 mg/dL.

**LINEARITY**

The reaction is linear up to 500 mg/dL. For samples with values above 500 mg/dL, dilute the sample with Sodium Chloride 0.85%, repeat the dosage and multiply the results obtained by the dilution factor.

**DIAGNOSTIC SIGNIFICANCE**

Glucose is the main source of energy for the body's tissues, which can be obtained through food, hepatic

metabolism of glycogen or by the production of glucose through amino acid precursors (gluconeogenesis). Blood glucose is directly controlled by the action of the hormones insulin and glucagon, but is also directly influenced by hormones such as catecholamines, growth hormones (GH) and corticosteroids.

Physiological hyperglycemia is usually associated with a persistent increase in blood glucose and its main causes are: Diabetes mellitus (type I - destruction of  $\beta$  cells in the pancreas, Type II diabetes mellitus - caused by insulin resistance), Hyperadrenocorticism, Acromegaly (excess growth hormone), Hyperglucagonemia, Hyperpituitarism/ Pituitary Pars Intermedia Dysfunction (PPID) (in horses, the Pituitary Pars Intermedia dysfunction that leads to increased cortisol and insulin resistance) and Pheochromocytoma (tumors that synthesize catecholamines).

Conditions such as hyperthyroidism in cats, sepsis and obstruction in the proximal portion of the duodenum (especially in cattle) are important causes of hyperglycemia, but it usually occurs transiently.

The increase in blood glucose can occur through physiological mechanisms such as: postprandial period, stress or pregnancy. The use of drugs such as xylasin, detomidine, propranolol, ketamine and megestrol acetate can cause an increase in blood glucose due to the pharmacological effect. In these situations, the increase in glucose is often transient or non-persistent.

Hypoglycemia can cause seizures, and is of great relevance in neonates, due to the immaturity of the hepatic pathways of gluconeogenesis. Glucose reduction can occur due to changes in insulin concentration (reduction in production or increase in consumption by the body) or reduction in obtaining glucose from the feed, but it can also be present in the occurrence of sepsis, ketosis in cattle, pregnancy toxemia in sheep and goats or in case of excessive physical exertion (more frequent in athletic dogs and horses).

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**QUALITY ASSURANCE**

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

**QUIBASA QUÍMICA BÁSICA Ltda**

Rua Teles de Menezes, 92 - Santa Branca

CEP 31565-130 - Belo Horizonte - MG - Brasil

Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 - Made in Brazil

**CUSTOMER SERVICE**

Customer Advisory Service

Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

**Registration Number:** Product exempt from registration with the Ministry of Agriculture, Livestock and Supply.

**Technical Responsible:** Dr. Camila Eckstein (CRMV/MG 20.611)

**Review:** June/2025

**UNIVERSAL SYMBOLOGY**

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N>- TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER