MAGNÉSIO MONORREAGENTE

REF K314

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FINALIDADE

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- 1 - Mann CK, Yoe J H. Ann. Chem. 28:202, 1956.
- 2 - Mann CK, Yoe JH. Anal. Chim. Acta. 16:155, 1957.
- 3 - Rice EW, Lapara CZ. Clin. Chim. Acta. 10:260, 1964.
- 4 - Weissmann N, Pileggi VJ. (1974) in Clinical Chemistry Principles and Technics 2nd. Ed. Henry R., Cannon, D. C. e Winkelman, J. W. p. 678. Haper and Row Publishers.
- 5 - Tonkks DB. Clin. Chem. 9:217, 1963.
- 6 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- 7 - ECLiNPath. Cornell University College of veterinary Medicine. 2013.
- 8 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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- 1 - Mann CK, Yoe J H. Ann. Chem. 28:202, 1956.
- 2 - Mann CK, Yoe JH. Anal. Chim. Acta. 16:155, 1957.
- 3 - Rice EW, Lapara CZ. Clin. Chim. Acta. 10:260, 1964.
- 4 - Weissmann N, Pileggi VJ. (1974) in Clinical Chemistry Principles and Technics 2nd. Ed. Henry R., Cannon, D. C. e Winkelman, J. W. p. 678. Haper and Row Publishers.
- 5 - Tonkks DB. Clin. Chem. 9:217, 1963.
- 6 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- 7- EclinPaath. Cornell University College of veterinary Medicine. 2013.
- 8 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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INSTRUCTIONS FOR USE

VETERINARY USE

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:
SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

FUNCTION

Method for determination of Magnesium in biological samples of serum, plasma (Heparin), liquor and urine. Colorimetric test, for *in vitro* diagnostics use only.

PRINCIPLE OF ACTION

Methodology: Mann and Yoe/ Magon Sulfonate/ Xilidil Blue

Mann and Yoe dye in alkaline pH and in the presence of Magnesium develop a red color. The intensity of the red color of the complex is proportional to the concentration of Magnesium. The method does not require deproteinization, being more sensitive than the Titan Yellow method, allowing therefore the use of only 10 µL of the sample to be analyzed, what turns the technique excellent for pediatric use. Besides its simplicity, the method does not suffer interference from Calcium Gluconate. The presence of surfactants rules out the interference of lipemic serum.

Reagents

Number 1 - Color Reagent - Store between 2 and 8°C. Contains: Sulphonated Magon < 1 mmol/L, Anhydrous Potassium Carbonate < 300 mmol/L, stabilizers, surfactant and preservative.

Number 2 - Standard* - Store between 2 and 8 °C. Contains: Magnesium 2 mg/dL and preservative.

PRESENTATION

Presentation	Reagent N° 1	Reagent N° 2*
K314-2-VET	2 x 100 mL	1 x 3 mL

* The R2 Standard is not present in all presentations.

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, watch or stopwatch, pipettes, test tubes, Biocal Kit, Biocontrol N and Biocontrol P Bioclin. Can be found in markets specialized in Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- Water used in the cleansing of the material should be recent and free of contaminants.
- 4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- 5- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 6- For biosafety information or in case of accidents with the product, consult the SDS (Safety Data Sheet)

available on the website www.bioclin.com.br or request it via Quibasa's Customer Support Service.

7- Do not use the product in case of damaged packaging.
8- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum obtained free of hemolysis, plasma collected with Heparin, liquor and 24h urine.
The serum and plasma sample is stable for 7 days at 2 to 8°C, and 1 year at -20 °C.⁶
The urine sample is stable for 3 days at 2 to 8 °C and 1 year at -20°C.⁶

PROCESS DESCRIPTION

The calibration stability of the Magnesium Single Reagent kit installed on refrigerated equipment is at least 03 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

TECHNIQUE

Bioclin recommends, use as calibrator the Biocal Bioclin Kit and control serum, Biocontrol N and P Bioclin Kits.

Mark 3 test tubes: B (Blank), A (Sample), S (Standard) and proceed as follows:

	Blank	Standard	Sample
Reagent N° 1	1.0 mL	1.0 mL	1.0 mL
Sample	----	----	10 µL
Reagent N° 2	----	10 µL	----

Homogenize well and read the absorbance from Standard and from the Sample in 500 nm, correcting the zero with the Blank. The immediate color is stable for 30 minutes.

24 hour urine – Sample preparation:

Homogenize urine and measure its volume. Add 1 drop of concentrated HCl in 20 mL of this urine. Dilute 1,0 mL of acidified urine with 4,0 mL of distilled water (1:5 dilution). Homogenize it well by agitating and proceed the determination of Magnesium the same way that was proposed to the serum. Multiply the results obtained by 5 (dilution factor).

CALCULATIONS

Serum

$$\text{Magnesium (mg/dL)} = \frac{\text{Sample Absorbance} \times 2}{\text{Standard Absorbance}}$$

Since the reaction follows the Lambert-Beer law, the Calibration Factor may be used.

$$\text{Calibration Factor} = \frac{\text{Standard Concentration (2 mg/dL)}}{\text{Standard Absorbance}}$$

$$\text{Magnesium (mg/dL)} = \text{Sample absorbance} \times \text{Calibration Factor}$$

Urine

$$\text{Magnesium (mg/24h)} = \frac{\text{mg/dL} \times \text{Urinary Volume (mL)}}{100}$$

Results are expressed in mg/dL for serum and mg/24h for urine.

INTERFERENCES

No interference was observed by Bilirubin until 8 mg/dL and Hemoglobin until 250 mg/dL.
Usage of hemolyzed samples may result in falsely elevated results, since the erythrocytes contain 3 times more Magnesium than serum or plasma.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 929A.

REFERENCE VALUES

Reference values were determined in serum samples of healthy male and female, adult dogs, cats, cattle and horses without predilection of breed.

Espécie	Valor de Referência
Dogs	1.435 - 2.11 mg/dL
Cats	1.623 - 2.483 mg/dL
Cattle	1.537 - 2.521 mg/dL
Horses	0.27 - 3 mg/dL

These values should be used as guidelines and each laboratory should establish its own reference value, according to the population served. The results provided by this kit must be interpreted by the veterinary medical professional, and are not the only criterion for determining the diagnosis and/or treatment of the patient.
In order to convert values of mg/dL into mmol/L (SI) multiply it by 0.41.

PRODUCT PERFORMANCE**ACCURACY****Comparison of Methods**

The Magnesium Single Reagent was compared with other commercially available methods for measurement of Magnesium. 42 analyzes were performed and the results were evaluated. The linear equation obtained was $Y = 1.0091X - 0.0683$ and the correlation coefficient 0.9923. With these results, we can conclude that the kit shows good methodological specificity.

PRECISION**Repeatability**

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	0.96	2.01	5.13
Standard Deviation (mg/dL)	0.05	0.03	0.10
Coefficient of Variation (%)	5.17	1.66	1.97

Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	0.96	2.00	5.04
Standard Deviation (mg/dL)	0.03	0.01	0.09
Coefficient of Variation (%)	3.41	0.66	1.79

SENSITIVITY

The sensitivity was calculated from 40 determinations of a sample free the Magnesium. The average found was 0.01 mg/dL and the standard deviation of 0.01 mg/dL. The sensitivity, that indicates the method detection limit, corresponds the average plus 3 times the standard deviation and is equal to 0.04 mg/dL.

LINEARITY

Reaction is linear up to the concentration of 5.0 mg/dL. For higher values, dilute the sample with NaCl 0.85%, repeat the determination and multiply the results by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

Magnesium is an important mineral that acts as a catalyst for several intracellular reactions. The balance of serum magnesium concentration occurs by the control between intestinal absorption and renal excretion (which in ruminants occurs mainly through saliva and in horses through sweat). The reduction in the concentration of serum magnesium is closely related to the reduction of serum albumin or the reduction in food. Cattle grazing in places with low magnesium content or fertilized with calcium, nitrates, ammonia or sulfates can lead to reduced magnesium absorption, especially in lactating or pregnant dairy cattle.

BIBLIOGRAPHIC REFERENCES

- 1 - Mann CK, Yoe J H. Ann. Chem. 28:202, 1956.
- 2 - Mann CK, Yoe JH. Anal. Chim. Acta. 16:155, 1957.
- 3 - Rice EW, Lapara CZ. Clin. Chim. Acta. 10:260, 1964.
- 4 - Weissmann N, Pileggi VJ. (1974) in Clinical Chemistry Principles and Techniques 2nd. Ed. Henry R., Cannon, D. C. e Winkelman, J. W. p. 678. Haper and Row Publishers.
- 5 - Tonkks DB. Clin. Chem. 9:217, 1963.
- 6 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- 7- EclinPaht. Cornell University College of veterinary Medicine. 2013.
- 8 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

**QUIBASA QUÍMICA BÁSICA Ltda**

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: + 55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

Registration Number: Product exempt from registration with the Ministry of Agriculture, Livestock and Supply.

Technical Responsible: Dr. Camila Eckstein (CRMV/MG 20.611)

Review: June/2025

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER