PROTEÍNAS TOTAIS MONORREAGENTE

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REF **K315**

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FINALIDADE

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REF **K315**

USAGE INSTRUCTIONS

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:
SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

FUNCTION

Method for determination of Proteins in biological serum, plasma (EDTA or Heparin), pleural, synovial and ascitic fluids samples. Colorimetric test. For *in vitro* diagnostic only.

PRINCIPLE OF ACTION

Methodology: Biuret

The peptidic bindings of the protein (-CONH-) reacts with the cupric ions, in an alkaline environment, forming a violet colored complex proportional to the protein content in the environment. The Sodium Tartrate and Potassium iodate establishes the reagent and the appropriate concentration of Potassium iodate prevines its self-reduction.

REAGENTS

Number 1 - Biuret - Store between 2 and 8°C. Contains: Sodium Hydroxide < 500 mmol/L, Sodium and Potassium Tartrate < 100 mmol/L, Potassium iodate < 20 mmol/L and Cooper Sulphate < 200 mmol/L.

Number 2 - Standard - Store between 2 and 8°C. Shake before use. Contains: Albumin 4g/dL, stabilizer and preservative.

PRESENTATION

Presentation	Reagent Nº 1	Reagent Nº 2
K315-2-VET	250 mL	2 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometers or colorimeter, watch and stopwatch, pipettes and test tubes, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized in clinical analysis laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature shall be between 2 and 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not Freeze.**

SPECIAL CARE

1- For *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- Reagent Nº 2 is caustic, handle carefully. Reagent Nº 1 contains sodium azide, irritating to skin and mucous membrane.

6- Do not use plasma.

7- Determine the factor periodically and in each product batch.

8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

9- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

10- Do not use the product in case of damaged packaging.

11- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum and plasma (EDTA or Heparin) obtained free of hemolysis and fluids (pleural, synovial and ascitic). The analyte is stable for 7 days between 2 and 8°C and 60 days at negative 10°C⁸, 9.

PROCESS DESCRIPTION

The calibration stability of the Total Protein Monoreagent kit installed on refrigerated equipment is up to 15 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

REAGENT PREPARATION

The reagents are ready for use.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Mark 3 test tubes: B (Blank), P (standard), A (Sample) and proceed as follows:

	Blank	Standard	Sample
Sample	---	---	50 µL
Standard	---	50 µL	---
Biuret	2.5 mL	2.5 mL	2.5 mL

Homogenize well and let it rest for 10 minutes. Read the sample and standard absorbances at 545 nm (510 - 550 nm), matching the zero with the blank. The color is stable for 30 minutes.

CALCULATIONS

Total Protein = $\frac{\text{Sample absorbance}}{\text{Standard Absorbance}} \times 4$ (g/dL)

As the reaction follows the Lambert-Beer law, the Calibration Factor can be used.

Calibration factor = $\frac{\text{Standard Concentration}}{\text{Standard Absorbance}}$ (4 g/dL)

Total Protein = Sample Absorbance x Calibration Factor (g/dL)

The results are express in g/dL

PROCEDURE LIMITATIONS

The described method is not applied to the Protein dosage in urine and liquor. For this samples use the Bioclin Bioprot U/LCR kit.

INTERFERENCES

Hemolyzed, icteric or lipemic samples interfere with the measurement by the colorimetric method, producing falsely increased results. There are no studies that prove the interference of drugs on the serum concentration of total protein, but it is important that the use of medications be considered when interpreting the test.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 927D.

REFERENCE VALUES

The reference values for were determined in healthy dogs, cats, cattle and horses, males and females, adults and with no predilection for breed.

Specie	Reference Value
Dogs	5.1 - 7.1 g/dL
Cats	5.5 - 8.7 g/dL
Cattle	5.77 - 9.03 g/dL
Horses	5.18 - 7.65 g/dL

These values should be used as a guideline and each laboratory should establish its own reference value, according to the population served. The results provided by this kit must be interpreted by the veterinary medical professional, and are not the only criterion for determining the diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

ACCURACY

Comparison of Methods

The Total Protein Monoreagent kit was compared wit other method for total protein dosage, commercially available. 42 analyses were performed and the results were evaluated. The linear equation obtained was $Y = 0.984X + 0.0315$ and the correlation coefficient 0.9978. With these results we can conclude that the kit shows good methodological specificity.

PRECISION

Repeatability

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average Concentration (g/dL)	4.99	8.87	6.94
Standard Deviation (g/dL)	0.04	0.03	0.02
Coefficient of Variation (%)	0.75	0.39	0.24

Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Average Concentration (g/dL)	4.95	8.84	6.74
Standard Deviation (g/dL)	0.04	0.06	0.35
Coefficient of Variation (%)	0.77	0.73	5.22

SENSITIVITY

The sensitivity was calculated from 40 determinations of a sample free of Total Protein. The average found was 0.036 g/dL with a standard deviation of 0.002 g/dL. The sensitivity, that indicates the method detection limit, corresponds to the average plus 3 times the standard deviation, and is equal to 0.043 g/dL.

LINEARITY

The reaction is linear up to 12 g/dL. For higher values dilute the serum with saline 0.85%, repeat the dosage and multiply the result obtained per the dilution factor.

DIAGNOSTIC SIGNIFICANCE

Proteins present in the serum of animals and humans perform essential functions for the body, such as: transport of substances, maintenance of osmotic pressure inside blood vessels, coagulation, etc. Serum total protein concentration is affected by factors such as age, race, sex, hormone concentration, lactation and pregnancy, nutritional status, stress, and fluid loss. Albumin (A) and globulin (G) constitute the largest fraction of total proteins. The unexpected profile of total serum protein concentration indicates the occurrence of a pathological process in the animal and must be preceded by the evaluation of the proportion between albumin and globulin, which complement the diagnosis. Therefore, it is important that the serum total protein profile is evaluated together with the A:G ratio.

The increase in total proteins can occur in cases of dehydration, while its reduction can occur due to excessive fluid therapy or fluid loss. In both cases the change in total protein concentration does not cause changes in the albumin:globulin ratio.

Hypoproteinemia is a common finding in disorders of the gastrointestinal tract, usually caused by a reduction in serum albumin. Chronic liver diseases can cause a reduction in total serum protein due to reduced albumin production. Elevation of globulins usually causes a reduction in the A:G ratio and occurs in inflammatory diseases, lymphomas, nephrotic syndrome and suppurative skin diseases. In horses, the elevation of β -globulins may result from infection by gastrointestinal parasites.

The reduction of globulins when the A:G ratio is increased is usually due to failures in colostrum consumption in cattle and horses. In foals, Arabian foal immunodeficiency, an autosomal recessive hereditary disease, should also be considered.

BIBLIOGRAPHIC REFERENCES

- 1 - Gornall AG, Bardawill CJ, David MMJ. Biol. Chem., 1977, 751.
- 2 - Weichselbaum TE, Amer J. Clin. Pathol., 1946, 16,40.
- 3 - Slater L, Carter P M, Hobbs JR. Ann. Clin.Biochem., 1975, 12,333.
- 4 - Batsakis JG, Arousohn RS, Walker WA, Barnes B, Amer J.Clin. Pathol., 1.976, 66,238.
- 5 - Hoel PG. Estatística Elementar, Ed. Fundo de Cultura S/A, 1969.
- 6 - Tonks DB. Clin. Chem., 1983, 9,217.
- 7 - Carl AB, Edward RA. Tietz Textbook of Clinical Chem. 2nd ed, 1994, 695-697.
- 8 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- 9 - Kaneko JJ, Harvey JW, Bruss ML. Clinical Biochemistry of Domestic Animals. San Diego: Academic Express, 6. ed., 1998. 936 p.
- 10 - EClinPahth. Cornell University College of veterinary Medicine. 2013.
- 11 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUIMICA BASICA Ltda.

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

Registration Number: Product exempt from registration with the Ministry of Agriculture, Livestock and Supply.

Technical Responsible: Dr. Camila Eckstein (CRMV/MG 20.611)

Review: June/2025

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER